










1/13/2022

**Worklist: 5512**

REVIEWED

By Brittany Wylie at 3:13 pm, Jan 13, 2022

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2021-2770	3	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-2790		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-2793		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-0006		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-0012		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-0032		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2022-0045		BCK	AM 27 Blood THC Quant by LC-QQQ	



# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 1/12/22  
Plate lot#: 211018

Analyst: Anne Nord  
Plate re-test: 4/18/22

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 21D52496 **Urine Blank:** 83121 **Column:** UCT Selectra DA 100 x 2.1mm 3um  
**LCMS-QQQ ID:** 69679

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.


## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *THC-OH not evaluated in this run, ratios indicate an interferant.*




	1	2	3	4	5	6
a	cal 1	Internal control (urine)	c2021-2790-1			
b	cal 2	negative blood	c2022-0032-1			
c	cal 3	c2021-2793-1				
d	cal 4	c2022-0006-1				
e	Cal 5	c2022-0012-1				
f	cal 6	c2022-0045-1				
g	cal 7	negative urine				
h	Internal control (blood)	c2021-2770-3				

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 12:35:21 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	46813	∞	1191.8	∞	2437701	5.355 ng/ml
THC-COOH	1.416	167772	466.5	34.0	238.3	751463	16.170 ng/ml
THC	3.137	418729	∞	24.6	50842 68000 1341.8	3564313	4.656 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-B2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 12:42:05 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.454	6801	∞	1600.3 <b>High</b>	∞	2714699	0.468 ng/ml <b>Low</b>

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj urine internal control.d
<b>Type</b>	Sample	<b>Sample</b>	mj urine internal control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 1:42:16 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	60185	∞	868.6	∞	3389819	4.931 ng/ml
THC-COOH	1.431	141706	446.7	32.8	378.7	578276	17.668 ng/ml
THC	3.137	75725	1695957797 3703.5	26.5	488.5	664892	4.524 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative urine.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative urine
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 1:48:58 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood end of run.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood end of run.d
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 3:08:58 PM		

**Sample Info.**


## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	48096	∞	1208.2 <b>High</b>	∞	2604917	5.138 ng/ml
THC-COOH	1.416	175926	206.3	35.4	244.8	832862	15.343 ng/ml
THC	3.137	542134	∞	24.3	3191.1	4548848	4.718 ng/ml




# Compound Calibration Report




**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Last Cal. Update** 1/12/2022 3:20 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.1	114.5
mj cal2	2	✓	3.0	2.9	95.2
mj cal 3	3	✓	5.0	4.7	94.0
mj cal 4	4	✓	10.0	9.6	96.1
mj cal 5	5	✓	25.0	24.8	99.1
mj cal 6	6	✓	50.0	50.2	100.4
mj cal 7	7	✓	100.0	100.7	100.7

# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Last Cal. Update** 1/12/2022 3:20 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.3	105.6
mj cal2	2	✓	10.0	9.9	99.3
mj cal 3	3	✓	20.0	19.7	98.7
mj cal 4	4	✓	50.0	48.1	96.2
mj cal 5	5	✓	75.0	74.0	98.6
mj cal 6	6	✓	100.0	100.5	100.5
mj cal 7	7	✓	250.0	252.5	101.0

# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Last Cal. Update** 1/12/2022 3:20 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.1	114.9
mj cal2	2	✓	3.0	3.0	101.3
mj cal 3	3	✓	5.0	4.8	95.2
mj cal 4	4	✓	10.0	9.0	90.1
mj cal 5	5	✓	25.0	24.2	96.7
mj cal 6	6	✓	50.0	49.9	99.8
mj cal 7	7	✓	100.0	102.0	102.0

Not evaluated


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 11:41:46 AM		

**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	1.409	11877	∞	2236.1 <b>High</b>	∞	2458138	1.149 ng/ml	<b>Low</b>
THC-COOH	1.431	46172	135.5	37.1	370.4	711232	5.282 ng/ml	
THC	3.167	73850	∞	25.1	294.1	3433497	1.145 ng/ml	

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

**Instrument** 69679 **Data File** mj cal2.d  
**Type** Cal **Sample** mj cal2  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-B1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 1/12/2022 11:48:31 AM  
**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.409	27587	∞	1439.3 High	∞	2442479	3.040 ng/ml
THC-COOH	1.431	93968	337.0	36.4	55.8	708809	9.934 ng/ml
THC	3.152	231009	36866.5	25.1	1076.7	3382812	2.856 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 11:55:13 AM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.409	41147	∞	1281.9 High	∞	2396011	4.761 ng/ml
THC-COOH	1.431	197446	169.4	34.7	151.2	717544	19.739 ng/ml
THC	3.152	407128	4618471638 55621.0	23.1	1375.3	3428688	4.702 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 12:01:55 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	79535	∞	1139.1	∞	2510238	9.008 ng/ml
THC-COOH	1.431	516754	719.4	37.9	2098.4	751425	48.106 ng/ml
THC	3.152	900426	∞	22.7	∞	3560556	9.610 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

**Instrument** 69679 **Data File** mj cal 5.d  
**Type** Cal **Sample** mj cal 5  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-E1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 1/12/2022 12:08:37 PM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	202323	∞	969.9	∞	2421822	24.184 ng/ml
THC-COOH	1.431	773516	2774.2	37.8	3371.2	727139	73.966 ng/ml
THC	3.152	2291536	∞	23.7	∞	3434828	24.766 ng/ml




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 12:15:19 PM		

**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	424473	∞	900.8	∞	2477158	49.883 ng/ml
THC-COOH	1.431	1043642	10625.2	37.8	302.3	719951	100.495 ng/ml
THC	3.152	5186164	∞	24.7	∞	3805085	50.223 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\011222\QuantResults\cann.batch.bin  
**Calibration Last Update** 1/12/2022 3:20:10 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/12/2022 12:22:03 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	827584	∞	858.5	∞	2368938	101.975 ng/ml
THC-COOH	1.431	2396980	1635.4	38.9	2037.5	654933	252.478 ng/ml
THC	3.136	10359047	∞	24.4	∞	3777114	100.698 ng/ml