



3/30/2022

Worklist: 5727

REVIEWED
By Brittany Wylie at 7:58 am, Mar 31, 2022

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>
C2022-0438	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-0553	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-0587	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-0614	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-0638	1	BCK	AM 27 Blood THC Quant by LC-QQQ





AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 3/30/22
Plate lot#: 211018

Analyst: Anne Nord
Plate re-test: 4/18/22

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 22B52020 **Urine Blank:** **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.


Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Blood only run*

Curve Limits: THC-OH 3-100

THC reported as less than 3 for values between 1 and 3.




	1	2	3	4	5	6
a	cal 1	Internal control				
b	cal 2	negative blood				
c	cal 3	553-1				
d	cal 4	587-1				
e	Cal 5	614-1				
f	cal 6	638-1				
g	cal 7	438-1				
h	Internal control (blood)					

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:46:20 AM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	322589	∞	12.6	∞	2375978	5.440 ng/ml
THC-COOH	1.431	151963	482.0	35.6	226.1	695576	15.331 ng/ml
THC	3.152	109046	∞	24.4	∞	920042	4.592 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj negative blood.d
Type	Sample	Sample	mj negative blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B2	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:53:04 AM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM


Instrument	69679	Data File	end of run blood control.d
Type	Sample	Sample	end of run blood control
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 1:06:35 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	321279	∞	13.8	∞	2464780	5.244 ng/ml
THC-COOH	1.431	163929	311.1	35.2	103.3	738621	15.562 ng/ml
THC	3.152	124003	∞	25.2	∞	1063115	4.526 ng/ml

Compound Calibration Report



Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Last Cal. Update 3/30/2022 1:39 PM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	120.6
mj cal2	2	✓	3.0	2.9	97.7
mj cal 3	3	✓	5.0	4.6	91.1
mj cal 4	4	✓	10.0	9.1	90.9
mj cal 5	5	✓	25.0	24.7	98.7
mj cal 6	6	✓	50.0	49.3	98.7
mj cal 7	7	✓	100.0	102.2	102.2

Compound Calibration Report




Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Last Cal. Update 3/30/2022 1:39 PM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.1	102.9
mj cal2	2	✓	10.0	10.2	102.0
mj cal 3	3	✓	20.0	19.3	96.3
mj cal 4	4	✓	50.0	48.5	97.0
mj cal 5	5	✓	75.0	78.0	104.0
mj cal 6	6	✓	100.0	97.0	97.0
mj cal 7	7	✓	250.0	251.9	100.7

Compound Calibration Report



Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Last Cal. Update 3/30/2022 1:39 PM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	x	1.0	1.3	133.1
mj cal2	2	✓	3.0	3.1	101.8
mj cal 3	3	✓	5.0	4.9	98.9
mj cal 4	4	✓	10.0	10.0	99.9
mj cal 5	5	✓	25.0	24.7	98.9
mj cal 6	6	✓	50.0	50.1	100.3
mj cal 7	7	✓	100.0	100.1	100.1


dropped cal 1 due to ratio being out of range.

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 1.d
Type	Cal	Sample	mj cal 1
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 10:52:41 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	1.406	46712	∞	21.4 High	139.0	2128417	1.331 ng/ml	Low
THC-COOH	1.431	40260	137.7	36.6	66.6	612252	5.144 ng/ml	
THC	3.182	16264	∞	28.5	26.2	740977	1.206 ng/ml	


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal2.d
Type	Cal	Sample	mj cal2
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 10:59:25 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	171260	∞	13.0	∞	2457894	3.054 ng/ml
THC-COOH	1.431	96026	571.1	32.9	194.1	678128	10.204 ng/ml
THC	3.182	59730	∞	29.1	1439.5	838951	2.932 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 3.d
Type	Cal	Sample	mj cal 3
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:06:09 AM		

Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	297207	∞	13.5	∞	2435048	4.945 ng/ml
THC-COOH	1.431	189824	12407.1	37.1	249.9	684518	19.256 ng/ml
THC	3.167	103666	∞	26.7	∞	882547	4.554 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 4.d
Type	Cal	Sample	mj cal 4
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:12:53 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	635724	∞	11.8	∞	2427110	9.995 ng/ml
THC-COOH	1.431	494253	1310.8	38.5	929.9	690272	48.520 ng/ml
THC	3.167	211899	∞	25.2	∞	858060	9.094 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 5.d
Type	Cal	Sample	mj cal 5
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:19:35 AM		

Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	1550741	∞	12.1	∞	2313942	24.735 ng/ml
THC-COOH	1.431	726474	1937.1	37.9	22707 9.6	627384	77.997 ng/ml
THC	3.167	604840	∞	24.2	∞	874361	24.686 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 6.d
Type	Cal	Sample	mj cal 6
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:26:18 AM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	3255713	∞	12.4	6164.1	2369472	50.147 ng/ml
THC-COOH	1.431	954171	354.7	38.9	2462.3	661220	97.014 ng/ml
THC	3.167	1266959	∞	24.0	∞	908027	49.349 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\033022\QuantResults\cann.batch.bin
Calibration Last Update 3/30/2022 1:39:43 PM

Instrument	69679	Data File	mj cal 7.d
Type	Cal	Sample	mj cal 7
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	3/30/2022 11:33:00 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.406	6239158	∞	12.5	∞	2262010	100.124 ng/ml
THC-COOH	1.431	2246332	8136.0	39.0	10862.2	596711	251.865 ng/ml
THC	3.152	2522578	∞	24.5	∞	869176	102.178 ng/ml