








REVIEWED
By Amber Gerheart at 9:39 am, Oct 21, 2022

10/19/2022

Worklist: 6138

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>	
C2022-2126		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-2212		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-2220		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-2295		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2022-2343		BCK	AM 27 Blood THC Quant by LC-QQQ	

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 10/19/22
Plate lot#: 220802

Analyst: Anne Nord
Plate re-test: 2/2/23

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 22B52016-1 **Urine Blank:** blood only **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.


Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL blood (calibrated pipette) Pipette ID: I41142J in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: THC-OH range 3-100 (cal 1 dropped due to ratio)



	1	2	3	4	5	6
a	cal 1	blood control internal				
b	cal 2	negative blood				
c	cal 3	2126-1				
d	cal 4	2212-1				
e	Cal 5	2220-1				
f	cal 6	2295-1				
g	cal 7	2343-1				
h	Internal control (blood)					

Plate position 3


c2022- ____ - _

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	10/19/2022 7:51:37 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	76707	∞	994.0	504.3	4875499	4.848 ng/ml
THC-COOH	1.403	100708	977.2	281.4	∞	1353158	13.775 ng/ml
THC	3.061	476440	2388.5	24.7	∞	3847868	4.673 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument	69679	Data File	mj negative blood.d
Type	Sample	Sample	mj negative blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B2	Comment	
Injection Volume	10		
Acq. Date-Time	10/19/2022 7:58:22 PM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument	69679	Data File	end of run blood control.d
Type	Sample	Sample	end of run blood control
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	10/19/2022 9:12:01 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	95468	∞	1014.3	∞	6094862	4.830 ng/ml
THC-COOH	1.403	137499	1044.7	262.6	∞	1682349	15.019 ng/ml
THC	3.061	736534	133210.2	23.4	1554.9	5780911	4.796 ng/ml

Compound Calibration Report

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Last Cal. Update 10/20/2022 9:52 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	119.8
mj cal 2	2	✓	3.0	2.9	96.0
mj cal 3	3	✓	5.0	4.5	91.0
mj cal 4	4	✓	10.0	9.4	94.0
mj cal 5	5	✓	25.0	24.4	97.8
mj cal 6	6	✓	50.0	50.0	99.9
mj cal 7	7	✓	100.0	101.6	101.6

Compound Calibration Report

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Last Cal. Update 10/20/2022 9:52 AM
Analyst Name ISP\datastor
Analyte THC-COOH

Internal Standard THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.5	109.1
mj cal 2	2	✓	10.0	9.9	99.4
mj cal 3	3	✓	20.0	18.9	94.6
mj cal 4	4	✓	50.0	48.4	96.8
mj cal 5	5	✓	75.0	74.4	99.2
mj cal 6	6	✓	100.0	99.6	99.6
mj cal 7	7	✓	250.0	253.3	101.3

Compound Calibration Report

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Last Cal. Update 10/20/2022 9:52 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	x	1.0	1.6	162.8
mj cal 2	2	✓	3.0	3.3	111.1
mj cal 3	3	✓	5.0	4.8	95.0
mj cal 4	4	✓	10.0	9.7	96.8
mj cal 5	5	✓	25.0	24.4	97.5
mj cal 6	6	✓	50.0	48.8	97.5
mj cal 7	7	✓	100.0	102.1	102.1

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679 **Data File** mj cal 1.d
Type Cal **Sample** mj cal 1
Acq. Method AM 27 THC quant.m **Operator** Anne Nord
Sample Position P3-A1 **Comment**
Injection Volume 10
Acq. Date-Time 10/19/2022 6:57:49 PM
Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	13219	109.3	1645.5 High	∞	3935637	1.628 ng/ml Low
THC-COOH	1.403	22558	∞	259.6	132.3	883097	5.456 ng/ml
THC	3.061	48764	423.4	30.1	64.1	2150716	1.198 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679 **Data File** mj cal 2.d
Type Cal **Sample** mj cal 2
Acq. Method AM 27 THC quant.m **Operator** Anne Nord
Sample Position P3-B1 **Comment**
Injection Volume 10
Acq. Date-Time 10/19/2022 7:04:33 PM
Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	43703	2315.7	1087.1	∞	4411522	3.332 ng/ml
THC-COOH	1.403	59768	423.4	257.8	∞	1151353	9.943 ng/ml
THC	3.077	207276	2442.0	23.9	372.7	2893487	2.880 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679 **Data File** mj cal 3.d
Type Cal **Sample** mj cal 3
Acq. Method AM 27 THC quant.m **Operator** Anne Nord
Sample Position P3-C1 **Comment**
Injection Volume 10
Acq. Date-Time 10/19/2022 7:11:18 PM
Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	70028	∞	1061.5	∞	4560041	4.750 ng/ml
THC-COOH	1.403	128722	1662.3	273.1	∞	1229917	18.921 ng/ml
THC	3.061	391957	2081.3	23.9	457.1	3261397	4.548 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679 **Data File** mj cal 4.d
Type Cal **Sample** mj cal 4
Acq. Method AM 27 THC quant.m **Operator** Anne Nord
Sample Position P3-D1 **Comment**
Injection Volume 10
Acq. Date-Time 10/19/2022 7:18:02 PM
Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	172280	∞	907.0	∞	5021383	9.682 ng/ml
THC-COOH	1.403	382140	∞	267.7	∞	1375853	48.383 ng/ml
THC	3.077	1058383	8058.2	23.8	1405.1	4047805	9.403 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679 **Data File** mj cal 5.d
Type Cal **Sample** mj cal 5
Acq. Method AM 27 THC quant.m **Operator** Anne Nord
Sample Position P3-E1 **Comment**
Injection Volume 10
Acq. Date-Time 10/19/2022 7:24:47 PM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	416728	∞	852.4	∞	4592058	24.368 ng/ml
THC-COOH	1.403	540385	∞	267.1	∞	1255400	74.374 ng/ml
THC	3.077	2530668	∞	24.7	3261.1	3619877	24.441 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-F1
Injection Volume 10
Acq. Date-Time 10/19/2022 7:31:31 PM
Sample Info.

Data File mj cal 6.d
Sample mj cal 6
Operator Anne Nord
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	813361	∞	833.1	5173.4	4408825	48.759 ng/ml
THC-COOH	1.403	682768	∞	266.4	∞	1180127	99.582 ng/ml
THC	3.061	4945440	∞	24.8	3476.0	3430117	49.960 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2022\am 27-28\101922\QuantResults\cann.batch.bin
Calibration Last Update 10/20/2022 9:52:17 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-G1
Injection Volume 10
Acq. Date-Time 10/19/2022 7:38:15 PM
Sample Info.

Data File mj cal 7.d
Sample mj cal 7
Operator Anne Nord
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	1608598	64427.9	811.8	∞	4129855	102.108 ng/ml
THC-COOH	1.403	1614582	∞	261.9	11335.3	1089529	253.341 ng/ml
THC	3.077	9129665	195053.7	24.8	9160.3	3101382	101.570 ng/ml