



**REVIEWED**  
By Tamara Salazar at 11:31 am, Dec 08, 2022

12/7/2022

**REVIEWED**  
By Brittany Wylie at 10:34 am, Dec 09, 2022

**Worklist: 6179**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>
C2022-2571	2	BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-2586		BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-2646		BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-2647		BCK	AM 27 Blood THC Quant by LC-QQQ
C2022-2650		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ
C2022-2664		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 12/6/22  
Plate lot#: 220802

Analyst: Anne Nord  
Plate re-test: 2/2/23

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE LCMS Methanol

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 22B52016-1 **Urine Blank:** 12522  
**LCMS-QQQ ID:** 69679

**Column:** UCT Selectra DA 100 x 2.1mm 3um

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.


## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette 1000µL blood (calibrated pipette) Pipette ID: I41142J in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: The negative urine, C2022-2264-1 and the end of run blood control evaporated and did not inject. The samples were reconstituted and re-injected. Those injections were evaluated.



	1	2	3	4	5	6
a	cal 1	Internal control urine	2650-1			
b	cal 2	negative blood				
c	cal 3	2571-2				
d	cal 4	2586-1				
e	Cal 5	2646-1				
f	cal 6	2647-1				
g	cal 7	negative urine				
h	Internal control (blood)	2664-1				

Plate position 3


c2022-\_\_\_\_-\_\_

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/6/2022 4:28:49 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	58311	∞	994.44	∞	3649948	4.442 ng/ml
THC-COOH	1.388	85715	∞	275.29	∞	1266078	13.515 ng/ml
THC	3.001	225327	35476.7	25.46	559.9	1821588	4.740 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-B2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/6/2022 4:35:34 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative urine g2.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative urine g2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/6/2022 5:35:50 PM		
<b>Sample Info.</b>			

Data not used

## Sample Chromatogram




Sample did not inject, sample was reconstituted and injected that injection was evaluated.

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative urine g2r.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative urine g2r
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/7/2022 11:07:49 AM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control urine end of run.d
<b>Type</b>	Sample	<b>Sample</b>	mj internal control urine end of run
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/6/2022 6:15:59 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	45084	∞	819.19	∞	2562733	4.857 ng/ml
THC-COOH	1.373	63257	∞	272.82	∞	872580	14.377 ng/ml
THC	3.001	77040	1027.5	24.63	∞	631877	4.679 ng/ml




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	end of run blood control blood.d
<b>Type</b>	Sample	<b>Sample</b>	end of run blood control blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/6/2022 6:22:44 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




sample did not inject, the sample was reconstituted and injected that injection was evaluated.

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

<b>Instrument</b>	69679	<b>Data File</b>	end of run blood control bloodr.d
<b>Type</b>	Sample	<b>Sample</b>	end of run blood control bloodr
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/7/2022 11:27:57 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	37478	∞	865.06	∞	2224162	4.666 ng/ml
THC-COOH	1.388	59091	144.2	263.24	∞	822320	14.262 ng/ml
THC	3.016	146526	3744.8	25.96	280.7	1229657	4.583 ng/ml








# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-A1  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 3:35:01 PM  
**Sample Info.**

**Data File** mj cal 1.d  
**Sample** mj cal 1  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	12050	$\infty$	967.44	$\infty$	3509925	1.221 ng/ml <b>Low</b>
THC-COOH	1.388	23253	$\infty$	261.69	93.5	1030063	5.401 ng/ml
THC	3.016	38547	533.4	24.65	165.7	1675538	1.255 ng/ml


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B1  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 3:41:45 PM  
**Sample Info.**

**Data File** mj cal 2.d  
**Sample** mj cal 2  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	1.364	33748	383.7	918.51	∞	3378506	2.905 ng/ml	Low
THC-COOH	1.388	51918	∞	255.41	385.2	1054384	10.195 ng/ml	
THC	3.016	124062	3758.4	23.97	247.1	1789472	2.858 ng/ml	

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679 **Data File** mj cal 3.d  
**Type** Cal **Sample** mj cal 3  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-C1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 3:48:29 PM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	59899	∞	978.27	∞	3646043	4.558 ng/ml
THC-COOH	1.388	120003	1432.3	263.22	1089.9	1207790	19.207 ng/ml
THC	3.016	223127	1777.6	24.68	498.0	1894267	4.536 ng/ml




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679 **Data File** mj cal 4.d  
**Type** Cal **Sample** mj cal 4  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-D1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 3:55:14 PM  
**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	124684	873.2	916.90	∞	3605891	9.219 ng/ml
THC-COOH	1.388	300980	186.6	269.28	∞	1191428	46.764 ng/ml
THC	3.016	470666	∞	24.36	709.1	1886356	9.096 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679 **Data File** mj cal 5.d  
**Type** Cal **Sample** mj cal 5  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-E1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 4:01:58 PM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	337741	∞	873.80	∞	3691825	23.832 ng/ml
THC-COOH	1.388	472632	∞	260.40	5045.3	1180201	73.347 ng/ml
THC	3.001	1322902	11855.8	23.87	882.1	1950453	23.938 ng/ml


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-F1  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 4:08:42 PM  
**Sample Info.**

**Data File** mj cal 6.d  
**Sample** mj cal 6  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	710804	15748.0	826.79	∞	3664944	50.144 ng/ml
THC-COOH	1.388	631581	∞	255.36	∞	1138680	101.072 ng/ml
THC	3.001	2644641	37376.7	24.15	1878.2	1866227	49.516 ng/ml


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2022\am 27-28\120622\QuantResults\cann.batch.bin  
**Calibration Last Update** 12/7/2022 12:09:45 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-G1  
**Injection Volume** 10  
**Acq. Date-Time** 12/6/2022 4:15:27 PM  
**Sample Info.**

**Data File** mj cal 7.d  
**Sample** mj cal 7  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.364	1447644	72604.8	816.69	∞	3652409	102.121 ng/ml
THC-COOH	1.388	1503454	386.3	251.93	∞	1069870	254.014 ng/ml
THC	3.016	5197568	1278794.5	24.67	5659.8	1758096	102.802 ng/ml