








REVIEWED

By Tamara Salazar at 10:53 am, Sep 16, 2022

9/12/2022



Worklist: 6100

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
P2022-2113	2	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2113	3	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2113	5	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2113	6	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2113	7	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2527	1	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	
P2022-2681	1	BCK	AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ	

AM# 28: Multi-Drug Quantitation by LC-MS/MS

Extraction Date: 09/09/2022

Analyst: Celena Shrum

Plate lot#: 220816

Plate Retest Date: 02/16/2023

Mobile phase A: 5mM Amm Form + 0.01% FA

Mobile phase B: 0.01% Formic Acid in MeOH

Blank Blood Lot: Lampire 22B52015-1

Column: Agilent 120 EC-C18 (2.1x 100-2.7um)

LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **250µL blood (calibrated pipette) or 250µL hydrolyzed urine** in wells of analytical (standards) plate.
Pipette ID: #42
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. (SKIPPED PER DEVIATION)
- 4. Pipette **250µL 0.5 M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **300µL of blood+base/urine+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right)
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate**.
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **900uL ethyl acetate**.
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate.
- 16. Add 50 ul 1% HCl in MeOH to all wells with samples.
- 17. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 18. Reconstitute in **100µL 20% MeOH** and heat seal plate with foil.

Post-Analytic

- 1. Create batch and process data.
- 2. Make necessary changes to integration limits
- 3. Integration linear and R² values ≥0.98 for each analyte.
- 4. For unknown samples and controls: response ratio within 20% of average of controls and standards, RT within +/- 5% (tramadol RT +/-2%), S/N for primary transition >10 and secondary transitions >5.
- 5. Did all QCs pass for each analyte? Yes, see comments Add Control data to QC tracking spreadsheet.
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports.

Compounds evaluated: 10-Hydroxycarbamazepine, alpha-hydroxymidazolam 5-500 (cal 8 removed due to accuracy), amitriptyline, carbamazepine, maprotiline 5-100 (cals 6, 7, and 8 removed due to accuracy), midazolam, nortriptyline, and topiramate 5-100 (cals 6, 7, and 8 removed due to accuracy).

P2022-2113-5 and P2022-2681-1 were reinjected due to carryover in the blank before the sample.

Maprotiline not evaluated in the reinjected samples due to the ratio being out in the QC100 reinject.

**Idaho State Police
Forensic Services**

Request for Departure from an Analytical Method or Quality Standard

Deviation Number (assigned by QM): TOX-22-01

Date of Request: **2/3/2022**

Requestor/Discipline: Celena Shrum/Toxicology

Analytical Method/Quality Standard, Revision #: AM #25, AM #28, AM #29, Revision 13

Temporary or Permanent Deviation: Permanent

Scope of Deviation (record specific information, e.g. affected programs, evidence types, expected end date; etc): Deviation will remain in place until the change is made in the next method revision.

Deviation Request (Describe detailed instructions of the changes being made; include reference to specific section number(s) in the method manual): 4.1.4 (Place plate on shaking incubator at approximately 900 rpm for approximately 15 minutes) of AM #25, AM # 28, and AM #29 is being removed. The removal of this step was tested in the validation “Addition of Compounds/Modifications for the MDS” (approved on 2/2/2022) and it was determined that that step is not necessary and can be removed.

Technical Justification for Analytical Method Deviations: Refer to validation “Addition of Compounds/Modifications for the MDS” (approved on 2/2/2022)

Technical Review

Departure approved
Comments:

Departure Not Approved
Comments:

Approver: Rachel Cutler
Title: Laboratory Manager



Date: 2/10/2022

Quality Review

Quality Approver: Jason Crowe
Title: Quality Manager
Date: 2/10/2022



	1	2	3	4	5	6	7	8	9	10	11	12
A										USED FOR ISTD BLANK	P2022-2113-3	IS + Cal. 8
B										USED FOR ISTD BLANK	P2022-2113-2	IS + Cal. 7
C										USED FOR ISTD BLANK	Blood NEG	IS + Cal. 6
D										P2022-2527-1	USED FOR ISTD BLANK	IS + Cal. 5
E										P2022-2681-1	IS + QC_4	IS + Cal. 4
F										P2022-2113-7	IS + QC_3	IS + Cal. 3
G										P2022-2113-6	IS + QC_2	IS + Cal. 2
H										P2022-2113-5	IS + QC_1	IS + Cal. 1

	1	2	3	4	5	6	7	8	9	10	11	12
A							IS + Cal. 8 (P2)	P2022-2113-3 (P2)	P2022-2747-1	P2022-2113-4	M2022-3195-1	IS + Cal. 8
B							IS + Cal. 7 (P2)	P2022-2113-2 (P2)	P2022-2721-1	P2022-2113-2	M2022-3141-1	IS + Cal. 7
C							IS + Cal. 6 (P2)	Blood NEG (P2)	P2022-2685-1	M2022-3631-3	Neg Ctrl	IS + Cal. 6
D						P2022-2527-1 (P2)	IS + Cal. 5 (P2)	IS + QC_4 (P2)	P2022-2681-1	M2022-3617-2	IS + QC_2	IS + Cal. 5
E						P2022-2681-1 (P2)	IS + Cal. 4 (P2)	IS + QC_3 (P2)	P2022-2673-1	M2022-3600-2	IS + QC_4	IS + Cal. 4
F						P2022-2113-7 (P2)	IS + Cal. 3 (P2)	IS + QC_2 (P2)	P2022-2667-1	M2022-3561-2	IS + QC_3	IS + Cal. 3
G						P2022-2113-6 (P2)	IS + Cal. 2 (P2)	IS + QC_1 (P2)	P2022-2644-1	M2022-3541-2	IS + QC_2	IS + Cal. 2
H						P2022-2113-5 (P2)	IS + Cal. 1 (P2)	P2022-2748-1	P2022-2113-5	M2022-3299-1	IS + QC_1	IS + Cal. 1

SLE Map


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument	Falco (069901)	Data File	p2 Negative Blood.d
Type	Sample	Sample	p2 Negative Blood
Acq. Method	AM 28 MDQ P2 Updated 081022 CS.m	Operator	Celena Shrum
Sample Position	P6-C8	Comment	
Injection Volume	5		
Acq. Date-Time	9/10/2022 4:50:19 AM		
Sample Info.			

Sample Chromatogram




AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 10.d
Type QC **Sample** p2 QC 10
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-G8 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 3:24:27 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.001	302764	518.36	90.3	11374.44	1220847	10.0344 ng/ml
alpha-hydroxymidazolam	5.807	70568	304.79	61.1	5244.34	548888	10.6154 ng/ml
Amitriptyline	5.701	100489	398.74	114.5	632.93	426891	10.6483 ng/ml
Carbamazepine	5.607	454054	2015.69	8.3	279.70	2392340	10.7463 ng/ml
Maprotiline	5.667	63137	392.05	204.1	646.62	426891	12.3736 ng/ml
Midazolam	5.842	46261	37055.29	91.3	5511.83	497152	10.8695 ng/ml
Nortriptyline	5.722	103975	2049.00	35.1	569.43	301606	10.3748 ng/ml
Topiramate	4.984	5388	10373.73	48.4	3844.62	39648	10.1959 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 100.d
Type QC **Sample** p2 QC 100
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-F8 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 7:52:43 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.008	2318152	1023.09	89.6	10302.46	906791	99.0988 ng/ml
alpha-hydroxymidazolam	5.814	605012	645.77	61.7	3036.83	489924	105.4338 ng/ml
Amitriptyline	5.701	715653	2839.49	110.3	6595.50	305794	107.2784 ng/ml
Carbamazepine	5.607	3631900	5082.38	8.0	1520.44	1841669	94.6975 ng/ml
Maprotiline	5.674	313960	5318.33	265.6	2158.60	305794	91.5158 ng/ml
Midazolam	5.849	367981	8244.38	95.0	26362.23	413771	98.6748 ng/ml
Nortriptyline	5.722	617063	473695.59	36.5	7600.36	179325	99.7570 ng/ml
Topiramate	4.984	31259	40527.65	45.3	9776.47	31455	92.7930 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 250.d
Type QC **Sample** p2 QC 250
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-E8 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 3:45:53 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.001	4153861	1100.19	87.7	7739.75	631157	254.3883 ng/ml
alpha-hydroxymidazolam	5.814	1146875	3653.89	63.3	681.18	409532	239.6075 ng/ml
Amitriptyline	5.701	1424858	3475.81	109.0	8919.76	243788	268.1525 ng/ml
Carbamazepine	5.607	6517055	4105.01	8.4	869.11	1279670	241.6928 ng/ml
Maprotiline	5.674	475058	787.29	335.6 High	6428.08	243788	174.5437 ng/ml*
Midazolam	5.849	717485	819.06	91.4	523.64	312703	253.6163 ng/ml
Nortriptyline	5.722	1147368	3927.90	36.2	4307.16	137201	241.8327 ng/ml
Topiramate	4.984	45639	39086.23	43.6	26013.52	17737	244.8506 ng/ml*

*Outside curve range.


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 1000.d
Type QC **Sample** p2 QC 1000
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
 CS.m
Sample Position P6-D8 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 4:07:22 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.008	7250650	5821.79	87.9	74181.13	284660	983.2000 ng/ml
alpha-hydroxymidazolam	5.807	1838863	249.24	62.7	26489.07	201377	782.2014 ng/ml *
Amitriptyline	5.701	3171287	11348.92	103.9	45811.46	141776	1026.7026 ng/ml
Carbamazepine	5.614	9341772	81230.77	8.5	888.43	435718	1011.7008 ng/ml
Maprotiline	5.674	475929	18040.64	686.5 High	2494.35	141776	301.3672 ng/ml *
Midazolam	5.849	1451816	3200.30	91.2	19646.12	149663	1070.2766 ng/ml
Nortriptyline	5.729	2402983	3026.71	36.2	363825.39	69131	1003.8579 ng/ml
Topiramate	4.991	38106	432.80	73.0 High	15.24	4723	773.9053 ng/ml *

*Outside curve range.




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 250_r.d
Type QC **Sample** p2 QC 250_r
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
CS.m
Sample Position P6-E8 **Comment**
Injection Volume 5
Acq. Date-Time 9/12/2022 12:08:36 PM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.022	5374220	1795.01	87.3	41298.36	814554	255.0210 ng/ml
alpha-hydroxymidazolam	5.814	1243578	3102.08	62.3	2527.58	430276	247.2980 ng/ml
Amitriptyline	5.708	1161724	1647.47	104.2	6796.70	195844	272.1564 ng/ml
Carbamazepine	5.614	8128303	11308.71	8.5	4139.56	1524277	252.9882 ng/ml
Maprotiline	5.680	312787	1192.77	395.3 High	2393.54	195844	142.8858 ng/ml *
Midazolam	5.849	807745	1695.18	93.3	548.75	371847	240.1404 ng/ml
Nortriptyline	5.729	732324	253995.37	36.4	7592.05	86434	245.0071 ng/ml
Topiramate	4.998	59551	613.17	44.4	25896.91	21123	268.5480 ng/ml *

*Outside curve range.




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 QC 100_r.d
Type QC **Sample** p2 QC 100_r
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
CS.m
Sample Position P6-F8 **Comment**
Injection Volume 5
Acq. Date-Time 9/12/2022 1:23:59 PM
Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.008	2775328	1185.09	85.9	6348.53	1043055	103.1242 ng/ml
alpha-hydroxymidazolam	5.814	687097	331.00	62.2	7422.52	564051	103.9971 ng/ml
Amitriptyline	5.701	693293	1472.19	110.9	2404.33	302096	105.1957 ng/ml
Carbamazepine	5.607	4274540	4520.48	8.4	687.68	2166998	94.7207 ng/ml
Maprotiline	5.674	253447	1162.97	315.2 High	10932.11	302096	74.6086 ng/ml
Midazolam	5.849	437243	3862.67	93.1	3091.58	497155	97.5891 ng/ml
Nortriptyline	5.722	516445	60967.03	36.3	6562.22	151107	99.0846 ng/ml
Topiramate	4.991	31284	43018.64	42.8	31760.18	28440	103.0189 ng/ml



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte 10-OH-Carbazepine **Internal Standard** 10-OH-Carbazepine-13-D6




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	5.0	100.0
p2 Cal 2-10ng	2	✓	10.0	10.4	104.1
p2 Cal 3 -25ng	3	✓	25.0	24.7	98.9
p2 Cal 4-50ng	4	✓	50.0	48.4	96.8
p2 Cal 5-100ng	5	✓	100.0	100.3	100.3
p2 Cal 6-250ng	6	✓	250.0	250.6	100.2
p2 Cal 7-500ng	7	✓	500.0	496.1	99.2
p2Cal 8-1000ng	8	✓	1000.0	1004.4	100.4



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte alpha-hydroxymidazolam **Internal Standard** alpha-hydroxymidazolam-D4




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	4.8	95.3
p2 Cal 2-10ng	2	✓	10.0	10.7	107.2
p2 Cal 3 -25ng	3	✓	25.0	26.0	104.1
p2 Cal 4-50ng	4	✓	50.0	51.2	102.5
p2 Cal 5-100ng	5	✓	100.0	103.6	103.6
p2 Cal 6-250ng	6	✓	250.0	248.6	99.5
p2 Cal 7-500ng	7	✓	500.0	439.5	87.9
p2Cal 8-1000ng	8	x	1000.0	725.9	72.6



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Amitriptyline **Internal Standard** Amitriptyline-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	4.8	95.6
p2 Cal 2-10ng	2	✓	10.0	10.5	105.4
p2 Cal 3 -25ng	3	✓	25.0	25.1	100.2
p2 Cal 4-50ng	4	✓	50.0	48.7	97.5
p2 Cal 5-100ng	5	✓	100.0	100.1	100.1
p2 Cal 6-250ng	6	✓	250.0	251.4	100.6
p2 Cal 7-500ng	7	✓	500.0	506.3	101.3
p2Cal 8-1000ng	8	✓	1000.0	993.0	99.3



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Carbamazepine **Internal Standard** Carbamazepine-13C6




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	5.9	118.4
p2 Cal 2-10ng	2	✓	10.0	10.8	108.0
p2 Cal 3 -25ng	3	✓	25.0	23.3	93.4
p2 Cal 4-50ng	4	✓	50.0	44.2	88.4
p2 Cal 5-100ng	5	✓	100.0	93.1	93.1
p2 Cal 6-250ng	6	✓	250.0	241.0	96.4
p2 Cal 7-500ng	7	✓	500.0	501.9	100.4
p2Cal 8-1000ng	8	✓	1000.0	1019.8	102.0



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Maprotiline **Internal Standard** Amitriptyline-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	4.7	93.4
p2 Cal 2-10ng	2	✓	10.0	11.1	111.2
p2 Cal 3 -25ng	3	✓	25.0	27.1	108.5
p2 Cal 4-50ng	4	✓	50.0	49.6	99.1
p2 Cal 5-100ng	5	✓	100.0	87.8	87.8
p2 Cal 6-250ng	6	✗	250.0	173.8	69.5
p2 Cal 7-500ng	7	✗	500.0	243.2	48.6
p2Cal 8-1000ng	8	✗	1000.0	317.5	31.7


CS



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Midazolam **Internal Standard** Midazolam-D4

Midazolam - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 6 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	5.0	100.5
p2 Cal 2-10ng	2	✓	10.0	10.7	107.1
p2 Cal 3 -25ng	3	✓	25.0	25.3	101.0
p2 Cal 4-50ng	4	✓	50.0	48.7	97.4
p2 Cal 5-100ng	5	✓	100.0	96.2	96.2
p2 Cal 6-250ng	6	✓	250.0	247.2	98.9
p2 Cal 7-500ng	7	✓	500.0	482.3	96.5
p2Cal 8-1000ng	8	✓	1000.0	1024.6	102.5



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Nortriptyline **Internal Standard** Nortriptyline-d3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	5.1	102.6
p2 Cal 2-10ng	2	✓	10.0	10.4	103.9
p2 Cal 3 -25ng	3	✓	25.0	25.1	100.4
p2 Cal 4-50ng	4	✓	50.0	48.4	96.8
p2 Cal 5-100ng	5	✓	100.0	98.5	98.5
p2 Cal 6-250ng	6	✓	250.0	243.0	97.2
p2 Cal 7-500ng	7	✓	500.0	496.6	99.3
p2Cal 8-1000ng	8	✓	1000.0	1012.9	101.3



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Last Cal. Update 9/16/2022 8:48 AM
Analyst Name ISP\datastor
Analyte Topiramate **Internal Standard** Topiramate-d12



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
p2 Cal 1-5ng	1	✓	5.0	3.8	75.9
p2 Cal 2-10ng	2	✓	10.0	10.9	108.6
p2 Cal 3 -25ng	3	✓	25.0	29.4	117.6
p2 Cal 4-50ng	4	✓	50.0	51.9	103.8
p2 Cal 5-100ng	5	✓	100.0	94.0	94.0
p2 Cal 6-250ng	6	✗	250.0	225.3	90.1
p2 Cal 7-500ng	7	✗	500.0	287.5	57.5
p2Cal 8-1000ng	8	✗	1000.0	716.3	71.6




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 1-5ng.d
Type Cal **Sample** p2 Cal 1-5ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
CS.m
Sample Position P6-H7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 1:26:30 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	4.995	154233	263.06	90.6	2967.66	1312809	4.9989 ng/ml
alpha-hydroxymidazolam	5.807	35062	128.91	60.4	6976.46	581563	4.7638 ng/ml
Amitriptyline	5.694	39303	467.77	115.9	550.74	365214	4.7822 ng/ml
Carbamazepine	5.601	227809	1247.91	8.4	117.04	2608413	5.9201 ng/ml
Maprotiline	5.667	22779	352.57	220.6	864.54	365214	4.6712 ng/ml
Midazolam	5.842	20676	2135.03	97.1	11105.14	516358	5.0243 ng/ml
Nortriptyline	5.716	40155	32484.19	35.1	462.88	246206	5.1312 ng/ml
Topiramate	4.977	2867	2321.10	43.3	197.03	41300	3.7955 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 2-10ng.d
Type Cal **Sample** p2 Cal 2-10ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-G7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 1:37:14 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	4.995	324672	228.41	87.8	1984.46	1260119	10.4070 ng/ml
alpha-hydroxymidazolam	5.807	72761	583.58	62.9	23023.21	560564	10.7211 ng/ml
Amitriptyline	5.694	99663	1511.90	112.0	2537.67	427806	10.5366 ng/ml
Carbamazepine	5.601	478778	3782.55	8.3	267.86	2506751	10.8028 ng/ml
Maprotiline	5.667	57300	365.05	214.2	575.85	427806	11.1163 ng/ml
Midazolam	5.842	46738	1739.38	93.8	40242.38	510152	10.7112 ng/ml
Nortriptyline	5.716	98953	40296.54	37.4	50131.39	286665	10.3877 ng/ml
Topiramate	4.984	6095	13006.63	38.3	6.69	42679	10.8615 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 3 -25ng.d
Type Cal **Sample** p2 Cal 3 -25ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-F7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 1:47:57 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	4.995	766091	2339.15	90.4	2037.47	1218543	24.7225 ng/ml
alpha-hydroxymidazolam	5.807	174110	1238.54	62.7	7406.91	564513	26.0301 ng/ml
Amitriptyline	5.694	220890	1568.79	112.7	2309.38	402143	25.0578 ng/ml
Carbamazepine	5.601	1137001	5939.83	8.3	433.66	2487357	23.3378 ng/ml
Maprotiline	5.667	125298	911.65	222.0	4791.37	402143	27.1138 ng/ml
Midazolam	5.842	112440	1702.30	93.0	93469.44	503120	25.2522 ng/ml
Nortriptyline	5.716	226259	91399.84	35.5	1599.68	264599	25.1077 ng/ml
Topiramate	4.977	12880	16061.09	43.9	3792.84	38397	29.4082 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 4-50ng.d
Type Cal **Sample** p2 Cal 4-50ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-E7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 1:58:41 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	4.995	1357571	925.25	89.5	7484.30	1092456	48.4113 ng/ml
alpha-hydroxymidazolam	5.807	346339	1331.04	62.5	112140.19	574844	51.2330 ng/ml
Amitriptyline	5.694	333364	2327.51	113.8	2343.42	312931	48.7463 ng/ml
Carbamazepine	5.601	2064318	5634.10	8.4	1258.44	2293992	44.1937 ng/ml
Maprotiline	5.667	175493	909.84	235.4	2222.04	312931	49.5584 ng/ml
Midazolam	5.842	200006	186775.75	92.8	58613.40	458562	48.7036 ng/ml
Nortriptyline	5.716	314579	3284.42	36.0	1510.53	189330	48.3875 ng/ml
Topiramate	4.977	20301	22274.78	40.8	4357.16	35670	51.9077 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 5-100ng.d
Type Cal **Sample** p2 Cal 5-100ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-D7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 2:09:24 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	4.995	2350466	382.08	88.8	6510.17	908052	100.3349 ng/ml
alpha-hydroxymidazolam	5.807	637949	1392.68	63.8	2873.99	525847	103.5718 ng/ml
Amitriptyline	5.694	481124	597.95	112.9	18367.88	220207	100.1429 ng/ml
Carbamazepine	5.601	3731404	8688.38	8.4	1899.60	1925683	93.0788 ng/ml
Maprotiline	5.667	217102	955.36	264.8	1338.18	220207	87.8413 ng/ml
Midazolam	5.842	347346	5183.82	94.8	3636.06	400714	96.1918 ng/ml
Nortriptyline	5.716	401124	27732.77	37.0	3133.90	118078	98.4892 ng/ml
Topiramate	4.977	30449	552.42	45.3	24.85	30249	94.0271 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 6-250ng.d
Type Cal **Sample** p2 Cal 6-250ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-C7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 2:20:07 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.001	4183659	1461.52	88.8	55949.66	645276	250.6142 ng/ml
alpha-hydroxymidazolam	5.807	1242337	5251.26	62.3	6899.38	427521	248.6453 ng/ml
Amitriptyline	5.701	1172758	2728.89	107.2	266.74	214019	251.3981 ng/ml
Carbamazepine	5.607	6754582	9317.55	8.4	2937.06	1330155	240.9993 ng/ml
Maprotiline	5.667	415287	10570.68	318.4 High	2889.75	214019	173.8026 ng/ml
Midazolam	5.842	713688	2563.79	92.2	3723.67	319191	247.1621 ng/ml
Nortriptyline	5.722	912691	14403.04	36.6	2965.23	108622	242.9818 ng/ml
Topiramate	4.977	45575	527.44	44.0	28991.34	19226	225.3363 ng/ml

CS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2 Cal 7-500ng.d
Type Cal **Sample** p2 Cal 7-500ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-B7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 2:30:49 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.001	5723172	4281.37	87.2	148319.31	445538	496.0733 ng/ml
alpha-hydroxymidazolam	5.807	1672564	14360.50	63.6	4843.48	325875	439.4767 ng/ml
Amitriptyline	5.701	2072052	3135.85	107.3	54.81	187816	506.3023 ng/ml
Carbamazepine	5.607	8820631	11411.27	8.4	1165.89	830880	501.8565 ng/ml
Maprotiline	5.674	509087	296.87	436.8 High	2191.36	187816	243.1584 ng/ml
Midazolam	5.849	1052457	4182.84	92.5	1233.71	240908	482.3407 ng/ml
Nortriptyline	5.722	1610400	47979.20	36.0	8129.22	93694	496.5953 ng/ml
Topiramate	4.984	34410	148465.05	51.3 High	22446.12	11409	287.5004 ng/ml


AM #28 Multi-Drug Quant. Results



Batch results D:\MassHunter\Data\2022\AM 27-28\090922 AM 28 P1 P2 CS\QuantResults\AM 28 P2.batch.bin
Calibration Last Update 9/16/2022 8:48:01 AM

Instrument Falco (069901) **Data File** p2Cal 8-1000ng.d
Type Cal **Sample** p2Cal 8-1000ng
Acq. Method AM 28 MDQ P2 Updated 081022 **Operator** Celena Shrum
Sample Position P6-A7 **Comment**
Injection Volume 5
Acq. Date-Time 9/10/2022 2:41:32 AM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
10-OH-Carbazepine	5.008	7389987	6095.93	86.6	328213.23	283993	1004.4379 ng/ml
alpha-hydroxymidazolam	5.807	1637251	2621.12	63.5	6559.52	193191	725.9223 ng/ml
Amitriptyline	5.701	3214230	1084.30	103.0	6945.54	148567	993.0338 ng/ml
Carbamazepine	5.614	9318220	38651.02	8.6	14632.18	431157	1019.8110 ng/ml
Maprotiline	5.674	525326	898.36	634.2 High	3494.08	148567	317.4909 ng/ml
Midazolam	5.849	1436137	6095.39	93.5	200.04	154649	1024.6141 ng/ml
Nortriptyline	5.722	2563869	41493.38	35.7	26485.63	73099	1012.9197 ng/ml
Topiramate	4.991	38788	196800.09	46.7	18607.81	5193	716.2619 ng/ml