


Worklist: 6244

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
P2023-0098	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0129	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0138	1	CBUK	AM 27 Blood THC Quant by LC-QQQ
P2023-0162	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0168	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0176	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0202	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0215	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2023-0226	2	BCK	AM 27 Blood THC Quant by LC-QQQ



AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 02/06/2023

Analyst: Tamara Salazar

Plate lot#: 230113

Plate Retest Date: 07/13/2023

Mobile phase A: 0.1% Formic Acid in LCMS Water

Mobile phase B: 0.1% Formic acid in Acetonitrile

Blank Blood Lot: Lampire 23A52593

Blank Urine Lot: N/A

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes.
- 3. Using a calibrated pipette, add **1000µl blood and urine (if applicable) (calibrated pipette)** into the appropriate wells of analytical (standards) plate. **Pipette ID: 42**
- 4. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 5. Add **500µL of 0.1% formic acid in water to blood samples, and 500µL of saturated phosphate buffer to urine samples** in the wells of the analytical plate.
- 6. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 7. Transfer **700-800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.
Amount transferred: 750µL
- 8. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right)
- 9. Wait 5 minutes.
- 10. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 11. Wait 5 minutes.
- 12. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 13. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 14. Wait 5 minutes.
- 15. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 16. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 17. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Did all QCs pass for each analyte? (if not, describe in comments section)
- 5. Enter QCs into control charting.
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: THC 3-100 – calibrator one dropped due to S/N

	1	2	3	4	5	6
A	IS + Cal. 1	IS + QC_1	P2023-0202-1	IS + Sample	IS + Sample	IS + QC_1
B	IS + Cal. 2	Neg Blood	P2023-0215-2	IS + Sample	IS + Sample	IS + Cal. 7
C	IS + Cal. 3	P2023-0098-1	P2023-0226-2	IS + Sample	IS + Sample	IS + Cal. 6
D	IS + Cal. 4	P2023-0129-2	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 5
E	IS + Cal. 5	P2023-0138-1	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 4
F	IS + Cal. 6	P2023-0162-2	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 3
G	IS + Cal. 7	P2023-0168-1	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 2
H	IS + QC_1	P2023-0176-1	IS + Sample	IS + Sample	IS + QC_1	IS + Cal. 1

All wells to contain 100 μ l of residual DMSO




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	MJ Negative Blood.d
Type	Sample	Sample	MJ Negative Blood
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-B2	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 2:05:12 PM		
Sample Info.			

Sample Chromatogram






AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	MJ QC Control Blood.d
Type	QC	Sample	MJ QC Control Blood
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-H1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:49:58 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	118512	∞	12.9	∞	1468159	4.5235 ng/ml
THC-COOH	1.614	30267	362.45	254.2	149.61	369232	14.9180 ng/ml
THC	3.510	95772	∞	24.5	119.85	2057324	5.3449 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	QC end MJ.d
Type	QC	Sample	QC end MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-A2	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 4:37:32 PM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	127485	∞	13.1	∞	1571290	4.5460 ng/ml
THC-COOH	1.614	34064	699.32	256.0	1547.45	424130	14.6365 ng/ml
THC	3.510	104167	∞	24.4	∞	2285228	5.2377 ng/ml

TS



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 2/6/2023 4:53 PM
Analyst Name ISP\tsalazar
Analyte THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1 MJ	1	x	1.0	1.1	114.2
Cal 2 MJ	2	✓	3.0	3.0	101.2
Cal 3 MJ	3	✓	5.0	5.0	100.4
Cal 4 MJ	4	✓	10.0	10.0	100.4
Cal 5 MJ	5	✓	25.0	23.7	94.6
Cal 6 MJ	6	✓	50.0	52.1	104.1
Cal 7 MJ	7	✓	100.0	99.2	99.2

TS



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 2/6/2023 4:53 PM
Analyst Name ISP\tsalazar
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1 MJ	1	✓	5.0	5.5	110.0
Cal 2 MJ	2	✓	10.0	9.6	96.3
Cal 3 MJ	3	✓	20.0	19.4	97.0
Cal 4 MJ	4	✓	50.0	47.6	95.2
Cal 5 MJ	5	✓	75.0	75.5	100.6
Cal 6 MJ	6	✓	100.0	100.0	100.0
Cal 7 MJ	7	✓	250.0	252.5	101.0

TS



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 2/6/2023 4:53 PM
Analyst Name ISP\tsalazar
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1 MJ	1	✓	1.0	1.1	107.7
Cal 2 MJ	2	✓	3.0	2.9	98.2
Cal 3 MJ	3	✓	5.0	4.8	96.3
Cal 4 MJ	4	✓	10.0	9.8	98.0
Cal 5 MJ	5	✓	25.0	24.9	99.5
Cal 6 MJ	6	✓	50.0	49.8	99.6
Cal 7 MJ	7	✓	100.0	100.7	100.7




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 1 MJ.d
Type	Cal	Sample	Cal 1 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-A1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 12:49:00 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.558	23458	∞	13.6	25.73	1339287	1.0768 ng/ml Low
THC-COOH	1.584	9349	159.33	256.1	∞	352674	5.5015 ng/ml
THC	3.495	19155	∞	22.9	6.64 Low	2230474	1.1421 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 2 MJ.d
Type	Cal	Sample	Cal 2 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-B1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 12:56:45 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	70145	∞	13.4	29.81	1353838	2.9470 ng/ml Low
THC-COOH	1.614	18538	95.16	255.9	∞	364878	9.6267 ng/ml
THC	3.510	64568	∞	25.3	40.48	2511583	3.0374 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 3 MJ.d
Type	Cal	Sample	Cal 3 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-C1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:04:22 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	118306	∞	13.0	176.43	1375248	4.8128 ng/ml
THC-COOH	1.599	39088	735.45	251.0	∞	360755	19.3966 ng/ml
THC	3.510	103338	∞	24.8	∞	2368945	5.0206 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 4 MJ.d
Type	Cal	Sample	Cal 4 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-D1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:11:58 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	242289	∞	13.0	264.25	1364476	9.8048 ng/ml
THC-COOH	1.599	97441	∞	252.0	1174.56	355131	47.5854 ng/ml
THC	3.510	209240	∞	23.7	116.15	2352428	10.0383 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 5 MJ.d
Type	Cal	Sample	Cal 5 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-E1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:19:33 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	649858	∞	13.1	2341.48	1431936	24.8700 ng/ml
THC-COOH	1.614	158892	1754.32	242.0	6108.74	362348	75.4507 ng/ml
THC	3.510	560176	∞	23.2	∞	2642470	23.6598 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 6 MJ.d
Type	Cal	Sample	Cal 6 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-F1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:27:10 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.588	1284908	∞	13.2	∞	1409998	49.8158 ng/ml
THC-COOH	1.599	199887	603.68	240.0	∞	342961	99.9540 ng/ml
THC	3.510	1067788	∞	22.5	∞	2278394	52.0744 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results G:\TOX\Pocatello\Falco\2023\AM 27 28\020623 AM 27 28 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 2/6/2023 4:53:02 PM

Instrument	Falco (069901)	Data File	Cal 7 MJ.d
Type	Cal	Sample	Cal 7 MJ
Acq. Method	AM 27 THCQ.m	Operator	Tamara Salazar
Sample Position	P1-G1	Comment	
Injection Volume	10		
Acq. Date-Time	2/6/2023 1:34:46 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.573	2567404	∞	13.4	∞	1392383	100.6728 ng/ml
THC-COOH	1.599	505807	13777.29	241.8	5538.05	341479	252.4852 ng/ml
THC	3.510	1766981	∞	22.4	∞	1976339	99.1695 ng/ml