


Bylee

Worklist: 6370

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2023-0949	1	BLOOD	AM 27 Blood THC Quant by LC-QQQ
C2023-0949	2	BLOOD	AM 27 Blood THC Quant by LC-QQQ
C2023-0949	3	BLOOD	AM 27 Blood THC Quant by LC-QQQ
C2023-0989	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-1001	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-1003	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-1049	1	BCK	AM 27 Blood THC Quant by LC-QQQ



AM# 27: Quantitation of THC and Metabolites in Blood and Urine *by Wylie*

LC-MS/MS

Extraction Date 5/11/23
Plate lot#: 230113

Analyst: Britany Wylie
Plate re-test: 7/13/23

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 22B52016-1 **Urine Blank:** 12522
LCMS-QQQ ID: 69679

Column: UCT Selectra DA 100 x 2.1mm 3um

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL (calibrated pipette) blood or 1000µL hydrolyzed urine Pipette ID: K52558G in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. (Add in 3 increments of 750uL)
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 12. Add 2.25mL Hexane. (Add in 3 increments of 750uL)
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators. SN > 10
- 4. Case sample response for THC 1ng/ml LOD ; 3ng/ml LOQ, OH-THC 3ng/mL LOD and LOQ, Carboxy-THC: 5 ng/mL (qualitative only). Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: .Blood only Run

end of run control- incorrect plate location was selected. QC was reconstituted and reinjected 5/12/23

BW

Bylee

	1	2	3	4	5	6
a	cal 1	Internal control urine	1049-1 (a2 on SLE and Ext)			
b	cal 2	negative blood				
c	cal 3	949-1				
d	cal 4	949-2				
e	Cal 5	949-3				
f	cal 6	989-1				
g	cal 7	1001-1				
h	Internal control (blood)	1003-1				

Plate position 3

c2023-0____-__


AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj negative blood.d
Type	Sample	Sample	mj negative blood
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-B2	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 11:18:45 PM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 11:12:10 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	47933	1643.8	938.64	∞	2640588	4.310 ng/ml
THC-COOH	1.358	57258	∞	252.22	∞	722013	15.653 ng/ml
THC	2.911	157332	∞	24.94	∞	1087455	4.881 ng/ml

AM #27 Cannabinoids


BWylee

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	end of run blood control blood A2.d
Type	Sample	Sample	end of run blood control blood A2
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	5/12/2023 12:57:34 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.358	17705	∞	303.14 High	∞	1010627	4.695 ng/ml Low
THC	2.911	37277	31629.9	33.20 High	∞	4484142	0.688 ng/ml Low

Incorrect plate location injected- sample was reinjected with the correct position.

This sample was not evaluated

BWylee


AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/15/2023 10:44:18 AM

Instrument	69679	Data File	end of run blood control blood h1.d
Type	Sample	Sample	end of run blood control blood h1
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	5/12/2023 9:03:33 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	52255	∞	936.75	∞	2801211	4.419 ng/ml
THC-COOH	1.358	55826	19.6	261.62	∞	837667	13.408 ng/ml
THC	2.911	159345	∞	25.23	∞	734618	7.103 ng/ml


Compound Calibration Report

Byylee

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Last Cal. Update 5/12/2023 9:17 AM
Analyst Name ISP\datastor
Analyte THC

Internal Standard THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	123.9
mj cal 2	2	✓	3.0	2.9	97.6
mj cal 3	3	✓	5.0	4.5	90.4
mj cal 4	4	✓	10.0	9.1	91.1
mj cal 5	5	✓	25.0	23.7	94.8
mj cal 6	6	✓	50.0	49.7	99.5
mj cal 7	7	✓	100.0	102.8	102.8


Compound Calibration Report

Byylee

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Last Cal. Update 5/12/2023 9:17 AM
Analyst Name ISP\datastor
Analyte THC-COOH

Internal Standard THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.5	110.2
mj cal 2	2	✓	10.0	9.7	97.4
mj cal 3	3	✓	20.0	19.5	97.3
mj cal 4	4	✓	50.0	47.8	95.7
mj cal 5	5	✓	75.0	73.6	98.1
mj cal 6	6	✓	100.0	99.4	99.4
mj cal 7	7	✓	250.0	254.4	101.8

Compound Calibration Report

Byylee

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Last Cal. Update 5/12/2023 9:17 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	120.3
mj cal 2	2	✓	3.0	3.0	98.4
mj cal 3	3	✓	5.0	4.6	92.8
mj cal 4	4	✓	10.0	9.1	90.9
mj cal 5	5	✓	25.0	23.9	95.8
mj cal 6	6	✓	50.0	49.6	99.2
mj cal 7	7	✓	100.0	102.5	102.5

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 1.d
Type	Cal	Sample	mj cal 1
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:25:58 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	10126	∞	895.94	∞	2592778	1.203 ng/ml Low
THC-COOH	1.358	14866	∞	284.68	∞	671892	5.512 ng/ml
THC	2.926	24487	∞	30.01	71.6	933754	1.239 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 2.d
Type	Cal	Sample	mj cal 2
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:32:42 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	31343	1254.8	913.07	∞	2627795	2.952 ng/ml Low
THC-COOH	1.358	31777	∞	277.40	∞	691047	9.743 ng/ml
THC	2.926	81474	925.0	25.81	∞	1004327	2.927 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 3.d
Type	Cal	Sample	mj cal 3
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:39:19 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	51266	∞	902.37	∞	2606092	4.641 ng/ml
THC-COOH	1.358	69402	∞	256.43	∞	688379	19.469 ng/ml
THC	2.926	135393	∞	26.66	574.3	1018324	4.521 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 4.d
Type	Cal	Sample	mj cal 4
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:45:53 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	115122	∞	918.70	∞	2872864	9.090 ng/ml
THC-COOH	1.358	198962	∞	254.90	∞	762835	47.848 ng/ml
THC	2.926	323500	2268.5	25.14	∞	1146595	9.109 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 5.d
Type	Cal	Sample	mj cal 5
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:52:28 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	306360	∞	870.68	122207.1	2831781	23.943 ng/ml
THC-COOH	1.358	306973	∞	255.89	∞	756209	73.587 ng/ml
THC	2.926	860065	19534.7	23.94	4478.1	1136383	23.707 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 6.d
Type	Cal	Sample	mj cal 6
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 10:59:02 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	632702	∞	856.66	∞	2800091	49.625 ng/ml
THC-COOH	1.358	403061	∞	250.39	∞	730936	99.393 ng/ml
THC	2.926	1813463	∞	23.94	∞	1131000	49.741 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2023\am 27-28\51123r\QuantResults\thcq.batch.bin
Calibration Last Update 5/12/2023 9:17:44 AM

Instrument	69679	Data File	mj cal 7.d
Type	Cal	Sample	mj cal 7
Acq. Method	AM 27 THC quant.m	Operator	Britany Wylie
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	5/11/2023 11:05:36 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	1246388	∞	835.59	∞	2659576	102.546 ng/ml
THC-COOH	1.358	952794	∞	250.65	∞	668336	254.447 ng/ml
THC	2.926	3697370	∞	24.25	17357.5	1111194	102.757 ng/ml