



data reviewed 6/16/16

B. Wylee

Worklist: 1157

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>	
C2016-0764	1	58598	3.10.3 Blood confirmation of co	
M2016-1315	1	58597	3.10.3 Blood confirmation of co	

**Idaho State Police
Forensic Services
Toxicology Discipline**

Request for Departure from an Analytical Method

Date of Request

05/16/2016

Forensic Scientist

Celena Shrum

Analytical Method

3.10 Manual Solid Phase Extraction (SPE) Methods

3.10.3 Extraction and Confirmation of Free (Unbound) Codeine and Morphine in Blood Employing the United Chemical Technologies (UCT) 200 mg CLEAN SCREEN® DAU Extraction Column

I request permission to use this method for full scan analysis only. A negative and positive control will still be run. All samples and controls will be run in full scan mode and any drugs in the samples will be confirmed qualitatively by comparing the full scan spectrum to that of authenticated reference material(s). This deviation is to remain in effect until permanent changes to the analytical method are completed.

Discipline Leader Review

Departure Approved

Comments:

Departure Not Approved

Comments:



Date: May 16, 2016

All vial positions verified.

simulate_sequence.log

Simulate Run Sequence Tue Jun 14 14:39:38 2016


Instrument Name: Major Mass Spec
 Sequence File: C:\Users\ISPuser\Desktop\Sequences\CS-0P061416.sequence.xml
 Comment: MassHunter sequence
 Operator: ISP\datastor
 Data Path: D:\DATA\CDS\2016\061416OP\
 Method Path: C:\Users\datastor\Desktop\OP Methods\

Line	Type	Vials	DataFile	Sample Name


Acquisition Method: AM3-10-3-FS.M				
1)	Sample	100	Blank 1	Pre-run Solvent Blank
2)	Sample	100	Blank 2	Pre-run Solvent Blank
3)	Sample	1	Negative Control	Negative Control -
...1013				
4)	Sample	100	C2016-0764-BLK	Lab No.: C2016-0764-1
5)	Sample	3	C2016-0764-FS	Lab No.: C2016-0764-1
6)	Sample	2	Spiked Positive Control	Positive Control
7)	Sample	99	Blank 3	Solvent Blank
8)	Sample	100	M2016-1315-BLK	Lab No.: M2016-1315-1
9)	Sample	4	M2016-1315-FS	Lab No.: M2016-1315-1
10)	Sample	99	Blank 4	Solvent Blank

megabytes Needed: 19 Space on drive D: 239346
 Sequence Verification Done!


File :I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Blank 1.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 14:46 using AcqMethod AM3-10-3-FS.M
Sample Name: Pre-run Solvent Blank
Misc Info : Chloroform





File :I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Blank 2.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 15:05 using AcqMethod AM3-10-3-FS.M
Sample Name: Pre-run Solvent Blank
Misc Info : Chloroform





File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Negative Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 15:23 using AcqMethod AM3-10-3-FS.M
Sample Name: Negative Control - Utak Lot B1013
Misc Info : UTAK B1013




File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Negative Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 15:23 using AcqMethod AM3-10-3-FS.M
Sample Name: Negative Control - Utak Lot B1013
Misc Info : UTAK B1013




File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Negative Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 15:23 using AcqMethod AM3-10-3-FS.M
Sample Name: Negative Control - Utak Lot B1013
Misc Info : UTAK B1013





File :I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Blank 3.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:36 using AcqMethod AM3-10-3-FS.M
Sample Name: Solvent Blank
Misc Info : Chloroform





File :I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Spiked Positive Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:18 using AcqMethod AM3-10-3-FS.M
Sample Name: Positive Control
Misc Info : UTAK B1013 + 061316 WS





File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Spiked Positive Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:18 using AcqMethod AM3-10-3-FS.M
Sample Name: Positive Control
Misc Info : UTAK B1013 + 061316 WS





File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Spiked Positive Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:18 using AcqMethod AM3-10-3-FS.M
Sample Name: Positive Control
Misc Info : UTAK B1013 + 061316 WS



File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Spiked Positive Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:18 using AcqMethod AM3-10-3-FS.M
Sample Name: Positive Control
Misc Info : UTAK B1013 + 061316 WS



File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Spiked Positive Control.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 16:18 using AcqMethod AM3-10-3-FS.M
Sample Name: Positive Control
Misc Info : UTAK B1013 + 061316 WS



File : I:\Instrument Data\Pocatello\Major Mass Spec\CDS\2016\061416
... OP\Blank 4.D
Operator : ISP\datastor
Instrument : Major Mass Spec
Acquired : 14 Jun 2016 17:31 using AcqMethod AM3-10-3-FS.M
Sample Name: Solvent Blank
Misc Info : Chloroform

