




Downloading performance verification data from the EasyCal calibration unit




Plug in a USB thumb drive into the EasyCal. This will automatically be recognized by the unit with a message "USB Drive plugged in"





Press Yes and a file will be downloaded to a USB key



For performance verification data from the EasyCal, the Tech log is not necessary so press "No"



Take the USB drive off and transfer the data into a computer/lap top



In a computer/laptop, open Excel. Open up the file that was in the USB drive. When finding the file choose "All Files" then select Print Log that was downloaded and select open. This file is a CSV or comma separated values file. The files are named at the time the file was downloaded in the USB drive. First choose "delimited" then next then check "comma" then finish. Check if this setting is working for your computer and if this is showing the data in a presentable manner. This will open an excel file.

SN	software	version	calDate	calTime	calStand	checkDate	checkTime	checkStan	checkRes	units	model	operator	expDate	lotCode	result	checkRes	checkRes	calibrateMode
90201895	6.5	#####	10:17	0.16	4/6/2016	11:15	0.08	0.077	BRAC	Sequence		Nov-17	29315080	Passed	0.077		2	
15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:20	0.08	0.08	BRAC	Sequence		Nov-17	29315080	Passed	0.079		2	
15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:24	0.2	0.195	BRAC	Sequence		Nov-17	29315200	Passed	0.193		2	
15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:34	0.04	0.04	BRAC	Sequence		Nov-17	29315040	Passed	0.04		2	
90205804	6.5	#####	9:53	0.2	#####	10:18	0.2	0.194	BRAC	Sequence		Nov-17	29315200	Passed	0.191		2	
90205741	6.5	#####	9:54	0.2	#####	10:24	0.2	0.196	BRAC	Sequence		Nov-17	29315200	Passed	0.192		2	
90205804	6.5	#####	9:53	0.2	#####	10:28	0.08	0.078	BRAC	Sequence		Nov-17	29315080	Passed	0.078		2	
90205741	6.5	#####	9:54	0.2	#####	10:32	0.08	0.079	BRAC	Sequence		Nov-17	29315080	Passed	0.078		2	
90205804	6.5	#####	9:53	0.2	#####	10:35	0.04	0.04	BRAC	Sequence		Nov-17	29315040	Passed	0.039		2	
90205741	6.5	#####	9:54	0.2	#####	10:40	0.04	0.04	BRAC	Sequence		Nov-17	29315040	Passed	0.04		2	
90205661	6.5	#####	10:32	0.2	#####	10:53	0.2	0.199	BRAC	Sequence		Nov-17	29315200	Passed	0.195		2	
15230014	8.7.21	#####	10:23	0.16	#####	10:57	0.08	0.081	BRAC	Sequence		Nov-17	29315080	Passed	0.08		2	
90205803	6.5	#####	10:10	0.2	#####	11:00	0.08	0.078	BRAC	Sequence		Nov-17	29315080	Passed	0.078		2	
90205663	6.5	#####	10:22	0.2	#####	11:04	0.08	0.079	BRAC	Sequence		Nov-17	29315080	Passed	0.078		2	
90205661	6.5	#####	10:32	0.2	#####	11:07	0.08	0.08	BRAC	Sequence		Nov-17	29315080	Passed	0.079		2	
15230014	8.7.21	#####	10:23	0.16	#####	11:11	0.04	0.041	BRAC	Sequence		Nov-17	29315040	Passed	0.04		2	
14200005	6.5	#####	10:43	0.2	#####	11:15	0.04	0.039	BRAC	Sequence		Nov-17	29315040	Passed	0.04		2	
90205803	6.5	#####	10:10	0.2	#####	11:18	0.04	0.039	BRAC	Sequence		Nov-17	29315040	Passed	0.039		2	
90205663	6.5	#####	10:22	0.2	#####	11:23	0.04	0.04	BRAC	Sequence		Nov-17	29315040	Passed	0.039		2	
90205661	6.5	#####	10:32	0.2	#####	11:26	0.04	0.039	BRAC	Sequence		Nov-17	29315040	Passed	0.04		2	
90205663	6.5	#####	13:43	0.16	#####	13:48	0.2	0.197	BRAC	Sequence		Nov-17	29315200	Passed	0.195		2	
90205803	6.5	#####	13:50	0.16	#####	13:52	0.2	0.199	BRAC	Sequence		Nov-17	29315200	Passed	0.197		2	

Adjust the width of the columns with "#####" to have enough space to show the information

PrintLog_041916_1419 - Excel

Bradley, Nikk

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	SN	softVersio	calDate	calTime	calStandai	checkDate	checkTime	checkStan	checkRes	units	model	operator	expDate	lotCode	result	checkRe
2	90201895	6.5	12/29/2015	10:17	0.16	4/6/2016	11:15	0.08	0.077	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
3	15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:20	0.08	0.08	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
4	15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:24	0.2	0.195	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
5	15360114	8.7.21	4/6/2016	11:12	0.16	4/6/2016	11:34	0.04	0.04	BrAC	Sequence		Nov-17 29315040A1	Passed	0.1	
6	90205804	6.5	4/19/2016	9:53	0.2	4/19/2016	10:18	0.2	0.194	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
7	90205741	6.5	4/19/2016	9:54	0.2	4/19/2016	10:24	0.2	0.196	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
8	90205804	6.5	4/19/2016	9:53	0.2	4/19/2016	10:28	0.08	0.078	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
9	90205741	6.5	4/19/2016	9:54	0.2	4/19/2016	10:32	0.08	0.079	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
10	90205804	6.5	4/19/2016	9:53	0.2	4/19/2016	10:35	0.04	0.04	BrAC	Sequence		Nov-17 29315040A1	Passed	0.0	
11	90205741	6.5	4/19/2016	9:54	0.2	4/19/2016	10:40	0.04	0.04	BrAC	Sequence		Nov-17 29315040A1	Passed	0.1	
12	90205661	6.5	12/4/2015	10:32	0.2	12/4/2015	10:53	0.2	0.199	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
13	15230014	8.7.21	12/4/2015	10:23	0.16	12/4/2015	10:57	0.08	0.081	BrAC	Sequence		Nov-17 29315080A2	Passed	0.1	
14	90205803	6.5	12/4/2015	10:10	0.2	12/4/2015	11:00	0.08	0.078	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
15	90205663	6.5	12/4/2015	10:22	0.2	12/4/2015	11:04	0.08	0.079	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
16	90205661	6.5	12/4/2015	10:32	0.2	12/4/2015	11:07	0.08	0.08	BrAC	Sequence		Nov-17 29315080A2	Passed	0.0	
17	15230014	8.7.21	12/4/2015	10:23	0.16	12/4/2015	11:11	0.04	0.041	BrAC	Sequence		Nov-17 29315040A1	Passed	0.1	
18	14200005	6.5	12/4/2015	10:43	0.2	12/4/2015	11:15	0.04	0.039	BrAC	Sequence		Nov-17 29315040A1	Passed	0.1	
19	90205803	6.5	12/4/2015	10:10	0.2	12/4/2015	11:18	0.04	0.039	BrAC	Sequence		Nov-17 29315040A1	Passed	0.0	
20	90205663	6.5	12/4/2015	10:22	0.2	12/4/2015	11:23	0.04	0.04	BrAC	Sequence		Nov-17 29315040A1	Passed	0.0	
21	90205661	6.5	12/4/2015	10:32	0.2	12/4/2015	11:26	0.04	0.039	BrAC	Sequence		Nov-17 29315040A1	Passed	0.1	
22	90205663	6.5	12/4/2015	13:43	0.16	12/4/2015	13:48	0.2	0.197	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
23	90205803	6.5	12/4/2015	13:50	0.16	12/4/2015	13:52	0.2	0.199	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	
24	90205661	6.5	12/4/2015	13:53	0.16	12/4/2015	13:55	0.2	0.195	BrAC	Sequence		Nov-17 29315200A4	Passed	0.1	

Save the file as either another CSV file, excel sheet or pdf depending on preference. Another option is to print the data to a regular printer. Page layout may need to be adjusted to allow based on presentation preferences.