


BWylee

Worklist: 1851

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2017-1469	2	92968	AM 27 Blood THC Quant by LC
C2017-1519	1	92965	AM 27 Blood THC Quant by LC
C2017-1656	1	92966	AM 27 Blood THC Quant by LC
M2017-3500	1	92967	AM 27 Blood THC Quant by LC
M2017-3613	1	92964	AM 27 Blood THC Quant by LC



[Handwritten signature]

Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 8-22-17

Analyst: Star Anne Nord

PRE-ANALYTIC

Plate Lot# Custom - 499102 Plate Exp. 1/29/2018 External QC Lot 21718 exp 2-17-18

- ✓ 1. Ensure all solutions are within expiration date.
 - Mobile Phase A: *0.1% Formic Acid in LCMS Water*
 - Mobile Phase B: *0.1% Formic Acid in LCMS Acetonitrile*
 - *LCMS Methanol*
 - **Blank/Negative Blood: Lot 321632-1**
 - *0.1% Formic Acid in water*
 - *MTBE*
 - *Hexane*Column: *UCT Selectra DA 100 x 2.1 mm 3um*
- ✓ 2. Check levels of mobile phases and needle wash and refill as necessary. Ensure waste is not full.
- ✓ 3. Purge Pump and Load appropriate Acq. Method, allow system to equilibrate for approx. 30 min.
- ✓ 4. Create worklist. Data path name: 82217 can screen+quant

ANALYTIC

- ✓ 1. Remove standards plate, blood, and samples from cold storage. Allow to reach room temperature.
- ✓ 2. Add **1000 µL blood** to wells of analytical (standards) plate. Place cover on Plate
 - Blank blood for locations containing standards/QCs and internal standards
 - Sample blood for locations containing only internal standards
- ✓ 3. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**. *Shaker ID 66759*
- ✓ 4. Pipette **500µL 0.1% formic acid** to all wells of standards plate.
- ✓ 5. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**.
- ✓ 6. Transfer **800µL of blood+acid mixture** to corresponding wells of SLE+ plate.
- ✓ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). Wait 5 min. (**Load blood samples at 85- 100 PSI- Selector to Right**)
- ✓ 8. Add **2.25mL MTBE** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 9. Apply positive pressure for approx. 15 seconds (**10-15 PSI- Selector to left -**).
- ✓ 10. Add **2.25mL Hexane** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 11. Apply positive pressure for approx. 15 seconds. (**10-15 PSI Selector to the left**)
- ✓ 12. Remove collection plate containing eluate.
- ✓ 13. Place collection plate on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID 66819*
- ✓ 14. Reconstitute in **100 µL MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

POST-ANALYTIC

- ✓ 1. Open quantitation software and create a new quantitation batch.
Batch name: 82217 can quant
- ✓ 2. Make any necessary integration changes. Limit curves based on validated linear ranges (3-50ng/mL).
- ✓ 3. Were all appropriate standards used in the curve for each analyte? Y / N
Are r^2 values ≥ 0.98 for each analyte? Y / N
- ✓ 4. Did all QCs pass for each analyte? (Y) / N Were QCs entered into QC charting? (Y) / N
- ✓ 5. Central File Packet to include: ✓ LIMS Worklist: ✓ Method Checklist ✓ Calibration and Control Reports

COMMENTS



ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 17:47 **Data File** Cal 1 - 3ng.d
Sample Type Calibration **Sample Name** Cal 1 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.416	6856	263226	0.0260	2.8415
THC-COOH	THC-COOH-d9	2.566	5826	86977	0.0670	2.9381
THC	THC-d3	6.813	3640	101956	0.0357	2.8331

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 17:59 **Data File** Cal 2 - 5ng.d
Sample Type Calibration **Sample Name** Cal 2 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.416	11924	264791	0.0450	4.7829
THC-COOH	THC-COOH-d9	2.586	10029	87899	0.1141	5.0523
THC	THC-d3	6.793	6073	99858	0.0608	4.8926

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 18:11 **Data File** Cal 3 - 10ng.d
Sample Type Calibration **Sample Name** Cal 3 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	22804	224255	0.1017	10.5753
THC-COOH	THC-COOH-d9	2.526	17223	77311	0.2228	9.9293
THC	THC-d3	6.733	9758	78433	0.1244	10.1073

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 18:23 **Data File** Cal 4 - 25ng.d
Sample Type Calibration **Sample Name** Cal 4 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	53183	213047	0.2496	25.7005
THC-COOH	THC-COOH-d9	2.486	46808	72772	0.6432	28.7960
THC	THC-d3	6.633	24856	80381	0.3092	25.2603

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 18:35 **Data File** Cal 5 - 50ng.d
Sample Type Calibration **Sample Name** Cal 5 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	116454	218360	0.5333	54.7038
THC-COOH	THC-COOH-d9	2.486	85644	73844	1.1598	51.9765
THC	THC-d3	6.593	54824	84238	0.6508	53.2670

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 18:46 **Data File** Cal 6 - 100ng.d
Sample Type Calibration **Sample Name** Cal 6 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	192630	219648	0.8770	89.8416
THC-COOH	THC-COOH-d9	2.486	138792	70165	1.9781	88.6957
THC	THC-d3	6.613	95626	77804	1.2291	100.6775

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 18:58 **Data File** Cal 7 - 250ng.d
Sample Type Calibration **Sample Name** Cal 7 - 250ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	521937	209777	2.4881	254.5544
THC-COOH	THC-COOH-d9	2.506	337919	64583	5.2323	234.7264
THC	THC-d3	6.613	228443	76121	3.0011	245.9622

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 19:22 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-COOH	THC-COOH-d9	2.406	14779	119770	0.1234	5.4697

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 19:34 **Data File** QC - 10ng.d
Sample Type QC **Sample Name** QC - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	22418	192586	0.1164	12.0800
THC-COOH	THC-COOH-d9	2.486	16169	63000	0.2567	11.4495
THC	THC-d3	6.653	11338	73029	0.1553	12.6352

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 19:46 **Data File** External Control lot 21717.d
Sample Type Sample **Sample Name** External Control lot 21717 8 A
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	15219	185762	0.0819	8.5552
THC-COOH	THC-COOH-d9	2.486	14097	61671	0.2286	10.1898
THC	THC-d3	6.633	7543	67271	0.1121	9.0999

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin
Analysis Time 8/24/2017 8:48 AM **Analyst Name** ISP Tox
Report Time 8/24/2017 8:50 AM **Reporter Name** ISP Tox
Last Calib Update 8/24/2017 8:48 AM **Batch State** Processed

Analysis Info

Acq Time 2017-08-23 21:56 **Data File** external control 91317 10 ng.d
Sample Type Sample **Sample Name** external control 91317 10 ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	29307	274107	0.1069	11.1101
THC-COOH	THC-COOH-d9	2.526	16604	86868	0.1911	8.5098
THC	THC-d3	6.713	11049	86275	0.1281	10.4065


ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin

Last Calib Update 8/24/2017 8:48 AM **Analyst Name** ISP TOX

Target Compound *THC-OH*

Internal Standard *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.8	94.7
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.8	95.7
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.6	105.8
QC - 10ng	3	<input checked="" type="checkbox"/>	10	12.1	120.8
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.7	102.8
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	54.7	109.4
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	89.8	89.8
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	254.6	101.8

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin

Last Calib Update 8/24/2017 8:48 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.9	97.9
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.1	101.0
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.9	99.3
QC - 10ng	3	<input checked="" type="checkbox"/>	10	11.4	114.5
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	28.8	115.2
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	52.0	104.0
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	88.7	88.7
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	234.7	93.9

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\82217 cann quant\QuantResults\82217 cann quant.batch.bin

Last Calib Update 8/24/2017 8:48 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.8	94.4
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.9	97.9
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.1	101.1
QC - 10ng	3	<input checked="" type="checkbox"/>	10	12.6	126.4
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.3	101.0
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	53.3	106.5
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	100.7	100.7
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	246.0	98.4