

Worklist: 2195

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
M2018-0220	1	107771	AM 27 Blood THC Quant by LC	
M2018-0229	1	107772	AM 27 Blood THC Quant by LC	
P2018-0245	2	107773	AM 27 Blood THC Quant by LC	
P2018-0332	1	107774	AM 27 Blood THC Quant by LC	
P2018-0372	1	107775	AM 27 Blood THC Quant by LC	
P2018-0381	1	107777	AM 27 Blood THC Quant by LC	
P2018-0462	1	107783	AM 27 Blood THC Quant by LC	

Quantitation of THC and Metabolites in Blood by LC-MS/MS

D

Extraction Date: 02/15/18

Analyst: Sarah Pickle

Plate lot#: 0515037

Plate Expiration: 09/28/18

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 361331-3

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 59740

Pre-Analytic:

1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
3. Create worklist:

Analytic:

1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
2. Pipette **1000µL blood (calibrated pipette) Pipette ID: #27** in wells of analytical (standards) plate.
3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
8. Wait 5 minutes.
9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
10. Wait 5 minutes.
11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
13. Wait 5 minutes.
14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 067103
16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

1. Create batch and process data.
Worklist path: *021518 THC Quant SP* Batch Name: *THC Quant SP 021518*
2. Make any necessary integration changes, r^2 values ≥ 0.98 for each analyte.
3. Did all QCs pass for each analyte? Y/N Enter QCs into control charting?
4. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:

Curve Range Limited: 5-250ng for THC-COOH

R

Toxicology AM method 27 external prep information
Working Solution 1 ug/ml in MeOH C-THC, THC-OH, THC
Stock Solution 1mg/ml 10 ul THC, 100ug/ml 100ul THC-OH, C-THC in 9790 ul MeOH (Fisher Lot # 172516)
Prepared: 02/07/18 Expires: 02/07/19 By: TS

Drug	Lot	Expiration
C-THC	497429	N/A
THC	FE04231406	4/30/2019
THC-OH	FE01121503	1/31/2020

AM 27 control 100 ul working solution lot (WS020718) in 9900ul blood (Hemostat Lot# 361331-2)
Prepared: 02/07/18 Expires: 02/07/19 By: TS

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


S

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:37 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 13:12	Data File	Negative.d
Sample Type	Sample	Sample Name	Negative
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-H11	Sample Info	
Inj Vol	-1	Comment	Hemostat 361331-3

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.025	716	299063	0.0024	0.5453 < 3
THC-COOH	THC-COOH-D9	2.299	3743	87486	0.0428	0.1554 < 10

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 12:49	Data File	QC Control.d
Sample Type	Sample	Sample Name	QC Control
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-A12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	33818	273974	0.1234	10.9096
THC-COOH	THC-COOH-D9	2.325	20547	79963	0.2570	10.5171
THC	THC-D3	5.919	13048	94663	0.1378	11.1901

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:37 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 13:36	Data File	External Control.d
Sample Type	Sample	Sample Name	External Control
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-G11	Sample Info	
Inj Vol	-1	Comment	Hemostat 361331-3 + WS 020718

Sample Chromatogram



Results


Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.212	28700	275771	0.1041	9.2516
THC-COOH	THC-COOH-D9	2.312	19011	77968	0.2438	9.8818
THC	THC-D3	5.812	11417	95101	0.1200	9.6547

Calibration Curve Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP
021518.batch.bin

Last Calib Update 2/16/2018 10:34 AM **Analyst Name** ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-D3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	3.0	99.4
Cal 2	2	<input checked="" type="checkbox"/>	5	5.2	104.6
Cal 3	3	<input checked="" type="checkbox"/>	10	10.0	100.4
Cal 4	4	<input checked="" type="checkbox"/>	25	26.4	105.6
Cal 5	5	<input checked="" type="checkbox"/>	50	45.9	91.9
Cal 6	6	<input checked="" type="checkbox"/>	100	95.1	95.1
Cal 7	7	<input checked="" type="checkbox"/>	250	257.3	102.9

Calibration Curve Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP
 021518.batch.bin

Last Calib Update 2/16/2018 10:34 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-D9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input type="checkbox"/>	3	2.5	81.8
Cal 2	2	<input checked="" type="checkbox"/>	5	5.6	112.2
Cal 3	3	<input checked="" type="checkbox"/>	10	10.0	99.7
Cal 4	4	<input checked="" type="checkbox"/>	25	25.1	100.2
Cal 5	5	<input checked="" type="checkbox"/>	50	44.5	89.0
Cal 6	6	<input checked="" type="checkbox"/>	100	95.0	95.0
Cal 7	7	<input checked="" type="checkbox"/>	250	259.9	104.0

Calibration Curve Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP
021518.batch.bin

Last Calib Update 2/16/2018 10:34 AM **Analyst Name** ISP TOX

Target Compound THC
Internal Standard THC-D3



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	2.5	82.2
Cal 2	2	<input checked="" type="checkbox"/>	5	5.2	104.7
Cal 3	3	<input checked="" type="checkbox"/>	10	10.5	105.1
Cal 4	4	<input checked="" type="checkbox"/>	25	27.9	111.6
Cal 5	5	<input checked="" type="checkbox"/>	50	48.0	95.9
Cal 6	6	<input checked="" type="checkbox"/>	100	101.5	101.5
Cal 7	7	<input checked="" type="checkbox"/>	250	247.4	99.0

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 11:14	Data File	Cal 1.d
Sample Type	Calibration	Sample Name	Cal 1
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-H12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.225	8838	286513	0.0308	2.9816
THC-COOH	THC-COOH-D9	2.325	7980	88380	0.0903	2.4538
THC	THC-D3	5.892	3815	103823	0.0367	2.4660

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


PS

Batch Data Path C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin
Analysis Time 2/16/2018 10:34 AM **Analyst Name** ISPUser
Report Time 2/16/2018 10:36 AM **Reporter Name** ISPUser
Last Calib Update 2/16/2018 10:34 AM **Batch State** Processed

Analysis Info

Acq Time 2018-02-15 11:26 **Data File** Cal 2.d
Sample Type Calibration **Sample Name** Cal 2
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-G12 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	15325	268266	0.0571	5.2319
THC-COOH	THC-COOH-D9	2.325	12285	79014	0.1555	5.6078
THC	THC-D3	5.932	6275	91148	0.0688	5.2356

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin
Analysis Time 2/16/2018 10:34 AM **Analyst Name** ISPUser
Report Time 2/16/2018 10:36 AM **Reporter Name** ISPUser
Last Calib Update 2/16/2018 10:34 AM **Batch State** Processed

Analysis Info

Acq Time 2018-02-15 11:38 **Data File** Cal 3.d
Sample Type Calibration **Sample Name** Cal 3
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-F12 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	28403	250651	0.1133	10.0433
THC-COOH	THC-COOH-D9	2.325	18307	74525	0.2457	9.9701
THC	THC-D3	5.892	11095	85384	0.1299	10.5089

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 11:50	Data File	Cal 4.d
Sample Type	Calibration	Sample Name	Cal 4
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-E12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	72577	238386	0.3045	26.4095
THC-COOH	THC-COOH-D9	2.339	39495	70854	0.5574	25.0536
THC	THC-D3	5.905	26798	80829	0.3315	27.9061

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 12:01	Data File	Cal 5.d
Sample Type	Calibration	Sample Name	Cal 5
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-D12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	154373	289820	0.5327	45.9496
THC-COOH	THC-COOH-D9	2.339	80338	83735	0.9594	44.5039
THC	THC-D3	5.919	55805	98960	0.5639	47.9603

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 12:13	Data File	Cal 6.d
Sample Type	Calibration	Sample Name	Cal 6
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-C12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	323069	291983	1.1065	95.0835
THC-COOH	THC-COOH-D9	2.339	159529	79672	2.0023	94.9598
THC	THC-D3	5.919	114025	96305	1.1840	101.4732

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\021518 THC Quant SP\QuantResults\THC Quant SP 021518.batch.bin		
Analysis Time	2/16/2018 10:34 AM	Analyst Name	ISPUser
Report Time	2/16/2018 10:36 AM	Reporter Name	ISPUser
Last Calib Update	2/16/2018 10:34 AM	Batch State	Processed

Analysis Info

Acq Time	2018-02-15 12:25	Data File	Cal 7.d
Sample Type	Calibration	Sample Name	Cal 7
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-B12	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.239	884341	294690	3.0009	257.3005
THC-COOH	THC-COOH-D9	2.339	433334	80075	5.4116	259.9049
THC	THC-D3	5.919	313326	108964	2.8755	247.4499