


Worklist: 2766

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2018-2060	1	130703	AM 27 Blood THC Quant by LC-QQQ
C2018-2089	1	130699	AM 27 Blood THC Quant by LC-QQQ
C2018-2101	1	130700	AM 27 Blood THC Quant by LC-QQQ
C2018-2113	1	130701	AM 27 Blood THC Quant by LC-QQQ
C2018-2163	1	130702	AM 27 Blood THC Quant by LC-QQQ
C2018-2172	1	130704	AM 27 Blood THC Quant by LC-QQQ



AN

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 10/31/18
Plate lot#: 0539904

Analyst: Britany Wylie
Plate Expiration: 09/10/2019

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile
LCMS Methanol
Hexane
Blank Blood Lot: 18G207D7
Column: UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 62340

Pre-Analytic:

1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
3. Create worklist:

Analytic:

1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
8. Wait 5 minutes.
9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
10. Wait 5 minutes.
11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
13. Wait 5 minutes.
14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

1. Create batch and process data.
Batch path: cann quant Batch Name: _____
2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
5. Did all QCs pass for each analyte? Y N
6. Enter QCs into control charting.
7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:

Toxicology AM method 27 external prep information
working solution 1 ug/ml in meoh C-THC, THC-OH, THC
Stock solution 1mg/ml 10 ul each THC, THC-OH 100 ug/ml 100 ul C-THC in 9890 ul meOH
Ppd 6/5/18 Exp: 4/1/19 lot 6518 by AMN

Drug	lot (certifiant)	expiration
C-THC	FE03121501	3/1/2020
THC-OH	FE01141502	1/1/2020
THC	FE04231406	4/1/2019

AM 27 control 50 ul working solution lot (6518) in 4950 ul blood lot (17120718)
ppd 6/5/18 Exp 4/1/19 lot 6518 Concentration 10 ng/ml each by AMN
ppd 9/20/18 Exp 4/1/19 lot 92018 Concentration 10 ng/ml each by AMN

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 11:29 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



ISP

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 11:41 **Data File** QC - 10ng.d
Sample Type QC **Sample Name** QC - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.256	4669	115678	0.0404	5.0459
THC-COOH	THC-COOH-d9	2.306	7340	46027	0.1595	9.2502
THC	THC-d3	5.833	2395	45457	0.0527	4.8300

BLW

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 11:53 **Data File** External control.d
Sample Type Sample **Sample Name** External control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	8131	100891	0.0806	9.8323
THC-COOH	THC-COOH-d9	2.326	6231	41071	0.1517	8.7483
THC	THC-d3	5.833	4844	38574	0.1256	12.0675

BW

ISP Forensics Calibration Curve Report

Batch Data Path

D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update

11/5/2018 8:30 AM

Analyst Name


ISP TOX

Target Compound

THC-OH

Internal Standard

THC-OH-d3




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 1ng	1	<input type="checkbox"/>	1	1.2	117.6
Cal 2 - 3ng	2	<input checked="" type="checkbox"/>	3	3.3	109.2
Cal 3 - 5ng	3	<input checked="" type="checkbox"/>	5	4.9	98.7
QC - 10ng	3	<input checked="" type="checkbox"/>	5	5.0	100.9
Cal 4 - 10ng	4	<input checked="" type="checkbox"/>	10	9.3	93.5
Cal 5 - 25ng	5	<input checked="" type="checkbox"/>	25	24.4	97.6
Cal 6 - 50ng	6	<input checked="" type="checkbox"/>	50	50.1	100.1
Cal 7 - 100ng	7	<input checked="" type="checkbox"/>	100	101.0	101.0

ISP Forensics Calibration Curve Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update 11/5/2018 8:30 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 1ng	1	<input type="checkbox"/>	1	1.2	122.0
Cal 2 - 3ng	2	<input checked="" type="checkbox"/>	3	3.2	106.8
Cal 3 - 5ng	3	<input checked="" type="checkbox"/>	5	4.8	95.9
QC - 10ng	3	<input checked="" type="checkbox"/>	5	9.3	185.0
Cal 4 - 10ng	4	<input checked="" type="checkbox"/>	10	9.9	99.1
Cal 5 - 25ng	5	<input checked="" type="checkbox"/>	25	24.7	98.9
Cal 6 - 50ng	6	<input checked="" type="checkbox"/>	50	48.9	97.9
Cal 7 - 100ng	7	<input checked="" type="checkbox"/>	100	101.4	101.4


BW

ISP Forensics Calibration Curve Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update 11/5/2018 8:30 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 1ng	1	<input type="checkbox"/>	1	0.9	93.0
Cal 2 - 3ng	2	<input checked="" type="checkbox"/>	3	3.0	100.6
Cal 3 - 5ng	3	<input checked="" type="checkbox"/>	5	5.1	101.3
QC - 10ng	3	<input checked="" type="checkbox"/>	5	4.8	96.6
Cal 4 - 10ng	4	<input checked="" type="checkbox"/>	10	9.9	98.8
Cal 5 - 25ng	5	<input checked="" type="checkbox"/>	25	24.8	99.4
Cal 6 - 50ng	6	<input checked="" type="checkbox"/>	50	49.8	99.5
Cal 7 - 100ng	7	<input checked="" type="checkbox"/>	100	100.4	100.4

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 09:54 **Data File** Cal 1 - 1ng.d
Sample Type Calibration **Sample Name** Cal 1 - 1ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	853	108735	0.0078	1.1762
THC-COOH	THC-COOH-d9	2.306	1519	43171	0.0352	1.2198
THC	THC-d3	5.773	543	40456	0.0134	0.9304

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 10:06 **Data File** Cal 2 - 3ng.d
Sample Type Calibration **Sample Name** Cal 2 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	3372	132315	0.0255	3.2752
THC-COOH	THC-COOH-d9	2.326	3482	52839	0.0659	3.2045
THC	THC-d3	5.833	1807	52462	0.0344	3.0177

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 10:18 **Data File** Cal 3 - 5ng.d
Sample Type Calibration **Sample Name** Cal 3 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.216	5287	134154	0.0394	4.9328
THC-COOH	THC-COOH-d9	2.306	4633	51162	0.0905	4.7971
THC	THC-d3	5.813	2823	51288	0.0550	5.0631

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 10:30 **Data File** Cal 4 - 10ng.d
Sample Type Calibration **Sample Name** Cal 4 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	9993	130607	0.0765	9.3467
THC-COOH	THC-COOH-d9	2.326	8614	50780	0.1696	9.9062
THC	THC-d3	5.833	5365	51799	0.1036	9.8820

BW

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 10:42 **Data File** Cal 5 - 25ng.d
Sample Type Calibration **Sample Name** Cal 5 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	25967	127904	0.2030	24.3994
THC-COOH	THC-COOH-d9	2.326	20064	50296	0.3989	24.7192
THC	THC-d3	5.853	12935	50880	0.2542	24.8423

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 10:54 **Data File** Cal 6 - 50ng.d
Sample Type Calibration **Sample Name** Cal 6 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	51633	123323	0.4187	50.0603
THC-COOH	THC-COOH-d9	2.326	36500	47177	0.7737	48.9298
THC	THC-d3	5.833	25332	50145	0.5052	49.7647

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\10312018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 11/5/2018 8:30 AM **Analyst Name** ISP Tox
Report Time 11/5/2018 8:31 AM **Reporter Name** ISP Tox
Last Calib Update 11/5/2018 8:30 AM **Batch State** Processed

Analysis Info

Acq Time 2018-11-02 11:06 **Data File** Cal 7 - 100ng.d
Sample Type Calibration **Sample Name** Cal 7 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.236	103226	121919	0.8467	100.9857
THC-COOH	THC-COOH-d9	2.326	73133	46096	1.5865	101.4433
THC	THC-d3	5.833	48477	47744	1.0154	100.4302