REVIEWED


By Sarah Pickle at 8:35 am, Aug 05, 2019



7/24/2019

Worklist: 3564

LAB CASE	<u>ITEM</u>	TASK ID	DESCRIPTION
M2019-2774	1	158081	AM 27 Blood THC Quant by LC-QQQ
M2019-3002	1	158082	AM 27 Blood THC Quant by LC-QQQ
P2019-1944	1	158083	AM 27 Blood THC Quant by LC-QQQ
P2019-2137	1	158084	AM 27 Blood THC Quant by LC-QQQ
P2019-2147	1	158085	AM 27 Blood THC Quant by LC-QQQ
P2019-2183	1	158086	AM 27 Blood THC Quant by LC-QQQ
			·
P2019-2188	1	158087	AM 27 Blood THC Quant by LC-QQQ





AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 07/25/2019

Analyst: Tamara Salazar

Plate lot#: 0539904

Plate Expiration: 09/10/19

MTBE

Mobile phase A: 0.1% Formic Acid in LCMS Water

Mobile phase B: 0.1% Formic acid in Acetonitrile

Hexane

Blank Blood Lot: Hemostat 445283-2

LCMS Methanol

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 069901

Pre-Analytic:

□ I. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.

⊠ 3. Create worklist:

Analytic:

☑ 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.

Σ 2. Pipette 1000μL blood/urine (calibrated pipette) Pipette ID: 3 in wells of analytical (standards) plate.

□ 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. Shaker ID: 067105

Δ 4. Pipette 500μL 0.1% formic acid in water in wells of analytical plate for blood samples.

(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104

⊠ 8. Wait 5 minutes.

□ Signature Street S

 \boxtimes 10. Wait 5 minutes.

□ 11. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).

□ 12. Add 2.25mL Hexane. (Add in 3 increments of 750uL)

 \boxtimes 13. Wait 5 minutes.

△ 14. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).

 □ 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. SPE Dry ID: 067103

Δ 16. Reconstitute in 100 μL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

Worklist path: D:\MassHunter\Data\2019\AM 27\072319 wklst 3556 MDS TS reinjects Batch Name THCQ TS

- \boxtimes 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- ☑ 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- ☑ 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- Did all QCs pass for each analyte? Y/N **⊠** 5.
- Enter OCs into control charting.
- \boxtimes 7 Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Curves limited: THC-COOH 10-100





Idaho State Police Forensic Services

AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

Methanol External Control Solution (Lot: WS041619)

10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH Approximate concentration 1ug/mL.

		T .	77 1 11
Component	Source	Source Lot Number	Expiration Date
Methanol (LCMS)	Fisher	184782	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2020
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	04/16/2019		
Prepared By:	Tamara Salaza	r	
Expires:	01/31/2020		

Blood External Control Solution (Lot: 072319)

100 ul of methanol external control solution was added to 9900 ul of blood.

Approximately 10ng/mL of each compound.

Component	Source	Source Lot Number
Negative Blood	Hemostat	445283-2
Methanol External Control Solution	*	WS041619
Prepared:	07/23/19	
Prepared by:	Tamara Salaza	ır
Expires:	01/31/2020	



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

7/24/2019 6:18:27 PM Sample P3-H1 Injection Volume Sample Position Acq. Date-Time Acq. Method Instrument Type

AM 27 THC quant.m

Comment

Data File

Sample

THC_QC Control.d THC_QC Control

Sample Chromatogram

Sample Info.

+ TIC MRM (** -> **) THC_QC Control.d (THC_QC Control)



Resp.	138003	52393	73305
RT	3.912	1.730	1.693

THC THC-COOH THC-OH

Name

































Acquisition Time (min)

3.8

3.6

3.4

3.2

2.8

2.6

2.4

2.2

2

Final Conc.

3684098 ISTD Resp.





D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

Sample

AM 27 THC quant.m

P3-A2

Injection Volume Sample Position

Acq. Method

Instrument

Type

Acq. Date-Time

Sample Info.

7/24/2019 6:33:37 PM

THC_Negative.d THC_Negative

Data File Sample

Comment

Sample Chromatogram

+ TIC MRM (** -> **) THC_Negative.d (THC_Negative)



Resp. 16242 **RT** 1.768 High

Name THC-OH

8

S/N

Ratio 3.3 Low

S/N 3.22 **Low**

ISTD Resp. 845587

Final Conc. 0.4328 ng/ml <

63

Acquisition Time (min)

3.8

3.4

3.2

က

2.8

2.6

2.4

2.2

2

<u>←</u>



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

7/24/2019 6:48:48 PM AM 27 THC quant.m Sample P3-B2 10 Injection Volume Acq. Method Sample Position Acq. Date-Time Sample Info. Instrument Type


Data File Sample

Comment

THC_External Control-1.d THC_External Control-1

Sample Chromatogram





Acquisition Time (min)

Final Conc. 9.2197 ng/ml 8.9869 ng/ml 8.7720 ng/ml

ISTD Resp.4326349
271688
962040

S/N 132.99 481.45 362.79

4.2

3.8

3.6

3.4

3.2

က

2.8

2.6

AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ

TS.batch.bin

Last Cal. Update

7/26/2019 11:40 AM

Analyst Name


ISP\datastor

Analyte

THC

Internal Standard

THC-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
THC_Cal 1-3ng	1	✓	3.0	3.0	100.2
THC_Cal 2- 5ng	2	✓	5.0	4.8	96.2
THC_Cal 3 -10ng	3	✓	10.0	9.5	95.4
THC_Cal 4-25ng	4	✓	25.0	28.0	112.0
THC Cal 5-50ng	5	√	50.0	48.5	97.0
THC_Cal 6-100ng	6	/	100.0	99.1	99.1

AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ

TS.batch.bin

Last Cal. Update

7/26/2019 11:40 AM

Analyst Name


ISP\datastor

Analyte

THC-COOH

Internal Standard

THC-COOH-D9



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
THC_Cal 1-3ng	1	×	3.0	3.2	105.4
THC_Cal 2- 5ng	2	×	5.0	5.0	99.3
THC_Cal 3 -10ng	3	✓	10.0	8.7	87.1
THC_Cal 4-25ng	4	✓	25.0	28.8	115.3
THC Cal 5-50ng	5	✓	50.0	50.1	100.2
THC_Cal 6-100ng	6	✓	100.0	97.3	97.3

AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ

TS.batch.bin

Last Cal. Update

7/26/2019 11:40 AM

Analyst Name


ISP\datastor

Analyte

THC-OH

Internal Standard

THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
THC_Cal 1-3ng	1	✓	3.0	3.0	99.8
THC_Cal 2- 5ng	2	✓	5.0	4.7	94.8
THC_Cal 3 -10ng	3	✓	10.0	9.7	96.9
THC_Cal 4-25ng	4	√	25.0	28.0	112.1
THC_Cal 5-50ng	5	✓	50.0	48.8	97.7
THC_Cal 6-100ng	6	1	100.0	98.7	98.7



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

AM 27 THC quant.m Sample Position Acq. Method Instrument

Type

P3-B1

Comment

THC_Cal 1-3ng_r.d THC_Cal 1-3ng

Data File


Sample

7/24/2019 8:58:19 PM

Acq. Date-Time Sample Info.

Injection Volume

Sample Chromatogram



178.00

Acquisition Time (min)

3.8



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Data File Sample Calibration Last Update Batch results


7/24/2019 9:05:54 PM AM 27 THC quant.m P3-C1 10 Injection Volume Sample Position Acq. Date-Time Acq. Method Instrument Гуре

THC_Cal 2- 5ng_r.d THC_Cal 2- 5ng

Comment

Sample Chromatogram

Sample Info.



Acquisition Time (min)

3.8

4.8106 ng/ml 4.9659 ng/ml 4.7410 ng/ml

Final Conc.



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

Comment Data File Sample 7/24/2019 5:47:59 PM AM 27 THC quant.m P3-D1 Injection Volume Sample Position Acq. Date-Time Acq. Method Instrument Type

THC_Cal 3 -10ng.d THC_Cal 3 -10ng

Sample Info.

Final Conc. 9.5407 ng/ml 8.7142 ng/ml 9.6944 ng/ml 3.8 **ISTD Resp.**3746092
252010
859623 3.4 3.2 **S/N** 495.15 370.41 269.50 က 2.8 5.6 Ratio 28.7 52.4 11.8 2.2 **S/N** 850.04 351.64 1298.63 2 6 **Resp.** 288271 49400 126790 9. + TIC MRM (** -> **) THC_Cal 3 -10ng.d (THC_Cal 3 -10ng) 1.4 1.2 **RT** 3.882 1.715 1.663 0.8 Sample Chromatogram 9.0 THC-COOH THC-OH 3.5 2.5 1.5 1.5 Counts × 105 3 Name

Acquisition Time (min)

4.2



D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

AM 27 THC quant.m

P3-E1 10 7/24/2019 5:55:43 PM

Injection Volume Sample Position

Acq. Method

Instrument

Type

Acq. Date-Time


Sample Info.

THC_Cal 4-25ng.d THC_Cal 4-25ng

Data File Sample

Comment

Sample Chromatogram



me RT Resp. S/N Ratio C 3.912 1056265 8187.51 27.6 C-COOH 1.730 163985 465.93 54.9 C-OH 1.693 419950 636.96 12.5	<u>+</u>	0.6	8:0	-	12	4.	1.6	8:	- 73	2:2	2.4	2.6	2
3.912 1056265 8187.51 1.730 163985 465.93 1.693 419950 636.96	пе			RT			tesp.		S/N		Ratio		
1.730 163985 465.93 1.693 419950 636.96			17)	3.912		105	6265	ω	187,51		27.6		
1.693 419950 636.96	-C00H		7-1	1.730		16	3385		465,93		54.9		
	등		•	1.693		41	.9950		96'989		12.5		

Acquisition Time (min)

3.8

28.0058 ng/ml 28.8170 ng/ml 28.0213 ng/ml

Final Conc.




D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

THC_Cal 5-50ng.d THC_Cal 5-50ng

Comment Data File Sample 7/24/2019 6:03:18 PM AM 27 THC quant.m P3-F1 Injection Volume Sample Position Acq. Date-Time Sample Info. Acq. Method Instrument Гуре

Sample Chromatogram





D:\MassHunter\Data\2019\AM 25\072319 wklst 3556 MDS TS_reinjects\QuantResults\THCQ TS.batch.bin 7/26/2019 11:40:49 AM Calibration Last Update Batch results

THC_Cal 6-100ng.d THC_Cal 6-100ng

 Instrument
 Falco
 Data File

 Type
 Cal
 Sample

 Acq. Method
 AM 27 THC quant.m
 Comment

 Sample Position
 P3-G1
 Comment

 Injection Volume
 10
 Comment

 Acq. Date-Time
 7/24/2019 6:10:52 PM
 Sample Info.

Sample Chromatogram

+ TIC MRM (** -> **) THC_Cal 6-100ng.d (THC_Cal 6-100ng)

