

Worklist: 4526

LAB CASE	<u>ITEM</u>	ITEM TYPE	DESCRIPTION	
P2020-2795	1.10	CBUK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
P2020-2795	1.3	CBUK	AM 27 Blood THC Quant by LC-QQQ	

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 09/17/20 Plate lot#: IDP-108-2-200723 Mobile phase A: 0.1% Formic Acid in LCMS Water Blank Blood Lot: Hemostat 445283-4 Blank Urine Lot: POCO31319 LCMS-QQQ ID: 069901 Analyst: <u>Sophia Jackson</u> Plate Expiration: 01/23/2021 **Mobile phase B:** 0.1% Formic acid in Acetonitrile **Column**: UCT Selectra DA 100 x 2.1mm 3um

Pre-Analytic:

- ☑ 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- ☑ 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- ☑ 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- ☑ 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes. Using a calibrated pipette, add 1000µl blood and urine (if applicable) (calibrated pipette) into the appropriate wells of analytical (standards) plate. Pipette ID: 3382167
- ☑ 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- A. Pipette 500μL 0.1% formic acid in water blood sample, 500 μL saturated phosphate buffer in urine in wells of analytical plate.
- \boxtimes 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- Δ 6. Transfer 700-800μL of blood+acid or urine+acid mixture to corresponding wells of SLE+ plate. Amount transferred: 800 uL
- ☑ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
 (Load at 85-100 PSI- Selector to the right)
- \boxtimes 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. (Add in 3 increments of 750uL)
- \boxtimes 10. Wait 5 minutes.
- ☑ 11. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- ☑ 12. Add 2.25mL Hexane. (Add in 3 increments of 750uL)
- \boxtimes 13. Wait 5 minutes.
- ☑ 14. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- ☑ 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- ☑ 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- \boxtimes 1. Create batch and process data.
- \boxtimes 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values \ge 0.98 for each analyte
- ☑ 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- ☑ 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 5ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- ☑ 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- \boxtimes 6 Enter QCs into control charting.
- 27 Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Did not evaluate THC-OH Curves limited: THC 3-100*




 Batch results
 D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin

 Calibration Last Update
 9/21/2020 9:59:51 AM

InstrumentFalcoTypeSampleAcq. MethodAM 27 TSample PositionP3-A6Injection Volume10Acq. Date-Time9/17/20Sample Info.Sample Info.

Falco Sample AM 27 THC quant.m P3-A6 10 9/17/2020 7:37:45 PM Data File Sample Operator Comment MJ_QC Control.d MJ_QC Control Sophia Jackson

Sample Chromatogram




* Compound not evaluated





Batch results D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Type Acq. Method Sample Position Injection Volume Acq. Date-Time Sample Info. Falco Sample AM 27 THC quant.m P3-H5 10 9/17/2020 7:52:57 PM Data File Sample Operator Comment MJ_Negative Blood.d MJ_Negative Blood Sophia Jackson





Idaho State Police Forensic Services

AM #26 Screening of THC and Metabolites and AM #27 Confirmation of THC and Metabolites Urine External Control Prep Sheet

Methanol External Control Solution (Lot: WS011620)

10 μL of 1mg/mL THC, 100 μL of 100 μg/mL THC-OH, C-THC in 9790 μL MeOH Approximate concentration 1ug/mL.

Component	Source	Source Lot Number	Expiration Date	
Methanol (LCMS)	Fisher	193941		
THC	Cerilliant	FE09101501	11/30/2020	
C-THC	Cerilliant	FE07171501	09/30/2020	
THC-OH	Cerilliant	FE07221601	07/31/2021	
Prepared:	01/16/2020			
Prepared By:	Tamara Salazar			
Expires:	09/30/2020			

Urine External Control Solution (Lot: 090120)

200 ul of methanol external control solution was added to 9800 ul of urine. Approximately 20ng/mL each

Component	Source	Source Lot Number	
Negative Urine	Pocatello Lab	POC031319	
Methanol External Control Solution	-	WS011620	
Prepared:	09/01/2020		
Prepared by:	Sarah Pickle		





Batch results


D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Type Acq. Method Sample Position **Injection Volume** Acq. Date-Time Sample Info.

Falco Sample AM 27 THC quant.m P3-G5 10 9/17/2020 8:23:22 PM

Data File Sample Operator Comment MJ external urine.d MJ external urine Sophia Jackson

Sample Chromatogram



* Compound not evaluated




 Batch results
 D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin

 Calibration Last Update
 9/21/2020 9:59:51 AM

Instrument Type Acq. Method Sample Position Injection Volume Acq. Date-Time Sample Info. Falco Sample AM 27 THC quant.m P3-F5 10 9/17/2020 8:08:10 PM


Data File Sample Operator Comment MJ_Negative urine.d MJ_Negative urine Sophia Jackson





AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results	D:\MassHunter\Data\2020\AM 27-28\091720 A 27.batch.bin	M 27 AM 28 P1 P2 SJ con	np test\QuantResults\AM
Last Cal. Update	9/21/2020 9:59 AM		
Analyst Name	ISP\datastor		
Analyte	THC	Internal Standard	THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	×	1.0	1.3	125.0
MJ_Cal 2	2	√	3.0	3.1	104.0
MJ_Cal 3	3	. 1	5.0	5.1	101.3
MJ_Cal 4	4	1	10.0	9.7	96.6
MJ_Cal 5	5	1	25.0	24.3	97.3
MJ_Cal 6	6	1	50.0	50.0	99.9
MJ_Cal 7	7	√	100.0	100.9	100.9



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results	D:\MassHunter\Data\2020\AM 27-28\091720 AI 27.batch.bin	/I 27 AM 28 P1 P2 SJ cor	np test\QuantResults\AM
Last Cal. Update	9/21/2020 9:59 AM		
Analyst Name	ISP\datastor		
Analyte	THC-COOH	Internal Standard	THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	1	5.0	5.3	105.8
MJ_Cal 2	2	1	10.0	10.4	103.7
MJ_Cal 3	3	1	20.0	19.6	98.0
MJ_Cal 4	4	√	50.0	47.9	95.7
MJ_Cal 5	5	1	75.0	72.4	96.5
MJ_Cal 6	6	1	100.0	97.5	97.5
MJ_Cal 7	7	√	250.0	257.0	102.8



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results	D:\MassHunter\Data\2020\AM 27-28\091720 A 27.batch.bin	M 27 AM 28 P1 P2 SJ con	np test\QuantResults\AM
Last Cal. Update	9/21/2020 9:59 AM		
Analyst Name	ISP\datastor		
Analyte	THC-OH	Internal Standard	THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	×	1.0	0.8	[.] 83.6
MJ_Cal 2	2	√	3.0	2.9	96.6
MJ_Cal 3	3	√	5.0	5.3	106.9
MJ_Cal 4	4	√	10.0	9.7	97.2
MJ_Cal 5	5	√	25.0	24.7	98.8
MJ_Cal 6	6	√	50.0	50.2	100.3
MJ_Cal 7	7	√	100.0	100.2	100.2

Did not evaluate






Batch results

D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Falco Type Cal Acq. Method **Sample Position** P3-H6 **Injection Volume** 10 Acq. Date-Time Sample Info.

AM 27 THC quant.m 9/17/2020 6:36:45 PM

Data File Sample Operator Comment MJ_Cal 1.d MJ_Cal 1 Sophia Jackson




Batch results

D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Falco Instrument Туре Cal Acq. Method Sample Position P3-G6 **Injection Volume** 10 Acq. Date-Time Sample Info.

AM 27 THC quant.m 9/17/2020 6:44:29 PM

Data File Sample Operator Comment MJ_Cal 2.d MJ Cal 2 Sophia Jackson






Batch results

D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Falco Type Cal Acq. Method Sample Position P3-F6 **Injection Volume** 10 Acq. Date-Time Sample Info.

AM 27 THC quant.m 9/17/2020 6:52:05 PM Data File Sample Operator Comment MJ_Cal 3.d MJ_Cal 3 Sophia Jackson






Batch results

D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Falco Type Cal Acq. Method AM 27 THC quant.m Sample Position P3-E6 **Injection Volume** 10 Acq. Date-Time 9/17/2020 6:59:41 PM Sample Info.

Data File Sample Operator Comment MJ Cal 4.d MJ_Cal 4 Sophia Jackson






 Batch results
 D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin

 Calibration Last Update
 9/21/2020 9:59:51 AM

Instrument	Falco
Туре	Cal
Acq. Method	AM 27 THC quant.m
Sample Position	P3-D6
Injection Volume	10
Acq. Date-Time	9/17/2020 7:07:18 PM
Sample Info.	

Data File Sample Operator Comment MJ_Cal 5.d MJ_Cal 5 Sophia Jackson








D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin **Batch results** Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Falco Type Cal Acq. Method **Sample Position** P3-C6 **Injection Volume** 10 Acq. Date-Time Sample Info.

AM 27 THC quant.m 9/17/2020 7:14:54 PM Data File Sample Operator Comment MJ_Cal 6.d MJ_Cal 6 Sophia Jackson





Batch results

D:\MassHunter\Data\2020\AM 27-28\091720 AM 27 AM 28 P1 P2 SJ comp test\QuantResults\AM 27.batch.bin Calibration Last Update 9/21/2020 9:59:51 AM

Instrument Falco Туре Cal Acq. Method AM 27 THC quant.m Sample Position P3-B6 **Injection Volume** 10 Acq. Date-Time 9/17/2020 7:22:31 PM Sample Info.

Data File Sample Operator Comment MJ_Cal 7.d MJ_Cal 7 Sophia Jackson

