

**REVIEWED**


By Anne Nord at 1:15 pm, Apr 30, 2021

4/29/2021

*Byylee*

**Worklist: 4939**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2021-0871	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0876	6	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0879	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0888	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0931	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0974	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0975	1	BCK	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood and Urine

## LC-MS/MS

Extraction Date: 04/29/21

Analyst: Britany Wylie

Plate lot#: 201412 210412 <sup>BW</sup>

Plate Expiration: ~~6-6-2021~~ 10-12-21 <sup>BW</sup>

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol

Hexane

**Blank Blood Lot:** 20K20702

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 69679

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Blood only batch; curves limitations: THC-OH 3-100

*BWylee*

	1	2	3	4	5	6
A	IS + Cal. 1	neg blood				IS + QC_1
B	IS + Cal. 2	871-1				IS + Cal. 7
C	IS + Cal. 3	876-6				IS + Cal. 6
D	IS + Cal. 4	879-1				IS + Cal. 5
E	IS + Cal. 5	888-1				IS + Cal. 4
F	IS + Cal. 6	931-1				IS + Cal. 3
G	IS + Cal. 7	974-1				IS + Cal. 2
H	IS + QC_1	975-1				IS + Cal. 1

All wells to contain 100 µl of residual DMSO

Case #:   C2021-0  -1


# AM #27 Cannabinoids

*BWylie*

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj blood neg.d
<b>Type</b>	Sample	<b>Sample</b>	mj blood neg
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:36:11 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

BWylie

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood 2nd.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood 2nd
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 9:16:17 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	317601	∞	11.0	∞	2275368	4.718 ng/ml
THC-COOH	1.476	163025	305.9	35.8	859.9	837437	14.308 ng/ml
THC	3.317	183243	∞	24.3	∞	1618508	4.295 ng/ml

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Last Cal. Update** 4/29/2021 9:26 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	1.0	1.2	121.9
mj cal2	2	✓	3.0	2.8	94.7
mj cal 3	3	✓	5.0	4.6	92.8
mj cal 4	4	✓	10.0	9.3	92.8
mj cal 5	5	✓	25.0	23.9	95.6
mj cal 6	6	✓	50.0	50.1	100.1
mj cal 7	7	✓	100.0	102.0	102.0


# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Last Cal. Update** 4/29/2021 9:26 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH

**Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	5.0	5.5	110.6
mj cal2	2	✓	10.0	9.8	97.9
mj cal 3	3	✓	20.0	19.1	95.4
mj cal 4	4	✓	50.0	48.4	96.9
mj cal 5	5	✓	75.0	74.6	99.4
mj cal 6	6	✓	100.0	97.9	97.9
mj cal 7	7	✓	250.0	254.7	101.9

# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Last Cal. Update** 4/29/2021 9:26 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	x	1.0	1.8	179.1
mj cal2	2	✓	3.0	2.8	94.1
mj cal 3	3	✓	5.0	5.3	106.4
mj cal 4	4	✓	10.0	10.1	101.3
mj cal 5	5	✓	25.0	24.7	98.6
mj cal 6	6	✓	50.0	49.5	99.0
mj cal 7	7	✓	100.0	100.6	100.6



# AM #27 Cannabinoids


BWylie

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal1.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 6:42:30 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	118549	∞	5.9 <b>Low</b>	∞	2214747	1.791 ng/ml <b>Low</b>
THC-COOH	1.476	51523	89630.9	33.1	81.5	746272	5.531 ng/ml
THC	3.317	31606	∞	25.7	1303986 460438. 0	1242743	1.219 ng/ml

# AM #27 Cannabinoids


BWylie

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal2.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 6:49:14 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	177345	∞	10.4	∞	2114097	2.824 ng/ml <b>Low</b>
THC-COOH	1.476	94701	521.3	32.0	27.2	728352	9.791 ng/ml
THC	3.317	91004	∞	25.7	4157842 3955399 .3	1268562	2.841 ng/ml


# AM #27 Cannabinoids

*BWylie*

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 6:55:57 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	329270	∞	9.6	∞	2093443	5.320 ng/ml
THC-COOH	1.476	192858	2147.5	35.7	121.4	733330	19.081 ng/ml
THC	3.317	150991	∞	25.1	∞	1226441	4.642 ng/ml

# AM #27 Cannabinoids


*BWylie*

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:02:39 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	615562	∞	10.4	∞	2061080	10.128 ng/ml
THC-COOH	1.476	495975	1091.9	35.9	393446.6	726099	48.430 ng/ml
THC	3.332	308170	∞	24.9	4079087.0874596.9	1206570	9.278 ng/ml

# AM #27 Cannabinoids


BWylie

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 5
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-E1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:09:21 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	1507312	∞	11.1	∞	2076830	24.655 ng/ml
THC-COOH	1.476	760398	1245184.3	36.0	782.7	719322	74.563 ng/ml
THC	3.317	859702	∞	23.8	∞	1278060	23.901 ng/ml

# AM #27 Cannabinoids


*BWylie*

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:16:03 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	3001292	∞	11.4	∞	2060597	49.508 ng/ml
THC-COOH	1.476	1006755	862.7	36.0	1738.5	723585	97.914 ng/ml
THC	3.317	1811336	∞	23.5	∞	1276122	50.070 ng/ml

# AM #27 Cannabinoids


*BWylie*

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:22:45 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	5866141	∞	11.9	∞	1983355	100.565 ng/ml
THC-COOH	1.476	2379178	8231.8	37.3	2326665	654467	254.690 ng/ml
THC	3.332	3658693	∞	23.8	∞	1260473	102.049 ng/ml

# AM #27 Cannabinoids


BWylie

**Batch results** D:\MassHunter\Data\2021\am 27-28\042921\QuantResults\cannq.batch.bin  
**Calibration Last Update** 4/29/2021 9:26:33 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/29/2021 7:29:29 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.451	331726	∞	9.9	∞	2160567	5.192 ng/ml
THC-COOH	1.476	153553	1605.7	34.6	532.1	793028	14.235 ng/ml
THC	3.317	160158	∞	25.8	∞	1412582	4.301 ng/ml