

## Section Three

### Blood Toxicology

#### 3.10 SPE Methods for Quantitative GC/MSD Confirmation

##### 3.10.2 Extraction and Quantitation of Methamphetamine and Amphetamine from Blood Employing the Bond Elut Certify™ Extraction Column

###### 3.10.2.1 BACKGROUND

The drug amphetamine dates back to 1887. It was used freely as a nasal decongestant, appetite suppressant, and to treat disorders such as narcolepsy in the early part of the 20th century until its potential for abuse was fully realized.<sup>4,5,6</sup> The use of amphetamine and methamphetamine to treat narcolepsy, attention deficit disorder and obesity continues in a more regulated environment. Amphetamine (figure 1) and Methamphetamine (figure 2)) are phenethylamines structurally related to norepinephrine and epinephrine, respectively.




figure 1




figure 2

The blood concentrations of methamphetamine and amphetamine should be considered in conjunction with all available information to determine the degree and nature of an individual's impairment.<sup>2,3</sup> Therapeutic levels for legitimate methamphetamine and amphetamine use are one to two orders of magnitude less than abuse and toxic levels.<sup>6</sup>

Consult provided references for additional information regarding the pharmacology of these compounds.

###### 3.10.2.2 PRINCIPLE

Methamphetamine and amphetamine are recovered through the application of the Varian Bond Elut Certify® solid phase extraction (SPE) cartridge. The Bond Elut Certify® SPE cartridge contains a sorbent which utilizes cation exchange and non-polar mechanisms to recover methamphetamine and amphetamine from blood. Following the addition of deuterated internal standard mixture, the blood proteins are precipitated with cold acetonitrile. Following centrifugation, the supernatant is decanted and made basic with a 100mM phosphate buffer (pH 6). The sample is loaded onto the SPE cartridge that has been conditioned with methanol and a 100mM phosphate buffer (pH 6). The methanol conditioning opens up the coiled hydrophobic portion of the sorbent so that it interacts with the polar, buffered blood

matrix. The addition of the buffer removes excess methanol and creates an environment similar to the matrix thus allowing for optimal interaction between the sorbent and the analytes of interest. The analyte is retained by ionic interaction of the cationic functional groups present on the drug and the anionic sulfonic acid exchanger on the sorbent.

The cartridge is subsequently washed with 100mM acetic acid followed by methanol, to selectively remove matrix components and interfering substances from the cartridge. The wash also disrupts the hydrophobic and adsorption interactions leaving behind the ionically bound material. Next, the sorbent is thoroughly dried to remove traces of aqueous and organic solvents which could adversely affect the analyte recovery. When the sorbent is dry, the analytes of interest are recovered from the cartridge with alkaline ethyl acetate. The alkaline environment serves to disrupt the ionic interactions of the analyte with the sorbent and the methanol disrupts the hydrophobic interactions. Following the elution from the SPE cartridge the evaporated extract is acylated for confirmation on the GC/MSD. The quantitation is accomplished through the use of a deuterated internal standard and a five-point calibration curve. This method is based on the method utilized by the Bioaeronautical Sciences Research laboratory.<sup>1</sup>

### 3.10.2.3 EQUIPMENT AND SUPPLIES

- 3.10.2.3.1 Varian Bond Elute Certify<sup>®</sup> SPE Cartridge  
Product No: 1210-2051 (Laboratory Robot Compatible (LRC)) or 1211-3050 (Straight barrel) or equivalent  
Sorbent type: Mixed mode octyl (C8) and benzenesulfonic acid (SCX), Sorbent mass: 130mg, Particle size: 40 µm
- 3.10.2.3.2 Drybath or laboratory oven capable of 70°C
- 3.10.2.3.3 Evaporative concentrator equipped with nitrogen tank.
- 3.10.2.3.4 Vacuum manifold/pump
- 3.10.2.3.5 Tube rocker
- 3.10.2.3.6 Vortex mixer
- 3.10.2.3.7 Laboratory centrifuge capable of 3400- 3500rpm
- 3.10.2.3.8 Fixed and adjustable volume single channel air displacement pipettors, and appropriate tips, capable of accurate and precise dispensing of volumes indicated.
- 3.10.2.3.9 16 x 100mm round bottom glass tube
- 3.10.2.3.10 Screw Cap for 16mm O.D. tube
- 3.10.2.3.11 GC/MS Automated Liquid Sample (ALS) vials
- 3.10.2.3.12 GC/MS Vial Microinsert
- 3.10.2.3.13 GC equipped with a mass selective detector and a nonpolar capillary column with a phase composition comparable to 95%-dimethyl-polysiloxane with 5%-diphenyl.

**3.10.2.4 REAGENTS**

*Refer to manual section 5.12 for solution preparation instructions.*

- 3.10.2.4.1 Deionized/distilled (DI) water
- 3.10.2.4.2 Methanol (Certified ACS grade or better)
- 3.10.2.4.3 Hexane (Certified ACS grade or better)
- 3.10.2.4.4 Ethyl Acetate (Certified ACS grade or better)
- 3.10.2.4.5 Acetonitrile (Certified ACS grade or better)
- 3.10.2.4.6 Ammonium Hydroxide (Certified ACS grade or better)
- 3.10.2.4.7 Concentrated HCl (Certified ACS grade or better)
- 3.10.2.4.8 1% HCl in Methanol
- 3.10.2.4.9 100mM Phosphate Buffer (pH 6.0)
- 3.10.2.4.10 1M Acetic Acid
- 3.10.2.4.11 Pentafluoropropionic acid anhydride (PFAA)

**3.10.2.5 QUALITY ASSURANCE MATERIAL****3.10.2.5.1 Calibrator and Control Solutions**

Corresponding calibrator and control reference materials must be obtained from different vendors, or be from different lot numbers if suitable second vendors are not available.

**3.10.2.5.1.1 Reference Material Stock Solutions**

Concentration: 1 mg/mL  
(±)-Methamphetamine  
(±)-Amphetamine

Store remaining stock solution in freezer.

**3.10.2.5.1.2 Reference Material Working Solutions**

Working solutions are stable for 6 months when stored under refrigeration.

Concentration: 10ng/μL

Add 100μL each 1mg/mL Amphetamine and Methamphetamine Stock Solution to ≈9mL Methanol in a 10mL volumetric class A flask. QS to 10mL.

Concentration: 1ng/μL

Add 1mL 10ng/μL working drug solution to ≈5mL Methanol in a 10mL volumetric class A flask. QS to 10mL.

3.10.2.5.2 Internal Standard Stock Solutions3.10.2.5.2.1 **Stock Solutions**

Concentration: 1mg/mL  
(±)-Methamphetamine-D<sub>8</sub>  
(±)-Amphetamine-D<sub>8</sub>

Store remaining stock solution in ALS vial in freezer.

3.10.2.5.2.2 **Working Internal Standard Solution**

Working internal standard solution is stable for 6 months when stored under refrigeration.

Concentration: 10ng/μL

Add 100μL each 1mg/mL Amphetamine-D<sub>8</sub> and Methamphetamine-D<sub>8</sub> Stock Solution to 9mL Methanol in a 10mL volumetric class A flask. QS to 10mL.

3.10.2.5.3 Vendor Obtained Whole Blood Controls3.10.2.5.5.1 **Negative Whole Blood**3.10.2.5.5.2 **Positive Whole Blood**

Control containing Amphetamine and Methamphetamine each at a specified target concentration. Refer to package insert for verified value and expected range.

**3.10.2.6 PROCEDURE**3.10.2.6.1 Initial set-up

Label extraction tubes, Bond Elut Certify<sup>®</sup> extraction columns, and GC/MSD vials with microinserts for calibrators, controls and case samples.

3.10.2.6.2 Calibrator Preparation

Use the same lot of negative blood used to prepare the negative control to prepare calibrators.

3.10.2.6.2.1 Add 2mL of negative whole blood to five screw-top extraction tubes. Use the same lot number of blood as negative control.

- 3.10.2.6.2.2 Add the volume of working 1ng/μL Amphetamine and Methamphetamine mixed reference material as indicated in the chart below.

Level	Desired ng/mL	μL Working Reference material
1	25	50
2	50	100

- 3.10.2.6.2.3 Add the volume of working 10ng/μL Amphetamine and Methamphetamine mixed reference material as indicated in the chart below.

Level	Desired ng/mL	μL Working Reference material
3	100	20
4	250	50
5	500	100

3.10.2.6.3 Positive Control Sample Preparation

Use the same lot of negative blood used to prepare the negative control for positive control preparation.

- 3.10.2.6.3.1 Add 2mL of negative whole blood to two screw top tubes.

- 3.10.2.6.3.2 Add indicated amount of working 10ng/μL mixed control solution.

Desired ng/mL	μL Working Control
75	15
300	60

3.10.2.6.4 Negative Control Sample Preparation

Add 2mL of negative whole blood to screw top tube.

3.10.2.6.5 Case Sample Preparation

- 3.10.2.6.5.1 Based on enzyme immunoassay screen results, samples may be diluted with negative whole blood prior to analysis.

- 3.10.2.6.5.2 Place sample container on tube rocker for a minimum of five minutes. If sample is clotted, homogenize as necessary.
- 3.10.2.6.5.3 Transfer 2mL neat or diluted sample to labeled screw top tube.
- 3.10.2.6.6 Internal Standard Addition
- 3.10.2.6.6.1 To calibrators, controls and case samples, add 20 $\mu$ L of internal standard mix (100ng/mL).
- 3.10.2.6.6.2 Vortex tube briefly and allow to stand 15 to 30 minutes for sample equilibration..
- 3.10.2.6.7 Protein Precipitation
- 3.10.2.6.7.1 While vortexing, add 5mL cold acetonitrile to case, calibrator and control samples.
- 3.10.2.6.7.2 Cap tubes and rock samples for approximately 15 minutes. Tubes should be at room temperature. Remove from rocker and place samples into centrifuge and let stand for 5 minutes.
- 3.10.2.6.7.3 Centrifuge at 3400 – 3500 rpm for 10 minutes.
- 3.10.2.6.7.4 Transfer organic supernatant into second labeled tapered bottom centrifuge tube.
- 3.10.2.6.7.5 Transfer tube to Evaporative Concentrator. Evaporate sample to approximately 1mL under nitrogen at approximately 40°C. *Do not allow extract to go to dryness.*
- 3.10.2.6.7.6 To evaporated extract add 2mL 100mM phosphate buffer (pH 6). Vortex to mix.
- 3.10.2.6.7.7 If needed, centrifuge an additional 5 minutes to remove blood fragments or foam.
- 3.10.2.6.8 SPE Column Preparation
- 3.10.2.6.8.1 Insert labeled Bond Elut Certify<sup>®</sup> Extraction column in the vacuum manifold.

- 3.10.2.6.8.2 Add 2mL methanol to the column. Aspirate at  $\leq 3$  in. Hg to prevent sorbent drying.
- 3.10.2.6.8.3 Add 2mL 100mM Phosphate buffer (pH 6.00) to the column. Aspirate at  $\leq 3$  in. Hg.
- 3.10.2.6.9 Blood Extract Loading  
Load buffered blood onto column and allow to gravity flow or apply minimal vacuum.
- 3.10.2.6.10 Column Clean-up
- 3.10.2.6.10.1 Add 1mL 1M Acetic Acid.  
Aspirate at  $\leq 3$  in. Hg.
- 3.10.2.6.10.2 Increase vacuum to  $\geq 10$  in. Hg ( $\geq 34$  kPa) for  $\geq 5$  minutes (disc should be dry).
- 3.10.2.6.10.3 Add 6mL methanol.
- 3.10.2.6.11 Pre-Elution Dry Disc  
Increase vacuum to  $\geq 10$  in. Hg ( $\geq 34$  kPa) for  $\cong 5$  minutes.
- 3.10.2.6.12 Compound Elution
- 3.10.2.6.12.1 Open vacuum manifold, wipe collection tips, and insert the collection rack containing the labeled tapered tip centrifuge tubes.
- 3.10.2.6.12.2 Add 4mL 2%  $\text{NH}_4\text{OH}$  in ethyl acetate elution solvent to the column.  
*Collect with gravity flow or apply minimal vacuum.*
- 3.10.2.6.12.3 Add 50 $\mu\text{L}$  1% HCl in Methanol into each tube to minimize analyte loss.
- 3.10.2.6.13 Eluate Evaporation  
Transfer centrifuge tube to evaporative concentrator ~~TurboVap~~. Take solvent to dryness, under a gentle stream of nitrogen at approximately 40°C.
- 3.10.2.6.14 Derivatization
- 3.10.2.6.14.1 In fume hood add 50 $\mu\text{L}$  ethyl acetate. Vortex for 15 seconds.

- 3.10.2.6.14.2 Add 50 $\mu$ L PFAA.
- 3.10.2.6.14.3 Cap tubes and vortex briefly.
- 3.10.2.6.14.4 Place tubes in 70°C dry bath or oven for 20 minutes.
- 3.10.2.6.14.5 Remove from heat and allow to cool to room temperature.
- 3.10.2.6.14.6 Return tubes to TurboVap and evaporate to dryness under nitrogen at approximately 40°C.
- 3.10.2.6.14.7 Reconstitute extract with 50 $\mu$ L ethyl acetate.
- 3.10.2.6.14.8 Transfer reconstituted extract to labeled GC/MSD ALS vial with microinsert.
- 3.10.2.6.15 Preparation for GC-MS Run
- 3.10.2.6.15.1 Perform an AUTOTUNE and TUNE EVALUATION.
- 3.10.2.6.15.2 When tune values are acceptable, program SEQUENCE TABLE with sample, calibrator and control information.
- 3.10.2.6.15.3 Load ALS vials into quadrant racks as indicated in the SEQUENCE TABLE.
- 3.10.2.6.16 GC-MS Calibration Curve
- 3.10.2.6.16.1 The calibration curve should be established with a minimum of four data points.
- 3.10.2.6.16.2 Calibrators should be analyzed in order of increasing concentration.
- 3.10.2.6.16.3 The least squares line resulting from the analysis of calibrators must have a coefficient of correlation of  $\geq 0.98$ .
- 3.10.2.6.16.4 If calibration reference materials are run in duplicate, it is not required that duplicate



calibration points be included as long as the linearity requirement is met.

### 3.10.2.7 GC and MSD ACQUISITION PARAMETERS

Critical parameters are specified below. Parameters not specified are at the discretion of the analyst and should be optimized for the particular GC-MSD instrument. Each laboratory should maintain a centrally stored printed or electronic copy of current and past GC-MSD methods. The data supporting the GC-MSD method should be stored centrally.

#### 3.10.2.7.1 GC Temperature Parameter

Injection Port: 250°C

#### 3.10.2.7.2 MSD Instrument Parameters

Detector/Transfer Line: 280°C

#### 3.10.2.7.3 ALS Parameters

Injection Volume: 1µL (1 stop)

Viscosity Delay: A minimum of 3 seconds

Solvent Washes (A & B): A minimum of 4 pre- and post-wash rinses.

#### 3.10.2.7.4 MS SIM Parameters

Analyte	Target Ion	Qualifier Ion 1	Qualifier Ion 2
Amphetamine	190	118	91
Amphetamine-D8	193	126	96
Methamphetamine	204	160	118
Methamphetamine-D8	211	163	123

### 3.10.2.8 REPORTING CRITERIA

#### 3.10.2.8.1 Qualitative Chromatographic Criteria

Acceptable retention time window established by calibrators is  $\pm 0.2$  minutes.

#### 3.10.2.8.2 Qualitative Mass Spectral SIM Criteria

Ion ratios for the analyte and its corresponding internal standard, established by calibrators for target and qualifier ions, must not differ by more than  $\pm 20\%$  (relative).

### 3.10.2.8.3 Quantitative Mass Spectral and Control Criteria

- 3.10.2.8.3.1 Quantitative results can be accepted if the calculated concentrations of all calibrator and control samples are within  $\pm 20\%$  of their respective concentrations.
- 3.10.2.8.3.2 Quantitation is achieved through the plotting of the target ion response ratio versus the concentration for each calibrator.
- 3.10.2.8.3.3 Quantitative values for case samples, calibrators and controls will be truncated for reporting purposes.
- 3.10.2.8.3.4 Administrative limit of detection (LOD) for Amphetamine and Methamphetamine is 25ng/mL. Results < this LOD should be reported as negative unless there are extenuating circumstances. The Toxicology Discipline Leader must be consulted to evaluate exceptions.
- 3.10.2.8.3.5 If the concentration exceeds the calibration range, the sample needs to be appropriately diluted with negative whole blood for reanalysis.

## 3.10.2.9 **REPORTING OF RESULTS**

### 3.10.2.9.1 Quantitative Value

Analysis results should be truncated and reported out without decimal places.

### 3.10.2.9.2 Uncertainty Value

Based on the current uncertainty assessment, the +/- range should be included on the analysis report. Refer to method quality monitoring spreadsheet for current uncertainty figure.

## 3.10.2.10 **QUALITY ASSURANCE REQUIREMENTS**

### 3.10.2.10.1 General

- 3.10.2.10.1.1 Blood samples are to be stored under refrigeration after aliquots are removed for analysis.

3.10.2.10.1.2 Refer to toxicology manual section 5.1 for pipette calibration and intermediate check options.

3.10.2.10.1.3 Refer to toxicology manual section 5.2 for balance calibration and intermediate check requirements.

3.10.2.10.1.4 Refer to toxicology manual section 5.8 for additional GC-MSD quality assurance requirements.

3.10.2.10.1.5 Refer to toxicology manual section 5.10 for reference material authentication requirements.

3.10.2.10.2 Per Analysis Run Quality Requirements

3.10.2.10.2.1 A solvent blank must follow the highest calibrator, as well as proceed each case sample.

3.10.2.10.2.2 A minimum of the spiked blood controls described in section 3.10.2.6.3 must be run per batch of samples.

3.10.2.10.2.3 If the number of case samples exceeds 10, in addition to the two spiked described in 3.10.2.6.3, one blood spiked or commercially obtained control must be run for each additional 10 case samples. Additional concentrations may be used.

3.10.2.10.3 Monitoring of Control Values

Upon the completion of analysis, input blood control values on spreadsheet used to assess uncertainty for this method.

**3.10.2.11 ANALYSIS DOCUMENTATION**

3.10.2.11.1 A packet containing original data for controls and calibrators will be prepared for each analysis run and stored centrally in the laboratory where the analysis was performed until archiving.

3.10.2.11.2 A copy of controls and calibrators need not be included in individual case files. When necessary, a copy of the control

and calibrator printouts can be prepared from the centrally stored document.

Property of Idaho State Police Forensic Services  
Uncontrolled Internet Copy  
OBSOLETE DOCUMENT 01/07/2013

**3.10.2.12 REFERENCES AND RECOMMENDED READING**

- 3.10.2.12.1 Chaturevidi, A.K., Cardona, P.S., Soper, J.W. and Canfield, D.V., *Distribution and Optical Purity of Methamphetamine Found in Toxic Concentration in a Civil Aviation Accident Pilot Fatality*, U.S. Department of Transportation Federal Aviation Administration Technical Report, December 2004.
- 3.10.2.12.2 Logan, B.K., *Methamphetamine - Effects on Human Performance and Behavior*, Forensic Science Rev. 14(1/2): 133-151, 2002.
- 3.10.2.12.3 Logan, B.K., *Methamphetamine and Driving Impairment*. J Forensic Sci, 1996, 41(3):457-464.
- 3.10.2.12.4 Drummer, O.H., *Stimulants*, pp. 49-96. *in*: The Forensic Pharmacology of Drugs of Abuse, Arnold: London, 2001.
- 3.10.2.12.5 Moore, K.A., *Amphetamine/Sympathomimetic Amines*. pp. 245-264. *in*: Principles of Forensic Toxicology. Levine, B. ed., AACC, 2003.
- 3.10.2.12.6 Baselt, R.C., d-Methamphetamine, pp. 683-685. and Amphetamine, pp. 66-69. *in*: Disposition of Toxic Drugs and Chemicals in Man, Seventh ed., 2004.

## *Revision History*

---

---

### Section Three

#### Blood Toxicology

---

#### **3.10 SPE Methods for Quantitative GC/MSD Confirmation**

#### **3.10.2 Extraction and Quantitation of Methamphetamine and Amphetamine from Blood Employing the Bond Elut Certify™ Extraction Column**

---

<b>Revision No.</b>	<b>Issue Date</b>	<b>Revision/Comments</b>
0	11-21-2006	Original Issue
1	07-28-2008	Clarified that negative blood used to prepare calibrators and positive controls is the same lot as used for negative control.
2	03-07-2011	Storage condition specifications updated, reformatted reference material section for clarity, emphasized need for sample homogeneity, updated nomenclature.