

APPROVED

By John Garner at 10:18 am, Sep 10, 2020

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB Liquid Processor/Dilutor Serial Number: ML600HC11378

Volatiles Quality Assurance Controls

Run Date(s): 9/3/20 - 9/4/20; calibration 9/3/20

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jul-23	1907006	0.0764	0.0688-0.0840	0.0755 g/100cc
					0.0736 g/100cc
					g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.2003 g/100cc
					g/100cc
					g/100cc
Multi-Component mixture:					acceptable
Curve Fit:		Column 1	Lot #	Column 2	0.99989
		0.99998	FN07101701		

Ethanol Calibration Reference Material						
Calibrator level	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
50	0.050	0.045 - 0.055	0.0509	0.0529	0.002	0.0519
100	0.100	0.090 - 0.110	0.0998	0.0997	0.0001	0.0997
200	0.200	0.180 - 0.220	0.1999	0.1982	0.0017	0.199
300	0.300	0.270 - 0.330	0.2984	0.2967	0.0017	0.2975
400	0.400	0.360 - 0.440			0	#DIV/0!
500	0.500	0.450 - 0.550	0.5009	0.5025	0.0016	0.5017

Aqueous Controls			
Control level	Target Value	Acceptable Range	Overall Results
80	0.080	0.076 - 0.084	0.081 g/100cc


On 9/3/20, after the run sequence has started, I realized that the 0.080 QA samples were accidentally swapped with the internal standard blank locations on the auto-sampler tray. Up to this point, the mixed volatile, QC-1-1, and M2020-2310-2 already completed with M2020-3318-1 equilibrating in the oven awaiting injection. Upon this discovery, I decided to abort the sequence and I re-extracted internal standard blanks, 0.080 QA samples, QC1-1, M2020-2310-2, and M2020-3318-1. I then recreated the run sequence and double checked all the vial positions in the batch and restarted the run.

MBradley
9/4/20

NB

Worklist: 4493

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>
M2020-2310	2	BCK	Alcohol Analysis
M2020-3318	1	BCK	Alcohol Analysis
M2020-3344	1	BCK	Alcohol Analysis
M2020-3351	1	BCK	Alcohol Analysis
M2020-3362	1	BCK	Alcohol Analysis
M2020-3366	1	BCK	Alcohol Analysis
M2020-3376	1	BCK	Alcohol Analysis
M2020-3382	2	BCK	Alcohol Analysis
M2020-3383	1	BCK	Alcohol Analysis
M2020-3387	1	BCK	Alcohol Analysis
M2020-3399	1	BCK	Alcohol Analysis
M2020-3400	1	BCK	Alcohol Analysis
P2020-2543	1	BCK	Alcohol Analysis
P2020-2544	1	BCK	Alcohol Analysis
P2020-2574	1	BCK	Alcohol Analysis
P2020-2575	1	BCK	Alcohol Analysis
P2020-2576	1	BCK	Alcohol Analysis
P2020-2590	1	BCK	Alcohol Analysis
P2020-2591	1	BCK	Alcohol Analysis



NB

=====
Calibration Table
=====

General Calibration Setting

Calib. Data Modified : Thursday, September 03, 2020 2:25:10 PM
Signals calculated separately : No

Rel. Reference Window : 0.000 %
Abs. Reference Window : 0.100 min
Rel. Non-ref. Window : 0.000 %
Abs. Non-ref. Window : 0.100 min
Uncalibrated Peaks : not reported
Partial Calibration : Yes, identified peaks are recalibrated
Correct All Ret. Times: No, only for identified peaks

Curve Type : Linear
Origin : Ignored
Weight : Equal

Recalibration Settings:
Average Response : Average all calibrations
Average Retention Time: Floating Average New 75%

Calibration Report Options :
Printout of recalibrations within a sequence:
 Calibration Table after Recalibration
 Normal Report after Recalibration
If the sequence is done with bracketing:
 Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-propanol
2	1.00000	n-propanol

Signal Details

Signal 1: FID1 A, Front Signal
Signal 2: FID2 B, Back Signal

Overview Table

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No 1	methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
3.075	1	1	5.00000e-2	4.35949	1.14692e-2	No	No 1	ethanol
		2	1.00000e-1	8.76935	1.14034e-2			
		3	2.00000e-1	17.73212	1.12790e-2			
		4	3.00000e-1	26.66915	1.12490e-2			
		5	5.00000e-1	43.90200	1.13890e-2			
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No 2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No 1	isopropyl alcohol
4.285	2	1	5.00000e-2	4.55851	1.09685e-2	No	No 2	ethanol
		2	1.00000e-1	9.04151	1.10601e-2			
		3	2.00000e-1	18.39215	1.08742e-2			
		4	3.00000e-1	27.80130	1.07909e-2			
		5	5.00000e-1	46.19029	1.08248e-2			
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No 1	acetone
4.620	1	1	1.00000	40.85903	2.44744e-2	No	Yes 1	n-propanol
		2	1.00000	40.98506	2.43991e-2			
		3	1.00000	40.87630	2.44641e-2			
		4	1.00000	41.02362	2.43762e-2			
		5	1.00000	40.09966	2.49379e-2			
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No 2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No 2	isopropyl alcohol
7.550	2	1	1.00000	42.51776	2.35196e-2	No	Yes 2	n-propanol
		2	1.00000	42.22186	2.36844e-2			
		3	1.00000	41.86564	2.38859e-2			
		4	1.00000	41.85574	2.38916e-2			
		5	1.00000	40.71259	2.45624e-2			


Peak Sum Table

No Entries in table

1 Warnings or Errors :


Warning : Curve requires more calibration points., (methanol)

Calibration Curves




methanol at exp. RT: 2.586
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 9.04744e-2
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio


NB




Acetaldehyde at exp. RT: 2.809
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 1.00217e-1
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio



Acetaldehyde at exp. RT: 2.977
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 1.00217e-1
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio




ethanol at exp. RT: 3.075
 FID1 A, Front Signal
 Correlation: 0.99998
 Residual Std. Dev.: 0.00265
 Formula: $y = mx + b$
 m: 2.19589
 b: -5.17961e-3
 x: Amount Ratio
 y: Area Ratio




methanol at exp. RT: 3.388
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 1.00208e-1
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio


NB




isopropyl alcohol at exp. RT: 3.628
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 2.38149e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio




ethanol at exp. RT: 4.285
FID2 B, Back Signal
Correlation: 0.99989
Residual Std. Dev.: 0.00710
Formula: $y = mx + b$
m: 2.28520
b: -1.37099e-2
x: Amount Ratio
y: Area Ratio




acetone at exp. RT: 4.308
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.59069e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio




n-propanol at exp. RT: 4.620
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio



acetone at exp. RT: 4.661
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.62121e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio



isopropyl alcohol at exp. RT: 4.969
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 2.51810e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio




n-propanol at exp. RT: 7.550
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio

=====

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.050 FN05211804
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014 -CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.35949	0.0509	g/100cc
2.	Ethanol	Column 2:	4.55851	0.0529	g/100cc
3.	n-Propanol	Column 1:	40.85903	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.51776	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN02271802
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.76935	0.0998	g/100cc
2.	Ethanol	Column 2:	9.04151	0.0997	g/100cc
3.	n-Propanol	Column 1:	40.98506	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.22186	1.0000	g/100cc

NS

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 FN06231704
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.73212	0.1999	g/100cc
2.	Ethanol	Column 2:	18.39215	0.1982	g/100cc
3.	n-Propanol	Column 1:	40.87630	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.86564	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 FN07311804
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	26.66915	0.2984	g/100cc
2.	Ethanol	Column 2:	27.80130	0.2967	g/100cc
3.	n-Propanol	Column 1:	41.02362	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.85574	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 FN08241801
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	43.90200	0.5009	g/100cc
2.	Ethanol	Column 2:	46.19029	0.5025	g/100cc
3.	n-Propanol	Column 1:	40.09966	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.71259	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.64092	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.48873	1.0000	g/100cc

NB

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\09-03-20_CAL\09-03-20_CAL 2020-09-03 13-19-55\09-03-20_CAL.S
 Data directory path: C:\Chem32\1\Data\09-03-20_CAL\09-03-20_CAL 2020-09-03 13-19-55\
 Logbook: C:\Chem32\1\Data\09-03-20_CAL\09-03-20_CAL 2020-09-03 13-19-55\09-03-20_CAL.LOG
 Sequence start: 9/3/2020 1:34:35 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM

Method file name: C:\Chem32\1\Data\09-03-20_CAL\09-03-20_CAL 2020-09-03 13-19-55\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.050 FN05211804	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN02271802	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN06231704	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN07311804	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN08241801	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-1

Analysis Date(s): 03 Sep 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.0750	0.0760	0.0010	0.0755	0.0000	0.0755
(g/100cc)	0.0752	0.0758	0.0006	0.0755		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.075	0.071	0.079	0.004


Reported Result	
0.075	

Calibration and control data are stored centrally.



ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-A
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.37433	0.0750	g/100cc
2.	Ethanol	Column 2:	6.57792	0.0760	g/100cc
3.	n-Propanol	Column 1:	39.95022	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.10364	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.44259	0.0752	g/100cc
2.	Ethanol	Column 2:	6.58231	0.0758	g/100cc
3.	n-Propanol	Column 1:	40.28199	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.27292	1.0000	g/100cc

RB

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-2

Analysis Date(s): 04 Sep 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.0735	0.0743	0.0008	0.0739	0.0005	0.0736
(g/100cc)	0.0728	0.0740	0.0012	0.0734		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.073	0.069	0.077	0.004


Reported Result	
0.073	

Calibration and control data are stored centrally.



ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-A
 Laboratory : Meridian
 Injection Date : Sep 4, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.53967	0.0735	g/100cc
2.	Ethanol	Column 2:	6.71160	0.0743	g/100cc
3.	n-Propanol	Column 1:	41.86920	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.99441	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-B
 Laboratory : Meridian
 Injection Date : Sep 4, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.48233	0.0728	g/100cc
2.	Ethanol	Column 2:	6.71334	0.0740	g/100cc
3.	n-Propanol	Column 1:	41.92204	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.20908	1.0000	g/100cc

NB

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1

Analysis Date(s): 03 Sep 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.2003	0.1994	0.0009	0.1998	0.0009	0.2003
(g/100cc)	0.2009	0.2006	0.0003	0.2007		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.200	0.190	0.210	0.010


Reported Result	
0.200	

Calibration and control data are stored centrally.



ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-A
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.25188	0.2003	g/100cc
2.	Ethanol	Column 2:	19.02520	0.1994	g/100cc
3.	n-Propanol	Column 1:	41.99362	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.05408	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.23156	0.2009	g/100cc
2.	Ethanol	Column 2:	19.05581	0.2006	g/100cc
3.	n-Propanol	Column 1:	41.81246	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.85102	1.0000	g/100cc

NB

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN04171701

Analysis Date(s): 03 Sep 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.0813	0.0820	0.0007	0.0816	0.0002	0.0815
(g/100cc)	0.0808	0.0820	0.0012	0.0814		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

Reported Result	
0.081	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.00751	0.0813	g/100cc
2.	Ethanol	Column 2:	7.19862	0.0820	g/100cc
3.	n-Propanol	Column 1:	40.43824	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.47297	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.93351	0.0808	g/100cc
2.	Ethanol	Column 2:	7.19094	0.0820	g/100cc
3.	n-Propanol	Column 1:	40.25801	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.41967	1.0000	g/100cc

MB

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK 1
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	40.55322	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.90674	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN07101701
 Laboratory : Meridian
 Injection Date : Sep 3, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.60192	0.1036	g/100cc
2.	Ethanol	Column 2:	8.88040	0.1039	g/100cc
3.	n-Propanol	Column 1:	38.69347	1.0000	g/100cc
4.	n-Propanol	Column 2:	39.69321	1.0000	g/100cc

NB

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK
 Laboratory : Meridian
 Injection Date : Sep 4, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.64141	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.74518	1.0000	g/100cc

NB

Sample Summary

Sequence table: C:\Chem32\1\Data\09-03-20b_SAMPLES\09-03-20b_SAMPLES 2020-09-03 17-19-51
 \09-03-20b_SAMPLES.S
 Data directory path: C:\Chem32\1\Data\09-03-20b_SAMPLES\09-03-20b_SAMPLES 2020-09-03 17-19-51\
 Logbook: C:\Chem32\1\Data\09-03-20b_SAMPLES\09-03-20b_SAMPLES 2020-09-03 17-19-51\
 \09-03-20b_SAMPLES.LOG
 Sequence start: 9/3/2020 5:34:34 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM
 Method file name: C:\Chem32\1\Data\09-03-20b_SAMPLES\09-03-20b_SAMPLES 2020-09-03 17-19-51\
 \ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D		2
2	2	1	MIX VOL FN071017	-	1.0000	002F0201.D		10
3	3	1	QC1-1-A	-	1.0000	003F0301.D		4
4	4	1	QC1-1-B	-	1.0000	004F0401.D		4
5	5	1	0.08 FN04171701-	-	1.0000	005F0501.D		4
6	6	1	0.08 FN04171701-	-	1.0000	006F0601.D		4
7	7	1	M2020-2310-2-A	-	1.0000	007F0701.D		2
8	8	1	M2020-2310-2-B	-	1.0000	008F0801.D		2
9	9	1	M2020-3318-1-A	-	1.0000	009F0901.D		4
10	10	1	M2020-3318-1-B	-	1.0000	010F1001.D		4
11	11	1	M2020-3344-1-A	-	1.0000	011F1101.D		2
12	12	1	M2020-3344-1-B	-	1.0000	012F1201.D		2
13	13	1	M2020-3351-1-A	-	1.0000	013F1301.D		4
14	14	1	M2020-3351-1-B	-	1.0000	014F1401.D		4
15	15	1	M2020-3362-1-A	-	1.0000	015F1501.D		4
16	16	1	M2020-3362-1-B	-	1.0000	016F1601.D		4
17	17	1	M2020-3366-1-A	-	1.0000	017F1701.D		4
18	18	1	M2020-3366-1-B	-	1.0000	018F1801.D		4
19	19	1	M2020-3376-1-A	-	1.0000	019F1901.D		4
20	20	1	M2020-3376-1-B	-	1.0000	020F2001.D		4
21	21	1	M2020-3382-2-A	-	1.0000	021F2101.D		2
22	22	1	M2020-3382-2-B	-	1.0000	022F2201.D		2
23	23	1	M2020-3383-1-A	-	1.0000	023F2301.D		5
24	24	1	M2020-3383-1-B	-	1.0000	024F2401.D		5
25	25	1	QC2-1-A	-	1.0000	025F2501.D		4
26	26	1	QC2-1-B	-	1.0000	026F2601.D		4
27	27	1	M2020-3387-1-A	-	1.0000	027F2701.D		4
28	28	1	M2020-3387-1-B	-	1.0000	028F2801.D		4
29	29	1	M2020-3399-1-A	-	1.0000	029F2901.D		4
30	30	1	M2020-3399-1-B	-	1.0000	030F3001.D		4
31	31	1	M2020-3400-1-A	-	1.0000	031F3101.D		4
32	32	1	M2020-3400-1-B	-	1.0000	032F3201.D		4
33	33	1	P2020-2543-1-A	-	1.0000	033F3301.D		4
34	34	1	P2020-2543-1-B	-	1.0000	034F3401.D		4
35	35	1	P2020-2544-1-A	-	1.0000	035F3501.D		4
36	36	1	P2020-2544-1-B	-	1.0000	036F3601.D		4
37	37	1	P2020-2574-1-A	-	1.0000	037F3701.D		4
38	38	1	P2020-2574-1-B	-	1.0000	038F3801.D		4
39	39	1	P2020-2575-1-A	-	1.0000	039F3901.D		4
40	40	1	P2020-2575-1-B	-	1.0000	040F4001.D		4
41	41	1	P2020-2576-1-A	-	1.0000	041F4101.D		4
42	42	1	P2020-2576-1-B	-	1.0000	042F4201.D		4
43	43	1	P2020-2590-1-A	-	1.0000	043F4301.D		4

NB

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
44	44	1	P2020-2590-1-B	-	1.0000	044F4401.D		4
45	45	1	P2020-2591-1-A	-	1.0000	045F4501.D		4
46	46	1	P2020-2591-1-B	-	1.0000	046F4601.D		4
47	47	1	QC1-2-A	-	1.0000	047F4701.D		4
48	48	1	QC1-2-B	-	1.0000	048F4801.D		4
49	49	1	INTERNAL STD BLK	-	1.0000	049F4901.D		2

Method file name: C:\Chem32\1\Data\09-03-20b_SAMPLES\09-03-20b_SAMPLES 2020-09-03 17-19-51
 \SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
50	50	1	EMPTY	-	1.0000	050F5001.D		0

NB