















Worklist: 1709

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
• C2017-0651	1	83774	AM 27 Blood THC Quant by LC	
- C2017-0695	1	83778	AM 27 Blood THC Quant by LC	
• C2017-0696	1	83776	AM 27 Blood THC Quant by LC	
, C2017-0708	1	83779	AM 27 Blood THC Quant by LC	
• C2017-0710	1	83773	AM 27 Blood THC Quant by LC	
• M2017-0629	1	79158	AM 27 Blood THC Quant by LC	
• M2017-1456	3	83775	AM 27 Blood THC Quant by LC	
, M2017-1640	2	83777	AM 27 Blood THC Quant by LC	
• M2017-1645	1	83772	AM 27 Blood THC Quant by LC	
• M2017-1763	3	83770	AM 27 Blood THC Quant by LC	
, P2017-0433	3	79159	AM 27 Blood THC Quant by LC	
• P2017-0910	1	83771	AM 27 Blood THC Quant by LC	



Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 5/3/17

Analyst: B. Wylie

PRE-ANALYTIC

0499102

1-28-2018

Plate Lot# Custom -0490364 * Plate Exp. 9-21-2017 * External QC Lot 61317, exp 6-13-17

- ✓ 1. Ensure all solutions are within expiration date.
 - Mobile Phase A: *0.1% Formic Acid in LCMS Water*
 - Mobile Phase B: *0.1% Formic Acid in LCMS Acetonitrile*
 - *LCMS Methanol*
 - **Blank/Negative Blood: Lot 321632-1**
 - *0.1% Formic Acid in water*
 - *MTBE*
 - *Hexane*
- Column: *UCT Selectra DA 100 x 2.1 mm 3um*
2. Check levels of mobile phases and needle wash and refill as necessary. Ensure waste is not full.
3. Purge Pump and Load appropriate Acq. Method, allow system to equilibrate for approx. 30 min.
4. Create worklist. Data path name: _____

ANALYTIC

- ✓ 1. Remove standards plate, blood, and samples from cold storage. Allow to reach room temperature.
- ✓ 2. Add **1000 µL blood** to wells of analytical (standards) plate. Mix via aspirate and dispense. Place cover on Plate
 - Blank blood for locations containing standards/QCs and internal standards
 - Sample blood for locations containing only internal standards
- ✓ 3. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**. *Shaker ID 66759*
- ✓ 4. Pipette **500µL 0.1% formic acid** to all wells of standards plate.
- ✓ 5. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**.
- ✓ 6. Transfer **800µL of blood+acid mixture** to corresponding wells of SLE+ plate.
- ✓ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). Wait 5 min. (**Load blood samples at 85- 100 PSI- Selector to Right**)
- ✓ 8. Add **2.25mL MTBE** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 9. Apply positive pressure for approx. 15 seconds (**10-15 PSI- Selector to left -**).
- ✓ 10. Add **2.25mL Hexane** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 11. Apply positive pressure for approx. 15 seconds. (**10-15 PSI Selector to the left**)
- ✓ 12. Remove collection plate containing eluate.
- ✓ 13. Place collection plate on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID 66819*
- ✓ 14. Reconstitute in **100 µL MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

POST-ANALYTIC

- ✓ 1. Open quantitation software and create a new quantitation batch.
 - Batch name: 5-3-17 THC Quant
- ✓ 2. Make any necessary integration changes. Limit curves based on validated linear ranges (3-50ng/mL).
- ✓ 3. Were all appropriate standards used in the curve for each analyte? (Y) / N
 - Are r^2 values ≥ 0.98 for each analyte? (Y) / N
- ✓ 4. Did all QCs pass for each analyte? (Y) / N Were QCs entered into QC charting? (Y) / N
- ✓ 5. Central File Packet to include: ✓ LIMS Worklist: ✓ Method Checklist ✓ Calibration and Control Reports

COMMENTS

BW

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 09:45 **Data File** Neg Control.d
Sample Type Sample **Sample Name** Neg Control
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-a2 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 09:57 **Data File** QC 10.d
Sample Type QC **Sample Name** QC 10
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	17193	179160	0.0960	10.1259
THC-COOH	THC-COOH-d9	2.446	12258	53618	0.2286	10.6597
THC	THC-d3	6.313	7064	55465	0.1274	10.6306

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 10:09 **Data File** External control lot 61317.d
Sample Type QC **Sample Name** External control lot 61317
Dilution 1 **Acq Method** Quant THC 2017.m
Position p2b2 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation 10 ng

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	19172	167558	0.1144	12.0192
THC-COOH	THC-COOH-d9	2.426	9071	47510	0.1909	8.8394
THC	THC-d3	6.293	6331	52086	0.1215	10.1387


ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin

Last Calib Update 5/5/2017 4:22 PM

Analyst Name ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1- 3ng/mL	1	<input checked="" type="checkbox"/>	3	3.1	101.8
cal 2- 5ng/mL	2	<input checked="" type="checkbox"/>	5	5.0	100.5
cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	9.9	98.9
QC 10	3	<input checked="" type="checkbox"/>	10	10.1	101.3
External control lot 61317	3	<input checked="" type="checkbox"/>	10	12.0	120.2
cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	25.0	100.1
cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	50.0	100.1
cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	97.7	97.7
cal 7 - 250ng/mL	7	<input checked="" type="checkbox"/>	250	252.3	100.9

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin

Last Calib Update 5/5/2017 4:22 PM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1- 3ng/mL	1	<input checked="" type="checkbox"/>	3	3.0	101.3
cal 2- 5ng/mL	2	<input checked="" type="checkbox"/>	5	4.9	97.8
cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	10.0	99.6
QC 10	3	<input checked="" type="checkbox"/>	10	10.7	106.6
External control lot 61317	3	<input checked="" type="checkbox"/>	10	8.8	88.4
cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	25.1	100.4
cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	50.9	101.8
cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	99.1	99.1
cal 7 - 250ng/mL	7	<input type="checkbox"/>	250	268.7	107.5


ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin

Last Calib Update 5/5/2017 4:22 PM

Analyst Name ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1- 3ng/mL	1	<input checked="" type="checkbox"/>	3	2.9	96.8
cal 2- 5ng/mL	2	<input checked="" type="checkbox"/>	5	5.3	105.2
cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	9.8	97.7
QC 10	3	<input checked="" type="checkbox"/>	10	10.6	106.3
External control lot 61317	3	<input checked="" type="checkbox"/>	10	10.1	101.4
cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	24.7	98.8
cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	51.1	102.3
cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	99.2	99.2
cal 7 - 250ng/mL	7	<input type="checkbox"/>	250	256.6	102.6

BW

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 09:33 **Data File** blank internal std.d
Sample Type Sample **Sample Name** blank internal std
Dilution 1 **Acq Method** Quant THC 2017.m
Position Vial 2 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



ISP

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 08:10 **Data File** Cal 1.d
Sample Type Calibration **Sample Name** Cal 1- 3ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.296	5600	207060	0.0270	3.0551
THC-COOH	THC-COOH-d9	2.426	4470	63066	0.0709	3.0398
THC	THC-d3	6.273	2249	62426	0.0360	2.9028

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 08:22 **Data File** cal 2.d
Sample Type Calibration **Sample Name** cal 2- 5ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.296	9437	204126	0.0462	5.0238
THC-COOH	THC-COOH-d9	2.386	6893	63149	0.1092	4.8893
THC	THC-d3	6.253	3700	57908	0.0639	5.2605

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 08:34 **Data File** cal 3.d
Sample Type Calibration **Sample Name** cal 3 - 10ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	19674	210069	0.0937	9.8888
THC-COOH	THC-COOH-d9	2.446	13557	63343	0.2140	9.9556
THC	THC-d3	6.373	7365	62823	0.1172	9.7748

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 08:46 **Data File** cal 4.d
Sample Type Calibration **Sample Name** cal 4 - 25ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	49897	206883	0.2412	25.0248
THC-COOH	THC-COOH-d9	2.446	32925	62408	0.5276	25.1020
THC	THC-d3	6.373	17950	61153	0.2935	24.6914

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 08:57 **Data File** cal 5.d
Sample Type Calibration **Sample Name** cal 5 - 50ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	97574	201154	0.4851	50.0463
THC-COOH	THC-COOH-d9	2.446	63597	59894	1.0618	50.9101
THC	THC-d3	6.353	36041	59457	0.6062	51.1449

ISP

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 09:09 **Data File** cal 6.d
Sample Type Calibration **Sample Name** cal 6 - 100ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.296	187826	197795	0.9496	97.7043
THC-COOH	THC-COOH-d9	2.406	120388	58457	2.0594	99.1032
THC	THC-d3	6.273	65847	56069	1.1744	99.2255

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\5-3-17 THC Quant\QuantResults\5-3-17 THC Quant.batch.bin
Analysis Time 5/5/2017 4:22 PM **Analyst Name** ISP Tox
Report Time 5/5/2017 4:23 PM **Reporter Name** ISP Tox
Last Calib Update 5/5/2017 4:22 PM **Batch State** Processed

Analysis Info

Acq Time 2017-05-05 09:21 **Data File** cal 7.d
Sample Type Calibration **Sample Name** cal 7 - 250ng/mL
Dilution 1 **Acq Method** Quant THC 2017.m
Position P2-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	502950	204780	2.4561	252.2569
THC-COOH	THC-COOH-d9	2.426	304736	54710	5.5701	268.6943
THC	THC-d3	6.293	178411	58795	3.0345	256.6134

BW

Run Date: 5-3-2017 extracted
 Worklist ID: _____

	1	2	3	4	5	6
A	Std 1 + ISTD ✓ (3ng/mL)	neg ✗	C651 ✓			QC 1 + ISTD (10ng/ml)
B	Std 2 + ISTD ✓ (5ng/ml)	ext ctrl ✗	m1456-3 ✓			Std 7 + ISTD (250ng/ml)
C	Std 3 + ISTD ✓ (10ng/ml)	P433 ✓	Purple for C696 (B) ✓			Std 6 + ISTD (100ng/ml)
D	Std 4 + ISTD ✓ (25ng/ml)	m629 ✓	m1640-2 ✓			Std 5 + ISTD (50ng/ml)
E	Std 5 + ISTD ✓ (50ng/ml)	m1763-3 ✓	C695 ✓			Std 4 + ISTD (25ng/ml)
F	Std 6 + ISTD ✓ (100ng/ml)	p910 ✓	C708 ✓			Std 3 + ISTD (10ng/ml)
G	Std 7 + ISTD ✓ (250ng/ml)	m1645-1 (B) ✓	p910R			Std 2 + ISTD (5ng/ml)
H	QC 1 + ISTD ✓ (10ng/ml)	C710 ✓				Std 1 + ISTD (3ng/mL)

