

*Bylee*

**Worklist: 1978**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2017-2094	1	98738	AM 27 Blood THC Quant by LC
C2017-2142	1	98739	AM 27 Blood THC Quant by LC
C2017-2143	1	98740	AM 27 Blood THC Quant by LC
C2017-2146	1	98741	AM 27 Blood THC Quant by LC



# Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 10-24-17

Analyst: Ame Nord

Plate lot#: 0499102

Plate Expiration: 1/29/2018

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol Hexane

**Blank Blood Lot:** 321632-1

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 62340

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 2609543** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
**(Load at 85-100 PSI- Selector to the right)** Manifold ID: 66729
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.  
Worklist path: 102417 cann quant Batch Name: 102417 cann quant
- 2. Make any necessary integration changes,  $r^2$  values  $\geq 0.98$  for each analyte.
- 3. Did all QCs pass for each analyte? Y/N Enter QCs into control charting? - *See comments*
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: [Click here to enter text.](#)

THC-OH Not evaluated. No samples screened positive  
for THC-OH.



# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 12:15 **Data File** Negative Control.d  
**Sample Type** Sample **Sample Name** Negative Control  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-A2 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
<del>THC-OH</del>	<del>THC-OH-d3</del>	<del>1.895</del>	<del>122453</del>	<del>131870</del>	<del>0.9286</del>	<del>92.8280</del>
THC-COOH	THC-COOH-d9	2.305	5088	40717	0.1249	0.0000

*Not evaluated  
Noise*

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 12:27 **Data File** QC - 10ng.d  
**Sample Type** QC **Sample Name** QC - 10ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-H1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	0	162237	0.0000	0.0000
THC-COOH	THC-COOH-d9	2.245	20242	50310	0.4024	11.2093
THC	THC-d3	5.552	7961	51684	0.1540	10.9249

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 12:39 **Data File** External Control lot 21718.d  
**Sample Type** Sample **Sample Name** External Control lot 21718  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B2 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.915	0	124491	0.0000	0.0000
THC-COOH	THC-COOH-d9	2.245	10958	37432	0.2927	5.7985
THC	THC-d3	5.532	4687	42012	0.1116	7.2823

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 10:40 **Data File** Cal 1 - 3ng.d  
**Sample Type** Calibration **Sample Name** Cal 1 - 3ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-A1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.035	38174	121261	0.3148	29.8301
THC-COOH	THC-COOH-d9	2.245	8561	35972	0.2380	3.0965
THC	THC-d3	5.552	1990	36333	0.0548	2.4102

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 10:52 **Data File** Cal 2 - 5ng.d  
**Sample Type** Calibration **Sample Name** Cal 2 - 5ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.995	50558	118639	0.4262	41.2582
THC-COOH	THC-COOH-d9	2.245	9571	35491	0.2697	4.6605
THC	THC-d3	5.552	2867	37739	0.0760	4.2303

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 11:04 **Data File** Cal 3 - 10ng.d  
**Sample Type** Calibration **Sample Name** Cal 3 - 10ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-C1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.095	15051	116680	0.1290	10.7583
THC-COOH	THC-COOH-d9	2.245	13260	34731	0.3818	10.1936
THC	THC-d3	5.552	6343	35843	0.1770	12.8912



# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 11:16 **Data File** Cal 4 - 25ng.d  
**Sample Type** Calibration **Sample Name** Cal 4 - 25ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-D1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	35803	142656	0.2510	23.2781
THC-COOH	THC-COOH-d9	2.245	28128	41282	0.6814	24.9809
THC	THC-d3	5.572	13616	42358	0.3214	25.2844

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 11:27 **Data File** Cal 5 - 50ng.d  
**Sample Type** Calibration **Sample Name** Cal 5 - 50ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-E1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	60705	132347	0.4587	44.5972
THC-COOH	THC-COOH-d9	2.245	47369	38005	1.2464	52.8696
THC	THC-d3	5.552	25540	39360	0.6489	53.3687

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 11:39 **Data File** Cal 6 - 100ng.d  
**Sample Type** Calibration **Sample Name** Cal 6 - 100ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-F1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	143138	127654	1.1213	112.6066
THC-COOH	THC-COOH-d9	2.225	82431	37637	2.1902	99.4502
THC	THC-d3	5.552	49060	40911	1.1992	100.5665

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Analysis Time** 10/26/2017 8:52 AM **Analyst Name** anord  
**Report Time** 10/26/2017 9:11 AM **Reporter Name** ISP Tox  
**Last Calib Update** 10/26/2017 8:52 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-10-25 11:51 **Data File** Cal 7 - 250ng.d  
**Sample Type** Calibration **Sample Name** Cal 7 - 250ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-G1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram




### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	293751	122442	2.3991	243.7599
THC-COOH	THC-COOH-d9	2.245	179017	35349	5.0643	241.3105
THC	THC-d3	5.572	110542	38458	2.8744	244.2487

# ISP Forensics Calibration Curve Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Last Calib Update** 10/26/2017 8:52 AM **Analyst Name** ISP TOX


**Target Compound** *THC*  
**Internal Standard** *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.4	80.3
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.2	84.6
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	12.9	128.9
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.9	109.2
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.3	101.1
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	53.4	106.7
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	100.6	100.6
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	244.2	97.7

# ISP Forensics Calibration Curve Report


**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Last Calib Update** 10/26/2017 8:52 AM **Analyst Name** ISP TOX  
**Target Compound** *THC-COOH*  
**Internal Standard** *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.1	103.2
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.7	93.2
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.2	101.9
QC - 10ng	3	<input checked="" type="checkbox"/>	10	11.2	112.1
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.0	99.9
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	52.9	105.7
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	99.5	99.5
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	241.3	96.5

# ISP Forensics Calibration Curve Report

**Batch Data Path** C:\MyFiles\anord\Documents\blood data\1024 cann quant\QuantResults\102417 cann quant.batch.bin  
**Last Calib Update** 10/26/2017 8:52 AM **Analyst Name** ISP TOX  
**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input type="checkbox"/>	3	29.8	994.3
Cal 2 - 5ng	2	<input type="checkbox"/>	5	41.3	825.2
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.8	107.6
QC - 10ng	3	<input checked="" type="checkbox"/>	10	0.0	0.0
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	23.3	93.1
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	44.6	89.2
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	112.6	112.6
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	243.8	97.5

Not evaluate too much noise.