


Worklist: 2585

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
C2018-1243	1	121980	AM 27 Blood THC Quant by LC-QQQ
C2018-1299	1	121981	AM 27 Blood THC Quant by LC-QQQ
C2018-1305	1	121982	AM 27 Blood THC Quant by LC-QQQ
C2018-1329	1	121983	AM 27 Blood THC Quant by LC-QQQ
C2018-1339	2	121984	AM 27 Blood THC Quant by LC-QQQ



A handwritten signature in black ink, consisting of a stylized, cursive letter 'A' followed by a horizontal line extending to the right.

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 7-19-18

Analyst: Ann Nord

Plate lot#: 0515037

Plate Expiration: 9/28/2018

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
LCMS Methanol Hexane

Blank Blood Lot: 17J0718

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 62340

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: 07192018 Cann quant Batch Name: Cann quant
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:

Toxicology AM method 27 external prep information

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 10 ul each THC, THC-OH 100 ug/ml 100 ul C-THC in 9890 ul meOH

Ppd 6/5/18 Exp: 4/1/19 lot 6518 by AMN

Drug	lot (cerilliant)	expiration
C-THC	FE03121501	3/1/2020
THC-OH	FE01141502	1/1/2020
THC	FE04231406	4/1/2019

AM 27 control 50 ul working solution lot (6518) in 4950 ul blood lot (17J20718)

ppd 6/5/18 Exp 4/1/19 lot 6518 Concentration 10 ng/ml each by AMN

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 18:14 **Data File** QC - 10ng.d
Sample Type QC **Sample Name** QC - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	14831	138468	0.1071	9.9669
THC-COOH	THC-COOH-d9	2.406	9534	47638	0.2001	9.9725
THC	THC-d3	6.093	4586	38796	0.1182	10.2144

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 18:03 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.416	1270	156370	0.0081	0.2373
THC-COOH	THC-COOH-d9	2.065	1083	53403	0.0203	0.6075

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 18:26 **Data File** External Control lot 6518.d
Sample Type Sample **Sample Name** External Control lot 6518
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results


Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	12840	152975	0.0839	7.6891
THC-COOH	THC-COOH-d9	2.406	9985	52669	0.1896	9.4233
THC	THC-d3	6.093	4467	38363	0.1164	10.0651

ISP Forensics Calibration Curve Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update 7/20/2018 7:02 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.2	106.1
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.9	97.5
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.7	96.9
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.2	102.1
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.8	99.3
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	49.3	98.7
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	101.7	101.7
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	249.4	99.8


ISP Forensics Calibration Curve Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update 7/20/2018 7:02 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*

Internal Standard *THC-COOH-d9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input type="checkbox"/>	3	3.0	98.4
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.1	102.1
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.6	105.7
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.0	99.7
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	23.5	94.0
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	49.0	98.1
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	99.1	99.1
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	252.7	101.1

ISP Forensics Calibration Curve Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin

Last Calib Update 7/20/2018 7:02 AM **Analyst Name** ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.1	102.7
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.9	98.3
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.1	101.2
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.0	99.7
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.3	97.3
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	49.8	99.5
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	101.3	101.3
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	249.5	99.8

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 16:39 **Data File** Cal 1 - 3ng.d
Sample Type Calibration **Sample Name** Cal 1 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	5199	140361	0.0370	3.0797
THC-COOH	THC-COOH-d9	2.406	3038	46538	0.0653	2.9512
THC	THC-d3	6.073	1313	37691	0.0348	3.1844

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 16:51 **Data File** Cal 2 - 5ng.d
Sample Type Calibration **Sample Name** Cal 2 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	7729	138799	0.0557	4.9126
THC-COOH	THC-COOH-d9	2.406	5020	47078	0.1066	5.1044
THC	THC-d3	6.093	2197	40049	0.0549	4.8740

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 17:03 **Data File** Cal 3 - 10ng.d
Sample Type Calibration **Sample Name** Cal 3 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	15979	147081	0.1086	10.1174
THC-COOH	THC-COOH-d9	2.406	10308	48705	0.2116	10.5722
THC	THC-d3	6.093	4517	40332	0.1120	9.6925

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 17:15 **Data File** Cal 4 - 25ng.d
Sample Type Calibration **Sample Name** Cal 4 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	33427	131974	0.2533	24.3342
THC-COOH	THC-COOH-d9	2.426	20490	44570	0.4597	23.4893
THC	THC-d3	6.113	11304	38764	0.2916	24.8370

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 17:27 **Data File** Cal 5 - 50ng.d
Sample Type Calibration **Sample Name** Cal 5 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	68348	133520	0.5119	49.7534
THC-COOH	THC-COOH-d9	2.406	42578	44798	0.9504	49.0408
THC	THC-d3	6.073	21548	37013	0.5822	49.3406

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 17:39 **Data File** Cal 6 - 100ng.d
Sample Type Calibration **Sample Name** Cal 6 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	140525	135659	1.0359	101.2554
THC-COOH	THC-COOH-d9	2.406	84401	44166	1.9110	99.0578
THC	THC-d3	6.093	44896	37336	1.2025	101.6511

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\07192018 cann quant\QuantResults\cann quant.batch.bin
Analysis Time 7/20/2018 7:02 AM **Analyst Name** ISP Tox
Report Time 7/20/2018 7:04 AM **Reporter Name** ISP Tox
Last Calib Update 7/20/2018 7:02 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-19 16:16 **Data File** Cal 7 - 250ng.d
Sample Type Calibration **Sample Name** Cal 7 - 250ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.316	347738	136658	2.5446	249.5472
THC-COOH	THC-COOH-d9	2.406	209879	43164	4.8624	252.7356
THC	THC-d3	6.073	109548	37074	2.9548	249.4204

**Idaho State Police
Forensic Services
Toxicology Discipline**

Request for Departure from an Analytical Method

Date of Request

06/28/18

Forensic Scientist

Anne Nord

Analytical Method

Method 26

Type of Deviation

Minor

Major

Departure

Allow the response from 5ng/ml for Carboxy-THC to be evaluated as positive at the analyst's discretion.

4.3 Evaluation of Results

4.3.1 Minimum Criteria:

4.3.1.1

Analyte Administrative Threshold

THC 3 ng/mL

Carboxy-THC 10 ng/mL

THC-OH 3 ng/mL

Calculated sample concentration of 3 ng/mL or greater for THC and THC-OH, a calculated sample concentration of 10 ng/mL or greater for Carboxy-THC.

Allow for the positive control to be prepared in house.

Rationale for Departure

AM 26 is a screening method for cannabinoids, the function of the method is to indicate if cannabinoids are present. There is a potential for a positive sample to give a response in this range in the screen.

Discipline Leader Review

Departure approved- this deviation will remain in effect until the next method revision at which time the change will be made in the analytical method.

Departure Not Approved

Comments:

Celena Shrum

Celena Shrum

Toxicology Program Discipline Leader

Date: 06/28/18

