


REVIEWED
By Tamara Salazar at 2:54 pm, Feb 11, 2019

2/11/2019

Worklist: 2941

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2019-0109	1	139399	AM 27 Blood THC Quant by LC-QQQ
C2019-0182	1	139400	AM 27 Blood THC Quant by LC-QQQ
C2019-0209	1	139401	AM 27 Blood THC Quant by LC-QQQ
C2019-0210	1	139402	AM 27 Blood THC Quant by LC-QQQ
C2019-0218	1	139403	AM 27 Blood THC Quant by LC-QQQ



A handwritten signature in black ink, consisting of a stylized, cursive letter 'A' with a long horizontal stroke extending to the right.

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 2/7/19
Plate lot#: 0539904

Analyst: Anne Nord
Plate Expiration: 09/10/19

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane
Blank Blood Lot: 18G207D7
LCMS-QQQ ID: 62340
Column: UCT Selectra DA 100 x 2.1mm 3um
LCMS Methanol

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: 2019 Data\AM 27\020719 Batch Name: cann quant
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:



Toxicology AMI method 27 external prep information

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 10 ul each THC, THC-OH 100 ug/ml 100 ul C-THC in 9890 ul meOH

Ppd 6/5/18 Exp: 4/1/19 lot 6518 by AMN

Drug	lot (certiliant)	expiration
C-THC	FE03121501	3/1/2020
THC-OH	FE01141502	1/1/2020
THC	FE04231406	4/1/2019

AM 27 control 50 ul working solution lot (6518) in 4950 ul blood lot (17J20718)

ppd 6/5/18 Exp 4/1/19	lot 6518	Concentration 10 ng/ml each	by AMN
ppd 9/20/18 Exp 4/1/19	lot 92018	Concentration 10 ng/ml each	by AMN
ppd 01/14/18 Exp 4/1/19	lot 011419	Concentration 10 ng/ml each	by AMN



ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 16:00 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.256	395	68828	0.0057	1.3601 $\left. \begin{array}{l} \leftarrow \\ \rightarrow \end{array} \right\}$

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 16:12 **Data File** QC - 10ng C-THC, 5 ng THC THC-OH.d
Sample Type QC **Sample Name** QC - 10ng C-THC, 5 ng THC THC-OH
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.035	2761	70473	0.0392	5.6513
THC-COOH	THC-COOH-d9	2.105	4622	30919	0.1495	9.9784
THC	THC-d3	5.032	1330	26954	0.0493	4.9038

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 16:24 **Data File** External control.d
Sample Type Sample **Sample Name** External control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.055	5836	77293	0.0755	10.3127
THC-COOH	THC-COOH-d9	2.085	4812	32679	0.1472	9.8182
THC	THC-d3	4.992	3420	30349	0.1127	11.4986


ISP Forensics Calibration Curve Report

Batch Data Path
Last Calib Update

D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
2/8/2019 10:24 AM

Analyst Name ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.1	103.5
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.3	105.7
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.1	90.7
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.6	98.6
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	50.6	101.2
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	100.3	100.3
QC - 10ng C-THC, 5 ng THC THC-OH	7	<input checked="" type="checkbox"/>	5	5.7	113.0


ISP Forensics Calibration Curve Report

Batch Data Path
Last Calib Update

D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
2/8/2019 10:24 AM

Analyst Name ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.1	104.7
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.0	100.2
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.4	94.0
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.3	101.3
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	49.8	99.5
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	100.4	100.4
QC - 10ng C-THC, 5 ng THC THC-OH	7	<input checked="" type="checkbox"/>	10	10.0	99.8


ISP Forensics Calibration Curve Report

Batch Data Path
Last Calib Update

D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
2/8/2019 10:24 AM

Analyst Name ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.2	107.9
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.7	95.0
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.4	93.9
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.3	101.2
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	51.7	103.3
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	98.7	98.7
QC - 10ng C-THC, 5 ng THC THC-OH	7	<input checked="" type="checkbox"/>	5	4.9	98.1

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 14:37 **Data File** Cal 1 - 3ng.d
Sample Type Calibration **Sample Name** Cal 1 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.035	1474	76173	0.0193	3.1060
THC-COOH	THC-COOH-d9	2.045	1790	33362	0.0537	3.1405
THC	THC-d3	4.952	910	27293	0.0333	3.2381

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 14:49 **Data File** Cal 2 - 5ng.d
Sample Type Calibration **Sample Name** Cal 2 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.015	2561	70526	0.0363	5.2836
THC-COOH	THC-COOH-d9	2.085	2502	31341	0.0798	5.0081
THC	THC-d3	4.952	1201	25107	0.0478	4.7491

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 15:01 **Data File** Cal 3 - 10ng.d
Sample Type Calibration **Sample Name** Cal 3 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.015	4261	64729	0.0658	9.0709
THC-COOH	THC-COOH-d9	2.085	4129	29221	0.1413	9.3958
THC	THC-d3	4.952	2281	24686	0.0924	9.3873

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 15:12 **Data File** Cal 4 - 25ng.d
Sample Type Calibration **Sample Name** Cal 4 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.015	12471	66614	0.1872	24.6468
THC-COOH	THC-COOH-d9	2.085	10235	28075	0.3646	25.3271
THC	THC-d3	4.972	6188	25215	0.2454	25.3087

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 15:24 **Data File** Cal 5 - 50ng.d
Sample Type Calibration **Sample Name** Cal 5 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.015	26826	68875	0.3895	50.6031
THC-COOH	THC-COOH-d9	2.085	21207	30000	0.7069	49.7591
THC	THC-d3	4.952	12592	25256	0.4986	51.6588

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2019 Data\AM 27\020719\QuantResults\cann quant.batch.bin
Analysis Time 2/8/2019 10:24 AM **Analyst Name** ISP Tox
Report Time 2/8/2019 10:25 AM **Reporter Name** ISP Tox
Last Calib Update 2/8/2019 10:24 AM **Batch State** Processed

Analysis Info

Acq Time 2019-02-07 15:36 **Data File** Cal 6 - 100ng.d
Sample Type Calibration **Sample Name** Cal 6 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.015	53053	68307	0.7767	100.2895
THC-COOH	THC-COOH-d9	2.105	41259	29136	1.4161	100.3693
THC	THC-d3	4.952	24094	25357	0.9502	98.6582