


**Worklist: 3318**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
M2019-1432	2	149340	AM 27 Blood THC Quant by LC-QQQ
P2019-1046	1	149341	AM 27 Blood THC Quant by LC-QQQ
P2019-1055	1	149342	AM 27 Blood THC Quant by LC-QQQ
P2019-1096	1	149343	AM 27 Blood THC Quant by LC-QQQ
P2019-1167	1	149344	AM 27 Blood THC Quant by LC-QQQ
P2019-1170	1	149345	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

2

Extraction Date: 4/17/19  
Plate lot#: 0539904

Analyst: Sarah Pickle  
Plate Expiration: 09/10/19

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE  
**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol Hexane  
**Blank Blood Lot:** 445283-1  
**Column:** UCT Selectra DA 100 x 2.1mm 3um  
**LCMS-QQQ ID:** 59740

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: Data Path: \_\_\_\_\_

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: #3** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data. C:\MassHunter\Data\2019\Am 27\041719 THCQ SP  
Worklist path: G:\TOX\Pocatello\Abby\2019\AM 27\041719 THCQ SP Batch Name: THCQ
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve range limited: Carboxy-THC 10-100*



# Idaho State Police Forensic Services

R

## AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

### Methanol External Control Solution (Lot: WS041619)

10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH  
Approximate concentration 1ug/mL.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	184782	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2019
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	04/16/2019		
Prepared By:	Tamara Salazar		
Expires:	01/31/2020		

### Blood External Control Solution (Lot: 041619)

100 ul of methanol external control solution was added to 9900 ul of blood.  
Approximately 10ng/mL of each compound.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-1
Methanol External Control Solution	-	WS041619
Prepared:	04/16/19	
Prepared by:	Tamara Salazar	
Expires:	01/31/2020	

S


**ISP FORENSICS**  
**Pocatello**  
**Benzodiazepine Analysis Report**

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM      **Analyst Name** datastor  
**Report Time** 4/22/2019 3:14 PM      **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM      **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-04-18 11:30      **Data File** Negative.d  
**Sample Type** Sample      **Sample Name** Negative  
**Dilution** 1      **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-H5      **Sample Info**  
**Inj Vol** -1      **Comment**

**Sample Chromatogram**



P


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:14 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-04-18 11:54 **Data File** External Control.d  
**Sample Type** Sample **Sample Name** External Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-G5 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.212	54425	594340	0.0916	8.0137
THC-COOH	THC-COOH-D9	2.299	39533	210250	0.1880	7.4009
THC	THC-D3	5.732	15930	224879	0.0708	8.5237

PS


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:14 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

## Analysis Info

**Acq Time** 2019-04-18 11:06 **Data File** QC Control.d  
**Sample Type** QC **Sample Name** QC Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-A6 **Sample Info**  
**Inj Vol** -1 **Comment**

## Sample Chromatogram



## Results


Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.212	37151	677073	0.0549	4.8844
THC-COOH	THC-COOH-D9	2.299	56117	237017	0.2368	10.5664
THC	THC-D3	5.759	10274	257740	0.0399	4.6692

# ISP Forensics Calibration Curve Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin

**Last Calib Update** 4/23/2019 11:35 AM      **Analyst Name** ISP TOX

**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input checked="" type="checkbox"/>	3	3.2	105.1
Cal 2-5ng	2	<input checked="" type="checkbox"/>	5	5.1	101.6
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	9.8	97.9
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	23.2	92.8
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	50.8	101.6
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	101.0	101.0
QC Control	8	<input checked="" type="checkbox"/>	5	4.9	97.7


# ISP Forensics Calibration Curve Report

P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin

**Last Calib Update** 4/23/2019 11:35 AM      **Analyst Name** ISP TOX

**Target Compound** *THC-COOH*  
**Internal Standard** *THC-COOH-D9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input type="checkbox"/>	3	0.6	21.0
Cal 2-5ng	2	<input type="checkbox"/>	5	5.1	101.4
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	9.6	96.3
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	23.8	95.3
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	56.9	113.8
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	94.7	94.7
QC Control	8	<input type="checkbox"/>	10	10.6	105.7




# ISP Forensics Calibration Curve Report

D

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin

**Last Calib Update** 4/23/2019 11:35 AM      **Analyst Name** ISP TOX

**Target Compound** *THC*  
**Internal Standard** *THC-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input checked="" type="checkbox"/>	3	2.9	97.9
Cal 2-5ng	2	<input checked="" type="checkbox"/>	5	4.8	95.7
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	10.3	102.5
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	25.5	102.1
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	52.3	104.6
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	97.2	97.2
QC Control	8	<input checked="" type="checkbox"/>	5	4.7	93.4

P


**ISP FORENSICS**  
**Pocatello**  
**Benzodiazepine Analysis Report**

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM      **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM      **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM      **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-04-18 09:43      **Data File** Cal 1-3ng.d  
**Sample Type** Calibration      **Sample Name** Cal 1-3ng  
**Dilution** 1      **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-G6      **Sample Info**  
**Inj Vol** -1      **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.199	19460	562982	0.0346	3.1532
THC-COOH	THC-COOH-D9	2.285	16711	199447	0.0838	0.6302
THC	THC-D3	5.772	6310	243291	0.0259	2.9361

P


**ISP FORENSICS**  
**Pocatello**  
**Benzodiazepine Analysis Report**

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-04-18 09:55 **Data File** Cal 2-5ng.d  
**Sample Type** Calibration **Sample Name** Cal 2-5ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-F6 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.199	35873	627552	0.0572	5.0800
THC-COOH	THC-COOH-D9	2.272	33241	218481	0.1521	5.0704
THC	THC-D3	5.692	9841	241230	0.0408	4.7850

R


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

## Analysis Info

**Acq Time** 2019-04-18 10:07 **Data File** Cal 3-10ng.d  
**Sample Type** Calibration **Sample Name** Cal 3-10ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-E6 **Sample Info**  
**Inj Vol** -1 **Comment**

## Sample Chromatogram



## Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.199	70936	630886	0.1124	9.7928
THC-COOH	THC-COOH-D9	2.285	48412	217739	0.2223	9.6294
THC	THC-D3	5.772	20413	240882	0.0847	10.2539

R


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

## Analysis Info

**Acq Time** 2019-04-18 10:19 **Data File** Cal 4-25ng.d  
**Sample Type** Calibration **Sample Name** Cal 4-25ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-D6 **Sample Info**  
**Inj Vol** -1 **Comment**

## Sample Chromatogram



## Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.185	188340	698209	0.2697	23.2050
THC-COOH	THC-COOH-D9	2.285	105021	238292	0.4407	23.8140
THC	THC-D3	5.759	54525	262848	0.2074	25.5223

D


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

## Analysis Info

**Acq Time** 2019-04-18 10:31 **Data File** Cal 5-50ng.d  
**Sample Type** Calibration **Sample Name** Cal 5-50ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-C6 **Sample Info**  
**Inj Vol** -1 **Comment**

## Sample Chromatogram



## Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.199	415638	700757	0.5931	50.7766
THC-COOH	THC-COOH-D9	2.285	226915	238855	0.9500	56.8931
THC	THC-D3	5.745	116272	275122	0.4226	52.3000

D


# ISP FORENSICS Pocatello Benzodiazepine Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\041719 THCQ SP\QuantResults\THCQ.batch.bin  
**Analysis Time** 4/22/2019 3:11 PM **Analyst Name** datastor  
**Report Time** 4/22/2019 3:13 PM **Reporter Name** datastor  
**Last Calib Update** 4/22/2019 3:11 PM **Batch State** Processed

## Analysis Info

**Acq Time** 2019-04-18 10:43 **Data File** Cal 6-100ng.d  
**Sample Type** Calibration **Sample Name** Cal 6-100ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-B6 **Sample Info**  
**Inj Vol** -1 **Comment**

## Sample Chromatogram



## Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.212	32896	27829	1.1821	100.9924
THC-COOH	THC-COOH-D9	2.312	14212	9280	1.5315	94.6635
THC	THC-D3	5.652	6720	8577	0.7835	97.2028