

**REVIEWED**


By Celena Shrum at 9:22 am, Jun 13, 2019



6/12/2019

**Worklist: 3471**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
C2019-0985	1	154412	AM 27 Blood THC Quant by LC-QQQ
C2019-1003	2	154413	AM 27 Blood THC Quant by LC-QQQ
C2019-1057	1	154414	AM 27 Blood THC Quant by LC-QQQ
C2019-1079	1	154415	AM 27 Blood THC Quant by LC-QQQ





## AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 06/11/19  
Plate lot#: 0539904

Analyst: Anne Nord  
Plate Expiration: 09/10/19

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 19A207P3

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 69679

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.  
Worklist path: Data\2019\AM 27\061119 Batch Name: cann quant
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:



## Acquisition Method Info

**Method Name** AM 27 THC quant.m  
**Method Path** D:\MassHunter\methods\AM 27 THC quant.m  
**Method Description** UCT Column  
**Device List**  
 Binary Pump

**Name: Binary Pump** **Module: G7120A**

**Flow** 0.600 mL/min  
**Use Solvent Types** Yes  
**Stroke Mode** Synchronized  
**Low Pressure Limit** 0.00 bar  
**High Pressure Limit** 1000.00 bar  
**Max. Flow Ramp Up** 100.000 mL/min<sup>2</sup>  
**Max. Flow Ramp Down** 100.000 mL/min<sup>2</sup>  
**Expected Mixer** No check  
**Stroke A**  
   **Automatic Stroke Calculation A** Yes  
**Stop Time**  
   **Stoptime Mode** Time set  
   **Stoptime** 7.00 min  
**Post Time**  
   **Posttime Mode** Off

**Solvent Composition**

	Channel	Ch. 1 Solv.	Name 1	Ch2 Solv.	Name 2	Selected	Used	Percent
1	A	100.0 % Water V.03		100.0 % Water V.03	0.1 Formic Acid	Ch. 2	Yes	40.00 %
2	B	100.0 % Acetonitrile V.03		100.0 % Acetonitrile V.03	0.1 Formic Acid	Ch. 2	Yes	60.00 %

**Timetable**

	Time	A	B	Flow
1	3.50 min	40.00 %	60.00 %	--- mL/min
2	3.51 min	5.00 %	95.00 %	--- mL/min
3	4.50 min	5.00 %	95.00 %	--- mL/min
4	4.51 min	45.00 %	55.00 %	--- mL/min



Toxicology AM method 27 external prep information

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 75 ul THC-OH in 9767.5 ul meOH

Ppd 3/13/19 Exp: 9/13/19 lot 91319 by AMN

Drug	lot	expiration
C-THC	FE01061702 cerillient	3/1/2022
THC-OH	318.1b18.1L1a	12/1/2021
THC	135.1b71.0L6	11/1/2021

AM 27 control 100 ul working solution lot (91319) in 9900 ul blood lot (19A207p3)

ppd 3/13/19 Exp 9/13/19 lot 31319 Concentration 7.5 ng/ml THC, THC-OH and 15 ng/ml C-TH by AMN

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
6/12/2019 10:22:38 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-H1  
**Injection Volume** 10  
**Acq. Date-Time** 6/11/2019 4:14:13 PM  
**Sample Info.**


**Data File** QC.d  
**Sample** QC  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.209	140574	∞	8.3	∞	2199810	4.553 ng/ml
THC-COOH	1.220	102360	128.1	35.1	1338.0	1092557	8.987 ng/ml <b>Low</b>
THC	2.326	36970	351.5	24.5	219.1	1141139	4.933 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	neg control.d
<b>Type</b>	Sample	<b>Sample</b>	neg control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 4:21:58 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.209	27689	$\infty$			2414957	0.910 ng/ml <b>Low</b>

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
6/12/2019 10:22:38 AM

**Instrument** 69679  
**Type** Sample  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B2  
**Injection Volume** 10  
**Acq. Date-Time** 6/11/2019 4:29:42 PM  
**Sample Info.**

**Data File** external control.d  
**Sample** external control  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	244036	∞	8.8	∞	2361134	7.294 ng/ml
THC-COOH	1.220	143620	516.5	35.4	283.7	1166812	11.771 ng/ml
THC	2.341	54416	3095.7	23.4	274.6	1077785	7.359 ng/ml








# Compound Calibration Report




**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Last Cal. Update** 6/12/2019 10:22 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 6 Levels, 6 Levels Used, 6 Points, 6 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
cal 2	2	✓	3.0	3.0	101.2
cal 3	3	✓	5.0	4.5	89.0
cal 4	4	✓	10.0	10.9	109.2
cal 5	5	✓	25.0	25.4	101.5
cal 6	6	✓	50.0	49.8	99.6
cal 7	7	✓	100.0	99.4	99.4


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	cal 2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:20:13 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	90717	∞	12.2	∞	2155886	3.037 ng/ml
THC-COOH	1.220	35565	59.4	30.7	117.8	1050189	3.320 ng/ml <b>Low</b>
THC	2.326	21516	1413.4	26.6	51.4	1072500	3.278 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:27:57 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	148383	∞	9.5	∞	2375946	4.452 ng/ml
THC-COOH	1.220	57433	179.6	34.8	681.4	1132761	4.915 ng/ml <b>Low</b>
THC	2.326	39834	405.6	23.7	45.0	1195105	5.058 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:35:41 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	328760	∞	9.2	∞	2113672	10.919 ng/ml
THC-COOH	1.220	101055	782.9	36.6	152.9	1038623	9.328 ng/ml <b>Low</b>
THC	2.326	68612	635.0	25.0	168.2	1052623	9.329 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	cal 5
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-E1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:43:25 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	729365	∞	10.4	∞	2005664	25.376 ng/ml
THC-COOH	1.220	249890	429.7	36.7	953.5	987414	24.085 ng/ml
THC	2.326	174653	10603.4	26.4	374.4	1004124	23.913 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:51:09 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	1539927	∞	11.3	∞	2152003	49.823 ng/ml
THC-COOH	1.220	555483	1373.4	37.7	3778.2	1053281	50.069 ng/ml
THC	2.326	389565	2057.7	24.2	532.2	1074754	49.197 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2019\AM 27\061119\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 6/12/2019 10:22:38 AM

<b>Instrument</b>	69679	<b>Data File</b>	cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/11/2019 3:58:53 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.224	3072826	∞	11.5	∞	2150118	99.393 ng/ml
THC-COOH	1.220	1078033	2049.5	38.1	2250.6	1009341	101.284 ng/ml
THC	2.326	832355	∞	24.2	3019.1	1098251	102.224 ng/ml