


cg

**Worklist: 3600**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
C2019-1460	2	159730	AM 3 Urine Carboxy-THC
M2019-3195	2	159735	AM 3 Urine Carboxy-THC
M2019-3289	2	159733	AM 3 Urine Carboxy-THC
P2019-2262	1	159736	AM 3 Urine Carboxy-THC



A method revision occurred after this batch had been done. The revision that applies to this batch is revision 8.

cg

## AM 3: Carboxy-THC Urine Extraction

**Extraction Date:** 08/09/2019  
**Negative Urine Lot:** POC031319  
**Positive Control Working Solution Lot:** WS012319  
**1N KOH Lot:** 091817  
**Ethyl Acetate Lot:** BDH 121615D  
**Potassium Phosphate Buffer Lot:** 020118

**Analyst:** Celena Shrum  
**GC/MS ID:** Major Mass Spec  
**BioRad C3 Control Lot:** 68460  
**Hexane Lot:** BDH 121015A  
**BSTFA + 1% TMCS Lot:** Cerilliant FN08231301

### Pre-Analytic:

- 1. *Positive Control Working Solution Preparation Instructions:*  
 Add 180uL of 100ug/mL 11-nor-9-carboxy- $\Delta$ 9-THC Stock Solution to 9.82mL Methanol. Other volumes may be prepared. Solution is stable for 1-year or the expiration of the stock reference material (whichever is sooner). Store under refrigeration.
- 2. Verify Tune and Tune evaluation completed within the previous 7 days. Tune and Tune evaluation reports initialed and filed.
- 3. Create GCMS sequence to include controls, case blanks and case samples.

### Analytic:

- 1. Remove working solutions, external control, negative urine and case samples from cold storage.
- 2. To each labeled round bottomed tube add 3mL sample, using negative urine sample for both negative and positive control. Positive control: spike negative urine with 100uL positive control working solution.
- 3. Add 500uL 1N KOH to all tubes. Check pH. *(If pH <12, add additional 500ul 1N KOH).* *(Note: put a mark on the tube or separate the tubes that have a pH<12 as you will need to know this in step 5).*
- 4. Place tubes in 40C water bath for ~15 min. Remove and allow to cool.
- 5. If original pH was >12, add 1.5mL pH 1.8 Saturated Phosphate Buffer and 3mL Hexane/Ethyl Acetate (87:13)  
 If original pH was <12, add 3mL pH 1.8 Saturated Phosphate Buffer and 4mL Hexane/Ethyl Acetate (87:13)
- 6. Rock at ambient temp for ~ 10 minutes.
- 7. Centrifuge for ~ 10 min at ~3500rpm.
- 8. Transfer solvent to tapered bottom tube and evaporate to dryness under nitrogen @ 37C.
- 9. Add 50uL Ethyl Acetate and 50uL silylating reagent, cap and vortex. Heat @ 95C for 15min, then allow to cool.
- 10. Transfer sample to labeled ALS vial with insert.
- 11. Place ALS Vials in appropriate location on GCMS rack and run using appropriate GCMS method.

### Post-Analytic

- 1. Complete Data analysis on all samples and corresponding sample blanks  
 GCMS Data path: D:\DATA\CDS\2019\am 3 worklist 3600
- 2. Did positive and negative control samples provide intended response? Y / N
- 3. **Criteria for ID:** RT +/- 0.1 min., Ion Ratio of 347:473 & 371:488 within +/- 20%.  
 Sample response greater than Min Corrected Area, Diluted samples
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, and Control sample GCMS data printouts

Comments:

Data Path : G:\TOX\Pocatello\MMS\CDS\2019\am 3 worklist 3600\  
Data File : THC-PC1.D  
Acq On : 09 Aug 2019 12:21  
Operator : ISP\Datastor  
Sample : Spiked Positive c-THC Control  
Misc : c-THC lot # 0497429 in Negative Lot # POC031319; Worklist 3600  
ALS Vial : 32 Sample Multiplier: 1

Integration Parameters: events.e  
Integrator: ChemStation

Method : C:\gcms\1\methods\TOX.M  
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): THC-PC1.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.200	307	388	406	BV	9831	172666	100.00%	100.000%

Sum of corrected areas: 172666

Signal : EIC Ion 473.00 (472.70 to 473.70): THC-PC1.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.201	373	388	413	BB	3260	56781	100.00%	100.000%


Sum of corrected areas: 56781

Signal : EIC Ion 488.00 (487.70 to 488.70): THC-PC1.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.202	372	388	409	BB	1865	31938	100.00%	100.000%

Sum of corrected areas: 31938

TOX.M Mon Aug 12 19:28:55 2019



Data Path : G:\TOX\Pocatello\MMS\CDS\2019\am 3 worklist 3600\  
Data File : THC-PC2.D  
Acq On : 09 Aug 2019 15:05  
Operator : ISP\Datastor  
Sample : Spiked Positive c-THC Control  
Misc : c-THC lot # 0497429 in Negative Lot # POC031319; Worklist 3600  
ALS Vial : 32 Sample Multiplier: 1

Integration Parameters: events.e  
Integrator: ChemStation

Method : C:\gcms\1\methods\TOX.M  
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): THC-PC2.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.200	315	387	434	BB 2	12670	237703	100.00%	100.000%

Sum of corrected areas: 237703

Signal : EIC Ion 473.00 (472.70 to 473.70): THC-PC2.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.201	372	388	411	BB 2	4262	72242	100.00%	100.000%


Sum of corrected areas: 72242

Signal : EIC Ion 488.00 (487.70 to 488.70): THC-PC2.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.202	373	388	410	BB 2	2396	40590	100.00%	100.000%

Sum of corrected areas: 40590

TOX.M Mon Aug 12 19:30:15 2019



Data Path : G:\TOX\Pocatello\MMS\CDS\2019\am 3 worklist 3600\  
Data File : THC-NC.D  
Acq On : 09 Aug 2019 12:11  
Operator : ISP\Datastor  
Sample : Negative Control  
Misc : Lot # POC031319; AM #3 Worklist 3600  
ALS Vial : 31 Sample Multiplier: 1

Integration Parameters: events.e  
Integrator: ChemStation

Method : C:\gcms\1\methods\TOX.M  
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): THC-NC.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
--------	----------	------------	----------	-----------	-------	-------------	------------	--------------	------------

No peaks were detected using the method integration parameters!  
Signal : EIC Ion 473.00 (472.70 to 473.70): THC-NC.D\data.ms


peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
--------	----------	------------	----------	-----------	-------	-------------	------------	--------------	------------

No peaks were detected using the method integration parameters!  
Signal : EIC Ion 488.00 (487.70 to 488.70): THC-NC.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
--------	----------	------------	----------	-----------	-------	-------------	------------	--------------	------------

No peaks were detected using the method integration parameters!

TOX.M Mon Aug 12 19:32:18 2019



Data Path : G:\TOX\Pocatello\MMS\CDS\2019\am 3 worklist 3600\  
Data File : THC-C3.D  
Acq On : 09 Aug 2019 12:43  
Operator : ISP\Datastor  
Sample : Biorad C3 Control  
Misc : C3 Lot # 68460; Worklist 3600  
ALS Vial : 33 Sample Multiplier: 1

Integration Parameters: events.e  
Integrator: ChemStation

Method : C:\gcms\1\methods\TOX.M  
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): THC-C3.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.198	338	387	408	BV	2921	50735	100.00%	100.000%

Sum of corrected areas: 50735

Signal : EIC Ion 473.00 (472.70 to 473.70): THC-C3.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.200	374	387	406	BB	979	16627	100.00%	100.000%


Sum of corrected areas: 16627

Signal : EIC Ion 488.00 (487.70 to 488.70): THC-C3.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	6.200	373	388	404	BB	556	9512	100.00%	100.000%

Sum of corrected areas: 9512

TOX.M Mon Aug 12 19:31:51 2019



AM 3 carboxy-THC Calculations Sheet

**Run Date:** 8/9/2019

**Worklist #:** 3600

**Laboratory Case #:** Bioroad C3 Lot 68460

**Case Calculations:**

EIC Ion 371 Corr. Area: 50735

EIC Ion 473 Corr. Area: 16627

EIC Ion 488 Corr. Area: 9512

473:371 Ratio: **0.328**

488:371 Ratio: **0.187**

**Quality Control Calculations:**

**1st Control**

**2nd Control**

EIC Ion 371 Corr. Area: 172666

237703

EIC Ion 473 Corr. Area: 56781

72242

EIC Ion 488 Corr. Area: 31938

40590

Ratio of 473:371: 0.329

0.304

Average: 0.316

Ratio of 488:371: 0.185

0.171

Average: 0.178

**Acceptable Retention Time Range:** 6.100 to 6.300

**Approximate Minimum Corrected Area of 371 Ion:** 34533

**Acceptable 473:371 Ratio Range:** 0.253 to 0.380

**Acceptable 488:371 Ratio Range:** 0.142 to 0.213