


TS

4/29/2020

REVIEWED
By Sarah Pickle at 3:57 pm, May 05, 2020

Worklist: 4205

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2020-1137	51	BCK	AM 27 Blood THC Quant by LC-QQQ
M2020-1303	3	BCK	AM 27 Blood THC Quant by LC-QQQ
M2020-1343	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0612	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0903	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1045	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1180	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1181	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1225	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1226	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1227	1	BCK	AM 27 Blood THC Quant by LC-QQQ



TS

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 04/30/2020

Analyst: Tamara Salazar

Plate lot#: IDP-108-200303

Plate Expiration: 09/03/2020

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: Hemostat 445283-3

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 42** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water for blood samples** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 067103
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: *_D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wklst 4205 4206 TS*
Batch Name: *AM 27*
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curves limited: THC: 1-100 3-100 05/05/2020*
THC-OH not evaluated. TS



AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM


Instrument Falco
Type Sample
Acq. Method AM 27 THC quant.m
Sample Position P3-H5
Injection Volume 10
Acq. Date-Time 5/1/2020 3:23:05 AM
Sample Info.

Data File MJ_Negative Blood.d
Sample MJ_Negative Blood
Operator Tamara Salazar
Comment

Sample failed to inject. The sample was reconstituted and re injected on 05/01/2020. Please refer to the reinjection data. TS

Sample Chromatogram

+ TIC MRM (** -> **) MJ_Negative Blood.d (MJ_Negative Blood)





AM #27 Cannabinoid Quant. Results


Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_QC Control.d
Type	Sample	Sample	MJ_QC Control
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-A6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 3:07:53 AM		
Sample Info.			

Sample failed to inject. The sample was reconstituted and re injected on 05/01/2020. Please refer to the reinjection data. TS

Sample Chromatogram

+ TIC MRM (** -> **) MJ_QC Control.d (MJ_QC Control)






AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wkst 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Negative Blood_r.d
Type	Sample	Sample	MJ_Negative Blood_r
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-H5	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 7:17:08 PM		
Sample Info.			

Sample Chromatogram



TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_QC Control_r.d
Type	Sample	Sample	MJ_QC Control_r
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-A6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 7:01:47 PM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
* THC-OH	1.468	85679	∞	6.8 Low	23.18	446559	5.6170 ng/ml
THC-COOH	1.474	39684	∞	50.5	∞	108389	16.4044 ng/ml
THC	3.345	41978	133.01	27.8	36.37	1119824	4.7138 ng/ml

* Compound not evaluated - 5/5/2020 TS



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 5/5/2020 10:18 AM
Analyst Name ISP\Datastor
Analyte THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	×	1.0	1.3	132.4
MJ Cal 2	2	✓	3.0	3.0	101.0
MJ Cal 3	3	✓	5.0	5.0	100.9
MJ Cal 4	4	✓	10.0	10.1	101.1
MJ Cal 5	5	✓	25.0	24.3	97.0
MJ Cal 6	6	✓	50.0	49.4	98.9
MJ Cal 7	7	✓	100.0	101.1	101.1



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk\st 4205 4206 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 5/5/2020 10:18 AM
Analyst Name ISP\Datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	5.0	5.0	100.9
MJ Cal 2	2	✓	10.0	10.5	104.9
MJ Cal 3	3	✓	20.0	20.2	100.9
MJ Cal 4	4	✓	50.0	46.3	92.5
MJ Cal 5	5	✓	75.0	75.1	100.2
MJ Cal 6	6	✓	100.0	99.2	99.2
MJ Cal 7	7	✓	250.0	253.7	101.5



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk\st 4205 4206 TS\QuantResults\AM 27.batch.bin
Last Cal. Update 5/5/2020 10:18 AM
Analyst Name ISP\Datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	×	1.0	0.4	38.4
MJ Cal 2	2	×	3.0	2.9	95.2
MJ Cal 3	3	✓	5.0	4.0	80.8
MJ Cal 4	4	✓	10.0	11.8	118.4
MJ Cal 5	5	✓	25.0	26.2	104.7
MJ Cal 6	6	✓	50.0	48.1	96.2
MJ Cal 7	7	✓	100.0	99.8	99.8

Compound not evaluated. Multiple ratios out of range.
 5/5/2020 TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
 Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 1.d
Type	Cal	Sample	MJ_Cal 1
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-H6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:07:03 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.513 High	81799	∞	5.0 Low	∞	835613	0.3838 ng/ml Low
THC-COOH	1.474	24264	∞	45.6	116.73	205276	5.0429 ng/ml Low
THC	3.345	15168	80.10	27.5	6.24 Low	1842046	1.3239 ng/ml Low

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
 Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 2.d
Type	Cal	Sample	MJ_Cal 2
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-G6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:14:47 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.528 High	104610	∞	6.4 Low	∞	735244	2.8557 ng/ml Low
THC-COOH	1.489	42134	∞	48.7	41.15	177772	10.4876 ng/ml
THC	3.375	35968	102.72	30.6	10.13	1567471	3.0288 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wklst 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 3.d
Type	Cal	Sample	MJ_Cal 3
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-F6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:22:21 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.483	113315	∞	8.4	103.59	692839	4.0404 ng/ml
THC-COOH	1.504	73445	∞	53.5	146.57	163778	20.1767 ng/ml
THC	3.360	51658	131.07	28.7	16.35	1279838	5.0472 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wk1st 4205 4206 TS\QuantResults\AM 27.batch.bin
 Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 4.d
Type	Cal	Sample	MJ_Cal 4
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-E6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:29:55 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	197420	∞	8.5	155.06	650176	11.8418 ng/ml
THC-COOH	1.504	158446	∞	57.5	1528.33	155687	46.2649 ng/ml
THC	3.360	100739	96.53	27.0	27.79	1198783	10.1080 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wklst 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 5.d
Type	Cal	Sample	MJ_Cal 5
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-D6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:37:30 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	407735	∞	10.6	403.59	726613	26.1819 ng/ml
THC-COOH	1.489	287060	∞	57.4	∞	174217	75.1353 ng/ml
THC	3.345	303277	1572.21	26.9	197.73	1471241	24.2580 ng/ml

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wklst 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 6.d
Type	Cal	Sample	MJ_Cal 6
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-C6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:45:06 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	698832	∞	12.7 High	∞	731817	48.1112 ng/ml
THC-COOH	1.489	379666	∞	58.6	∞	174743	99.1943 ng/ml
THC	3.360	591472	1734.37	27.4	∞	1397231	49.4260 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\043020 AM 27 28 wklst 4205 4206 TS\QuantResults\AM 27.batch.bin
Calibration Last Update 5/5/2020 10:18:03 AM

Instrument	Falco	Data File	MJ_Cal 7.d
Type	Cal	Sample	MJ_Cal 7
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-B6	Comment	
Injection Volume	10		
Acq. Date-Time	5/1/2020 2:52:41 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	1327538	∞	13.1 High	∞	704810	99.8247 ng/ml
THC-COOH	1.504	901283	∞	58.8	∞	162563	253.6984 ng/ml
THC	3.360	1158257	3578.90	26.1	∞	1332097	101.1320 ng/ml