

REVIEWED

By Tamara Salazar at 7:38 pm, Apr 19, 2021

4/19/2021

CS

Worklist: 4911

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2021-0429	4	BCK	AM 27 Blood THC Quant by LC-QQQ



M2021-0777-1, M2021-0921-1, P2021-0673-1, P2021-0762-1, P2021-0779-1, and P2021-0782-1 were also included in this run.

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 04/14/2021

Analyst: Celena Shrum

Plate lot#: IDP-108-2-201206

Plate Expiration: 06/06/2021

Mobile phase A: 0.1% Formic Acid in LCMS Water

Mobile phase B: 0.1% Formic acid in Acetonitrile

Blank Blood Lot: Lampire 20L20724

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Using a calibrated pipette, add **1000µl blood and urine (if applicable) (calibrated pipette)** into the appropriate wells of analytical (standards) plate. **Pipette ID: #42**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample** of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). **(Load at 85-100 PSI- Selector to the right)**
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/mL and OH-THC 3ng/mL (quantitative), Carboxy-THC: 5ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Curve Limitations: THC 3-100, c-THC 5-250, THC-OH 3-100.

A possible manufacturing error with some of the analytical plates in this lot was discovered and as such this run was done to ensure that the error did not affect the plate(s) that had been used for Worklist 4863. The results of the samples from this plate were consistent with the previous run and the concentrations were within the expected uncertainty ranges. Since both runs were valid, the lower value of the runs will be reported. With further investigation, it was determined that the possible manufacturing error had no impact. Note: M2021-0429-4 was run with Worklist 4820 (not Worklist 4863).

**Idaho State Police
Forensic Services
Toxicology Discipline**

Request for Departure from an Analytical Method

Date of Request

02/24/2021

Forensic Scientist

Anne Nord

Analytical Method

Toxicology AM #27: Quantitative Analysis of THC and Metabolites in Blood and Urine by LCMS-QQQ

Request

The method currently reads:

4.3.2.5 If any points are dropped from the approved quantitative range of the curve, the compound will be reported qualitatively. For calibrators and controls 10 ng and below, the accuracy must be within 30%, for calibrators and controls greater than 10 ng/mL the accuracy must be within 20%. If a control falls outside the accuracy range, at the analyst's discretion, the compound may be reported qualitatively.

I would like to add in the following exception:

If the 1ng/ml point is dropped for THC, the quantitative range will be 3-50 ng/ml.

Discipline Leader Review

Departure approved

Comments: This deviation is approved and will remain in effect until it is changed in the actual method.

Departure Not Approved

Comments:

Celena Shrum

Toxicology Discipline Lead

Date: 02/24/2021

	1	2	3	4	5	6
a	cal 1ng	Blood NC				
b	cal 3 ng	M2021-0429-4				
c	cal 5 ng	M2021-0777-1				
d	cal 10ng	M2021-0921-1				
e	cal 25 ng	P2021-0673-1				
f	cal 50 ng	P2021-0762-1				
g	cal 100 ng	P2021-0779-1				
h	QC 1	P2021-0782-1				

Samples were moved to columns 5 and 6 on the SPE plate (A1 became A5, A2 became A6, etc.).

Cg


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_Negative Blood.d
Type	Sample	Sample	MJQ_Negative Blood
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-A6	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 1:16:56 PM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_QC Start.d
Type	Sample	Sample	MJQ_QC Start
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-H5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 1:32:09 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	171165	∞	8.7	∞	2202796	3.9472 ng/ml
THC-COOH	1.489	169372	∞	54.4	625.52	456151	14.4573 ng/ml
THC	3.330	100363	∞	34.2	∞	2660507	4.2734 ng/ml


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_QC End.d
Type	Sample	Sample	MJQ_QC End
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-H5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 3:34:57 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	178967	∞	8.8	∞	2268005	4.0380 ng/ml
THC-COOH	1.474	179126	1052.22	53.6	498.59	492041	14.1709 ng/ml
THC	3.315	107803	∞	35.4	175.34	2967237	4.1262 ng/ml



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Last Cal. Update 4/19/2021 8:38 AM
Analyst Name ISP\Datastor
Analyte THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	1.3	128.5
MJQ_Cal 2	2	✓	3.0	3.1	103.9
MJQ_Cal 3	3	✓	5.0	5.1	102.1
MJQ_Cal 4	4	✓	10.0	9.8	97.5
MJQ_Cal 5	5	✓	25.0	23.9	95.5
MJQ_Cal 6	6	✓	50.0	49.8	99.6
MJQ_Cal 7	7	✓	100.0	101.4	101.4

Calibrator 1 dropped due to poor secondary peak shape/response and ratio.



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Last Cal. Update 4/19/2021 8:38 AM
Analyst Name ISP\Datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	✓	5.0	5.3	106.1
MJQ_Cal 2	2	✓	10.0	9.7	97.5
MJQ_Cal 3	3	✓	20.0	19.7	98.7
MJQ_Cal 4	4	✓	50.0	47.7	95.5
MJQ_Cal 5	5	✓	75.0	75.6	100.8
MJQ_Cal 6	6	✓	100.0	101.2	101.2
MJQ_Cal 7	7	✓	250.0	250.6	100.3



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Last Cal. Update 4/19/2021 8:38 AM
Analyst Name ISP\Datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	0.9	86.3
MJQ_Cal 2	2	✓	3.0	2.9	97.0
MJQ_Cal 3	3	✓	5.0	5.1	102.3
MJQ_Cal 4	4	✓	10.0	10.2	101.5
MJQ_Cal 5	5	✓	25.0	24.8	99.0
MJQ_Cal 6	6	✓	50.0	50.1	100.1
MJQ_Cal 7	7	✓	100.0	100.0	100.0

Calibrator 1 dropped due to not meeting ratio requirement.


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_Cal 1.d
Type	Cal	Sample	MJQ_Cal 1
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-A5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 12:15:58 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.513 High	93737	∞	5.4 Low	24.42	2549172	0.8627 ng/ml Low
THC-COOH	1.474	73028	∞	48.7	118.03	524503	5.3026 ng/ml
THC	3.315	28565	∞	38.4 High	∞	3022496	1.2851 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument Instrument 1
Type Cal
Acq. Method AM 27 THCQ.m
Sample Position P1-B5
Injection Volume 10
Acq. Date-Time 4/14/2021 12:23:44 PM
Sample Info.

Data File MJQ_Cal 2.d
Sample MJQ_Cal 2
Operator Celena Shrum
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	154701	∞	8.0	∞	2419479	2.9100 ng/ml Low
THC-COOH	1.489	128130	∞	52.2	723.08	508732	9.7455 ng/ml
THC	3.330	75573	∞	33.4	31.68	2820655	3.1180 ng/ml


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_Cal 3.d
Type	Cal	Sample	MJQ_Cal 3
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-C5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 12:31:20 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	224290	∞	8.7	∞	2406033	5.1163 ng/ml
THC-COOH	1.489	254281	∞	55.2	4725.29	503227	19.7430 ng/ml
THC	3.345	127480	∞	35.9	∞	2797015	5.1035 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument Instrument 1
Type Cal
Acq. Method AM 27 THCQ.m
Sample Position P1-D5
Injection Volume 10
Acq. Date-Time 4/14/2021 12:38:56 PM
Sample Info.

Data File MJQ_Cal 4.d
Sample MJQ_Cal 4
Operator Celena Shrum
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	378470	∞	9.8	∞	2364899	10.1513 ng/ml
THC-COOH	1.489	600363	∞	57.6	∞	494144	47.7370 ng/ml
THC	3.345	244364	∞	30.1	∞	2728098	9.7536 ng/ml


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_Cal 5.d
Type	Cal	Sample	MJQ_Cal 5
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-E5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 12:46:32 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	811623	∞	10.7	278.79	2293503	24.7581 ng/ml
THC-COOH	1.489	890257	∞	56.7	∞	463167	75.6323 ng/ml
THC	3.330	605327	∞	27.3	121.21	2711517	23.8819 ng/ml




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument Instrument 1
Type Cal
Acq. Method AM 27 THCQ.m
Sample Position P1-F5
Injection Volume 10
Acq. Date-Time 4/14/2021 12:54:08 PM
Sample Info.

Data File MJQ_Cal 6.d
Sample MJQ_Cal 6
Operator Celena Shrum
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	1565013	∞	11.1	∞	2269098	50.0640 ng/ml
THC-COOH	1.489	1166439	∞	58.1	12117.5	453787	101.2077 ng/ml
THC	3.345	1284586	∞	25.7	∞	2742760	49.7890 ng/ml


AM #27 Cannabinoid Quant. Results



Batch results D:\MassHunter\Data\2021\AM 27-28\041421 AM 27 CS\QuantResults\AM 27 THCQ.batch.bin
Calibration Last Update 4/19/2021 8:38:16 AM

Instrument	Instrument 1	Data File	MJQ_Cal 7.d
Type	Cal	Sample	MJQ_Cal 7
Acq. Method	AM 27 THCQ.m	Operator	Celena Shrum
Sample Position	P1-G5	Comment	
Injection Volume	10		
Acq. Date-Time	4/14/2021 1:01:43 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	2870788	∞	11.5	∞	2122738	100.0003 ng/ml
THC-COOH	1.489	2643653	∞	58.5	∞	415774	250.6318 ng/ml
THC	3.330	2511517	∞	25.3	∞	2626504	101.3540 ng/ml