


**REVIEWED**  
By Sarah Collins at 7:22 am, Feb 16, 2023

2/14/2023

**Worklist: 6248**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>
C2023-0084	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0168	2	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0212	2	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0268	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0293	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0317	3	BCK	AM 27 Blood THC Quant by LC-QQQ
C2023-0352	1	BCK	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 2/14/23  
Plate lot#: 230113

Analyst: Anne Nord  
Plate re-test: 7/13/23

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 22B52016-1 **Urine Blank:** blood only run  
3um

**Column:** UCT Selectra DA 100 x 2.1mm

**LCMS-QQQ ID:** 69679

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette **1000µL blood (calibrated pipette) Pipette ID: I41142J** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators. SN > 10
- 4. Case sample response for THC 1ng/ml LOD 3ng/ml LOQ, OH-THC 3ng/mL LOD and LOQ, Carboxy-THC: 5 ng/mL (qualitative only). Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *autosampler failed prior to injecting any samples for this run. It was reset and samples were injected 2/15/23.*

*Curve limits THC-OH 3-100 dropped cal 1 due to poor response for qualifier ion.*

	1	2	3	4	5	6
a	cal 1	Internal blood control	0352-1 mixing plate			
b	cal 2	negative blood	0352-1 SLE and Injection plate			
c	cal 3	0084-1				
d	cal 4	0168-2				
e	Cal 5	0212-2				
f	cal 6	0268-1				
g	cal 7	0293-1				
h	Internal control (blood)	0317-3				

Plate position 3


c2023-\_\_\_\_-

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control blood.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 10:04:29 AM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	70074	∞	1005.2 <sub>5</sub>	∞	3620141	4.635 ng/ml
THC-COOH	1.358	97745	1847.4	246.03	∞	1266291	15.332 ng/ml
THC	2.956	231143	∞	24.39	470.6	1793212	4.904 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-B2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 10:11:05 AM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	end of run blood control blood A2.d
<b>Type</b>	Sample	<b>Sample</b>	end of run blood control blood A2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 11:49:59 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.334	92432	440.7	996.36	∞	4646959	4.737 ng/ml
THC-COOH	1.373	129906	854.0	263.01	∞	1680758	15.350 ng/ml
THC	2.956	373878	3654.0	24.01	11804. 2	2887125	4.925 ng/ml






# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Last Cal. Update** 2/15/2023 12:17 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	×	1.0	1.7	171.1
mj cal 2	2	✓	3.0	3.3	110.0
mj cal 3	3	✓	5.0	5.0	99.2
mj cal 4	4	✓	10.0	9.4	93.8
mj cal 5	5	✓	25.0	24.2	96.8
mj cal 6	6	✓	50.0	49.0	98.1
mj cal 7	7	✓	100.0	102.1	102.1

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 9:18:14 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	1.364	17583	575.9	857.30	∞	4277666	1.711 ng/ml	Low
THC-COOH	1.373	31527	∞	270.89	∞	1456206	5.377 ng/ml	
THC	2.986	51869	736.5	23.87	47.5	2378987	1.210 ng/ml	


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B1  
**Injection Volume** 10  
**Acq. Date-Time** 2/15/2023 9:24:58 AM  
**Sample Info.**

**Data File** mj cal 2.d  
**Sample** mj cal 2  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	54257	717.9	1035.5 0	∞	4377907	3.299 ng/ml
THC-COOH	1.373	68510	∞	261.21	702.3	1510496	9.626 ng/ml
THC	2.971	177180	3471.7	24.31	99.9	2518068	2.885 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 9:31:34 AM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	76576	149104.6	996.55	∞	3637777	4.959 ng/ml
THC-COOH	1.373	134017	859.6	258.18	2882.1	1346197	19.340 ng/ml
THC	2.956	254385	3461.4	25.24	141.8	2116680	4.604 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

**Instrument** 69679 **Data File** mj cal 4.d  
**Type** Cal **Sample** mj cal 4  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-D1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 2/15/2023 9:38:08 AM  
**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	167287	2023.8	932.66	∞	3792243	9.382 ng/ml
THC-COOH	1.373	359796	∞	249.43	∞	1336548	49.745 ng/ml
THC	2.956	572710	2049.5	24.27	1029.4	2250792	9.235 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

**Instrument** 69679 **Data File** mj cal 5.d  
**Type** Cal **Sample** mj cal 5  
**Acq. Method** AM 27 THC quant.m **Operator** Anne Nord  
**Sample Position** P3-E1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 2/15/2023 9:44:44 AM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	467480	∞	844.86	∞	3850022	24.207 ng/ml
THC-COOH	1.373	528635	∞	254.12	∞	1289620	74.965 ng/ml
THC	2.956	1558558	9615.2	23.89	3486.3	2272017	24.121 ng/ml


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-F1  
**Injection Volume** 10  
**Acq. Date-Time** 2/15/2023 9:51:18 AM  
**Sample Info.**

**Data File** mj cal 6.d  
**Sample** mj cal 6  
**Operator** Anne Nord  
**Comment**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	1002690	3717.5	832.05	∞	3996696	49.032 ng/ml
THC-COOH	1.373	723525	∞	250.70	14214.7	1323802	99.454 ng/ml
THC	2.956	3338557	3101017656 989630.0	24.04	∞	2327010	49.948 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2023\am 27-28\021423\QuantResults\cann.batch.bin  
**Calibration Last Update** 2/15/2023 12:17:38 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	2/15/2023 9:57:54 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.349	1908758	87863.1	806.67	∞	3616979	102.120 ng/ml
THC-COOH	1.373	1595889	25629.9	241.82	∞	1144117	251.494 ng/ml
THC	2.971	5791492	∞	24.30	12883.8	1967523	101.996 ng/ml