



9/12/2024

**Worklist: 6926****REVIEWED***By Britany Wylie at 2:16 pm, Sep 13, 2024*

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2024-1658	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-1666	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-1808	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-1809	1	BCK	AM 27 Blood THC Quant by LC-QQQ
M2024-3128	1	BCK	AM 27 Blood THC Quant by LC-QQQ





## AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 9/11/24

Plate lot#: 240513

Mobile phase A: 0.1% Formic Acid in LCMS Water

Blank Blood Lot: 24C52043

Column: UCT Selectra DA 100 x 2.1mm 3um

Analyst: Anne Nord

Plate Retest Date: 11/13/2024

Mobile phase B: 0.1% Formic acid in Acetonitrile

Blank Urine Lot: blood run only

LCMS-QQQ ID: 69679

### Pre-Analytic:

- ☒ 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- ☒ 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- ☒ 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- ☒ 2. Using a calibrated pipette, pipette 1000µL blood or 1000µL urine in wells of analytical (standards) plate. **Pipette ID: K52558G**
- ☐ 3. **Urine hydrolysis add 100 ul BG turbo, and 200 ul BG turbo buffer to the urine samples in wells of the analytical plate.**
- ☒ 4. Add **500µL of 0.1% formic acid in water** in the wells of the analytical plate.
- ☒ 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- ☒ 6. Transfer **700-800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.  
Amount transferred: **750 µL**
- ☒ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
**(Load at 85-100 PSI- Selector to the right)**
- ☒ 8. Wait 5 minutes.
- ☒ 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- ☒ 10. Wait 5 minutes.
- ☒ 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- ☒ 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- ☒ 13. Wait 5 minutes.
- ☒ 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- ☒ 15. Remove plate containing eluate. Place on SPE Dry 75401 and evaporate to dryness at approx. 35°C.
- ☒ 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- ☒ 1. Create batch and process data.
- ☒ 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- ☒ 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- ☒ 4. Did all QCs pass for each analyte? (if not, describe in comments section)
- ☒ 5. Enter QCs into control charting.
- ☒ 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: C2024-1666-1 did not inject, it was reconstituted along with the control and re-injected on 9/12/24.



	1	2	3	4	5	6
a	cal 1	blood end of run control				
b	cal 2	negative blood				
c	cal 3	M2024-3128-1				
d	cal 4	1658-1				
e	cal 5	1666-1				
f	cal 6	1808-1				
g	cal 7	1809-1				
h	Internal control (blood)					

Plate position 3

c2024-\_\_\_\_-\_\_

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-H1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 7:25:56 PM  
**Sample Info.**

**Data File** mj internal control blood.d  
**Sample** mj internal control blood  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	50059	234.2	831.85	∞	3395234	4.613 ng/ml
THC-COOH	2.077	78123	431.9	278.99	182576.5	1419018	13.694 ng/ml
THC	3.528	89997	1625.2	28.53	∞	762811	4.488 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Sample  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-B2  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 7:32:30 PM  
**Sample Info.**

**Data File** mj negative blood b2.d  
**Sample** mj negative blood b2  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal blood control end of run.d
<b>Type</b>	QC	<b>Sample</b>	mj internal blood control end of run
<b>Acq. Method</b>	thc quant 50 50.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	9/11/2024 8:45:02 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	43488	∞	883.88	∞	3040753	4.486 ng/ml
THC-COOH	2.077	78598	1684.0	287.95	1663.7	1259495	15.362 ng/ml
THC	3.528	60704	145.5	27.38	314.4	547771	4.240 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-A2  
**Injection Volume** 10  
**Acq. Date-Time** 9/12/2024 9:21:38 AM  
**Sample Info.**

**Data File** mj internal blood control end of runr.d  
**Sample** mj internal blood control end of run re-inject  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.036	22685	∞	791.24	∞	1456489	4.851 ng/ml
THC-COOH	2.107	43933	506.2	270.05	12508.3	773247	14.093 ng/ml
THC	3.543	45595	∞	29.55	∞	369691	4.673 ng/ml

this sample was reconstituted and injected on 9-12-24

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Last Cal. Update** 9/12/2024 9:52 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 3 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	117.4
mj cal 2	2	✓	3.0	2.8	94.7
mj cal 3	3	✓	5.0	4.8	95.2
mj cal 4	4	✓	10.0	9.3	93.1
mj cal 5	5	✓	25.0	24.4	97.6
mj cal 6	6	✓	50.0	50.5	101.0
mj cal 7	7	✓	100.0	101.0	101.0



# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Last Cal. Update** 9/12/2024 9:52 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 3 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.5	111.0
mj cal 2	2	✓	10.0	9.3	93.0
mj cal 3	3	✓	20.0	19.3	96.3
mj cal 4	4	✓	50.0	49.8	99.6
mj cal 5	5	✓	75.0	77.0	102.7
mj cal 6	6	✓	100.0	96.4	96.4
mj cal 7	7	✓	250.0	252.7	101.1

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Last Cal. Update** 9/12/2024 9:52 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 3 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	117.6
mj cal 2	2	✓	3.0	3.0	98.4
mj cal 3	3	✓	5.0	4.6	92.5
mj cal 4	4	✓	10.0	9.4	93.6
mj cal 5	5	✓	25.0	23.9	95.6
mj cal 6	6	✓	50.0	50.3	100.5
mj cal 7	7	✓	100.0	101.7	101.7

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-A1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 6:39:40 PM  
**Sample Info.**

**Data File** mj cal 1.d  
**Sample** mj cal 1  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	2.021	10945	∞	737.65	∞	3963394	1.176 ng/ml	Low
THC-COOH	2.092	28369	826.0	283.82	38539.3	1479380	5.549 ng/ml	
THC	3.528	20603	∞	24.20	∞	926654	1.174 ng/ml	

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-B1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 6:46:24 PM  
**Sample Info.**

**Data File** mj cal 2.d  
**Sample** mj cal 2  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	2.021	37330	∞	873.70	∞	4168692	2.952 ng/ml	Low
THC-COOH	2.077	56171	163.6	290.48	243588.3	1574192	9.296 ng/ml	
THC	3.528	87281	∞	23.70	∞	1239697	2.841 ng/ml	

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-C1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 6:53:00 PM  
**Sample Info.**

**Data File** mj cal 3.d  
**Sample** mj cal 3  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	51836	∞	862.70	∞	3505205	4.625 ng/ml
THC-COOH	2.077	119192	116119.8	274.30	238.5	1498012	19.258 ng/ml
THC	3.528	130546	397.5	23.77	∞	1036861	4.762 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-D1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 6:59:36 PM  
**Sample Info.**

**Data File** mj cal 4.d  
**Sample** mj cal 4  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	118998	796.2	858.14	5876.6	3802408	9.360 ng/ml
THC-COOH	2.077	306850	10396.4	276.45	38433.9	1432934	49.806 ng/ml
THC	3.528	272823	∞	23.62	1227.6	1060387	9.310 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-E1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 7:06:12 PM  
**Sample Info.**

**Data File** mj cal 5.d  
**Sample** mj cal 5  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	352962	1712.1	842.02	∞	4304111	23.903 ng/ml
THC-COOH	2.077	507809	717686.8	279.22	1447.1	1519885	77.039 ng/ml
THC	3.528	835065	∞	24.18	202.0	1204774	24.398 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-F1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 7:12:48 PM  
**Sample Info.**

**Data File** mj cal 6.d  
**Sample** mj cal 6  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	645339	∞	814.02	∞	3710310	50.268 ng/ml
THC-COOH	2.077	602947	1168230.6	276.07	967.0	1438300	96.356 ng/ml
THC	3.528	1510556	∞	25.10	∞	1044288	50.477 ng/ml



# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\091124\QuantResults\am 27.batch.bin  
**Calibration Last Update** 9/12/2024 9:52:40 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-G1  
**Injection Volume** 10  
**Acq. Date-Time** 9/11/2024 7:19:22 PM  
**Sample Info.**

**Data File** mj cal 7.d  
**Sample** mj cal 7  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.021	1291783	14034.1	809.61	∞	3656179	101.715 ng/ml
THC-COOH	2.077	1481654	2317.7	274.69	489416 0.7	1337301	252.697 ng/ml
THC	3.528	2358172	∞	27.28	∞	811186	101.038 ng/ml