












10/29/2024

**Worklist: 6962****REVIEWED***By Brittany Wylie at 1:53 pm, Oct 30, 2024*

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2024-1882	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-1975	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2002	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2021	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2029	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2024-2043	1	AVK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2062	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2069	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2024-2099	1	BCK	AM 27 Blood THC Quant by LC-QQQ	



## AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 10/29/24

Plate lot#: 240513

Mobile phase A: 0.1% Formic Acid in LCMS Water

Blank Blood Lot: 24C52043

Column: UCT Selectra DA 100 x 2.1mm 3um

Analyst: Anne Nord

Plate Retest Date: 11/13/2024

Mobile phase B: 0.1% Formic acid in Acetonitrile

Blank Urine Lot: 6524

LCMS-QQQ ID: 69679

### Pre-Analytic:

- ☒ 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- ☒ 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- ☒ 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- ☒ 2. Using a calibrated pipette, pipette 1000µL blood or 1000µL urine in wells of analytical (standards) plate. **Pipette ID: K52558G**
- ☒ 3. **Urine hydrolysis add 100 ul BG turbo, and 200 ul BG turbo buffer to the urine samples in wells of the analytical plate.**
- ☒ 4. Add **500µL of 0.1% formic acid in water** in the wells of the analytical plate.
- ☒ 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- ☒ 6. Transfer **700-800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.  
Amount transferred: **750 µL**
- ☒ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
**(Load at 85-100 PSI- Selector to the right)**
- ☒ 8. Wait 5 minutes.
- ☒ 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- ☒ 10. Wait 5 minutes.
- ☒ 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- ☒ 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- ☒ 13. Wait 5 minutes.
- ☒ 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- ☒ 15. Remove plate containing eluate. Place on SPE Dry 75401 and evaporate to dryness at approx. 35°C.
- ☒ 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- ☒ 1. Create batch and process data.
- ☒ 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- ☒ 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- ☒ 4. Did all QCs pass for each analyte? (if not, describe in comments section)
- ☒ 5. Enter QCs into control charting.
- ☒ 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *The end of run blood control did not inject it was re-constituted and injected 10/30/24*

Curve limitations: THC 3-100 cal 1 dropped due to poor response, qualifier ion indistinguishable from baseline.

	1	2	3	4	5	6
a	cal 1	Internal control urine	2069-1			
b	cal 2	negative blood	2099-1			
c	cal 3	1882-1	negative urine			
d	cal 4	1975-1	2029-1			
e	cal 5	2002-1				
f	cal 6	2021-1				
g	cal 7	2043-1				
h	Internal control (blood)	2062-1				

Plate position 3

c2024-\_\_\_\_-\_\_

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-H1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 6:10:12 PM  
**Sample Info.**

**Data File** mj internal control blood.d  
**Sample** mj internal control blood  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.991	26029	90.1	900.28	278.4	1705872	4.506 ng/ml
THC-COOH	2.062	46825	531.8	281.79	270.4	725142	15.149 ng/ml
THC	3.483	30363	219.6	23.87	204.7	233437	5.217 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Sample  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-B2  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 6:16:46 PM  
**Sample Info.**

**Data File** mj negative blood b2.d  
**Sample** mj negative blood b2  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative urine c3.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative urine c3
<b>Acq. Method</b>	thc quant 50 50.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-C3	<b>Comment</b>	Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	10/29/2024 8:15:35 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.825 <b>Low</b>	28100	73.1	837.77	∞	2142166	3.916 ng/ml

Negative for THC-OH peak is outside accepted retention time window.

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-A2  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 8:35:23 PM  
**Sample Info.**

**Data File** mj internal urine control a2.d  
**Sample** mj internal urine control a2  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.991	74880	∞	900.96	2069.8	2692805	7.966 ng/ml
THC-COOH	2.062	67851	236.7	261.68	84.1	1032319	15.378 ng/ml
THC	3.498	17586	135.7	25.39	16.9	159638	4.529 ng/ml


*[Handwritten signature]*

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal blood control end of run.d
<b>Type</b>	QC	<b>Sample</b>	mj internal blood control end of run
<b>Acq. Method</b>	thc quant 50 50.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-H1	<b>Comment</b>	Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	10/29/2024 8:41:59 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



sample did not inject, it was reconstituted and injected the next morning.



# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-H1  
**Injection Volume** 10  
**Acq. Date-Time** 10/30/2024 9:17:29 AM  
**Sample Info.**

**Data File** mj internal blood control end of runr.d  
**Sample** mj internal blood control end of run r  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	10838	59.4	741.01	∞	634134	5.011 ng/ml
THC-COOH	2.092	21677	335.4	276.29	90.2	349048	14.659 ng/ml
THC	3.543	20319	367.1	26.18	228.2	145569	5.546 ng/ml

re-constituted and injected 10/30/24

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Last Cal. Update** 10/30/2024 10:07 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 3 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	x	1.0	1.5	150.8
mj cal 2	2	✓	3.0	3.3	110.6
mj cal 3	3	✓	5.0	5.1	102.4
mj cal 4	4	✓	10.0	9.1	90.9
mj cal 5	5	✓	25.0	24.1	96.3
mj cal 6	6	✓	50.0	48.4	96.9
mj cal 7	7	✓	100.0	103.0	103.0

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Last Cal. Update** 10/30/2024 10:07 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 3 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.8	115.4
mj cal 2	2	✓	10.0	10.0	100.4
mj cal 3	3	✓	20.0	18.7	93.4
mj cal 4	4	✓	50.0	45.6	91.1
mj cal 5	5	✓	75.0	73.6	98.2
mj cal 6	6	✓	100.0	98.4	98.4
mj cal 7	7	✓	250.0	257.9	103.2

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Last Cal. Update** 10/30/2024 10:07 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 3 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	117.9
mj cal 2	2	✓	3.0	3.0	99.3
mj cal 3	3	✓	5.0	4.6	91.7
mj cal 4	4	✓	10.0	9.4	94.1
mj cal 5	5	✓	25.0	23.9	95.5
mj cal 6	6	✓	50.0	49.6	99.1
mj cal 7	7	✓	100.0	102.4	102.4

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-A1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:23:59 PM  
**Sample Info.**

**Data File** mj cal 1.d  
**Sample** mj cal 1  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	2.006	5394	71.5	647.41	28.1	1692364	1.179 ng/ml	Low
THC-COOH	2.062	12034	69.8	252.75	30.7	695065	5.769 ng/ml	
THC	3.498	5295	18.3	106.95 High	52.3	233318	1.508 ng/ml	

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-B1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:30:43 PM  
**Sample Info.**

**Data File** mj cal 2.d  
**Sample** mj cal 2  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	2.006	16769	158.3	741.92	∞	1726488	2.978 ng/ml	Low
THC-COOH	2.062	25036	44291.3	269.27	255.8	645114	10.035 ng/ml	
THC	3.513	16515	348.9	30.33	71.2	219991	3.317 ng/ml	

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-C1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:37:19 PM  
**Sample Info.**

**Data File** mj cal 3.d  
**Sample** mj cal 3  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	35079	227.1	877.74	∞	2255980	4.587 ng/ml
THC-COOH	2.062	70929	267.6	288.84	749.4	861304	18.677 ng/ml
THC	3.528	45782	478.5	26.41	21.2	359912	5.118 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-D1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:43:54 PM  
**Sample Info.**

**Data File** mj cal 4.d  
**Sample** mj cal 4  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	65807	2203.3	851.31	277.4	1990386	9.415 ng/ml
THC-COOH	2.062	179344	275.3	276.15	2172.7	823613	45.551 ng/ml
THC	3.528	65937	458.5	24.77	46.4	272289	9.090 ng/ml



# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-E1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:50:30 PM  
**Sample Info.**

**Data File** mj cal 5.d  
**Sample** mj cal 5  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	140199	361.3	868.76	∞	1640232	23.863 ng/ml
THC-COOH	2.062	255475	392914.1	271.94	86897.9	711138	73.633 ng/ml
THC	3.498	159361	1422.4	24.26	87.5	235660	24.085 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-F1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 5:57:04 PM  
**Sample Info.**

**Data File** mj cal 6.d  
**Sample** mj cal 6  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	309651	4250.5	836.87	∞	1733059	49.555 ng/ml
THC-COOH	2.062	354984	1416.3	265.53	14848.8	733069	98.441 ng/ml
THC	3.483	329848	1432.9	24.29	1062.7	238819	48.439 ng/ml

# AM #27 Cannabinoids


**Batch results** D:\MassHunter\Data\2024\am 27-28\102924\QuantResults\am 27.batch.bin  
**Calibration Last Update** 10/30/2024 10:07:08 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** thc quant 50 50.m  
**Sample Position** P3-G1  
**Injection Volume** 10  
**Acq. Date-Time** 10/29/2024 6:03:38 PM  
**Sample Info.**

**Data File** mj cal 7.d  
**Sample** mj cal 7  
**Operator** Anne Nord  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.006	549492	2933.9	829.78	∞	1483292	102.424 ng/ml
THC-COOH	2.062	838802	1663616.3	258.12	243325 2.3	651423	257.893 ng/ml
THC	3.483	631204	18410.7	23.70	17971. 7	213311	102.951 ng/ml