



**Worklist: 6727**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
M2024-0015	2	UCK	AM 25/AM 26 Urine MultiDrug/THC Screen by LC-QQQ	
M2024-0020	2	UCK	AM 25/AM 26 Urine MultiDrug/THC Screen by LC-QQQ	
M2024-0433	3	UCK	AM 25/AM 26 Urine MultiDrug/THC Screen by LC-QQQ	
M2024-0741	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2024-0742	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2024-0755	5	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2024-0929	6	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2024-0930	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2024-0931	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0113	2	UCK	AM 25/AM 26 Urine MultiDrug/THC Screen by LC-QQQ	
P2024-0343	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0457	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0480	1	UCK	AM 25/AM 26 Urine MultiDrug/THC Screen by LC-QQQ	
P2024-0495	1	BLOOD	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0637	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0677	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0678	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0679	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0680	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0699	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0709	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	

**Worklist: 6727**

TS

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
P2024-0730	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2024-0755	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	

## AM# 25: Multi-Drug Screen in Blood and Urine by LC-MS/MS

Extraction Date: 03/20/2024

Analyst: Tamara Salazar

Plate lot#: 231213

Plate Retest Date: 06/13/2024

Mobile phase A: 10mM Amm Form in LCMS Water

Mobile phase B: 0.1% Formic Acid in MeOH

Blank Blood Lot: Lampire 23A52595

Blank Urine Lot: POC021022

LCMS-QQQ ID: 069901

Column: Phenomenex Phenyl Hexyl (4.6x50mm, 2.6um)

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. **Urine Hydrolysis: In blank well, add 250µL urine, 40µL BG Turbo, and 100µL Instant Buffer I. Place on plate shaker for 5 minutes.**
- 3. Using a calibrated pipette, pipette 250µL blood or 250µL hydrolyzed urine in wells of analytical (standards) plate. Pipette ID: 42
- 4. Pipette **250µL 0.5 M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **200-450µL of blood+base and urine+base (if applicable)** mixture to corresponding wells of SLE+ plate. Amount transferred: 250uL
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). **(Load at 85-100 PSI- Selector to the right).**
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate.**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **900uL ethyl acetate.**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Add 50µL 1% HCl in MeOH to wells and place plate cover on plate before drying. This step is required for urine samples, but optional for blood samples.
- 17. Reconstitute in **100µL 20% LC MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Open quantitation software and create a new quantitation batch.
- 2. Make necessary changes to integration limits
- 3. Evaluate samples, S/N of primary transition >5 and S/N of secondary transition >3 or evaluation of peak symmetry and resolution. Within +/- 2% or 0.1 min RT of administrative control. Calculated concentration of 5 or greater or 2-5 for discretionary range.
- 4. Did all QCs pass for each analyte? If no, describe issue in comments (below).
- 5. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Benzoylcegonine note evaluated for M2024-0433-3 due to poor internal standard response.

TS

	1	2	3	4	5	6	7	8	9	10	11	12	
A	IS + Cal. 1	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	Neg Urine	P2024-0677-1	M2024-0755-5	
B	IS + Cal. 1	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	P2024-0755-2	P2024-0637-1	M2024-0742-1	
C	IS + Control 1	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	P2024-0730-1	P2024-0495-1	M2024-0741-1	
D	IS + Control 1	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	P2024-0480-1	P2024-0709-1	P2024-0457-1	Neg Blood
E	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	M2024-0433-3	P2024-0699-1	P2024-0343-2	IS + Control 1
F	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	P2024-0113-2	P2024-0680-1	M2024-0931-2	IS + Control 1
G	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	M2024-0015-2	P2024-0679-1	M2024-0930-2	IS + Cal. 1
H	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Sample	M2024-0020-2	P2024-0678-1	M2024-0929-6	IS + Cal. 1

All wells to contain 60 µl of residual DMSO

TS



# AM #25 Multi-Drug Screen. Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 25.batch.bin  
**Calibration Last Update** 3/27/2024 10:00:30 AM

**Instrument** Falco (069901)  
**Type** Sample  
**Acq. Method** AM 25 MDS.m  
**Sample Position** P2-D12  
**Injection Volume** 5  
**Acq. Date-Time** 3/20/2024 6:54:32 PM  
**Sample Info.**

**Data File** Negative Blood.d  
**Sample** Negative Blood  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

## Sample Chromatogram



TS




# AM #25 Multi-Drug Screen. Results

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 25.batch.bin  
**Calibration Last Update** 3/27/2024 10:00:30 AM

**Instrument** Falco (069901) **Data File** Cal.d  
**Type** Cal **Sample** Cal  
**Acq. Method** AM 25 MDS.m **Operator** Tamara Salazar  
**Sample Position** P2-H12 **Comment**  
**Injection Volume** 5  
**Acq. Date-Time** 3/20/2024 6:45:55 PM  
**Sample Info.**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

## Sample Chromatogram



Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
10-OH-Carbamazepine	3.795	2106212	129.74	247.23	20721299	10.0000
6-MAM	3.341	50672	21073.89	19682.43	2323025	10.0000
7-aminoclonazepam	3.592	1324344	186.45	482.35	7978088	10.0000
7-aminoflunitrazepam	3.792	2878534	5363.85	209.12	7978088	10.0000
9-Hydroxyrisperidone	4.062	9595031	724.06	91039.55	36423001	10.0000
Acetyl Fentanyl	4.129	417804	14.06	199209.7	35650437	10.0000
Acetyl Norfentanyl	3.028	259187	682.75	175.47	35650437	10.0000
a-hydroxyalprazolam	4.529	170077	20055.79	77.61	7978088	10.0000
alpha-hydroxymidazolam	4.605	2009152	780.31	308.36	7978088	10.0000
Alpha-PHP	4.027	3288416	5968.29	694.71	35650437	10.0000
alpha-PVP	3.782	4280484	15479.27	493.16	16940785	10.0000
Alprazolam	4.624	2155490	242.74	575.52	18530261	10.0000
Amitriptyline	4.566	1641852	145.59	299.53	10256289	10.0000
Amphetamine	3.047	3581332	672.84	344.37	16940785	10.0000
Benzoylecgonine	3.407	91018	109945.62	17.47	492717	10.0000
Bromazolam	4.696	938760	4141.18	19533.33	18530261	10.0000
Brompheniramine	4.160	100848	260.84	547.75	51590091	10.0000
Buprenorphine	4.998	104223	8596.75	9617.16	5593388	10.0000
Bupropion	4.012	5698643	14820.33	294.94	23783070	10.0000
Carbamazepine	4.261	10021748	∞	505.12	224431	10.0000
Carisoprodol	4.244	1215545	476.24	61.72	10432009	10.0000
Chlordiazepoxide	4.764	808905	346.48	115.30	18530261	10.0000
Chlorpheniramine	4.072	8186539	202.26	24.35	16205525	10.0000
Chlorpromazine	4.761	1976581	434851.30	2382.88	11749809	10.0000
Citalopram	4.159	3227377	534.26	3058254.73	51590091	10.0000
Clomipramine	4.762	2750083	15476.78	4994.52	51590091	10.0000
Clonazepam	4.453	904177	422.37	130.57	224431	10.0000
Clonazolam	4.373	1068784	2522.04	217616.96	18530261	10.0000
Clozapine	4.559	4022740	407.64	128289.59	22521529	10.0000

TS



# AM #25 Multi-Drug Screen. Results

Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
Cocaethylene	3.928	4219025	1046383.87	740.12	34561084	10.0000
Cocaine	3.760	5040726	29340.86	690.07	34561084	10.0000
Codeine	3.344	310425	378.30	171357.10	11913260	10.0000
Cyclobenzaprine	4.458	2780529	423.93	147.18	10256289	10.0000
Desipramine	4.459	4043880	453.49	132.14	10256289	10.0000
Dextromethorphan	4.196	2171742	13301.90	406453.81	16205525	10.0000
Dextrorphan	3.502	2482149	2925305.72	824.63	16205525	10.0000
Diazepam	4.872	1063930	1642.27	894.02	18530261	10.0000
Dihydrocodeine	3.128	880284	9719.55	439.03	11913260	10.0000
Diphenhydramine	4.151	10664824	713.32	2487.87	51590091	10.0000
DMT	3.168	551318	848.59	1734.34	16205525	10.0000
Doxepin	4.272	2506561	187.91	15.24	32774156	10.0000
Doxylamine	3.763	8617103	9544.11	6002063.18	16205525	10.0000
Duloxetine	4.410	106145	57005.82	16329.91	1591683	10.0000
EDDP	4.149	221225	368.70	42058.10	1541462	10.0000
Estazolam	4.549	3818737	1447.60	587692.21	18530261	10.0000
Etizolam	4.635	313682	195188.56	311856.35	18530261	10.0000
Fentanyl	4.328	413463	249.41	707.67	27436507	10.0000
Flualprazolam	4.498	814509	308690.01	71698.94	18530261	10.0000
Flunitrazepam	4.577	1984473	10189.66	∞	18530261	10.0000
Fluorofentanyl	4.373	692755	∞	650.81	27436507	10.0000
Fluoxetine	4.408	2877816	911.74	59.06	4729233	10.0000
Flurazepam	4.387	3628735	1671408.36	183892.02	18530261	10.0000
Hydrocodone	3.468	1780306	367.62	439.71	11913260	10.0000
Hydromorphone	3.054	1000913	336.26	389.13	225042	10.0000
Hydroxyzine	4.664	3519464	1020.49	1145807.46	22521529	10.0000
Imipramine	4.519	7242770	1014.28	568.04	10256289	10.0000
Ketamine	3.966	4273960	3321.40	103.81	15103514	10.0000
Lamotrigine	3.687	3446790	4714.52	11479.73	51590091	10.0000
Levamisole	3.445	3387890	524.64	168.28	34561084	10.0000
Levetiracetam	2.694	2996884	566.47	470.88	51590091	10.0000
Lorazepam	4.453	316088	95.19	56.05	18530261	10.0000
Maprotiline	4.474	445339	34.79	72.81	10256289	10.0000
MDA	3.137	2773706	625.18	44.44	32931834	10.0000
MDEA	3.351	5456029	4593.76	136.42	32931834	10.0000
MDMA	3.229	6514132	539.34	169.58	32931834	10.0000
Meperidine	3.780	3250031	1029.85	1470.00	16205525	10.0000
Meprobamate	3.705	967617	353.22	109.59	10432009	10.0000
Methadone	4.485	7609184	40587.78	1917.42	1541462	10.0000
Methamphetamine	3.154	6728134	887.04	1262.28	32931834	10.0000
Methocarbamol	3.595	294745	95.56	171368.35	1541462	10.0000
Methylphenidate	3.658	11023642	578.25	585.90	22092007	10.0000
Metoprolol	3.531	867135	198.49	526232.73	16205525	10.0000
Midazolam	4.775	774944	56739.79	147891.22	18530261	10.0000
Mirtazapine	4.352	3964824	1966.16	847544.82	16205525	10.0000
Mitragynine	4.386	655491	255590.69	513956.86	16205525	10.0000
Morphine	2.917	257995	135.15	592.99	225042	10.0000
Norbuprenorphine	3.924	108803	32542.80	49215.26	5593388	10.0000
Nordiazepam	4.720	1378126	11204.25	444.15	18530261	10.0000
Norfentanyl	3.443	7230364	11260.54	2427.97	35650437	10.0000

Cal

TS



# AM #25 Multi-Drug Screen. Results

Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
Norhydrocodone	3.085	213724	112.93	76.52	225042	10.0000
Norketamine	4.029	757706	728.95	2459.14	15103514	10.0000
Normeperidine	3.704	3131856	175.26	191.11	51590091	10.0000
Noroxycodone	3.006	1408759	∞	889.98	15103514	10.0000
Nortriptyline	4.506	1892852	630.37	308.34	10256289	10.0000
O-desmethyl-tramadol	3.056	6928011	5127.23	263.47	51590091	10.0000
O-desmethylvenlafaxine	3.362	1772978	453.53	50185.98	10306283	10.0000
Olanzapine	4.068	1263517	567648.69	725.38	224431	10.0000
Oxazepam	4.534	1384448	315.71	125.93	11239047	10.0000
Oxycodone	3.296	2874773	353.08	281.35	15103514	10.0000
Oxymorphone	2.823	1879643	∞	6761.29	225042	10.0000
Paroxetine	4.420	539379	207.35	220703.3	4729233	10.0000
Phenazepam	4.650	2172165	489.09	2530.55	18530261	10.0000
Phencyclidine	4.012	6827520	276.85	235.65	16205525	10.0000
Phentermine	3.277	1992751	266.32	40.24	22092007	10.0000
Phenytoin	4.152	375678	768.46	26.73	224431	10.0000
Primidone	3.505	15392140	5274.58	1149.41	224431	10.0000
Promethazine	4.503	6716324	45433.61	866.59	51590091	10.0000
Pseudoephedrine	2.862	59808303	5257.00	22890.71	32931834	10.0000
Quetiapine	4.694	5661988	389.34	1316400.69	38521642	10.0000
Risperidone	4.263	7998157	15524.50	341.21	36423001	10.0000
Sertraline	4.671	905847	435701.41	975.96	4729233	10.0000
Sufentanil	4.694	347314	175159.60	146.46	35650437	10.0000
Tapentadol	3.551	5666415	1089.35	3135.41	15103514	10.0000
Temazepam	4.687	3730004	485.89	173.95	18530261	10.0000
Topiramate	3.880	774451	322.34	374270.85	373996	10.0000
Tramadol	3.547	13409513	162.85	53.00	51590091	10.0000
Trazodone	4.863	7215458	9336322.03	1342797.42	32774156	10.0000
Venlafaxine	3.917	7939411	917.47	210.38	10306283	10.0000
Xylazine	3.551	2649779	611.25	28.26	15103514	10.0000
Zaleplon	4.348	2150507	705.21	750.48	38521642	10.0000
Zolpidem	4.486	10337772	3647280.16	1557236.77	38521642	10.0000
Zopiclone	4.433	731692	196352.82	266.10	3840864	10.0000



TS



# AM #25 Multi-Drug Screen. Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 25.batch.bin  
**Calibration Last Update** 3/27/2024 10:00:30 AM

**Instrument** Falco (069901)  
**Type** Sample  
**Acq. Method** AM 25 MDS.m  
**Sample Position** P2-A10  
**Injection Volume** 5  
**Acq. Date-Time** 3/20/2024 9:43:00 PM  
**Sample Info.**

**Data File** Negative Urine.d  
**Sample** Negative Urine  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

## Sample Chromatogram



TS




# AM #25 Multi-Drug Screen. Results

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 25.batch.bin  
**Calibration Last Update** 3/27/2024 10:00:30 AM

**Instrument** Falco (069901) **Data File** Internal Urine Control.d  
**Type** Sample **Sample** Internal Urine Control  
**Acq. Method** AM 25 MDS.m **Operator** Tamara Salazar  
**Sample Position** P2-F12 **Comment**  
**Injection Volume** 5  
**Acq. Date-Time** 3/20/2024 9:34:34 PM  
**Sample Info.**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

## Sample Chromatogram



Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
Alprazolam	4.624	24334550	5956.13	952.55	15474555	135.1888
Amphetamine	3.047	34489729	1176.29	549.77	12348388	132.1200
Codeine	3.344	3653178	1941540.37	6516.91	9321724	150.4001
Diphenhydramine	4.136	57041775	26507.71	2718.23	45567161	60.5555
Zolpidem	4.486	50817093	885436.73	4274.96	34775073	54.4527

# AM# 26: Screening of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 03/20/2024  
Plate lot#: 231212  
Mobile phase A: 10mM Amm Form in LCMS Water  
Blank Blood Lot: Lampire 23A52595  
LCMS-QQQ ID: 069901

Analyst: Tamara Salazar  
Plate Retest Date: 06/12/2024  
Mobile phase B: 0.1% Formic acid in MeOH  
Blank Urine Lot: POC021022  
Column: Phenomenex Phenyl Hexyl (4.6x50mm, 2.6um)

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250ul 1N KOH. Shake and incubate at 40 degrees for 15 minutes.
- 3. Using a calibrated pipette, pipette 1000µL blood or 1000µL hydrolyzed urine in wells of analytical (standards) plate. Pipette ID: 42
- 4. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 5. Add 500µL of 0.1% formic acid in water to blood samples, and 500µL of saturated phosphate buffer to urine samples in the wells of the analytical plate.
- 6. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 7. Transfer 700-800µL of blood+acid or urine+acid mixture to corresponding wells of SLE+ plate. Amount transferred: 750 µL
- 8. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). (Load at 85-100 PSI- Selector to the right)
- 9. Wait 5 minutes.
- 10. Add 2.25mL MTBE. (Add in 3 increments of 750uL)
- 11. Wait 5 minutes.
- 12. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 13. Add 2.25mL Hexane. (Add in 3 increments of 750uL)
- 14. Wait 5 minutes.
- 15. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 16. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 17. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, R<sup>2</sup> values ≥0.98 for each analyte
- 3. RT +/- 2% or 0.100 min, whichever is greater
- 4. Confirmation testing on case samples with a response for THC and OH-THC of 3ng/mL or greater and/or Carboxy-THC at 10ng/mL or greater (analyst discretion between 5-10ng/mL) may be pursued.
- 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports


COMMENTS: THC 3-100 – calibrator 1 dropped due to retention time shifts and poor chromatography.

TS

Analytical Plate Map

	1	2	3	4	5	6
A	IS + Cal. 1	IS + QC_1	P2024-0343-2	P2024-0699-1	M2024-0433-3	IS + QC_1
B	IS + Cal. 2	Neg Blood	P2024-0457-1	P2024-0709-1	P2024-0480-1	IS + Cal. 7
C	IS + Cal. 3	M2024-0741-1	P2024-0495-1	P2024-0730-1		IS + Cal. 6
D	IS + Cal. 4	M2024-0742-1	P2024-0637-1	P2024-0755-2		IS + Cal. 5
E	IS + Cal. 5	M2024-0755-5	P2024-0677-1	Neg Urine		IS + Cal. 4
F	IS + Cal. 6	M2024-0929-6	P2024-0678-1	M2024-0020-2		IS + Cal. 3
G	IS + Cal. 7	M2024-0930-2	P2024-0679-1	M2024-0015-2		IS + Cal. 2
H	IS + QC_1	M2024-0931-2	P2024-0680-1	P2024-0113-2	IS + QC_1	IS + Cal. 1

All wells to contain 100 µl of residual DMSO




SLE Plate Map

TS

	1	2	3	4	5	6
A	IS + Cal. 1	IS + QC_1	P2024-0343-2	P2024-0699-1	M2024-0433-3	
B	IS + Cal. 2	Neg Blood	P2024-0457-1	P2024-0709-1	P2024-0480-1	
C	IS + Cal. 3	M2024-0741-1	P2024-0495-1	P2024-0730-1		
D	IS + Cal. 4	M2024-0742-1	P2024-0637-1	P2024-0755-2		
E	IS + Cal. 5	M2024-0755-5	P2024-0677-1	Neg Urine		
F	IS + Cal. 6	M2024-0929-6	P2024-0678-1	M2024-0020-2		
G	IS + Cal. 7	M2024-0930-2	P2024-0679-1	M2024-0015-2		
H	IS + QC_1	M2024-0931-2	P2024-0680-1*	P2024-0113-2	P2024-0680-1*	

\*Sample moved during step 7 of the extraction due to clotting.



TS




# AM #26 Cannabinoids Screen Results

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJ Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	MJ Negative Blood
<b>Acq. Method</b>	AM 26 THC.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P1-B2	<b>Comment</b>	Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/20/2024 2:16:33 PM		
<b>Sample Info.</b>			

### Sample Chromatogram



TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** QC  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-H1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 2:01:21 PM  
**Sample Info.**

**Data File** MJ QC Control.d  
**Sample** MJ QC Control  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	9339	279293	6.30 ng/ml
THC-COOH	4.596	103301	682508	12.79 ng/ml
THC-OH	4.422	34948	3737336	4.83 ng/ml

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Sample  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-E4  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 4:55:43 PM  
**Sample Info.**

**Data File** MJ Negative Urine.d  
**Sample** MJ Negative Urine  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

## Sample Chromatogram





TS




# AM #26 Cannabinoids Screen Results

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901) **Data File** MJ QC Control\_Urine.d  
**Type** Sample **Sample** MJ QC Control\_Urine  
**Acq. Method** AM 26 THC.m **Operator** Tamara Salazar  
**Sample Position** P1-A2 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 4:40:32 PM  
**Sample Info.**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram




Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	2316	64262	6.73 ng/ml
THC-COOH	4.576	53519	474647	9.59 ng/ml
THC-OH	4.422	50469	5364270	4.86 ng/ml

TS



# AM #26 Cannabinoids Screen Calibration Curve Report

Batch results D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
 Last Cal. Update 3/21/2024 8:07 AM  
 Analyst Name ISP\datastor  
 Analyte THC Internal Standard THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	x	1.0	1.1	113.4
MJ Cal 2	2	✓	3.0	3.1	103.8
MJ Cal 3	3	✓	5.0	4.6	91.3
MJ Cal 4	4	✓	10.0	9.2	91.6
MJ Cal 5	5	✓	25.0	27.3	109.2
MJ Cal 6	6	✓	50.0	55.2	110.5
MJ Cal 7	7	✓	100.0	93.6	93.6

TS



# AM #26 Cannabinoids Screen Calibration Curve Report

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Last Cal. Update** 3/21/2024 8:07 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	5.0	5.9	118.6
MJ Cal 2	2	✓	10.0	8.6	85.6
MJ Cal 3	3	✓	20.0	19.1	95.7
MJ Cal 4	4	✓	50.0	49.9	99.8
MJ Cal 5	5	✓	75.0	77.5	103.4
MJ Cal 6	6	✓	100.0	95.5	95.5
MJ Cal 7	7	✓	250.0	253.4	101.4

TS



# AM #26 Cannabinoids Screen Calibration Curve Report

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Last Cal. Update** 3/21/2024 8:07 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	1.0	1.2	119.6
MJ Cal 2	2	✓	3.0	2.6	87.5
MJ Cal 3	3	✓	5.0	4.7	94.4
MJ Cal 4	4	✓	10.0	10.0	99.6
MJ Cal 5	5	✓	25.0	24.3	97.4
MJ Cal 6	6	✓	50.0	50.3	100.7
MJ Cal 7	7	✓	100.0	100.8	100.8

TS




# AM #26 Cannabinoids Screen Results

**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJ Cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 1
<b>Acq. Method</b>	AM 26 THC.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P1-A1	<b>Comment</b>	Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/20/2024 1:08:08 PM		
<b>Sample Info.</b>			

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.710 <b>High</b>	872	393297	1.13 ng/ml <b>Low</b>
THC-COOH	4.716 <b>High</b>	190833	2780159	5.93 ng/ml
THC-OH	4.543 <b>High</b>	23634	11001491	1.20 ng/ml <b>Low</b>

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-B1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:15:52 PM  
**Sample Info.**

**Data File** MJ Cal 2.d  
**Sample** MJ Cal 2  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	9107	641740	3.11 ng/ml
THC-COOH	4.616	177820	1772233	8.56 ng/ml
THC-OH	4.442	43374	8707897	2.63 ng/ml <b>Low</b>

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-C1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:23:29 PM  
**Sample Info.**

**Data File** MJ Cal 3.d  
**Sample** MJ Cal 3  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.570	18564	808929	4.56 ng/ml
THC-COOH	4.596	237075	1040204	19.15 ng/ml
THC-OH	4.422	43325	4748300	4.72 ng/ml

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-D1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:31:03 PM  
**Sample Info.**

**Data File** MJ Cal 4.d  
**Sample** MJ Cal 4  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.570	40182	791871	9.16 ng/ml
THC-COOH	4.596	515275	860570	49.92 ng/ml
THC-OH	4.422	85713	4394160	9.96 ng/ml



TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-E1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:38:38 PM  
**Sample Info.**

**Data File** MJ Cal 5.d  
**Sample** MJ Cal 5  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	74303	463129	27.31 ng/ml
THC-COOH	4.596	637102	683912	77.53 ng/ml
THC-OH	4.422	158281	3298572	24.34 ng/ml

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-F1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:46:13 PM  
**Sample Info.**

**Data File** MJ Cal 6.d  
**Sample** MJ Cal 6  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	155206	471394	55.23 ng/ml
THC-COOH	4.596	996733	868491	95.46 ng/ml
THC-OH	4.422	402499	4047485	50.33 ng/ml

TS



# AM #26 Cannabinoids Screen Results


**Batch results** D:\MassHunter\Data\2024\AM 25 26\031924 AM 25 26 TS\QuantResults\AM 26.batch.bin  
**Calibration Last Update** 3/21/2024 8:07:51 AM

**Instrument** Falco (069901)  
**Type** Cal  
**Acq. Method** AM 26 THC.m  
**Sample Position** P1-G1  
**Injection Volume** 10  
**Acq. Date-Time** 3/20/2024 1:53:47 PM  
**Sample Info.**

**Data File** MJ Cal 7.d  
**Sample** MJ Cal 7  
**Operator** Tamara Salazar  
**Comment**

Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods.

### Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	5.590	152764	272161	93.62 ng/ml
THC-COOH	4.596	2477502	811872	253.45 ng/ml
THC-OH	4.422	809183	4057486	100.82 ng/ml