

GA

10/1/2020

REVIEWED
By Brittany Wylie at 7:37 am, Oct 05, 2020

Worklist: 4552

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2020-1805	1	UCK	AM 2 Urine Toxi A
C2020-1823	2	UCK	AM 2 Urine Toxi A



AM 2: De-Tox Tube A Urine Extraction

Extraction Date: 10/1/20

Analyst: Anne Nord

Negative Urine Lot: 10120

GC/MS ID: 65198

(Optional Hydrolysis) Glucuronidase Lot: _____

Pre-Analytic:

- 1. *Positive Control Working Solution Preparation Instructions:*
Tube A positive control may be commercially obtained or prepared in-house. At a minimum, the control must contain at least one phenethylamine at an approximate concentration between 500 and 3000 ng/mL, and one opiate at an approximate concentration between 300 and 3000 ng/mL.
- 2. Verify Tune and Tune evaluation completed within the previous 7 days. Tune and Tune evaluation reports initialed and filed.
- 3. Create GCMS sequence to include controls, case blanks and case samples.

Analytic:

- 1. Remove working solutions, controls, and samples from cold storage.

(Optional Steps for Enzyme Hydrolysis- completed in addition to General extraction without Hydrolysis)

- 2a. In labeled round bottom Extraction tubes: add 4.5mL of case samples, and controls.
- 2b. Add 150uL of 2M acetate buffer, vortex.
- 2c. Add 100uL glucuronidase, cap and rock gently.
- 2d. Heat at 60C for 2 hours. Allow to cool before proceeding to step 3.
- 3. To each labeled De-Tox Tube add 5mL sample, Positive control: spike positive control working solution.
- 4. Place on tube rocker at ambient temp for approx. 10 minutes.
- 5. Centrifuge for approx. 10 min at ~2500-3000rpm.
- 6. Transfer solvent (upper layer) to new tube, and evaporate to ~100-300uL.
- 7. Transfer to labeled ALS vial with insert.
- 8. Place ALS Vials in appropriate location on GCMS rack and run using appropriate GCMS method.

Post-Analytic

- 1. Complete Data analysis on all samples and corresponding sample blanks
- 2. Did positive and negative control samples provide intended response? Yes
- 3. Sample Criteria for ID: RT +/- 0.2 min. (or 0.1 min. for phenethylamines)
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Working solution prep sheet(s), Positive control GCMS data printouts,

COMMENTS: Samples were extracted 9/25/20 the morphine in the positive control was not detected, the liner and septa on the GCMS was changed; the morphine still did not come out the samples were not evaluated in this run. The gold seal was changed and the front of the column was clipped. The sample^s were re-extracted on 9/29/20 the morphine still did not come out the samples were not evaluated in this run. The liner, septa, and gold seal were changed; the injection port was cleaned. I did a test 10/1/20 the salts have settled in the bottom of the de-tox tubes. I extracted the positive control using a new lot of tubes and ensuring the salts moved from the bottom of the tube. Morphine was detected. The samples were re-extracted 10/1/20 and evaluated. (the negative control used for the 9/25 and 9/29 extraction was lot 73020)

am 2 checklist

Toxicology AM method 2 control prep info 10000 ng/ml



working solution ~~20000~~ ng/ml in meoh methamphetamine, morphine
Stock solution 1mg/ml 100 ul each in 9800ul meOH

ppd 9/25/20: Exp: 9/25/21 lot 92520 by amn



Drug	lot	expiration
Methamphetamine	FE08101708	10/1/2022
Morphine	FE08221801	1/1/2024

AM 2 control add 500 ul working solution to 4500 ul negative urine and extract.
approximate concentration ~~2000~~ ng/ml
1000

File :D:\DATA\2020\am 2\100120 am2\00201001.D
Operator : Instrument 65198
Acquired : 01 Oct 2020 11:40 using AcqMethod TOXI-A 10115.M
Instrument : Instrument 65198 GCMS CdA
Sample Name: negative control
Misc Info : am 2
Vial Number: 2



File :D:\DATA\2020\am 2\100120 am2\00101002.D
Operator : Instrument 65198
Acquired : 01 Oct 2020 12:02 using AcqMethod TOXI-A 10115.M
Instrument : Instrument 65198 GCMS CdA
Sample Name: positive control
Misc Info : am 2
Vial Number: 1



File :D:\DATA\2020\am 2\100120 am2\00101002.D
Operator : Instrument 65198
Acquired : 01 Oct 2020 12:02 using AcqMethod TOXI-A 10115.M
Instrument : Instrument 65198 GCMS CdA
Sample Name: positive control
Misc Info : am 2
Vial Number: 1

