

Worklist: 1947

| <u>LAB CASE</u> | <u>ITEM</u> | <u>TASK ID</u> | <u>DESCRIPTION</u> | |
|-----------------|-------------|----------------|-------------------------------|--|
| M2017-3775 | 1 | 97419 | AM 25 Blood Multi-Drug Screer | |
| M2017-3896 | 1 | 97420 | AM 25 Blood Multi-Drug Screer | |
| M2017-3943 | 1 | 97421 | AM 25 Blood Multi-Drug Screer | |
| M2017-4078 | 1 | 97422 | AM 25 Blood Multi-Drug Screer | |
| M2017-4091 | 6 | 97423 | AM 25 Blood Multi-Drug Screer | |
| M2017-4285 | 1 | 97424 | AM 25 Blood Multi-Drug Screer | |
| P2017-2087 | 1 | 97425 | AM 25 Blood Multi-Drug Screer | |
| P2017-2088 | 1 | 97426 | AM 25 Blood Multi-Drug Screer | |
| P2017-2089 | 1 | 97427 | AM 25 Blood Multi-Drug Screer | |
| P2017-2090 | 1 | 97428 | AM 25 Blood Multi-Drug Screer | |
| P2017-2091 | 1 | 97429 | AM 25 Blood Multi-Drug Screer | |
| P2017-2092 | 1 | 97430 | AM 25 Blood Multi-Drug Screer | |
| P2017-2093 | 4 | 97431 | AM 25 Blood Multi-Drug Screer | |
| P2017-2118 | 1 | 97432 | AM 25 Blood Multi-Drug Screer | |
| P2017-2137 | 2 | 97433 | AM 25 Blood Multi-Drug Screer | |
| P2017-2161 | 1 | 97434 | AM 25 Blood Multi-Drug Screer | |
| P2017-2172 | 1 | 97435 | AM 25 Blood Multi-Drug Screer | |
| P2017-2227 | 1 | 97436 | AM 25 Blood Multi-Drug Screer | |
| P2017-2228 | 1 | 97437 | AM 25 Blood Multi-Drug Screer | |
| P2017-2231 | 1 | 97438 | AM 25 Blood Multi-Drug Screer | |
| P2017-2253 | 2 | 97439 | AM 25 Blood Multi-Drug Screer | |
| P2017-2254 | 1 | 97440 | AM 25 Blood Multi-Drug Screer | |
| P2017-2265 | 1 | 97441 | AM 25 Blood Multi-Drug Screer | |

Worklist: 1947

| <u>LAB CASE</u> | <u>ITEM</u> | <u>TASK ID</u> | <u>DESCRIPTION</u> |
|-----------------|-------------|----------------|-------------------------------|
| P2017-2266 | 1 | 97442 | AM 25 Blood Multi-Drug Screer |



P

SP

Multi-Drug Screen in Blood by LC-MS/MS

Extraction Date: 10/16/17

Analyst: Sarah Pickle

Plate lot#: 0495940

Plate Expiration: 12/12/17

Mobile phase A: 10mM Ammonium Formate

Mobile phase B 0.1% Formic acid in MeOH

0.5M Ammonium Hydroxide

Ethyl Acetate

Blank Blood Lot: 361331-3

Column: Phenomenex Phenyl Hexyl (4.6x50mm: 2.6 um)

LCMS-QQQ ID: 59740

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature. Each removed own samples.
- 2. Pipette **250µL blood (calibrated pipette)** Pipette ID: ~~067106~~ ⁸³⁵ in wells of analytical (standards) plate. Each pipetted own samples.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105* Performed by SP.
- 4. Pipette **250µL 0.5M ammonium hydroxide** buffer in wells of analytical plate. Performed by SP.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. Performed by SP.
- 6. Transfer **300µL of blood+base** mixture to corresponding wells of SLE+ plate. Each transferred own samples.
- 7. Apply positive pressure for approx. 4 seconds (or until no liquid remains on top of sorbent). Performed by SP.
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
Wait 5 minutes
- 8. Add **900µL ethyl acetate** and allow to flow for 5 minutes under gravity. Performed by SP.
- 9. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left) Performed by SP.
- 10. Add **900µL ethyl acetate** and allow to flow for 5 minutes under gravity. Performed by TS. §
- 11. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left) Performed by TS. §
- 12. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. Performed by TS. §
SPE Dry ID 067103
- 13. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist. Each reconstituted own samples.

POST- ANALYTIC

- 1. Create batch and process data.
Worklist path: C:\MassHunter\Data\2017\MD Screen\101617 SP and TS MD Screen Batch Name: 101617 MD Screen Panel 1 SP
101617 MD Screen Panel 2 SP
- 2. Evaluate samples, S/N of primary transition >5 and S/N of secondary transition >3 or evaluation of peak symmetry and resolution. Within +/- 2% or 0.1 min RT of administrative control. > 1/5 the response of administrative control.
- 3. Did all QCs pass for each analyte? Y (N) Did not evaluate Carisoprodol
Did not evaluate Norpropoxyphene §
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: and extraction

The analysis[^] was performed in conjunction with Tamara Salazar. Steps requiring pipetting or transferring of case samples was performed by the analyst with custody of the sample. The parts of the method not requiring direct interaction with case samples were split among Tamara and I with the other analyst observing.

J

Multidrug Screen External Control Preparation Information

Lot 053117

| Drug | Lot | Expiration date |
|------------|-----------------------|-----------------|
| Temazepam | Cerilliant FE08191401 | October 2019 |
| Quetiapine | Cerilliant FN02021502 | February 2020 |
| Zolpidem | Cerilliant FE12181402 | February 2020 |
| Zopiclone | Cerilliant FE12031501 | December 2020 |

Ext. control

100ul of 1mg/mL stock was added for each drug to 9960ul of LC MeOH (BDH Lot 15H250514). 100ul of the working solution was added to 9900ul of blood (Hemostat Lot 361331-3).

↳ W/S 092517

Created: 09/25/17

By: TS

Expires: 03/25/18

| Name | Transition | Acq. Date-Time | Pos | RT | Resp | S/N | Transition | S/N | Name | Transition | RT | Resp |
|-------------|-----------------|---------------------|--------|-------|--------|----------|-----------------|----------|-----------------|-----------------|-------|----------|
| Tramadol | 264.2 -> 58.3 | 10/17/2017 12:11 AM | P2-H12 | 5.154 | 209954 | Infinity | 264.2 -> 43.3 | 22.91 | Tramadol-13C-D3 | 268.2 -> 58.3 | 5.153 | 11735... |
| Tazodone | 372.2 -> 176... | 10/17/2017 12:11 AM | P2-H12 | 5.955 | 193194 | 1534.75 | 372.2 -> 148... | Infinity | Sufentanil-D5 | 392.2 -> 238... | 5.893 | 21442 |
| Venlafaxine | 278.2 -> 58.3 | 10/17/2017 12:11 AM | P2-H12 | 5.479 | 224816 | 31942... | 278.2 -> 260... | 245.21 | Venlafaxine-D6 | 284.2 -> 64.3 | 5.457 | 114256 |
| Zolpidem | 308.2 -> 235... | 10/17/2017 12:11 AM | P2-H12 | 5.738 | 161076 | Infinity | 308.2 -> 263... | 2279.94 | Zolpidem-D6 | 314.2 -> 235... | 5.716 | 736234 |

Did not evaluate carisoprodol.

2

| P1 negative-S1 | | P1 negative-S1 results | | Quarmer 1 M... | | Quarmer 2 M... | | IS1U Results | | | |
|----------------------|-----------------|------------------------|-------|----------------|-------|----------------|-----------------|--------------|--------------------------|-----------------|-------|
| Name | Transition | Acq. Date-Time | Pos | RT | Resp | S/N | Transition | S/N | Name | Transition | Resp |
| 6-MAM | 328.2 -> 165... | 10/16/2017 2:47 PM | P2-A5 | 4.635 | 271 | 105.05 | 328.2 -> 211... | | 6-MAM-D6 | 334.2 -> 165... | 4.635 |
| 7-Aminoclonazepam | 286.1 -> 121... | 10/16/2017 2:47 PM | P2-A5 | 5.313 | 33 | 0.88 | 286.1 -> 222... | | 7-Aminoclonazepam-D4 | 290.1 -> 121... | 5.231 |
| Acetyl Fentanyl | 323.2 -> 105... | 10/16/2017 2:47 PM | P2-A5 | | | | 323.2 -> 188... | | Acetyl Fentanyl-D5 | 328.2 -> 105... | 5.490 |
| Acetyl Norfentanyl | 219.1 -> 84.2 | 10/16/2017 2:47 PM | P2-A5 | | | | 219.1 -> 56.3 | | Acetyl Norfentanyl-D5 | 224.2 -> 84.2 | 4.609 |
| a-hydroxyalprazolam | 325.1 -> 297... | 10/16/2017 2:47 PM | P2-A5 | | | | 325.1 -> 215... | | a-hydroxyalprazolam-D5 | 330.1 -> 302... | 6.037 |
| alpha-PVP | 232.2 -> 91.0 | 10/16/2017 2:47 PM | P2-A5 | | | | 232.2 -> 77.1 | | alpha-PVP-d8 | 240.2 -> 91.1 | 5.238 |
| Alprazolam | 309.1 -> 281... | 10/16/2017 2:47 PM | P2-A5 | 6.165 | 223 | 0.96 | 309.1 -> 205... | | Alprazolam-D5 | 314.1 -> 286... | 6.102 |
| Amphetamine | 136.1 -> 91.1 | 10/16/2017 2:47 PM | P2-A5 | 4.714 | 24030 | 12.17 | 136.1 -> 119... | 2.32 | Amphetamine-D11 | 147.2 -> 130... | 4.531 |
| Bupropion | 240.1 -> 184... | 10/16/2017 2:47 PM | P2-A5 | | | | 240.1 -> 131... | | Ketamine-D4 | 242.1 -> 129... | 5.156 |
| Carisoprodol | 261.2 -> 176... | 10/16/2017 2:47 PM | P2-A5 | | | | 261.2 -> 55.3 | | Carisoprodol-D7 | 268.2 -> 183... | 5.765 |
| Citalopram | 325.2 -> 109... | 10/16/2017 2:47 PM | P2-A5 | 5.612 | 50 | 2.14 | 325.2 -> 262... | | Citalopram-D6 | 331.2 -> 109... | 5.692 |
| Clonazepam | 316.1 -> 269... | 10/16/2017 2:47 PM | P2-A5 | | | | 316.1 -> 213... | | Clonazepam-D4 | 320.1 -> 217... | 5.979 |
| Codeine | 300.2 -> 215... | 10/16/2017 2:47 PM | P2-A5 | | | | 300.2 -> 165... | | Codeine-D6 | 306.2 -> 218... | 4.515 |
| Cyclobenzaprine | 276.2 -> 214... | 10/16/2017 2:47 PM | P2-A5 | | | | 276.2 -> 231... | | Cyclobenzaprine-D3 | 279.2 -> 215... | 5.947 |
| Dextromethorphan | 272.2 -> 171... | 10/16/2017 2:47 PM | P2-A5 | | | | 272.2 -> 147... | | Dextromethorphan-D3 | 275.2 -> 171... | 5.723 |
| Dextroproprian | 258.2 -> 157... | 10/16/2017 2:47 PM | P2-A5 | | | | 258.2 -> 133... | | Dextroproprian-D3 | 261.2 -> 157... | 5.113 |
| Diazepam | 285.1 -> 193... | 10/16/2017 2:47 PM | P2-A5 | | | | 285.1 -> 153... | | Diazepam-D5 | 290.1 -> 198... | 6.329 |
| Dihydrocodeine | 302.2 -> 198... | 10/16/2017 2:47 PM | P2-A5 | | | | 302.2 -> 128... | | Dihydrocodeine-D6 | 308.2 -> 202... | 4.475 |
| Diphenhydramine | 256.2 -> 167... | 10/16/2017 2:47 PM | P2-A5 | | | | 256.2 -> 152... | | Diphenhydramine-D3 | 259.2 -> 167... | 5.665 |
| Doxylamine | 271.2 -> 182... | 10/16/2017 2:47 PM | P2-A5 | | | | 271.2 -> 167... | | Doxylamine-D5 | 276.2 -> 187... | 5.335 |
| EDDP | 279.2 -> 235... | 10/16/2017 2:47 PM | P2-A5 | | | | 279.2 -> 250... | | EDDP-D3 | 282.2 -> 235... | 5.721 |
| Fentanyl | 337.2 -> 105... | 10/16/2017 2:47 PM | P2-A5 | | | | 337.2 -> 188... | | Fentanyl-D5 | 342.3 -> 105... | 5.671 |
| Fluoxetine | 310.1 -> 44.3 | 10/16/2017 2:47 PM | P2-A5 | | | | 310.1 -> 148... | | Fluoxetine-D6 | 316.2 -> 44.3 | 5.899 |
| Hydrocodone | 300.2 -> 198... | 10/16/2017 2:47 PM | P2-A5 | | | | 300.2 -> 128... | | Hydrocodone-D6 | 306.2 -> 202... | 4.740 |
| Hydromorphone | 286.2 -> 184... | 10/16/2017 2:47 PM | P2-A5 | | | | 286.2 -> 157... | | Hydromorphone-D6 | 292.2 -> 185... | 4.157 |
| Ketamine | 238.1 -> 125... | 10/16/2017 2:47 PM | P2-A5 | 5.300 | 172 | 0.36 | 238.1 -> 220... | | Ketamine-D4 | 242.1 -> 129... | 5.156 |
| Meprobamate | 219.1 -> 158... | 10/16/2017 2:47 PM | P2-A5 | 5.301 | 48 | 12.34 | 219.1 -> 97.1 | | Carisoprodol-D7 | 268.2 -> 183... | 5.765 |
| Methadone | 310.2 -> 265... | 10/16/2017 2:47 PM | P2-A5 | 6.164 | 152 | 3.13 | 310.2 -> 105... | 1.16 | Methadone-D9 | 319.3 -> 268... | 5.980 |
| Methamphetamine | 150.1 -> 91.1 | 10/16/2017 2:47 PM | P2-A5 | | | | 150.1 -> 119... | | Methamphetamine-D11 | 161.2 -> 97.1 | 4.672 |
| Morphine | 286.2 -> 201... | 10/16/2017 2:47 PM | P2-A5 | | | | 286.2 -> 165... | | Morphine-D6 | 292.2 -> 152... | 3.894 |
| Naloxone | 328.2 -> 309... | 10/16/2017 2:47 PM | P2-A5 | 4.594 | 105 | 0.45 | 328.2 -> 211... | | Naltrexol-D3 | 347.2 -> 329... | 4.654 |
| Naltrexol | 344.2 -> 326... | 10/16/2017 2:47 PM | P2-A5 | 4.634 | 5244 | 6.64 | 344.2 -> 308... | | Naltrexol-D3 | 347.2 -> 329... | 4.654 |
| Naltrexone | 342.2 -> 324... | 10/16/2017 2:47 PM | P2-A5 | | | | 342.2 -> 55.3 | | Naltrexol-D3 | 347.2 -> 329... | 4.654 |
| Norbuprenorphine | 414.3 -> 101... | 10/16/2017 2:47 PM | P2-A5 | | | | 414.3 -> 57.3 | | Norbuprenorphine-D3 | 417.3 -> 101... | 5.506 |
| Nordiazepam | 271.1 -> 140... | 10/16/2017 2:47 PM | P2-A5 | | | | 271.1 -> 208... | | Nordiazepam-D5 | 276.1 -> 140... | 6.171 |
| Norfentanyl | 233.2 -> 84.2 | 10/16/2017 2:47 PM | P2-A5 | | | | 233.2 -> 55.3 | | Norfentanyl-D5 | 238.2 -> 84.2 | 5.055 |
| Norhydrocodone | 286.2 -> 199... | 10/16/2017 2:47 PM | P2-A5 | 4.662 | 133 | 1.02 | 286.2 -> 171... | | Norhydrocodone-D3 | 289.2 -> 202... | 4.681 |
| Noroxycodone | 302.1 -> 284... | 10/16/2017 2:47 PM | P2-A5 | 4.659 | 839 | 0.84 | 302.1 -> 186... | 60.07 | Noroxycodone-D3 | 305.2 -> 287... | 4.617 |
| O-desmethyl-tramadol | 250.2 -> 58.3 | 10/16/2017 2:47 PM | P2-A5 | 4.647 | 1203 | 0.87 | 250.2 -> 42.3 | 0.68 | O-desmethyl-tramadol-... | 256.2 -> 64.3 | 4.647 |
| Oxazepam | 287.1 -> 240... | 10/16/2017 2:47 PM | P2-A5 | | | | 287.1 -> 268... | | Oxazepam-D5 | 292.1 -> 246... | 6.025 |
| Oxycodone | 316.2 -> 298... | 10/16/2017 2:47 PM | P2-A5 | 4.656 | 1071 | 0.58 | 316.2 -> 241... | | Oxycodone-D6 | 322.2 -> 304... | 4.636 |
| Oxymorphone | 302.1 -> 284... | 10/16/2017 2:47 PM | P2-A5 | 3.952 | 92 | 0.64 | 302.1 -> 227... | 4.50 | Oxymorphone-D3 | 305.2 -> 287... | 3.971 |
| Phentermine | 150.1 -> 65.2 | 10/16/2017 2:47 PM | P2-A5 | 4.895 | 4400 | Infinity | 150.1 -> 133... | | Phentermine-D5 | 155.2 -> 96.2 | 4.854 |
| Promethazine | 285.1 -> 86.2 | 10/16/2017 2:47 PM | P2-A5 | | | | 285.1 -> 71.3 | | Promethazine-D3 | 288.2 -> 89.2 | 5.924 |
| Quetiapine | 384.2 -> 253... | 10/16/2017 2:47 PM | P2-A5 | | | | 384.2 -> 221... | | Quetiapine-D8 | 392.2 -> 226... | 5.913 |
| Serrraline | 306.1 -> 158... | 10/16/2017 2:47 PM | P2-A5 | | | | 306.1 -> 275... | | Serrraline-D3 | 309.1 -> 275... | 6.124 |
| Sufentanil | 387.2 -> 238... | 10/16/2017 2:47 PM | P2-A5 | 5.954 | 53 | 0.31 | 387.2 -> 111... | 202.01 | Sufentanil-D5 | 392.2 -> 238... | 5.933 |
| Temazepam | 301.1 -> 255... | 10/16/2017 2:47 PM | P2-A5 | 6.044 | 539 | 0.75 | 301.1 -> 283... | | Temazepam-D5 | 306.1 -> 260... | 6.166 |

| p1 negative-S1 | | | p1 negative-S1 results | | | Qualifier 1 M | | | Qualim | | | ISILV results | | |
|----------------|-----------------|--------------------|------------------------|-------|------|---------------|-----------------|------|-----------------|-----------------|-------|---------------|--|--|
| Name | Transition | Acq. Date-Time | Pos | RT | Resp | S/N | Transition | S/N | Name | Transition | RT | Resp | | |
| Tramadol | 264.2 -> 58.3 | 10/16/2017 2:47 PM | P2-A5 | 5.154 | 6640 | Infinity | 264.2 -> 43.3 | 1.30 | Tramadol-13C-D3 | 268.2 -> 58.3 | 5.153 | 12822... | | |
| Trazodone | 372.2 -> 176... | 10/16/2017 2:47 PM | P2-A5 | 5.813 | 137 | 0.21 | 372.2 -> 148... | 2.95 | Sufentanil-D5 | 392.2 -> 238... | 5.933 | 35068 | | |
| Venlafaxine | 278.2 -> 58.3 | 10/16/2017 2:47 PM | P2-A5 | | | | 278.2 -> 260... | | Venlafaxine-D6 | 284.2 -> 64.3 | 5.478 | 165066 | | |
| Zolpidem | 308.2 -> 235... | 10/16/2017 2:47 PM | P2-A5 | 5.881 | 105 | 0.89 | 308.2 -> 263... | | Zolpidem-D6 | 314.2 -> 235... | 5.737 | 868771 | | |

| Name | Transition | Acq. Date-Time | Pos. | RT | Resp. | S/N | Transition | S/N | Name | Transition | RT | Resp. |
|----------------------|-----------------|--------------------|-------|-------|--------|----------|-----------------|-----------|--------------------------|-----------------|-------|----------|
| 6-MAM | 328.2 -> 165... | 10/16/2017 3:07 PM | P2-B5 | | | | 328.2 -> 211... | | 6-MAM-D6 | 334.2 -> 165... | 4.655 | 38752 |
| 7-Aminoclonazepam | 286.1 -> 121... | 10/16/2017 3:07 PM | P2-B5 | | | | 286.1 -> 222... | | 7-Aminoclonazepam-D4 | 290.1 -> 121... | 5.292 | 3208 |
| Acetyl Fentanyl | 323.2 -> 105... | 10/16/2017 3:07 PM | P2-B5 | 5.512 | 40 | 0.40 | 323.2 -> 188... | | Acetyl Fentanyl-D5 | 328.2 -> 105... | 5.490 | 12361... |
| Acetyl Norfentanyl | 219.1 -> 84.2 | 10/16/2017 3:07 PM | P2-B5 | 4.630 | 83 | 0.64 | 219.1 -> 56.3 | | Acetyl Norfentanyl-D5 | 224.2 -> 84.2 | 4.630 | 343073 |
| a-hydroxyalprazolam | 325.1 -> 297... | 10/16/2017 3:07 PM | P2-B5 | 5.957 | 297 | 4.30 | 325.1 -> 215... | | a-hydroxyalprazolam-D5 | 330.1 -> 302... | 6.037 | 26649 |
| alpha-PVP | 232.2 -> 91.0 | 10/16/2017 3:07 PM | P2-B5 | | | | 232.2 -> 77.1 | | alpha-PVP-d8 | 240.2 -> 91.1 | 5.258 | 415940 |
| Alprazolam | 309.1 -> 281... | 10/16/2017 3:07 PM | P2-B5 | 6.185 | 88 | 0.60 | 309.1 -> 205... | 0.84 | Alprazolam-D5 | 314.1 -> 286... | 6.102 | 214247 |
| Amphetamine | 136.1 -> 91.1 | 10/16/2017 3:07 PM | P2-B5 | 4.734 | 26390 | 4.45 | 136.1 -> 119... | | Amphetamine-D11 | 147.2 -> 130... | 4.531 | 8152 |
| Bupropion | 240.1 -> 184... | 10/16/2017 3:07 PM | P2-B5 | | | | 240.1 -> 131... | | Ketamine-D4 | 242.1 -> 129... | 5.156 | 437930 |
| Carisoprodol | 261.2 -> 176... | 10/16/2017 3:07 PM | P2-B5 | | | | 261.2 -> 55.3 | | Carisoprodol-D7 | 268.2 -> 183... | 5.765 | 7917 |
| Citalopram | 325.2 -> 109... | 10/16/2017 3:07 PM | P2-B5 | 5.734 | 73 | 0.41 | 325.2 -> 262... | | Citalopram-D6 | 331.2 -> 109... | 5.692 | 51247 |
| Clonazepam | 316.1 -> 269... | 10/16/2017 3:07 PM | P2-B5 | | | | 316.1 -> 213... | | Clonazepam-D4 | 320.1 -> 217... | 5.959 | 10831 |
| Codeine | 300.2 -> 215... | 10/16/2017 3:07 PM | P2-B5 | | | | 300.2 -> 165... | | Codeine-D6 | 306.2 -> 218... | 4.536 | 28126 |
| Cyclobenzaprine | 276.2 -> 214... | 10/16/2017 3:07 PM | P2-B5 | 5.948 | 277 | 72.02 | 276.2 -> 231... | | Cyclobenzaprine-D3 | 279.2 -> 215... | 5.926 | 326693 |
| Dextromethorphan | 272.2 -> 171... | 10/16/2017 3:07 PM | P2-B5 | | | | 272.2 -> 147... | | Dextromethorphan-D3 | 275.2 -> 171... | 5.723 | 198288 |
| Dextrophan | 258.2 -> 157... | 10/16/2017 3:07 PM | P2-B5 | | | | 258.2 -> 133... | | Dextrophan-D3 | 261.2 -> 157... | 5.113 | 199394 |
| Diazepam | 285.1 -> 193... | 10/16/2017 3:07 PM | P2-B5 | 6.290 | 424 | 0.48 | 285.1 -> 153... | 0.75 | Diazepam-D5 | 290.1 -> 198... | 6.309 | 160120 |
| Dihydrocodeine | 302.2 -> 198... | 10/16/2017 3:07 PM | P2-B5 | | | | 302.2 -> 128... | | Dihydrocodeine-D6 | 308.2 -> 202... | 4.475 | 125090 |
| Diphenhydramine | 256.2 -> 167... | 10/16/2017 3:07 PM | P2-B5 | | | | 256.2 -> 152... | | Diphenhydramine-D3 | 259.2 -> 167... | 5.665 | 360133 |
| Doxylamine | 271.2 -> 182... | 10/16/2017 3:07 PM | P2-B5 | | | | 271.2 -> 167... | | Doxylamine-D5 | 276.2 -> 187... | 5.356 | 19792 |
| EDDP | 279.2 -> 235... | 10/16/2017 3:07 PM | P2-B5 | | | | 279.2 -> 250... | | EDDP-D3 | 282.2 -> 235... | 5.721 | 86419 |
| Fentanyl | 337.2 -> 105... | 10/16/2017 3:07 PM | P2-B5 | | | | 337.2 -> 188... | | Fentanyl-D5 | 342.3 -> 105... | 5.671 | 877032 |
| Fluoxetine | 310.1 -> 44.3 | 10/16/2017 3:07 PM | P2-B5 | 5.880 | 321 | 1.69 | 310.1 -> 148... | | Fluoxetine-D6 | 316.2 -> 44.3 | 5.899 | 11211... |
| Hydrocodone | 300.2 -> 198... | 10/16/2017 3:07 PM | P2-B5 | | | | 300.2 -> 128... | | Hydrocodone-D6 | 306.2 -> 202... | 4.740 | 182790 |
| Hydromorphone | 286.2 -> 184... | 10/16/2017 3:07 PM | P2-B5 | | | | 286.2 -> 157... | | Hydromorphone-D6 | 292.2 -> 185... | 4.157 | 47710 |
| Ketamine | 238.1 -> 125... | 10/16/2017 3:07 PM | P2-B5 | 5.055 | 1716 | 88.04 | 238.1 -> 220... | | Ketamine-D4 | 242.1 -> 129... | 5.156 | 437930 |
| Meprobamate | 219.1 -> 158... | 10/16/2017 3:07 PM | P2-B5 | 5.240 | 28 | 0.23 | 219.1 -> 97.1 | | Carisoprodol-D7 | 268.2 -> 183... | 5.765 | 7917 |
| Methadone | 310.2 -> 265... | 10/16/2017 3:07 PM | P2-B5 | 6.042 | 291 | 2.78 | 310.2 -> 105... | | Methadone-D9 | 319.3 -> 268... | 5.980 | 280697 |
| Methamphetamine | 150.1 -> 91.1 | 10/16/2017 3:07 PM | P2-B5 | | | | 150.1 -> 119... | | Methamphetamine-D11 | 161.2 -> 97.1 | 4.692 | 171103 |
| Morphine | 286.2 -> 201... | 10/16/2017 3:07 PM | P2-B5 | 3.854 | 31 | 0.49 | 286.2 -> 165... | 0.68 | Morphine-D6 | 292.2 -> 152... | 3.914 | 15514 |
| Naloxone | 328.2 -> 309... | 10/16/2017 3:07 PM | P2-B5 | 4.493 | 126 | 0.46 | 328.2 -> 211... | | Naloxone-D3 | 347.2 -> 329... | 4.674 | 57113 |
| Naltrexol | 344.2 -> 326... | 10/16/2017 3:07 PM | P2-B5 | 4.634 | 5014 | 11.82 | 344.2 -> 308... | | Naltrexol-D3 | 347.2 -> 329... | 4.674 | 57113 |
| Naltrexone | 342.2 -> 324... | 10/16/2017 3:07 PM | P2-B5 | | | | 342.2 -> 55.3 | | Naltrexone-D3 | 347.2 -> 329... | 4.674 | 57113 |
| Norbuprenorphine | 414.3 -> 101... | 10/16/2017 3:07 PM | P2-B5 | | | | 414.3 -> 57.3 | | Norbuprenorphine-D3 | 417.3 -> 101... | 5.506 | 4743 |
| Nordiazepam | 271.1 -> 140... | 10/16/2017 3:07 PM | P2-B5 | | | | 271.1 -> 208... | | Nordiazepam-D5 | 276.1 -> 140... | 6.171 | 22162 |
| Norfentanyl | 233.2 -> 84.2 | 10/16/2017 3:07 PM | P2-B5 | | | | 233.2 -> 55.3 | | Norfentanyl-D5 | 238.2 -> 84.2 | 5.055 | 293194 |
| Norhydrocodone | 286.2 -> 199... | 10/16/2017 3:07 PM | P2-B5 | 4.662 | 205 | 0.47 | 286.2 -> 171... | | Norhydrocodone-D3 | 289.2 -> 202... | 4.681 | 31815 |
| Noroxycodone | 302.1 -> 284... | 10/16/2017 3:07 PM | P2-B5 | | | | 302.1 -> 186... | | Noroxycodone-D3 | 305.2 -> 287... | 4.115 | 342 |
| O-desmethyl-tramadol | 250.2 -> 58.3 | 10/16/2017 3:07 PM | P2-B5 | 4.647 | 910 | 3.57 | 250.2 -> 42.3 | | O-desmethyl-tramadol-... | 256.2 -> 64.3 | 4.667 | 306164 |
| Oxazepam | 287.1 -> 240... | 10/16/2017 3:07 PM | P2-B5 | | | | 287.1 -> 268... | | Oxazepam-D5 | 292.1 -> 246... | 6.025 | 22186 |
| Oxycodone | 316.2 -> 298... | 10/16/2017 3:07 PM | P2-B5 | 4.656 | 1320 | 4.93 | 316.2 -> 241... | | Oxycodone-D6 | 322.2 -> 304... | 4.656 | 21206 |
| Oxymorphone | 302.1 -> 284... | 10/16/2017 3:07 PM | P2-B5 | | | | 302.1 -> 227... | | Oxymorphone-D3 | 305.2 -> 287... | 3.991 | 193881 |
| Phentermine | 150.1 -> 65.2 | 10/16/2017 3:07 PM | P2-B5 | 4.895 | 5662 | Infinity | 150.1 -> 133... | | Phentermine-D5 | 155.2 -> 96.2 | 4.854 | 427492 |
| Promethazine | 285.1 -> 86.2 | 10/16/2017 3:07 PM | P2-B5 | | | | 285.1 -> 71.3 | | Promethazine-D3 | 288.2 -> 89.2 | 5.924 | 368498 |
| Quetiapine | 384.2 -> 253... | 10/16/2017 3:07 PM | P2-B5 | 5.954 | 566930 | 1685.95 | 384.2 -> 221... | 376641... | Quetiapine-D8 | 392.2 -> 226... | 5.913 | 182876 |
| Serrraline | 306.1 -> 158... | 10/16/2017 3:07 PM | P2-B5 | 6.105 | 108 | 0.29 | 306.1 -> 275... | | Serrraline-D3 | 309.1 -> 275... | 6.124 | 225528 |
| Sufentanil | 387.2 -> 238... | 10/16/2017 3:07 PM | P2-B5 | | | | 387.2 -> 111... | | Sufentanil-D5 | 392.2 -> 238... | 5.933 | 32676 |
| Temazepam | 301.1 -> 255... | 10/16/2017 3:07 PM | P2-B5 | 6.167 | 502105 | Infinity | 301.1 -> 283... | Infinity | Temazepam-D5 | 306.1 -> 260... | 6.146 | 279789 |

12

| Compound Method | | P1 external control-SR | | | P1 external control-SR | | | P1 external control-SR | | | IS1U Method | | | IS1U Results | | |
|-----------------|-----------------|------------------------|-------|-------|------------------------|----------|-----------------|------------------------|-----------------|-----------------|-------------|----------|--|--------------|--|--|
| Name | Transition | Acq. Date-Time | Pos. | RT | Resp. | S/N | Transition | S/N | Name | Transition | RT | Resp. | | | | |
| Tramadol | 264.2 -> 58.3 | 10/16/2017 3:07 PM | P2-B5 | 5.174 | 6349 | Infinity | 264.2 -> 43.3 | 3.71 | Tramadol-13C-D3 | 268.2 -> 58.3 | 5.153 | 11753... | | | | |
| Trazodone | 372.2 -> 176... | 10/16/2017 3:07 PM | P2-B5 | 5.975 | 941 | 0.98 | 372.2 -> 148... | 1.06 | Sufentanil-D5 | 392.2 -> 238... | 5.933 | 32676 | | | | |
| Venlafaxine | 278.2 -> 58.3 | 10/16/2017 3:07 PM | P2-B5 | | | | 278.2 -> 260... | | Venlafaxine-D6 | 284.2 -> 64.3 | 5.478 | 155491 | | | | |
| Zolpidem | 308.2 -> 235... | 10/16/2017 3:07 PM | P2-B5 | 5.759 | 19624... | 12359... | 308.2 -> 263... | 1151.09 | Zolpidem-D6 | 314.2 -> 235... | 5.737 | 839092 | | | | |

| p2 cai 10-15 | | | | p2 cai 10-15 Results | | | | IS11D Method | | | | | | |
|-----------------|------------------------|-------------------------|----------------------------|----------------------|--------------|--------------|-----------------|---------------------------|----------------|----------------------------|-------------------------|--------------|--------------|------|
| Compound Method | Name | Transition | Acq. Date-Time | Pos | RT | Resp | S/N | Quamier 1 M... | S/N | Quam... | Name | Transition | RT | Resp |
| | 7-aminoflunitrazepam | 284.1 -> 135... | 10/17/2017 12:20 AM | P2-H12 | 5.404 | 30040 | 173.90 | 284.1 -> 93.2 | 25.57 | 7-aminoflunitrazepam-D7 | 291.2 -> 138... | 5.403 | 474576 | |
| | alpha-hydroxymidazolam | 342.1 -> 324... | 10/17/2017 12:20 AM | P2-H12 | 6.048 | 47340 | Infinity | 342.1 -> 203... | 53.28 | alpha-hydroxymidazolam-... | 346.1 -> 328... | 6.048 | 112201 | |
| | Amirtipryline | 278.2 -> 91.2 | 10/17/2017 12:20 AM | P2-H12 | 5.997 | 156621 | 22.73 | 278.2 -> 105... | 43.68 | Amirtipryline-D3 | 281.2 -> 91.2 | 5.996 | 427006 | |
| | Amoxapine | 314.1 -> 271... | 10/17/2017 12:20 AM | P2-H12 | 5.949 | 10890 | 35.57 | 314.1 -> 193... | 174.41 | Amirtipryline-D3 | 281.2 -> 91.2 | 5.996 | 427006 | |
| | Benzoylcegonine | 290.1 -> 168... | 10/17/2017 12:20 AM | P2-H12 | 5.058 | 7115 | 698.79 | 290.1 -> 105... | 77.36 | Benzoylcegonine-D8 | 298.2 -> 171... | 5.057 | 3162 | |
| | Carbamazepine | 237.1 -> 193... | 10/17/2017 12:20 AM | P2-H12 | 5.759 | 21190 | 26.78 | 237.1 -> 194... | 125.87 | Carbamazepine-13C6 | 243.1 -> 200... | 5.779 | 806281 | |
| | Chlordiazepoxide | 300.1 -> 227... | 10/17/2017 12:20 AM | P2-H12 | 6.177 | 8010 | 17.58 | 300.1 -> 282... | 19.28 | Chlordiazepoxide-D5 | 305.1 -> 286... | 6.176 | 210933 | |
| | Clomipramine | 315.2 -> 86.2 | 10/17/2017 12:20 AM | P2-H12 | 6.174 | 267713 | 541.47 | 315.2 -> 58.3 | 237.07 | Clomipramine-D3 | 318.2 -> 89.2 | 6.173 | 334238 | |
| | Cocaine | 304.2 -> 182... | 10/17/2017 12:20 AM | P2-H12 | 5.279 | 85475 | 373.02 | 304.2 -> 105... | 104.47 | Cocaine-D3 | 307.2 -> 185... | 5.299 | 211146 | |
| | Desipramine | 267.2 -> 72.3 | 10/17/2017 12:20 AM | P2-H12 | 5.937 | 154950 | 335.75 | 267.2 -> 44.3 | Infinity | Desipramine-D3 | 270.2 -> 75.3 | 5.937 | 70646 | |
| | Doxepin | 280.2 -> 107... | 10/17/2017 12:20 AM | P2-H12 | 5.753 | 12535 | 48.27 | 280.2 -> 77.2 | Infinity | Doxepin-D3 | 283.2 -> 107... | 5.752 | 2968 | |
| | Estazolam | 295.1 -> 267... | 10/17/2017 12:20 AM | P2-H12 | 6.014 | 17366 | 125.31 | 295.1 -> 205... | 122.05 | Estazolam-D5 | 300.1 -> 272... | 6.013 | 353711 | |
| | Flunitrazepam | 314.1 -> 268... | 10/17/2017 12:20 AM | P2-H12 | 6.051 | 73782 | 8941... | 314.1 -> 239... | 62.17 | Flunitrazepam-D7 | 321.1 -> 245... | 6.030 | 23456 | |
| | Flurazepam | 388.2 -> 315... | 10/17/2017 12:20 AM | P2-H12 | 5.742 | 150199 | Infinity | 388.2 -> 317... | 17.96 | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 10140 | |
| | Imipramine | 281.2 -> 86.2 | 10/17/2017 12:20 AM | P2-H12 | 5.956 | 173136 | 273.79 | 281.2 -> 58.3 | 346.22 | Imipramine-D3 | 284.2 -> 61.3 | 5.955 | 440825 | |
| | Levamisole | 205.1 -> 91.2 | 10/17/2017 12:20 AM | P2-H12 | 4.704 | 38888 | 99.56 | 205.1 -> 178... | 97.35 | Cocaine-D3 | 307.2 -> 185... | 5.299 | 211146 | |
| | Maprotiline | 278.2 -> 91.2 | 10/17/2017 12:20 AM | P2-H12 | 5.692 | 106965 | 206.08 | 278.2 -> 117... | Infinity | Meperidine-D4 | 252.2 -> 224... | 5.287 | 36586 | |
| | MDA | 180.1 -> 163... | 10/17/2017 12:20 AM | P2-H12 | 4.727 | 52229 | 36.69 | 180.1 -> 105... | 25.05 | MDA-D5 | 185.1 -> 168... | 4.706 | 88263 | |
| | MDEA | 208.1 -> 163... | 10/17/2017 12:20 AM | P2-H12 | 4.968 | 92567 | 209.84 | 208.1 -> 105... | 121.46 | MDEA-D6 | 214.2 -> 166... | 4.967 | 210837 | |
| | MDMA | 194.1 -> 163... | 10/17/2017 12:20 AM | P2-H12 | 4.827 | 101626 | 126.59 | 194.1 -> 105... | 139.66 | MDMA-D6 | 200.2 -> 166... | 4.806 | 10548 | |
| | Meperidine | 248.2 -> 220... | 10/17/2017 12:20 AM | P2-H12 | 5.288 | 39383 | 106.74 | 248.2 -> 174... | 46.94 | Meperidine-D4 | 252.2 -> 224... | 5.287 | 36586 | |
| | Methyphenidate | 234.2 -> 84.2 | 10/17/2017 12:20 AM | P2-H12 | 5.229 | 150247 | 58.67 | 234.2 -> 56.3 | 252.60 | Meperidine-D4 | 252.2 -> 224... | 5.287 | 36586 | |
| | Metoprolol | 268.2 -> 116... | 10/17/2017 12:20 AM | P2-H12 | 5.162 | 22622 | 102.45 | 268.2 -> 56.3 | 78.75 | Meperidine-D4 | 252.2 -> 224... | 5.287 | 36586 | |
| | Midazolam | 326.1 -> 223... | 10/17/2017 12:20 AM | P2-H12 | 6.192 | 24214 | 46.51 | 326.1 -> 249... | 178.62 | Midazolam-D4 | 330.1 -> 295... | 6.191 | 343675 | |
| | Mirtazapine | 266.2 -> 195... | 10/17/2017 12:20 AM | P2-H12 | 5.490 | 46291 | 118.69 | 266.2 -> 72.3 | 62.85 | Midazolam-D4 | 330.1 -> 295... | 6.191 | 343675 | |
| | Mirtazapine | 399.2 -> 174... | 10/17/2017 12:20 AM | P2-H12 | 5.761 | 59091 | 516.24 | 399.2 -> 238... | 62.34 | Carbamazepine-13C6 | 243.1 -> 200... | 5.779 | 806281 | |
| | Normeperidine | 234.1 -> 42.3 | 10/17/2017 12:20 AM | P2-H12 | 5.290 | 45356 | 143.25 | 234.1 -> 160... | 52.68 | Normeperidine-D4 | 238.2 -> 164... | 5.289 | 241168 | |
| | Norpropoxyphene | 326.2 -> 44.3 | 10/17/2017 12:20 AM | P2-H12 | 5.846 | 13507 | Infinity | 326.2 -> 252... | 1116... | Norpropoxyphene-D5 | 331.2 -> 44.3 | 5.825 | 53687 | |
| | Nortriptyline | 264.2 -> 91.2 | 10/17/2017 12:20 AM | P2-H12 | 5.897 | 331629 | 348.32 | 264.2 -> 232... | 1116... | Nortriptyline-D3 | 267.2 -> 233... | 5.978 | 154746 | |
| | Ondansetron | 294.2 -> 184... | 10/17/2017 12:20 AM | P2-H12 | 5.443 | 71794 | 579.61 | 294.2 -> 123... | 33.17 | Pentazocine-13C3 | 289.2 -> 72.3 | 5.362 | 506170 | |
| | Pentazocine | 286.2 -> 218... | 10/17/2017 12:20 AM | P2-H12 | 5.362 | 43291 | 117.08 | 286.2 -> 41.3 | 115.17 | Pentazocine-13C3 | 289.2 -> 72.3 | 5.362 | 506170 | |
| | Phenazepam | 349.0 -> 206... | 10/17/2017 12:20 AM | P2-H12 | 6.108 | 7815 | 103.63 | 349.0 -> 183... | 31.82 | Phenazepam-D4 | 353.0 -> 183... | 6.087 | 5658 | |
| | Phencyclidine | 244.2 -> 86.2 | 10/17/2017 12:20 AM | P2-H12 | 5.575 | 13711 | 28.51 | 244.2 -> 91.1 | 210.06 | Phencyclidine-D5 | 249.2 -> 86.2 | 5.554 | 216754 | |
| | Phenylethylamine | 253.1 -> 104... | 10/17/2017 12:20 AM | P2-H12 | 5.675 | 836 | 224.14 | 253.1 -> 182... | 64.36 | Phenylethylamine-D10 | 263.2 -> 192... | 5.655 | 2900 | |
| | Prazepam | 325.1 -> 271... | 10/17/2017 12:20 AM | P2-H12 | 6.577 | 151997 | 539.73 | 325.1 -> 140... | 510.02 | Prazepam-D5 | 330.1 -> 276... | 6.576 | 12664... | |
| | Primidone | 219.1 -> 91.2 | 10/17/2017 12:20 AM | P2-H12 | 5.129 | 3603 | 6.68 | 219.1 -> 162... | 11.72 | Prazepam-D5 | 330.1 -> 276... | 6.576 | 12664... | |
| | Propoxyphene | 340.2 -> 58.3 | 10/17/2017 12:20 AM | P2-H12 | 5.845 | 202384 | 417.90 | 340.2 -> 266... | 235.01 | Propoxyphene-D11 | 351.3 -> 277... | 5.824 | 3742 | |
| | Protriptyline | 264.2 -> 155... | 10/17/2017 12:20 AM | P2-H12 | 5.897 | 630854 | Infinity | 264.2 -> 161... | Infinity | Nortriptyline-D3 | 267.2 -> 233... | 5.978 | 154746 | |
| | Pseudoephedrine | 166.1 -> 148... | 10/17/2017 12:20 AM | P2-H12 | 4.425 | 1792205 | 59.92 | 166.1 -> 133... | Infinity | Pseudoephedrine-D3 | 169.1 -> 151... | 4.424 | 666427 | |
| | Tapentadol | 222.2 -> 107... | 10/17/2017 12:20 AM | P2-H12 | 5.149 | 96664 | 120.51 | 222.2 -> 121... | 68.47 | Tapentadol-D3 | 225.2 -> 107... | 5.149 | 253158 | |
| | Trimipramine | 295.2 -> 100... | 10/17/2017 12:20 AM | P2-H12 | 6.034 | 190641 | 833.66 | 295.2 -> 58.3 | 124.91 | Trimipramine-D3 | 298.2 -> 103... | 6.034 | 12526... | |
| | Verapamil | 455.3 -> 165... | 10/17/2017 12:20 AM | P2-H12 | 5.801 | 204171 | Infinity | 455.3 -> 150... | 461.60 | Trimipramine-D3 | 298.2 -> 103... | 6.034 | 12526... | |
| | Zaleplon | 306.1 -> 236... | 10/17/2017 12:20 AM | P2-H12 | 5.849 | 44544 | 313.98 | 306.1 -> 264... | 246.38 | Zaleplon-D4 | 310.2 -> 240... | 5.848 | 34010 | |
| | Zopiclone | 389.1 -> 244... | 10/17/2017 12:20 AM | P2-H12 | 5.620 | 7802 | 25.53 | 389.1 -> 216... | 701.48 | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 10140 | |

Did not evaluate norpropoxyphene.

| p2 negative-Sr | | | | p2 negative-Sr Results | | | | ISUD Method | | | | ISUD Results | | | |
|------------------------|-----------------|--------------------|-------|------------------------|-------|------|-----------------|-------------|-------|----------------------------|-----------------|--------------|----------|--|--|
| Name | Transition | Acq. Date-Time | Pos. | RT | Resp. | S/N | Transition | Qualit. | S/N | Name | Transition | RT | Resp. | | |
| 7-aminoflunitrazepam | 284.1 -> 135... | 10/16/2017 2:57 PM | P2-A5 | 5.282 | 82 | 0.79 | 284.1 -> 93.2 | | | 7-aminoflunitrazepam-D7 | 291.2 -> 138... | 5.423 | 538831 | | |
| alpha-hydroxymidazolam | 342.1 -> 324... | 10/16/2017 2:57 PM | P2-A5 | 5.946 | 238 | 1.58 | 342.1 -> 203... | 0.40 | 0.40 | alpha-hydroxymidazolam-... | 346.1 -> 328... | 6.088 | 134063 | | |
| Amiripiline | 278.2 -> 91.2 | 10/16/2017 2:57 PM | P2-A5 | 6.120 | 3306 | 1.94 | 278.2 -> 105... | 0.54 | 0.54 | Amiripiline-D3 | 281.2 -> 91.2 | 5.996 | 504023 | | |
| Amoxapine | 314.1 -> 271... | 10/16/2017 2:57 PM | P2-A5 | 6.153 | 111 | 0.45 | 314.1 -> 193... | 15.00 | 15.00 | Amiripiline-D3 | 281.2 -> 91.2 | 5.996 | 504023 | | |
| Benzylecgonine | 290.1 -> 168... | 10/16/2017 2:57 PM | P2-A5 | 5.139 | 200 | 0.91 | 290.1 -> 105... | 2.66 | 2.66 | Benzylecgonine-D8 | 298.2 -> 171... | 5.077 | 3912 | | |
| Carbamazepine | 237.1 -> 193... | 10/16/2017 2:57 PM | P2-A5 | 5.963 | 432 | 2.19 | 237.1 -> 194... | 0.62 | 0.62 | Carbamazepine-D8 | 243.1 -> 200... | 5.799 | 975096 | | |
| Chlordiazepoxide | 300.1 -> 227... | 10/16/2017 2:57 PM | P2-A5 | 6.034 | 141 | 5.36 | 300.1 -> 282... | Infini... | | Chlordiazepoxide-D5 | 305.1 -> 286... | 6.196 | 328872 | | |
| Clomipramine | 315.2 -> 86.2 | 10/16/2017 2:57 PM | P2-A5 | 6.498 | 21 | 1.15 | 315.2 -> 58.3 | 0.40 | 0.40 | Clomipramine-D3 | 318.2 -> 89.2 | 6.173 | 413119 | | |
| Cocaine | 304.2 -> 182... | 10/16/2017 2:57 PM | P2-A5 | 5.401 | 364 | 3.72 | 304.2 -> 105... | | | Cocaine-D3 | 307.2 -> 185... | 5.299 | 296938 | | |
| Desipramine | 267.2 -> 72.3 | 10/16/2017 2:57 PM | P2-A5 | 5.978 | 1100 | 6.10 | 267.2 -> 44.3 | 13.70 | 13.70 | Desipramine-D3 | 270.2 -> 75.3 | 5.957 | 111241 | | |
| Doxepin | 280.2 -> 107... | 10/16/2017 2:57 PM | P2-A5 | 5.956 | 14 | 0.27 | 280.2 -> 77.2 | 8.11 | 8.11 | Doxepin-D3 | 283.2 -> 107... | 5.772 | 6196 | | |
| Estazolam | 295.1 -> 267... | 10/16/2017 2:57 PM | P2-A5 | | | | 295.1 -> 205... | | | Estazolam-D5 | 300.1 -> 272... | 6.034 | 417997 | | |
| Flunitrazepam | 314.1 -> 268... | 10/16/2017 2:57 PM | P2-A5 | 6.113 | 1479 | 2.89 | 314.1 -> 239... | 1.00 | 1.00 | Flunitrazepam-D7 | 321.1 -> 245... | 6.071 | 26943 | | |
| Flurazepam | 388.2 -> 315... | 10/16/2017 2:57 PM | P2-A5 | 5.762 | 181 | 0.65 | 388.2 -> 317... | | | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 22053 | | |
| Imipramine | 281.2 -> 86.2 | 10/16/2017 2:57 PM | P2-A5 | 6.017 | 1746 | 7.35 | 281.2 -> 58.3 | 4.51 | 4.51 | Imipramine-D3 | 284.2 -> 61.3 | 5.975 | 519005 | | |
| Levamisole | 205.1 -> 91.2 | 10/16/2017 2:57 PM | P2-A5 | 4.745 | 100 | 0.68 | 205.1 -> 178... | | | Cocaine-D3 | 307.2 -> 185... | 5.299 | 296938 | | |
| Maprotiline | 278.2 -> 91.2 | 10/16/2017 2:57 PM | P2-A5 | 5.733 | 35 | 0.79 | 278.2 -> 117... | | | Meperidine-D4 | 282.2 -> 224... | 5.287 | 42823 | | |
| MDA | 180.1 -> 163... | 10/16/2017 2:57 PM | P2-A5 | 5.112 | 545 | 0.71 | 180.1 -> 105... | 1.27 | 1.27 | MDA-D5 | 185.1 -> 168... | 4.706 | 100415 | | |
| MDEA | 208.1 -> 163... | 10/16/2017 2:57 PM | P2-A5 | 4.886 | 511 | 0.44 | 208.1 -> 105... | 0.15 | 0.15 | MDEA-D6 | 214.2 -> 166... | 4.967 | 267106 | | |
| MDMA | 194.1 -> 163... | 10/16/2017 2:57 PM | P2-A5 | 4.887 | 428 | 0.84 | 194.1 -> 105... | 0.43 | 0.43 | MDMA-D6 | 200.2 -> 166... | 4.827 | 12650 | | |
| Meperidine | 248.2 -> 220... | 10/16/2017 2:57 PM | P2-A5 | 5.615 | 1884 | 4.33 | 248.2 -> 174... | 0.30 | 0.30 | Meperidine-D4 | 282.2 -> 224... | 5.287 | 42823 | | |
| Methylphenidate | 234.2 -> 84.2 | 10/16/2017 2:57 PM | P2-A5 | 5.637 | 232 | 0.73 | 234.2 -> 56.3 | 0.64 | 0.64 | Meperidine-D4 | 282.2 -> 224... | 5.287 | 42823 | | |
| Metoprolol | 268.2 -> 116... | 10/16/2017 2:57 PM | P2-A5 | 4.798 | 74 | 0.71 | 268.2 -> 56.3 | 4.01 | 4.01 | Meperidine-D4 | 282.2 -> 224... | 5.287 | 42823 | | |
| Midazolam | 326.1 -> 223... | 10/16/2017 2:57 PM | P2-A5 | 6.395 | 90 | 1.52 | 326.1 -> 249... | 0.82 | 0.82 | Midazolam-D4 | 330.1 -> 295... | 6.212 | 422168 | | |
| Mirtazapine | 266.2 -> 195... | 10/16/2017 2:57 PM | P2-A5 | 5.428 | 56 | 0.51 | 266.2 -> 72.3 | 0.68 | 0.68 | Midazolam-D4 | 330.1 -> 295... | 6.212 | 422168 | | |
| Mirtagynine | 399.2 -> 174... | 10/16/2017 2:57 PM | P2-A5 | | | | 399.2 -> 238... | | | Carbamazepine-13C6 | 243.1 -> 200... | 5.799 | 975096 | | |
| Norpropoxyphene | 234.1 -> 42.3 | 10/16/2017 2:57 PM | P2-A5 | 5.168 | 56 | 0.53 | 234.1 -> 160... | 0.73 | 0.73 | Norpropoxyphene-D4 | 238.2 -> 164... | 5.310 | 278797 | | |
| Nortriptyline | 326.2 -> 44.3 | 10/16/2017 2:57 PM | P2-A5 | 5.887 | 33 | 0.31 | 326.2 -> 252... | | | Norpropoxyphene-D5 | 331.2 -> 44.3 | 5.866 | 92573 | | |
| Ondansetron | 294.2 -> 91.2 | 10/16/2017 2:57 PM | P2-A5 | 5.836 | 82 | 0.26 | 294.2 -> 233... | | | Nortriptyline-D3 | 287.2 -> 233... | 5.998 | 189720 | | |
| Pentazocine | 286.2 -> 184... | 10/16/2017 2:57 PM | P2-A5 | 5.709 | 218 | 0.28 | 294.2 -> 212... | 0.60 | 0.60 | Pentazocine-13C3 | 289.2 -> 72.3 | 5.362 | 590615 | | |
| Phenazepam | 349.0 -> 206... | 10/16/2017 2:57 PM | P2-A5 | 6.068 | 36 | 0.64 | 286.2 -> 41.3 | 3.96 | 3.96 | Pentazocine-13C3 | 289.2 -> 72.3 | 5.362 | 590615 | | |
| Phencyclidine | 244.2 -> 86.2 | 10/16/2017 2:57 PM | P2-A5 | | 601 | 3.02 | 349.0 -> 183... | | | Phenazepam-D4 | 353.0 -> 183... | 6.128 | 5557 | | |
| Phenylethanolamine | 253.1 -> 104... | 10/16/2017 2:57 PM | P2-A5 | | | | 244.2 -> 91.1 | | | Phencyclidine-D5 | 249.2 -> 86.2 | 5.574 | 246764 | | |
| Prazeepam | 325.1 -> 271... | 10/16/2017 2:57 PM | P2-A5 | 6.638 | 112 | 0.79 | 325.1 -> 182... | | | Phenylethanolamine-D5 | 263.2 -> 192... | 5.696 | 4620 | | |
| Primidone | 219.1 -> 91.2 | 10/16/2017 2:57 PM | P2-A5 | 5.271 | 368 | 1.76 | 219.1 -> 162... | | | Prazeepam-D5 | 330.1 -> 276... | 6.555 | 14942... | | |
| Propoxyphene | 340.2 -> 58.3 | 10/16/2017 2:57 PM | P2-A5 | 5.825 | 160 | 0.41 | 340.2 -> 266... | | | Prazeepam-D5 | 330.1 -> 276... | 6.555 | 14942... | | |
| Protriptyline | 264.2 -> 155... | 10/16/2017 2:57 PM | P2-A5 | | | | 264.2 -> 161... | | | Propoxyphene-D11 | 351.3 -> 277... | 5.864 | 5229 | | |
| Pseudoephedrine | 166.1 -> 148... | 10/16/2017 2:57 PM | P2-A5 | 4.707 | 723 | 0.94 | 166.1 -> 133... | | | Nortriptyline-D3 | 267.2 -> 233... | 5.998 | 189720 | | |
| Tapentadol | 222.2 -> 107... | 10/16/2017 2:57 PM | P2-A5 | 5.618 | 753 | 0.73 | 222.2 -> 121... | 0.96 | 0.96 | Pseudoephedrine-D3 | 169.1 -> 151... | 4.404 | 682432 | | |
| Trimipramine | 295.2 -> 100... | 10/16/2017 2:57 PM | P2-A5 | | | | 295.2 -> 58.3 | | | Tapentadol-D3 | 225.2 -> 107... | 5.169 | 308957 | | |
| Verapamil | 455.3 -> 165... | 10/16/2017 2:57 PM | P2-A5 | 5.923 | 182 | 0.90 | 455.3 -> 150... | | | Trimipramine-D3 | 298.2 -> 103... | 6.034 | 15239... | | |
| Zaleplon | 306.1 -> 236... | 10/16/2017 2:57 PM | P2-A5 | 5.686 | 198 | 0.72 | 306.1 -> 264... | | | Trimipramine-D3 | 298.2 -> 103... | 6.034 | 15239... | | |
| Zopiclone | 389.1 -> 244... | 10/16/2017 2:57 PM | P2-A5 | | | | 389.1 -> 216... | | | Zaleplon-D4 | 310.2 -> 240... | 5.889 | 39183 | | |
| | | | | | | | | | | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 22053 | | |

Compound Method

P2 external control-S1

P2 external control-S1 H...

Qualifier 1 M...

Qualim

IS11U Method

IS11U Results

| Name | Transition | Acq. Date-Time | Pos. | RT | Resp. | S/N | Transition | S/N | Name | Transition | RT | Resp. |
|------------------------|-----------------|--------------------|-------|-------|-------|--------|-----------------|--------|----------------------------|-----------------|-------|----------|
| 7-aminoflunitrazepam | 284.1 -> 135... | 10/16/2017 3:16 PM | P2-B5 | 5.445 | 490 | 1.97 | 284.1 -> 93.2 | | 7-aminoflunitrazepam-D7 | 291.2 -> 138... | 5.423 | 561590 |
| alpha-hydroxymidazolam | 342.1 -> 324... | 10/16/2017 3:16 PM | P2-B5 | 6.170 | 44 | 0.42 | 342.1 -> 203... | 3.75 | alpha-hydroxymidazolam-... | 346.1 -> 328... | 6.068 | 141036 |
| Amiripryline | 278.2 -> 91.2 | 10/16/2017 3:16 PM | P2-B5 | 6.100 | 2333 | 2.84 | 278.2 -> 105... | 1.14 | Amiripryline-D3 | 281.2 -> 91.2 | 5.976 | 492101 |
| Amoxapine | 314.1 -> 271... | 10/16/2017 3:16 PM | P2-B5 | 6.092 | 101 | 0.72 | 314.1 -> 193... | 161.04 | Amiripryline-D3 | 281.2 -> 91.2 | 5.976 | 492101 |
| Benzoylcegonine | 290.1 -> 168... | 10/16/2017 3:16 PM | P2-B5 | 5.058 | 278 | 3.25 | 290.1 -> 105... | 1.16 | Benzoylcegonine-D8 | 298.2 -> 171... | 5.057 | 5035 |
| Carbamazepine | 237.1 -> 193... | 10/16/2017 3:16 PM | P2-B5 | 5.310 | 1782 | 3.92 | 237.1 -> 194... | 1.46 | Carbamazepine-13C6 | 243.1 -> 200... | 5.799 | 930164 |
| Chlordiazepoxide | 300.1 -> 227... | 10/16/2017 3:16 PM | P2-B5 | 6.013 | 1148 | 2.55 | 300.1 -> 282... | 10.92 | Chlordiazepoxide-D5 | 305.1 -> 286... | 6.176 | 305981 |
| Clomipramine | 315.2 -> 86.2 | 10/16/2017 3:16 PM | P2-B5 | 6.194 | 54 | 0.30 | 315.2 -> 58.3 | 0.74 | Clomipramine-D3 | 318.2 -> 89.2 | 6.153 | 389736 |
| Cocaine | 304.2 -> 182... | 10/16/2017 3:16 PM | P2-B5 | 5.320 | 895 | 2.89 | 304.2 -> 105... | 0.62 | Cocaine-D3 | 307.2 -> 185... | 5.299 | 264869 |
| Desipramine | 267.2 -> 72.3 | 10/16/2017 3:16 PM | P2-B5 | 5.978 | 1062 | 4.14 | 267.2 -> 44.3 | 5.75 | Desipramine-D3 | 270.2 -> 75.3 | 5.937 | 105329 |
| Doxepin | 280.2 -> 107... | 10/16/2017 3:16 PM | P2-B5 | 5.487 | 88 | 0.61 | 280.2 -> 77.2 | 22.19 | Doxepin-D3 | 283.2 -> 107... | 5.732 | 3349 |
| Estazolam | 295.1 -> 267... | 10/16/2017 3:16 PM | P2-B5 | 5.627 | 212 | 3.78 | 295.1 -> 205... | 54.59 | Estazolam-D5 | 300.1 -> 272... | 6.013 | 414007 |
| Flunitrazepam | 314.1 -> 268... | 10/16/2017 3:16 PM | P2-B5 | 6.092 | 1824 | 3.68 | 314.1 -> 239... | 1.53 | Flunitrazepam-D7 | 321.1 -> 245... | 6.051 | 27241 |
| Flurazepam | 388.2 -> 315... | 10/16/2017 3:16 PM | P2-B5 | 5.925 | 47 | 0.45 | 388.2 -> 317... | | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 16782 |
| Imipramine | 281.2 -> 86.2 | 10/16/2017 3:16 PM | P2-B5 | 5.996 | 1036 | 6.35 | 281.2 -> 58.3 | 2.58 | Imipramine-D3 | 284.2 -> 61.3 | 5.955 | 488860 |
| Levamisole | 205.1 -> 91.2 | 10/16/2017 3:16 PM | P2-B5 | 5.110 | 1236 | 0.76 | 205.1 -> 178... | 0.74 | Cocaine-D3 | 307.2 -> 185... | 5.299 | 264869 |
| Maprotiline | 278.2 -> 91.2 | 10/16/2017 3:16 PM | P2-B5 | 5.651 | 372 | 2.68 | 278.2 -> 117... | | Meperidine-D4 | 252.2 -> 224... | 5.287 | 42002 |
| MDA | 180.1 -> 163... | 10/16/2017 3:16 PM | P2-B5 | | | | 180.1 -> 105... | | MDA-D5 | 185.1 -> 168... | 4.706 | 96810 |
| MDEA | 208.1 -> 163... | 10/16/2017 3:16 PM | P2-B5 | | | | 208.1 -> 105... | | MDEA-D6 | 214.2 -> 166... | 4.967 | 260051 |
| MDMA | 194.1 -> 163... | 10/16/2017 3:16 PM | P2-B5 | 4.827 | 686 | 0.81 | 194.1 -> 105... | | MDMA-D6 | 200.2 -> 166... | 4.806 | 11570 |
| Meperidine | 248.2 -> 220... | 10/16/2017 3:16 PM | P2-B5 | 5.635 | 3131 | 20.52 | 248.2 -> 174... | 0.55 | Meperidine-D4 | 252.2 -> 224... | 5.287 | 42002 |
| Methylphenidate | 234.2 -> 84.2 | 10/16/2017 3:16 PM | P2-B5 | 5.392 | 173 | 2.59 | 234.2 -> 56.3 | | Meperidine-D4 | 252.2 -> 224... | 5.287 | 42002 |
| Metoprolol | 268.2 -> 116... | 10/16/2017 3:16 PM | P2-B5 | 5.101 | 19 | 0.54 | 268.2 -> 56.3 | 4.51 | Meperidine-D4 | 252.2 -> 224... | 5.287 | 42002 |
| Midazolam | 326.1 -> 223... | 10/16/2017 3:16 PM | P2-B5 | 6.375 | 40 | 0.75 | 326.1 -> 249... | | Midazolam-D4 | 330.1 -> 295... | 6.191 | 426706 |
| Mirtazapine | 266.2 -> 195... | 10/16/2017 3:16 PM | P2-B5 | 5.388 | 69 | 0.48 | 266.2 -> 72.3 | 0.49 | Midazolam-D4 | 330.1 -> 295... | 6.191 | 426706 |
| Mirtazapine | 399.2 -> 174... | 10/16/2017 3:16 PM | P2-B5 | 6.126 | 39 | 0.74 | 399.2 -> 238... | | Carbamazepine-13C6 | 243.1 -> 200... | 5.799 | 930164 |
| Normeperidine | 234.1 -> 42.3 | 10/16/2017 3:16 PM | P2-B5 | 5.229 | 109 | 0.50 | 234.1 -> 160... | 0.42 | Normeperidine-D4 | 238.2 -> 164... | 5.289 | 269319 |
| Normeperidine | 326.2 -> 44.3 | 10/16/2017 3:16 PM | P2-B5 | 5.765 | 75 | 0.41 | 326.2 -> 252... | | Normeperidine-D4 | 238.2 -> 164... | 5.289 | 269319 |
| Normpropoxyphene | 264.2 -> 91.2 | 10/16/2017 3:16 PM | P2-B5 | 6.041 | 158 | 1.35 | 264.2 -> 233... | 3.77 | Normpropoxyphene-D5 | 331.2 -> 44.3 | 5.846 | 79874 |
| Nortriptyline | 294.2 -> 184... | 10/16/2017 3:16 PM | P2-B5 | 5.525 | 483 | 0.74 | 294.2 -> 212... | | Nortriptyline-D3 | 267.2 -> 233... | 5.978 | 180997 |
| Ondansetron | 286.2 -> 218... | 10/16/2017 3:16 PM | P2-B5 | 6.068 | 713 | 190.95 | 286.2 -> 41.3 | | Pentazocine-13C3 | 289.2 -> 72.3 | 5.362 | 558166 |
| Pentazocine | 349.0 -> 206... | 10/16/2017 3:16 PM | P2-B5 | 5.595 | 20 | 0.22 | 349.0 -> 183... | 0.34 | Pentazocine-13C3 | 353.0 -> 183... | 6.108 | 6678 |
| Phenazepam | 244.2 -> 86.2 | 10/16/2017 3:16 PM | P2-B5 | 5.595 | 20 | 0.22 | 244.2 -> 91.1 | | Phenazepam-D4 | 249.2 -> 86.2 | 5.574 | 248762 |
| Phencyclidine | 253.1 -> 104... | 10/16/2017 3:16 PM | P2-B5 | 6.497 | 70 | 0.31 | 253.1 -> 182... | 0.07 | Phencyclidine-D5 | 263.2 -> 192... | 5.696 | 4202 |
| Phenylethanolamine | 325.1 -> 271... | 10/16/2017 3:16 PM | P2-B5 | 6.497 | 70 | 0.31 | 325.1 -> 140... | | Phenylethanolamine-D10 | 263.2 -> 192... | 5.696 | 4202 |
| Prazepam | 219.1 -> 91.2 | 10/16/2017 3:16 PM | P2-B5 | 5.557 | 235 | 0.79 | 219.1 -> 162... | | Prazepam-D5 | 330.1 -> 276... | 6.536 | 14469... |
| Primidone | 340.2 -> 58.3 | 10/16/2017 3:16 PM | P2-B5 | 5.805 | 85 | 0.51 | 340.2 -> 266... | | Prazepam-D5 | 330.1 -> 276... | 6.536 | 14469... |
| Propoxyphene | 264.2 -> 155... | 10/16/2017 3:16 PM | P2-B5 | | | | 264.2 -> 161... | | Propoxyphene-D11 | 351.3 -> 277... | 5.844 | 3143 |
| Protriptyline | 166.1 -> 148... | 10/16/2017 3:16 PM | P2-B5 | 4.606 | 154 | 0.31 | 166.1 -> 133... | 0.72 | Nortriptyline-D3 | 267.2 -> 233... | 5.978 | 180997 |
| Pseudoephedrine | 222.2 -> 107... | 10/16/2017 3:16 PM | P2-B5 | | | | 222.2 -> 121... | | Pseudoephedrine-D3 | 169.1 -> 151... | 4.404 | 701939 |
| Tapentadol | 295.2 -> 100... | 10/16/2017 3:16 PM | P2-B5 | | | | 295.2 -> 58.3 | | Tapentadol-D3 | 225.2 -> 107... | 5.149 | 301907 |
| Trimipramine | 455.3 -> 165... | 10/16/2017 3:16 PM | P2-B5 | 6.308 | 86 | 0.73 | 455.3 -> 150... | | Trimipramine-D3 | 298.2 -> 103... | 6.014 | 14784... |
| Verapamil | 306.1 -> 236... | 10/16/2017 3:16 PM | P2-B5 | 5.707 | 72 | 0.78 | 306.1 -> 264... | 1.29 | Trimipramine-D3 | 298.2 -> 103... | 6.014 | 14784... |
| Zaleplon | 389.1 -> 244... | 10/16/2017 3:16 PM | P2-B5 | 5.620 | 71094 | 568.00 | 389.1 -> 216... | 45.50 | Zaleplon-D4 | 310.2 -> 240... | 5.869 | 38502 |
| Zopiclone | | | P2-B5 | 5.620 | 71094 | 568.00 | 389.1 -> 216... | 45.50 | Zopiclone-D4 | 393.1 -> 244... | 5.620 | 16782 |