

ISP FORENSIC BIOLOGY CODIS METHODS MANUAL



Forensic Biology CODIS Methods Manual

Revision #12



APPROVED
December 7, 2011

Cynthia Cunnington Forensic Scientist IV Forensic Biology Supervisor/Technical Leader

Forensic Biology CODIS Methods Manual

REVISION RECORD

The following table must be filled out when revisions to the Biology Quality/Procedure Manual are made.

Date: The effective date of the revision(s).

Revision #: The manual revision number.

Description: A brief description of the changes made to the manual.

Addition: This column is checked if the revision reflects an addition (e.g. new SOP

or form) to the manual.

Deletion: This column is checked if the revision reflects a deletion (e.g. SOP or

form no longer in use) from the manual.

Initials: Initials of the Technical Leader making the revisions.

Date	Revision #	Description	Addition	Deletion	Initials
12/7/11	12	Separated CODIS and database methods into	X		CRC
	(See database manual for previous	two separate manuals. Renumbered methods.			
	revisions)	Cito Co Illa			
		9 110 60			
		0,0			
		7.0.			
		0, 0, 1, 4			
	X	, 70 1			
	<u> </u>				

Forensic Biology CODIS Methods Manual Table of Contents

i.	Signature/Approval page
ii.	Revision History
iii.	Table of Contents
1-5	CODIS Sample Data Entry and Upload
6-9	CODIS Database Hit Verification
10-14	CODIS Sample Removal
	A CONTROL OF THE POLITICE REIL A ON A STATE POLITICE REIL A STATE POLITIC



CODIS SAMPLE DATA ENTRY AND UPLOAD

1.0 BACKGROUND:

The implementation of the Combined DNA Index System (CODIS) in forensic DNA laboratories has provided an additional tool in assisting law enforcement agencies in solving or linking crimes that otherwise may not have resulted in the identification of a suspect. Accurate data entry for upload to NDIS (National DNA Index System) is essential.

It is the responsibility of the Forensic Scientist to generate profiles from convicted offender and/or forensic samples, determine CODIS eligibility, review other scientist's profiles, enter their profiles into CODIS and notify the CODIS Administrator of any potential problems associated with their CODIS DNA profiles. The CODIS Manager is responsible, in part, to serve as the laboratory's chief point of contact for all CODIS communications, maintain CODIS software updates and security, update/add/remove CODIS users, and upload profiles to NDIS.

DNA profiles entered into CODIS will be done so according to CODIS and NDIS Operational Procedures. Idaho currently allows for data entry into the following Indices: Offender (convicted offender); Unidentified Human Remains (deduced missing person, unidentified person); Missing Persons (alleged father/mother, biological child/father/mother, sibling, deduced missing person, missing person, maternal/paternal relative); Forensic Mixture; and Forensic (known and unknown). Profiles entered into the Forensic Index include both solved and unsolved cases in which the profile is associated with a crime and believed to be attributable to the putative perpetrator. Suspect reference samples and profiles matching the victim(s) and/or any elimination samples (e.g. consensual partner samples) may not be entered.

CODIS Methods: BI-500 Revision 12 Page 1 of 14

12/7/11 Issuing Authority: Quality Manager

2.0 SCOPE:

To provide a mechanism to ensure accurate data entry for all offender and forensic database samples

3.0 EQUIPMENT/REAGENTS:

CODIS Computer Workstation

Allele Tables and/or CMF files generated from sample analyses.

CODIS Training Manual

4.0 PROCEDURE:

4.1 SAMPLE DATA VERIFICATION.

- 4.1.1 Genotypic data (allele calls) are checked for accuracy and verified during the CODIS/Casework Review process and documented on the appropriate form (Form 306-BI and Form 214-BI, respectively).
- 4.1.2 Forensic and Missing Person/Unidentified Human profiles are verified for CODIS eligibility and Index and documented on the CODIS entry form (Form 218-BI).
- 4.1.3 Additionally, when 'STR Data Entry' is used to enter individual sample data (generally forensic samples) verification of allele calls and specimen category is performed prior to an SDIS search and NDIS upload.

4.2 SAMPLE DATA ENTRY IN CODIS:

4.2.1 Once sample data accuracy (forensic and offender) has been verified, data may be entered into the CODIS database either by use of the 'STR Data Entry' module to enter specimen ID/info and allele calls for individual samples or by using the 'Import' function for the entry of batches of data contained in a .cmf (or equivalent) file. Refer to the CODIS Training Manual and course

CODIS Methods: BI-500 Revision 12 Page 2 of 14 12/7/11

documentation for specific steps in accomplishing these tasks.

- 4.2.2 A second 'read' must also be entered for the individual samples through 'STR Data Entry'. This may be done by the analyst entering the sample data, or by another analyst logging on and entering the allelic data in the second read box. If a single individual enters both reads, the entire profile (i.e. all loci) for the first read must be entered prior to entering any data into the second read. A 'check' indicates agreement between readers at individual loci and discrepancies in entry must be rectified before saving to SDIS.
- 4.2.3 When using the 'Import' function, the scientist should validate the .cmf file prior to executing. This will ensure that any typos or inconsistencies (i.e., variant allele definitions/equivalencies) will be identified and corrected prior to importing the batch file data into SDIS.

4.3 CODIS DATA UPLOAD:

- 4.3.1 Prior to uploading eligible samples to NDIS, a search of SDIS is to be performed to identify any potential matches at the state level. This is typically done by running the 'Autosearcher' function in the software. See BI-303 for verification of any candidate matches identified during the autosearch.
- 4.3.2 There are various reasons that some samples present at SDIS should not be uploaded to NDIS. Prior to NDIS upload, these samples will be selected in Specimen Manager and 'unmarked for upload'. Generally speaking, an incremental upload will be performed. In Specimen Manager, 'incremental upload' is checked on the 'upload' pull-down menu and 'send upload' is selected. The upload is sent to NDIS as a message attachment via DNACOMM. If any 'candidate matches' are identified at NDIS, a match message will appear in DNACOMM and they will also be

CODIS Methods: BI-500 Page 3 of 14

reflected in Match Manager. For hit verification see BI-303.

5.0 COMMENTS:

- **5.1** Refer to CODIS Training Manual and course documentation for more specifics if necessary.
- 5.2 The CODIS software is redundant and there is generally more than one way to accomplish many tasks. Using a mechanism other than that listed here is acceptable.
- 5.3 The CODIS software is updated periodically and any necessary changes in procedure provided with new updates supercedes those in procedures written prior to update, if appropriate.
- 5.4 The sample history for convicted offender samples in DNA Tracker will be updated to reflect dates of sample analysis (in progress, complete, failed as appropriate), date sample was outsourced for testing, and/or date profile was entered into CODIS. This is accomplished for individual samples under the 'Samples' tab for the Offender by clicking on the appropriate sample and choosing 'Add History Event'. Multiple samples may be updated using the 'Batch Samples' function at the top of the screen. The current disposition box must be checked as appropriate for the history(s) added.
- 5.5 Eligible forensic profiles of unknown origin (i.e. profiles not matching any of the known reference samples) identified in the course of casework analysis must first be compared against the 'staff' index before entering the profile into the database. The 'searcher' function is used by hand entering the profile and choosing the 'staff' index. This profile will not be saved upon completion of the search and will need to be entered as outlined above for storage in the database.
- 5.6 'Searcher' is also used for performing one-time keyboard searches of forensic unknown profiles. A one-time search will only be performed under exigent circumstances and only on single source samples with a minimum of 8 loci present. The DNA profile must have been generated by or on behalf of a criminal justice

CODIS Methods: BI-500 Page 4 of 14

agency, by an accredited laboratory, in compliance with the FBI Director's Quality Assurance Standards (QAS), and must be maintained by a Federal, State, or Local NDIS participating laboratory. Additionally, requests for one-time searches of non-ISPFS profiles must be received from the CODIS Administrator of the laboratory that maintains the profile. When performing a one-time search, the requesting agency will be notified either verbally, or in writing in the event of a negative

cl.rmin, notific of a neg rification o. ring the one-ti

CODIS Methods: BI-500 Page 5 of 14

12/7/11 Issuing Authority: Quality Manager

Revision 12



CODIS DATABASE HIT VERIFICATION

1.0 BACKGROUND:

The implementation of the Combined DNA Index System (CODIS) in forensic DNA laboratories has provided an additional tool in assisting law enforcement agencies in solving or linking crimes that otherwise may not have resulted in the identification of a suspect. This is accomplished by the electronic storage and maintenance of DNA profiles at the local, state and national levels. Hits are obtained when a candidate match(es) is identified through a database search at any level. Hit verification involves evaluating the candidate match to determine if it is a true match and verification of CODIS offender sample data where necessary and possible.

NDIS CODIS Hit Disposition Reporting & Confirm an Interstate Candidate Match - Operational Procedures

2.0 SCOPE:

To provide a method of sample verification to be performed prior to law enforcement agency notification of a database hit to ensure reporting of only true, confirmed matches.

3.0 EQUIPMENT/REAGENTS:

CODIS

NDIS Procedure Manual; CODIS Training Manual/documentation DNA Tracker

AFIS

Offender Database Sample(s)
Equipment/Reagents for STR Analysis

4.0 PROCEDURE:

4.1 MATCH VERIFICATION (Forensic):

4.1.1 For 'hits'/matches involving an ISP Forensic Biology evidentiary sample (either case-to-case or case-to-offender and including one-time keyboard searches) the primary responsibility for

Revision 12 CODIS Methods: BI-501 12/7/11 Page 6 of 14 Issuing Authority: Quality Manager

match verification follow-up and disposition lies with the Idaho CODIS Administrator.

- 4.1.2 The CODIS Administrator or designee (typically the case analyst) will first evaluate the 'Candidate Match' in Match Manager to determine if there is a basis for exclusion and, therefore, disposition as 'No Match'. A 'No Match' disposition doesn't require further verification or, where applicable, notification of the other laboratory involved with the match.
- 4.1.3 If evaluation in Match Manager demonstrates that the candidate match consists of potential high stringency (e.g., exclusions attributable to different typing systems, Promega v. Applied Biosystems), or possibly moderate stringency in the event of a forensic mixture or degraded sample, the disposition is changed from 'Candidate Match' to 'Pending' until the verification process is complete. In general, for case-to-case matches, the verification will consist of communication between scientists regarding the data and case status, while case-to-offender matches typically necessitate sample verification at the 'offender lab'.
- 4.1.4 Once the status of the 'candidate match' has been
 resolved, the disposition is set accordingly
 (e.g., 'No Match', 'Offender Hit', 'Forensic
 Hit', 'Conviction Match', 'Investigative
 Information', etc.) and 'Investigations Aided'
 filled in as appropriate and as outlined in the
 NDIS 'CODIS Hit Disposition Reporting' procedure.
- 4.1.5 If verification results in a 'hit', a hit report is issued by the case analyst. A copy of the hit report, along with the CODIS match report, is filed in the CODIS file. The original hit report and a copy of the CODIS match report are placed in the associated case file. The appropriate law enforcement agency is notified of the 'hit'. If the law enforcement agency submits a sample from the identified offender, appropriate analysis and issuance of a supplemental report will be performed as in 4.2.5.

CODIS Methods: BI-501 Revision 12
Page 7 of 14

Revision 12
12/7/11
Issuing Authority: Quality Manager

4.2 MATCH VERIFICATION (Offender):

- 4.2.1 For 'hits'/matches involving an ISP Forensic Biology convicted offender sample the primary responsibility for match verification follow-up lies with the CODIS Administrator for the laboratory with the forensic (evidentiary) sample. However, the initial evaluation in Match Manager (see 4.1.2-4) and AFIS sample verification (see 4.2.2) will be initiated as soon as feasible once a verification request has been received from the forensic laboratory.
- 4.2.2 Once a potential match has been confirmed and a verification request received, the associated offender sample will be retrieved from the laboratory. If fingerprint verification was not already done, the sample card will be taken to the latent fingerprint section for an AFIS search of the thumbprint to verify identification of the offender.
- 4.2.3 Following AFIS verification of the thumbprint, re-analysis of the offender sample will be performed as appropriate (i.e., if duplicate analysis has already been performed either as a QC function or as the result of a duplicate sample, analysis will not be repeated) prior to agency notification. In situations where a thumbprint was not received with the DNA sample, or is of insufficient quality for verification, a notation will be made and re-analysis for confirmation may proceed. The forensic laboratory (or law enforcement agency for Idaho cases) will be notified that the Offender could not be verified through thumbprint confirmation.
- 4.2.4 Following sample verification (AFIS and reanalysis as appropriate) the forensic case laboratory, in the case of an interstate hit, or submitting law enforcement agency will be notified of the confirmed hit. Laboratory notification may be made verbally and relevant documentation will be provided to the forensic

CODIS Methods: BI-501 Revision 12
Page 8 of 14

Revision 12
12/7/11
Issuing Authority: Quality Manager

case laboratory as requested. In Idaho, initial notification as well as the request for a new DNA sample from the identified offender, may be made verbally. However, written notification and a formal request for a new DNA sample, in the form of a hit report, will be sent to the appropriate law enforcement agency.

- 4.2.5 For intrastate offender hits, where possible, a newly obtained DNA sample from the offender will be analyzed with all deliberate speed. The analysis of the forensic sample may also be repeated, though this is not required. A supplemental report will be issued delineating the match in the usual manner (i.e., same treatment as for matching suspect sample submitted with case evidence).
- **4.2.6** The sample history in DNA Tracker will be updated to reflect the date the hit was confirmed for the offender sample.

4.0 COMMENTS:

5.1 A good faith effort will be made to resolve matches within 30 business days from the date of the match or from the date of the verification request for interstate matches involving an Idaho offender.

CODIS Methods: BI-501 Revision 12
Page 9 of 14

Revision 12
12/7/11
Issuing Authority: Quality Manager

BI-502



CODIS SAMPLE REMOVAL

1.0 BACKGROUND:

Participation in the National DNA database, in accordance with the DNA Analysis Backlog Elimination Act of 2000, necessitates provisions for DNA profile expungement in the event that a qualifying offender's conviction is overturned. Additionally, the Idaho DNA Database Act of 1996 addresses court-granted expungement requests (I.C. §19-5513). Removal of DNA profile data and/or destruction of biological samples obtained from Convicted Offenders may be necessary as a result of conviction reversal or sample collection/submission errors.

Expungement is defined as the removal of DNA profile data from local (LDIS), state (SDIS) and national (NDIS) databases in response to a court order overturning the offender's conviction of a qualifying offense. Expungement will include the removal of identifying information from other laboratory documentation and destruction of the biological sample from which the offender database DNA profile was generated.

Administrative removal is defined as the destruction of a DNA sample and removal of any records relating to that sample. Examples for which administrative removal may be warranted include, but are not limited to, the collection of a sample from a non-qualifying offender, or the notification by the collection agency that removal is warranted. Generally, the determination that a DNA sample does not qualify for inclusion in the database occurs prior to entry of the profile into DNA Tracker and subsequently CODIS (see BI-301); however, there may be instances when the collection agency provides notification of an error after the profile has been generated. In these circumstances, the profile will be removed from the local, state, and national databases as part of the administrative removal.

NDIS Expunge a DNA Profile - Operational Procedures

Idaho Statutes: Title 19, Criminal Procedure, Chapter 55 "The Idaho DNA Database Act of 1996"

2.0 SCOPE:

To provide a protocol for CODIS sample expungement and administrative removal that protects an individual's rights to privacy and maintains the integrity of the Idaho DNA database program. All procedures will be documented on the CODIS Sample Removal Checklist (Form 310-BI).

3.0 EQUIPMENT/REAGENTS:

CODIS Workstation DNA Tracker Database

4.0 PROCEDURE:

4.1 EXPUNGEMENT

Forensic Servins

A pri Prior to removal of any DNA profile data, source identification, or biological sample destruction, the CODIS Administrator or designee will verify: 1) the request for expungement is accompanied by a certified court order that meets the requirements outlined in I.C. §19-5513, 2) that the offender has no other qualifying offense(s) and 3) the identification of the original DNA database sample where possible.

- 4.1.1 A search of DNA Tracker is performed to establish whether or not the specified sample has been received in the laboratory and if so, whether or not a DNA profile has been generated and/or entered into CODIS.
- **4.1.2** If the sample has been received in the laboratory, a criminal history check in ILETS will be performed to ensure that the offense for which the expungement is requested is the only qualifying offense. If the Offender has a separate qualifying offense, the sample will not be expunged and a copy of the criminal history check will be retained for documentation.
- **4.1.3** The associated offender sample folder will be retrieved from the laboratory and taken to the latent fingerprint section for an AFIS search of the thumbprint to verify identification of the offender, if not already done. The absence of a

Revision 12 CODIS Methods : BI-502 12/7/11 Page 11 of 14

Issuing Authority: Quality Manager

- thumbprint or a poor quality thumbprint does not preclude the sample expundement from proceeding.
- 4.1.4 The offender record will be removed from DNA Tracker and a sample removal report generated.
- **4.1.5** The DNA profile, if applicable, will be deleted from CODIS, followed by an upload to NDIS for removal at that level. A written notification of the expungement will be made to the NDIS Custodian and a request for verification of deletion at the national level.
- 4.1.6 The original DNA Collection Report received with the sample (or a copy of the sample folder when no report is received will be marked expunged, along with the scientist's initials and date of expungement.
- 4.1.7 The DNA sample will be destroyed and witnessed by a second scientist
- 4.1.8 Administrative paperwork documenting the event will be retained in the CODIS file. The paperwork may contain some administrative information about the sample/individual it was collected from and will include the following, as applicable: the court order/request for expungement, criminal history check and thumbprint verification, DNA Tracker sample removal report, CODIS deletion report(s), NDIS correspondence, CODIS sample removal checklist, and copies of any correspondence with the requesting party.
- .1.9 An expungement notification letter (or letter indicating why expungement did not occur), the DNA collection report or folder copy, and a copy of the CODIS sample removal checklist will be sent to the requesting party.

4.2 ADMINSTRATIVE REMOVAL

Prior to removal of any DNA profile data, source identification, or biological sample destruction, the CODIS Administrator or designee will verify: 1) the request for administrative removal is in writing from

Revision 12 CODIS Methods : BI-502 12/7/11 Issuing Authority: Quality Manager

the collecting agency and includes a description of the error resulting in the removal request and 2) that the offender has no other qualifying offense(s).

Removal resulting from an entry error on the part of ISP Forensic Biology staff does not require a written request; although, the nature of the error should be documented. Requests for administrative removal from a party other than the collecting agency or due to biology staff error will be referred to ISP legal staff for a determination of the appropriate action to be taken.

- 4.2.1 A search of DNA Tracker is performed to establish whether or not the specified sample has been received in the laboratory and if so, whether or not a DNA profile has been generated and/or entered into CODIS.
- 4.2.2 If the sample has been received in the laboratory, a criminal history check in ILETS will be performed to ensure that the offense for which the administrative removal is requested is the only qualifying offense. If the Offender has a separate qualifying offense, the sample will not be removed and a copy of the criminal history check will be retained for documentation.
- 4.2.3 The offender record will be removed from DNA Tracker and a sample removal report generated.
- 4.2.4 The DNA profile, if applicable, will be deleted from CODIS, followed by an upload to NDIS for removal at that level.
- 2.5 The original DNA Collection Report received with the sample (or a copy of the sample folder when no report is received) will be marked expunged, along with the scientist's initials and date of removal.
- 4.2.6 The DNA sample will be destroyed and witnessed by a second scientist.
- **4.2.7** Administrative paperwork documenting the event will be retained in the CODIS file. The paperwork may contain some administrative

Revision 12 12/7/11

CODIS Methods : BI-502 Page 13 of 14

information about the sample/individual it was collected from and will include the following, as applicable: the request for administrative removal, criminal history check, DNA Tracker sample removal report, CODIS deletion report(s), CODIS sample removal checklist, and copies of any correspondence with the collection agency.

4.2.8 A sample removal notification letter (or letter indicating why removal did not occur), the DNA collection report or folder copy, and a copy of the CODIS sample removal checklist will be sent to the requesting party, or submitting agency in the event of an internal entry error.

5.0 COMMENTS:

- 5.1 An Offender cannot be deleted from DNA Tracker if there are any samples and/or court orders associated with the offender.
- 5.2 A sample cannot be expunded or deleted from DNA Tracker if there are qualifying offenses associated with the Offender.
- 5.3 The 'Delete Sample' function will be used to remove a sample from DNA Tracker for an Administrative Removal.
- 5.4 The 'Expunse Sample' function will be used to remove a sample from DNA Tracker for an Expunsement.

Revision 12 CODIS Methods: BI-502 12/7/11 Page 14 of 14 Issuing Authority: Quality Manager