TRAINING MANUAL LATENT PRINT EXAMINER

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 1 of 109 Issu All printed copies are uncontrolled

Table of Contents

Table	e of Contents	2
Revis	sion History	4
1.0	Introduction	6
2.0	Roles and Responsibilities	10
3.0	Module 1: History and Background of Friction Ridge Identification	12
4.0	Module 2: Other Scientific Personal Identification Methods	15
5.0	Module 3: Safety Training	17
6.0	Module 4: Case Management and Reporting for Processing	18
7.0	Module 5: Digital Preservation of Latent Prints	21
8.0	Module 6: General Latent Print Processing	25
9.0	Module 6: General Latent Print Processing	27
10.0	Module 8: Processing Technique – Amido Black	29
11.0	Module 9: Processing Technique – 1, 8-Diazafluoren-9-One (DFO), 1, 2 – Indanedione, and 2 Indanedione TP	
12.0	Module 10: Processing Technique – Dye Stains – Rhodamine 6G and RAM	34
13.0	Module 11: Processing Technique – Gentian Violet/Crystal Violet	37
14.0	Module 12: Processing Technique – Iodine	39
15.0	Module 13: Processing Technique - Leuco Crystal Violet (LCV)	41
16.0	Module 14: Processing Technique - Ninhydrin and ThermaNin	43
17.0	Module 15: Processing Technique - Powder Development of Latent Prints	46
18.0	Module 16: Processing Technique – Physical Developer (PD)	48
19.0	Module 17: Processing Technique – RECOVER LFT	50
20.0	Module 18: Processing Technique – Small Particle Reagent (SPR)	52
21.0	Module 19: Processing Technique – Sticky Side Powder/Sticky Side Powder Equivalent	54
22.0	Module 20: Processing Technique – Sudan Black	56
23.0	Module 21: Processing Technique – Cyanoacrylate Ester (Super Glue®)	58
24.0	Module 22: Digital Imaging	60
25.0	Module 23: Biology and Physiology of Friction Ridge Skin	62
26.0	Module 24: Recording Inked Fingerprints, Palm Prints, and Footprints	64
27.0	Module 25: Friction Ridge Pattern Recognition and Interpretation	66
28.0	Module 26: Introduction to Latent Prints and the State of the Science	68

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 2 of 109

29.0	Module 27: Human Factors	70
30.0	Module 28: Analysis, Comparison, Evaluation, and Verification (ACE-V)	72
31.0	Module 29: Case Management and Reporting for Comparison and/or MBIS	76
32.0	Module 30: Court Procedures, Related Laws, Expert Testimony, Criminal and Civil Procedures Applicable to Latent Prints	78
33.0	Module 31: Multimodal Biometric Identification System (MBIS)	81
34.0	Module 32: DNA Database Fingerprint Comparison	84
Apper	ndix I – Reading Lists	86

Property of Idaho State Police Forensic Services

Revision History

Revision #	Description of Changes				
1	Ready for Qualtrax – no content changes				
2	Updated introduction to include requirements for DNA Database Card Comparisons; added Module 33: DNA Database Fingerprint Comparisons, added practical exercises for ThermaNin and 1,2, Indanedione TP and associated reading in appendix I				
3	Break out modules for Latent Print Field Service Response and ABIS; further define general grading policy and applicability to individual assignments; slight wording and grammatical changes throughout.				
4	Numbered practical exercises, modified introduction, removed Introduction to Crime Scenes unit, combined Taking Post Mortem Exemplars with unit on processing bodies for latent prints into new module - Advanced Latent Print Field Service Response, updated numbering, added readings to modules: 6, 14, 20, & 31, removed one reading from module 4, slight wording and grammatical changes throughout.				
5	Convert to pdf following automated conversion system error - no other changes were made				
6	Corrected info on Vucetich in Module 1, added written test for module 31, added/modified readings in modules 4, 28, & 29, slight wording and grammatical changes throughout.				
7	Minor wording changes throughout, added sign offs for exercises, updated background in module 7 updated objectives in modules 1, 5, 7, & 24 updated practical exercises for modules 1-27, & 29, updated readings for modules 1, 3, 5-21, 25, & 32. Removed KSI from ALS module.				
8	Minor wording changes throughout; changed ABIS to MBIS throughout, added column for written test grades; updated sections 6.4.2, 7.4.5, 16.3.1, 16.3.2, 18.4.3, 24.2.2, 32.4.2; updated title 11.0, 16.0, & 20.0; added sections 11.3.7, 11.3.8, 16.3.5-16.3.7, 31.4.13; removed Module 30.0 Advanced latent Print Field Service Response and associated readings and associated references; and updated readings in modules 5, 7, 9, 14, 16, 18, 20, 29 & 30.				

9	Document reformatted to fix issues with Qualtrax PDF conversion. No technical content was changed.
10	Minor wording changes and updated numbering throughout; updated title 1.0; updated section 6.4.2; added Module 17- RECOVER LFT; added section 27.4.2; and updated readings in modules 4, 7, 17, 22, & 28.
11	Minor wording changes and updated numbering throughout; updated sections 6.2.6, 7.4.3.2, 8.1, 15.4.4, 28.2.3, 34.2.2; added sections 6.4.3-6.4.5, 30.4.7, 30.4.8, and updated readings in modules 4-6, 26, 28, and 31.

Property of Idaho State Police Copy

1.0 Introduction

The purpose of this manual is to provide an in-house training program that will result in a competent and qualified expert Latent Print Examiner. This expert shall possess specialized knowledge, skills and training in the sub-disciplines of Latent Print Processing and Latent Print Comparison. In addition to establishing a minimum standard of professional competency, completion of this manual shall aid in maintaining quality and consistency among examiners within the section.

The training program, in its entirety, is designed for the Trainee who has little to no prior background or experience in the subject matter. The training program consists of two main segments: Latent Print Processing and Latent Print Comparison and two supplementary modules: Multimodal Biometric Identification System and DNA Database Card Comparison that may be used depending on work duties. Each segment is composed of a series of modules on specific topics. These modules consist of reading materials, observation and demonstration, and/or practical exercises. Each module has an associated test. Module tests shall evaluate the ability of the trainee to properly perform examinations and may be written, oral, hands-on or a combination thereof. They shall not be reviewed or verified prior to submission to the Trainer.

The modules outlined are the minimum requirements for completion of training. Additional exercises or readings may be assigned at the discretion of the Technical Lead, if necessary. The training may be abbreviated for trainees with prior experience and training or for those individuals who perform only limited duties. The background and experience of each individual will be assessed by the Technical Lead prior to the trainee beginning the training program. Training modules do not need to be completed in sequence. The order of completion may vary depending on the Trainee and/or operational needs.

All cases processed and examinations performed during training will be with the Trainee working as "the hands of the Trainer" as defined by the ISPFS Quality/Procedure Manual.

External training is used to supplement and/or meet certain portions of the training program. Trainees should attend workshops and/or training classes in the areas of latent print processing, latent print comparison, courtroom testimony, digital imaging, and photography. Attendance of outside training courses/workshops is subject to course availability and budget constraints. Requests for training shall be approved through the chain of command.

Progress is monitored by the Trainer, who reports to the Technical Lead and/or Supervisor. The Trainee must pass each written test with a minimum score of 80%. All tests are closed book unless otherwise noted. Many practical exercises require that the

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 6 of 109

Issuing Authority: Quality Manager

All printed copies are uncontrolled

Trainee search out or participate in a particular activity. These exercises are not graded and the Trainee or Trainer need only to document the date of occurrence. Other practical exercises will be graded "pass" or "fail" as noted. In order to receive a passing mark, the Trainee must demonstrate comprehension of the subject and demonstrate to the Trainer that they are able to complete the assignment with satisfactory results. If a practical exercise is assessed as "fail" the Trainee will be given additional training and/or additional exercises until competency is achieved. The Trainee must pass a final competency test and mock court in each of the sub-disciplines: Latent Print Processing and Latent Print Comparison. Competency tests and mock courts are also "pass" or "fail". Should the Trainee provide incorrect results or inaccurate testimony during these exercises additional training or testing will be necessary and mock courts may be repeated. Training is considered complete upon formal approval by the Quality Manager. This training program is estimated to last 18-24 months. The actual pace of instruction is dictated by agency resources and needs, as well as the Trainee's progress and demonstrated proficiency.

As de sionstra single police p

Modules for Latent Print Processing Sign Off Module 1: History and Background of Friction Ridge Identification Date Trainer Module 2: Other Scientific Personal Identification Methods Date Trainer Module 3: Safety Training Trainer Date Module 4: Case Management and Reporting for Processing Trainer Date Module 5: Digital Preservation of Latent Prints Trainer Date Module 6: General Latent Print Processing Trainer Date Module 7: Processing Technique - Alternate Light Sources Trainer Date Module 8: Processing Technique - Amido Black Trainer Date Module 9: Processing Techniques – 1, 8, Diazafluoren-9-One (DFO) and 1, 2, Indanedione Trainer Date Module 10: Processing Technique - Dye Stains - Rhodamine 6G and RAM Trainer Date Module 11: Processing Technique – Gentian Violet/Crystal Violet Date Trainer Module 12: Processing Technique - Iodine Trainer Date Module 13: Processing Technique - Leuco Crystal Violet (LCV Trainer Date Module 14: Processing Technique - Ninhydrin/Thermanin Trainer Date Module 15: Processing Technique – Powder Development of Latent Prints Trainer Date Module 16: Processing Technique - Physical Developer (PD) Trainer Date Module 17: Processing Technique - RECOVER LFT Trainer Date Module 18: Processing Technique – Small Particle Reagent (SPR) Trainer Date Module 19: Processing Technique – Sticky Side Powder Date Trainer Module 20: Processing Technique – Sudan Black Trainer Date Module 21: Processing Technique – Cyanoacrylate Ester (Super Glue®) Trainer Date Module 26: Introduction to Latent Prints and the State of the Science Trainer Date Module: 30: Court Procedures, Related Laws, Expert Testimony, Criminal and Civil Procedures Applicable to Latent Prints (reading & processing portions only) Trainer Date Latent Print Examiner Training Manual Revision 11 Issue Date: 12/15/2023

Page 8 of 109

Modules for Latent Print Comparison Sign Off		
Module 1: History and Background of Friction Ridge Identification		
	Trainer	Date
Module 2: Other Scientific Personal Identification Methods	Trainer	Date
Module 5: Digital Preservation of Latent Prints	Trainer	Date
. 104410 012-9.441 1 10001 1441011 01 2440110 1 11410	Trainer	Date
Module 22: Digital Imaging		
Maria and the state of the stat	Trainer	Date
Module 23: Biology and Physiology of Friction Ridge Skin	 Trainer	Date
Module 24: Recording Inked Fingerprints, Palm Prints and Footprints	Trailier	Daile S
	Trainer	Date
Module 25: Friction Ridge Pattern Recognition and Interpretation		
Madula 26. Introduction to Latent Drints and the State of the Science	Trainer	Date
Module 26: Introduction to Latent Prints and the State of the Science	Trainer	 Date
Module 27: Human Factors	S	
	Trainer	Date
Module 28: Analysis, Comparison, Evaluation, and Verification (ACE-V)		
Module 29: Case Management and Reporting for Comparison and/or MBIS	Trainer	Date
Trouble 27 date Numagement and Reporting for comparison analysis 1750	Trainer	Date
Module 30: Court Procedures, Related Laws, Expert Testimony, Criminal an		es
Applicable to Latent Prints (reading & comparison and/or MBIS portions or	nly)	
	Trainer	 Date
Module for Multimodal Biometric Identification System Sig	en Off	
Module 31: Multimodal Biometric Identification System – NOTE completion		Comparison is a
pre-requisite for MBIS.		
* /O.	Trainer	Date
Module for DNA Database Card Comparison Sign Off		
Module 32: DNA Database Card Comparison Module 32: DNA Database Fingerprint Comparison		
Arounie 32. Divi Buttabase i ingerprint comparison	Trainer	Date

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 9 of 109 Issu All printed copies are uncontrolled

2.0 Roles and Responsibilities

2.1 Supervisor

The Supervisor shall maintain an employee training file with all associated authorizations and shall evaluate mock court testimony.

2.2 Technical Lead

The Technical Lead shall assess any prior applicable training, review and/or modify the current training plan to reflect the trainee's prior training, assign the appropriate modules, and organize the training. The Technical Lead should regularly monitor the Trainee's progress and review their training record for completeness and accuracy, procure final competency tests, and schedule mock courts. The Technical Lead shall provide input regarding mock court performance to the Supervisor and/or other members of management. At the completion of Latent Print Processing and/or Latent Print Comparison training, the Technical Lead shall review all documentation regarding training to determine if the Trainee performed all required training and is competent to perform analysis. If the Trainee is competent to perform analysis, the Technical Lead shall forward all required documentation to the Quality Manager. The Technical Lead may designate an onsite Trainer.

2.3 Trainer

The Trainer shall provide a copy of the training plan to the Trainee with an anticipated timeline for completion. The Trainer is responsible for coordination of practical exercises, demonstrating techniques, reviewing assignments, providing feedback, and administration of module tests. The Trainer should monitor for comprehension and competency in theoretical knowledge and basic practical skills. The Trainer shall communicate progress, delays, or the need for supplemental activities to the Technical Lead and/or Supervisor. Deficiencies should be openly discussed among the Trainee, Trainer, Technical Lead and/or Supervisor in an attempt to rectify them.

2.4 Trainee

The Trainee shall maintain a record of training. This record shall include, but is not limited to: daily training received, observed events, activities performed by the Trainee, court testimony observed or performed, field cases observed or performed, completed assignments, and checklists. All steps in training shall be documented as they are completed. The record shall include a list of training samples that are utilized for hands-on processing exercises as well as the methods used to process them. With regards to comparison and/or MBIS training, the

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 10 of 109

record will include a list of cases utilized as practical comparison exercises and associated statistics (number of latent prints examined, number of comparisons performed, and number of identifications). The ILIMS training program may be utilized to record events or specific conclusions during training.

The Trainee should provide a weekly report to the Technical Lead and/or Trainer to include activities accomplished during the week (readings/exercises completed, casework observed, classes attended, etc.). They should keep the Technical Lead and/or their Trainer informed of any problems or questions that may arise.

At the completion of the Latent Print Processing or Latent Print Comparison segment, the Trainee will advance to supervised case work. Supervised case work will not commence until approval has been granted by the Quality Manager. At such time, a record of all cases, associated statistics, and the identity of the supervising examiner will be kept for all Latent Print Processing or Latent Print Comparison supervised case work.

The Trainee shall ensure that all training records for outside classes are ately refle forwarded to the Quality Manger for inclusion in his/her training file and shall ensure that their curriculum vitae accurately reflects successfully completed

3.0 Module 1: History and Background of Friction Ridge Identification

3.1 Background and Theory

Friction ridge identification has been relied upon for over 100 years to provide accurate identification. Fingerprints were originally used as signatures when signing business transactions and official government documents. In 1686, Professor Malpighi at the University of Bologna in Italy made observation of spirals, loops and ridges in fingerprints using the newly invented microscope. In 1858, Sir William Herschel was using fingerprints to "sign" documents. It was during this time that he noticed that no two prints were exactly alike and realized that they could be used for personal identification purposes. In the 1880's Henry Faulds was studying the permanency of friction ridge skin and was the first to publicly suggest that fingerprints could be used to identify criminals.

In 1888, Sir Francis Galton became the first person to provide evidence that no two fingerprints were exactly the same and that the prints remain the same throughout a person's lifetime. He calculated that the odds of finding two identical fingerprints were 1 in 64 billion. He went on to publish the first book on the subject titled "Finger Prints" in 1892, in which he detailed the first classification system for fingerprints. In his book, he identified three pattern types (loop, whorl, and arch).

In South America, Juan Vucetich developed his own system of classification by 1891 and published a book "Comparative Fingerprinting" (Dactiloscopia Comparada) in 1904. The first criminal fingerprint identification in a murder investigation came in 1892 by Police Inspector Alvarez, an Argentine police official trained by Vucetich.

In 1896, Sir Edward Richard Henry created a fingerprint classification system of his own in British India, which later spread to England. The Henry Classification system was used to establish a Fingerprint Bureau at Scotland Yard.

In 1902, New York was the first state in the United States to start implementing the new fingerprint technology. Within the next year, law enforcement agencies and military branches all over the United States started implementing their own identification departments.

Between 1911 and 1914, Edmund Locard established the first set of rules for fingerprint identification. Locard claimed that if there were 12 points of

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 12 of 109

agreement between prints with no disagreements, the identity was confirmed beyond doubt. This standard was formally adopted in many countries except for the United States who moved away from a standard based on counting points.

By the 1990's, Automated Fingerprint Identification Systems (AFIS) were being widely used. Currently, tens of thousands of individuals are added to repositories daily. These fingerprint collections provide the basis for criminal history records maintained by local, state, and federal law enforcement agencies.

The basic methodology for friction ridge identification has remained relatively unchanged. As other disciplines of forensic science continue to develop accurate statistics for their results, friction ridge identification seeks to quantify their own results. While still in its infancy, studies are beginning to surface based around this type of research.

- 3.2 Objectives, Principles, and Knowledge
 - 3.2.1 Understand the purpose of early methods of personal identification (Bertillon system, photography, scars, tattoos, sight recognition, marks, and mutilations).
 - 3.2.2 Knowledge of the earliest recorded awareness of fingerprints (cliff dwellers-Chinese) and be able to recall the earliest known uses of friction ridge impressions as a means of identification in China, Japan, and India.
 - 3.2.3 Knowledge of early anatomical observations (Grew, Malpighi, Purkinje, et. al.) and understand the biological significance of friction skin ridge patterns and their formation.
 - 3.2.4 Understand the scientific observations and use of fingerprints leading to modern friction ridge identification. Be able to recall the contributions of notable friction ridge pioneers to include: Locard, Herschel, Faulds, Galton, Vucetich, Henry, Holland, Cummins, and Ashbaugh.
 - 3.2.5 Knowledge of the historical events that led to the introduction and use of fingerprints in England (Belper Committee, Troup Committee) and in the United States (Thompson, Twain, DeForest, Ferrier, NY Prison System, Will/William West, establishment of the FBI Identification Division).
 - 3.2.6 Knowledge of the current criminal and civil applications of fingerprints, palm prints, and footprints and how these applications developed in the United States.
 - 3.2.7 Knowledge of the existence and development of various criminal and civil fingerprint files (FBI, U.S. military medical records, state and local fingerprint and palm print repositories).
- 3.3 Health and Safety Hazards
 - 3.3.1 N/A

3.4	Reading a	ınd Practical Exe	ercises			
	3.4.1	Complete Modul	e 1 reading list			
			Trainee	Trainer	Date	
	3.4.2	Practical Exercise following figures Pass/Fail.				
			Trainee	Trainer	Date	
	3.4.3	Practical Exercis	se II - visit <u>http:</u>	//onin.com to fa	amiliarize yours	self with this
		web site; with re	egards to this m ———— Trainee	nodule visit: <u>httr</u> Trainer	o://onin.com/fp	o/fphistory.html
	3.4.4	Practical Exercise of the names and Principles, and K reading that you trainer prior to it fun activity that tool, think Jeopa this activity will to it with each me	d historical eve Knowledge Sect feel are partic mplementing/ will allow you a rdy, Trivial Pur be carried thro	nts that are detation above as we ularly notable. It creating the gan and others to us rout, crossword oughout your tra	niled in the Object of the second sec	ectives, from your ea with your his should be a on as a study g etc. NOTE:
3.5	Written	Test – Module 1	Trainee Trainee	Trainer Trainer	Date	 Grade

4.0 Module 2: Other Scientific Personal Identification Methods

4.1 Background and Theory

Great strides have been made with regards to personal identification methods. In the late 1800's to early 1900's, agencies relied upon various methods of personal identification, including photography and anthropometry. The most common of these was the Bertillon method that utilized a person's physical measurements to prove identity. Those systems were replaced in the early 1900's by fingerprint identification. While fingerprint identification is still the most widely used system for personal identification, there are a number of other current personal identification methods of which a practitioner should be aware. These include DNA, odontology, handwriting and voice analysis, as well as various biometric techniques. Biometric verification is becoming increasingly popular in corporate and public security systems due to the rise in security breaches and transaction fraud. Biometrics use distinctive, measurable, physical, and behavioral characteristics to differentiate individuals. The physical characteristics used for biometric authentication include fingerprints, palm veins, facial recognition, DNA, palm print, hand geometry as well as iris or retina recognition. This information is often interpreted by a computer system that confirms identity.

- 4.2 Objectives, Principles, and Knowledge
 - 4.2.1 Awareness of personal identification methods other than friction ridge skin to include biometrics, iris recognition, face recognition, vascular pattern recognition, hand geometry question document analysis, voice analysis, odontology, and DNA.
 - 4.2.2 Awareness of the advantages/disadvantages of each.
- 4.3 Health and Safety Hazards
 - 4.3.1 N/A

4.4	Reading and	Practical	Exercises
-----	-------------	-----------	-----------

4.4.1	Complete Mod	dule 2 Reading	List		
		Trainee	Trainer	Date	
4.4.2	developed in l Objectives, Pr	Module 1. Incor inciples, and Kı	rporate each of nowledge Secti	game or other activity you the relevant terms located on above as well as any otl arly notable. Pass/Fail.	d in the
		Trainee	Trainer	 Date	

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 15 of 109

	Trainee	Trainer	Date	-
5 Written Test - Mod	lule 2			
	Trainee	Trainer	Date	Grade
. OPERTY OF 105				dices
			:0	3
		2,0	SUS,	
		dice of	S	
	, de	ej.		
	the luffer			
0/10				

5.0 Module 3: Safety Training

5.1 Background and Theory

Safety in the laboratory is an essential part of the job of a Forensic Scientist. The Occupational Safety & Health Administration (OSHA) was created in 1970 to protect workers. It mandates that each laboratory worker be knowledgeable about blood borne pathogens, chemical hygiene, universal precautions, biohazard disposal, decontamination, and vaccinations. It requires that all of the applicable information for the lab is given to the employee so that they may maintain safety in the workplace. It is also imperative that employees are able to access the Safety Data Sheets (SDS) in their laboratory in order to maintain safety around applicable chemicals.

5.2	Objectives.	Principles	, and Know	ledge
0.2		I I III CIPICO	, alla ittiovv.	10450

- 5.2.1 Understand safety hazards associated with the latent prints laboratory.
- 5.2.2 Knowledge of spill procedures/equipment and the use of personal protective equipment.
- 5.2.3 Knowledge of the potential explosion, fire, and contamination safety hazards associated with latent print development powders, solvents and chemicals.
- 5.2.4 Proper disposal of chemicals

5.3	Health	and Sa	fety Haz	ards
-----	--------	--------	----------	------

5.3.1 N/A

5.4	Dooding	and Practical	Evergices
5.4	Reading	anu Practica	i cxei@ises

J.4.1	Complete Modu	ne 5 Keaunig	List	
	7.0.	Trainee	Trainer	Date

Practical Exercise I - Trainer led session on section safety equipment (location of Safety Data Sheets, spill kits, eye washes, fire extinguishers); chemical storage and disposal; and forms and labeling requirements (bottle labels, hazard labels, reagent logs, equipment maintenance logs, control test logs, image deletion logs, etc.). The trainee shall demonstrate this knowledge by guiding the Discipline Lead or designee on a tour of the above listed items and showing them how/where to access the items. Pass/Fail.

		Trainee	Trainer	Date	
5.5	Written Test - Module 3				
		Trainee	Trainer	Date	Grade

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 17 of 109

6.0 Module 4: Case Management and Reporting for Processing

6.1 Background and Theory

In forensic science, it is imperative that procedures are accurately followed and documented appropriately. All documentation done for a case is subject to scrutiny by peers, the laboratory system, the courts, and accrediting bodies. Documentation should be as precise and error-free as possible.

It is important that measures are taken to prevent loss, deleterious change or tampering of evidence. Evidence should be tracked both internally (within the lab) and externally, as it transitions from agency to agency or person to person. This is done through chain of custody. When in the custody of an examiner, evidence integrity shall be ensured by properly securing, processing, marking, documenting, and re-sealing the evidence.

The system that is used to track information regarding a case is the Idaho Laboratory Information Management System (ILIMS). This system includes the internal chain of custody, information given to ISPFS by the submitting agency regarding the case, case correspondence, examiner generated notes and/or photographs, and all reports generated in relation to the evidence. The ILIMS system was implemented in 2013 to make all evidence processing paperless, efficient, and to afford timely access of records to submitting agencies and officers of the court. Comparison quality images are maintained in the Foray Digital Workplace database.

6.2 Objectives, Principles, and Knowledge

- 6.2.1 Knowledge of, and the ability to demonstrate, proper procedures for maintaining chain of custody (documentation and physical control).
- Knowledge of, and the ability to demonstrate, proper procedures for handling and marking physical evidence received for examination.
- 6.2.3 Ability to navigate and query ILIMS for latent print processing cases.
- Ability to demonstrate proper procedures for documentation of latent print processing casework. Documentation shall be such that another qualified Latent Print Examiner could evaluate what was done and why.
- 6.2.5 Understand how to prevent contamination.
- 6.2.6 Ability to demonstrate proper procedures for the collection and preservation of potential DNA from firearms and drink containers.
- 6.2.7 Knowledge of, and the ability to demonstrate, proper procedures for reporting latent print processing examination findings in an accurate, concise, and clear manner.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 18 of 109 Issuing Authority: Quality Manager

		results may be §	given to custon	ners.		
6.3	Health an	nd Safety Hazard N/A	ds			
6.4	Reading a	and Practical Ex	ercises			
	6.4.1	Complete Modu	ıle 4 Reading L	ist		
			Trainee	Trainer	Date	
	6.4.2				entation - shadow	
				=	least two process	
					out and writing la	
		processing repo	orts in ILIMS –	Trainer led di	scussion and dem	onstration.
			Trainee	Trainer	Date	 Case
			Transec	Tramer	: C)	dase
			Trainee	Trainer	Date	Case
			————— Trainee	Trainer	Date	 Case
	6.4.3	Practical Everci			what types of evic	
	0.4.5				nt examiners and	=
			=		iology section. T	
			-		ning/workstation	•
		· · · · · · · · · · · · · · · · · · ·	_ (sample packaging	
		creation in ILIM		eutheritation,	sample packaging	, sub item
			X	2		
			Trainee	Trainer	Date	
	6.4.4	Practical Exerci	se III – Trainee	shall demons	strate to trainer	
		cleaning/works	station set-up, _l	proper use of	PPE, swab collecti	ion, and
		documentation	for one firearn	n and one drir	ık container.	
		100				
			Trainee	Trainer	Date	
	6.4.5				ncy Exercise" Tra	inee will assess
	<i>X</i>	100 prints for c	omparison util	ity. Passing s	core is 90%.	
	0,					Grade
	616	Practical Exerci	Trainee	Trainer The Trainer	Date	Grade
0	0.4.0				modules, the trair	nga will process
X					servation by the tr	
		=			the items. The tra	
				-	The case examine	
					notes indicating a	
					rvision of the case	=
		report will be is		=		cammer. The
		report will be is	saca by the qu	aiiiicu taaiiiii	ner / transer.	
ont [Print Evamin	er Training Manus	al	R	evision 11	

Understand release of information policies, i.e. with whom, when, and how

Latent Print Examiner Training Manual

6.2.8

Issue Date: 12/15/2023 Page 19 of 109

		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
6.4.7		rcise VI – Traine se reports. This	=		ce three latent print
		Trainee	Trainer	Date	•
6.4.8	Practical Exe	rcise VII – Techr			essing cases -
		iscussion and/o			ajice
		Trainee	Trainer	Date	
6.4.9	Practical Exe	rcise VIII – Traiı	nee shall perfor	m administra	tive and technical
	review on at l	least ten proces	sing case repor	ts, authored b	y examiners other
	than their Tra	ainer(s). The Tra	ainer will be the	e reviewer of	record and
	ultimately res	sponsible for the	e review on the	se cases. Pass	/Fail.
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
	6/9.0	Trainee	Trainer	Date	Case
	,0	Trainee	Trainer	Date	Case
O'	3	Trainee	Trainer	Date	Case
06		Trainee	Trainer	Date	Case
Written	Test - Module	4			

Date

Grade

Trainer

Trainee

7.0 Module 5: Digital Preservation of Latent Prints

7.1 Background and Theory

Photography is widely used in Forensic Science. It dates back to the 1800s, when collections of photographs of criminals would hang in police stations for identification purposes. Today, we use digital photography for documentation of crime scenes, victim injuries and/or death, retrieval of evidence, and preservation of evidence. Digital cameras contain a sensor that records color and brightness values. These values are stored electronically and interpreted by computers. In general, the higher the resolution, the more information captured.

As with other evidence related to a case, evidentiary photographs should be properly captured, stored, and tracked to ensure their admissibility in court. Photography may be utilized at any point in the processing of evidence for latent prints, e.g. overall documentation of the evidence item, photographs of particular latent prints, to show orientation on an object, or final condition of an item. When photographing latent print evidence for comparison purposes, it is important to include both the impression and a scale. A variety of photographic techniques may be employed and will depend largely on the substrate as well as the particular development technique utilized on the item. Some of these techniques will require the use of an alternate light source (ALS) and specialized camera filters.

7.2 Objectives, Principles, and Knowledge

- 7.2.1 Understand the proper procedures for camera capture and digital scanning of latent and inked print images.
- 7.2.2 Familiarization with common digital photography terminology to include camera parts (body, lens, shutter diaphragm and shutter release) and function, file types (JPEG, RAW, TIF), compression, resolution, depth of field, bracketing, f-stop, shutter speed, aperture, exposure, etc.
- 7.2.3 Understand the different types of cameras and their suitability for latent print photography.
- 7.2.4 Understand the interplay between aperture and depth of field, aperture and shutter speed, and ISO. They shall know how to change these settings and why it may be applicable to do so.
- 7.2.5 Knowledge of and ability to apply special requirements for category 1 vs. category 2 images.
- 7.2.6 Understand the properties of light and how those properties relate to the use of filters and lighting techniques (oblique lighting, diffuse lighting, co-axial lighting, ALS lighting with appropriate filters, bounce lighting, etc.)

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 21 of 109 Issuing Authority: Quality Manager

- 7.2.7 Ability to photograph chemically treated and powder developed latent prints of various colors.
- 7.2.8 Ability to photograph three dimensional impressions (plastic prints).
- 7.2.9 Use and maintenance of cameras and other equipment.

7.3 Health and Safety Hazards

- 7.3.1 As with all electrical appliances, guard against electrical shock. This can be accomplished by ensuring that all connections are proper and that no loose, damaged, or frayed wires exist. Make sure the camera, scanner, and/or ALS is unplugged before attempting any maintenance and do not use outdoors if wet conditions exist.
- 7.3.2 The eyes are generally more vulnerable than the skin, and appropriate eye protection must be used to protect them. Permanent eye damage can occur from reflected, refracted, or direct illumination to the eye. Most of the light emitted by an ALS is not absorbed but is reflected and scattered off the surface being examined. Extreme care should be taken around highly reflective surfaces. Never look directly into the light or allow beams to bounce off the surface into your eyes or the eyes of another person in the vicinity. Filtered goggles or shields shall be utilized when using this equipment as they provide protection from potentially harmful rays and provide additional enhancement for viewing latent prints.
- 7.3.3 The nature and extent of all potential hazards are not yet known because indepth assessments have not been made on most of the high intensity light sources used in forensic identification work.

				· (2)	X O		
7.4	Readi	ing and l	Practical Exe	ercises			
	7.4.1	Cor	nplete Modu	le 5 Reading	List		
			0	Trainee	 Trainer	 Date	
	7.4.2	Pra	ctical Exercis	se I – continu	ue adding to the	game or other activi	ty you
		dev	veloped in Mo	dule 1. Inco	rporate each of	the relevant terms in	the
		Obj	ectives, Prin	ciples, and K	nowledge Secti	on above as well as a	ny others
		fro	m your readi	ng that you f	eel are particu	arly notable. Pass/Fa	ıil.
	~			Trainee	Trainer	Date	
	7.4.3	Pra	ctical Exercis	se II – Photog	graphy - Traine	e will need to familia	rize
0	(0)	the	mselves with	the camera	equipment (ca	meras, lenses, copy st	and) and
X		can	nera softwar	e utilized in 1	the laboratory.		
		7.4.3.1	Trainee wil	l be able to i	dentify basic ca	amera components (b	ody, lens,
			shutter dia	phragm and	shutter release	and demonstrate ho	w to change
			lenses, shu	tter speed ar	nd aperture bot	h manually and via co	mputer
			camera sof	tware. Pass/	Fail.		
				Trainee	Trainer	Date	

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 22 of 109 Issuing Authority: Quality Manager

	7.4.3.2	interplay speed, and known f-s	between apert d ISO. This sho	ould be done to tops)/shutter	of field, apert	cure and shutter es of photographs with
	7.4.3.3	formats (J lossless, f-	PEG, TIFF, & F	RAW), and be a	ible to define o	to the different file compression, lossy vs., bit vs. Byte, SLR,
	7.4.3.4	technique	s to include: o	blique lighting	g, diffuse lighti	various lighting ng, co-axial lighting, g, etc. Pass/Fail.
	7.4.3.5	photograp garnered	ohy." Practice	taking macro p . Present three		nt is macro g the information must be evidentiary in
7.4.4	acq utili prir	uisition dev ize these de nts, and pri	vices to include evices on train nts developed training samp	e flatbed scaniing samples to with a variety les will be eva	include pater of processing luated by the	on on digital ras. The Trainee will nt prints, plastic techniques. Images Trainer. Pass/Fail.
7.4.5	digi sett	tal imaging ings and lo	system to inc cations, etc. Tl	lude navigatio ne Trainee wil	n, features, ho	ed lesson on the low to upload, storage ling images from lice. Pass/Fail.
7.4.6	spe		ould attend a	weeklong basi	ic photography	y course or a more lable (attach copy of
			Trainee	Trainer	Date	_

Trainee	Trainer	Date	Grade

Property of Idaho State Police Forensic Services

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 24 of 109 Issuing Authority: Quality Manager

All printed copies are uncontrolled

8.0 Module 6: General Latent Print Processing

8.1 Background and Theory

Latent print visualization may be achieved using various visual, physical, or chemical processes, most of which have evolved during the past century. There are three types of friction ridge impressions: latent, patent, and plastic. Latent prints are generally hidden until a physical or chemical process makes them visible. Although latent means hidden, it has become synonymous will all types of crime scene and evidence impressions. A patent print is a visible print. Examples of patent prints may be those left in blood, paint, dust, etc. A plastic print is a three-dimensional print, for example, those left in clay, wax, melted plastic, or tacky paint.

Prior to any latent print processing, a thorough visual inspection of the evidence shall be conducted, using a strong light source.

Deciding what technique(s) to use to develop latent print evidence depends on several factors including: type of latent print residue, type of substrate, texture of substrate, condition of substrate (clean, dirty, sticky), known environmental conditions during or following latent print deposition, length of time since deposition, consequences of destructive processing methods, subsequent forensic examinations, and sequential ordering of reagents/development techniques.

8.2 Objectives, Principles, and Knowledge

- 8.2.1 Knowledge of the generally accepted techniques for the detection and visualization of friction ridge impressions.
- 8.2.2 Knowledge of latent print residue components targeted by different chemical development procedures.
- 8.2.3 Ability to assess the effectiveness and results of applied processing techniques.
- 8.2.4 Understand generally accepted preservation methods for friction ridge impressions.
- 8.2.5 Knowledge of surface and environmental factors affecting selection and sequencing of chemical development procedures.
- 8.2.6 Knowledge of the effects of various solvents on evidence surfaces (inks, plastics, varnishes, etc.).
- 8.2.7 Knowledge of equipment maintenance relative to chemical development of latent prints.

8.3 Health and Safety Hazards

8.3.1 N/A

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 25 of 109 Issu All printed copies are uncontrolled

8.4	Reading a	ınd Practical Exe	ercises						
	8.4.1	Complete Modu	le 6 Reading Lis	st					
			Trainee	Trainer	Date				
	8.4.2	The Trainee sho	ould attend a La	tent Fingerprin	t Processing/Cher	nical course			
		(36 hour minim	um - attach cer	tificate when co	mpleted).				
			Trainee	Trainer	 Date				
	8.4.3	Practical Exercis			ding the water cor	ntent of			
				_	sider how the info				
		regarding this to	opic has change	d. Prepare a tal	lk, power point, oi	r poster that			
		will allow you to	o present your f	indings to the s	ection. Pass/Fail	65			
						.0			
		_	Trainee	Trainer	Date	_			
	8.4.4			•	essing plan on ho	•			
		•	_		river rock, a dark				
	magazine page, a tree branch with bloody impressions. Preser								
	your trainer and explain your reasoning for why you chose that particu								
		sequence. Pass,	/Faii.	(0)					
			Trainee	Trainer	Date				
0.5	Maitten T	Soat Module (-8)					
8.5	written i	'est – Module 6		10 4					
			Trainee	Trainer	Date C	Grade			
8.6	Processin	g Competency T	Test-Trainee y	vill independe	ntly process a m	ock case. A			
0.0		Processing Competency Test-Trainee will independently process a mock case. A minimum of two item types will be processed using sequential processing. This							
	competency test will be entered into ILIMS, and as such, Trainee will need to complete all appropriate documentation and attachments and issue a report.								
	complete	an appropriate	uocumentatio	ii and attaciiii	erits and issue a	тероге.			
		19.0.	Trainee	Trainer	Date				
		4	11411160		Date				
8.7	Supervise	ed Cases – Comp	lete 20 Super	vised Processii	ng Cases. Traine	e shall record			
	all case n	umbers, associa	ted stats, and	the identity of	the supervising	examiner.			
	0								
	608		Trainee	Trainer	Date				

9.0 Module 7: Processing Technique – Alternate Light Sources

9.1 Background and Theory

Visible light consists of electromagnetic radiation of differing colors and wavelengths. Wavelengths at approximately 700 nm are viewed as red light while wavelengths approximate to 400nm are viewed as violet light. To visualize latent prints via fluorescence, a specific wavelength of radiation is absorbed by either an untreated latent print or one treated with a fluorescent chemical or powder and then re-emitted at a differing wavelength. The wavelengths chosen on the Alternate Light Source (ALS) may be determined by the inherent luminescent nature of the print, the specific chemical or powder utilized for processing, or the luminescent nature of the substrate. Evidence is viewed and photographed with various filters dependent upon the specific wavelength used.

9.2 Objectives, Principles, and Knowledge

- 9.2.1 Knowledge of luminescence, fluorescence, inherent luminescence, light wavelengths, band-pass filters, and light delivery systems as they relate to ALS detection of latent prints.
- 9.2.2 Knowledge of dye stain procedures used post-cyanoacrylate and the need for ALS processing.
- 9.2.3 Knowledge of 1, 8-Diazafluoren-9-One (DFO), 1, 2 Indanedione, and the need for ALS processing.
- 9.2.4 Knowledge of equipment maintenance relative to ALS detection of latent prints.

9.3 Health and Safety Hazards

- 9.3.1 As with other electrical appliances, guard against electrical shock. This can be accomplished by ensuring that all connections are proper and that no loose, damaged, or frayed wires exist. Make sure the ALS is unplugged before attempting any maintenance and do not use outdoors if wet conditions exist.
- 9.3.2 The eyes are generally more vulnerable than the skin, and appropriate eye protection must be used to protect them. Permanent eye damage can occur from reflected, refracted, or direct illumination to the eye. Most of the light emitted by an ALS is not absorbed but is reflected and scattered off the surface being examined. Extreme care should be taken around highly reflective surfaces. Never look directly into the light or allow beams to bounce off the surface into your eyes or the eyes of another person in the vicinity. Filtered goggles or shields shall be utilized when using this equipment as they provide protection from potentially harmful rays and provide additional enhancement for viewing latent prints.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

		sources used in	forensic ident	ification work.		
9.4	Reading a	and Practical Exe	ercises			
	9.4.1	Complete Modul		ist		
			Trainee	Trainer	Date	
	9.4.2	Practical Exercis				= =
		developed in Mo		•		
		Objectives, Prince from your reading	•	O		
		iroini your reaun	ng mat you ie	ei ai e pai ticulai	ly ilotable. Pas	SS/ Fall:
			Trainee	Trainer	Date	1110
	9.4.3	Practical Exercis			/	
		preservation of a followed by hand				
		•				rainer the process,
		0 1			•	ed in a processing
		sequence. Pass/	Fail.	COIL		
			Trainee	Trainer	Date	
0.5	XA7 - 44 7	5 M. J. 1. 7		iico V		
9.5	written i	Cest – Module 7	0	0,06,		
			Trainee	Trainer	Date	Grade
		of Idaho	X O.	0,		
			C			
		~0	1010			
		701				
		\$ 10.				
		0)				
		3				
	0					
al	,00					
0						
*	OPERT					

The nature and extent of all potential hazards are not yet known because indepth assessments have not been made on most of the high intensity light

9.3.3

10.0 Module 8: Processing Technique – Amido Black

10.1Background and Theory

Blood is composed of red blood cells, white blood cells and platelets, suspended in plasma. Red blood cells contain hemoglobin, a protein that carries oxygen from the respiratory organs to the remainder of the body. This protein is made up of four heme groups. There are two types of blood enhancement methods used in forensics: ones that react with the heme group to imply that blood is present and ones that react with proteins and their breakdown products. The methods that react with proteins are not specific to blood, but still tend to be sensitive methods due to the quantity of protein and protein breakdown products available in blood. Amido Black is a stain used in the latent print section to enhance the protein component of bloody prints. If blood is suspected, other presumptive blood testing techniques may need to be utilized.

10.2 Objectives, Principles, and Knowledge

- 10.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 10.2.2 Demonstrate proper chemical application and preservation of developed prints.
- Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

10.3 Health and Safety Hazards

- 10.3.1 Gloves, lab coats, goggles, and respirators (if there is a chance of the reagents becoming airborne) are worn when mixing or using Amido Black.
- Glacial acetic acid is corrosive and extremely irritating to the eyes and respiratory system. Avoid breathing the vapors and use in a fume hood, with a respirator, or with adequate ventilation. Glacial Acetic Acid will cause burns if it comes in contact with skin.
- Methanol is flammable. It needs to be handled carefully with gloves during the mixing and use of Amido Black. Methanol is toxic in quantities as small as 30 ml and should not be allowed to come in contact with the skin, eyes, or mouth. It is possible for methanol to be absorbed through the skin. If methanol comes into contact with the eyes or mouth, the area should be flushed with generous amounts of water and a doctor may be consulted. Inhalation of methanol vapors should be kept at a minimum and the solution should be used in a hood or well-ventilated area.
- In addition, trainees must be aware of the biological hazards associated with blood and other body fluids and take extra precautions to protect themselves.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 29 of 109

10.4.1	Complete Module	e 8 Reading	List		
		Trainee	Trainer	Date	
10.4.2	Practical Exercise				activity you
101112	developed in Mod		_	_	= =
	that you feel are		=	=	om your rough
		pui 010 uiui 1)	110 000 101 1 0.00	,	
		Trainee	Trainer	Date	
10.4.3	Practical Exercise	e II – locate	and read Safet	y Data Sheet – <i>A</i>	Amido Black an
	carrier solvents.				
					65
		Trainee	Trainer	Date	: (0)
10.4.4	Practical Exercise	e III – Train	er led lesson o	n the mixing of	Amido Black.
10.15		Trainee	Trainer	Date	
10.4.5	Practical Exercise				· -
	preservation of A		-		
	by the Trainee, u	tilizing traii	ning samples. 'I	'he trainee will	be able to expla
				· (/)	
	the trainer the pr	ocess, what	t it may be read	_	where it is gene
		ocess, what	t it may be read	_	where it is gene
	the trainer the pr	ocess, what essing sequ 	t it may be read ence. Pass/Fai	l. - —————	where it is gene
	the trainer the pr	ocess, what	t it may be read	_	where it is gene
10.5 Written	the trainer the pr	ocess, what essing sequ 	t it may be read ence. Pass/Fai	l. - —————	where it is gene
10.5 Written	the trainer the prutilized in a process	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. - —————	where it is gene
10.5 Written	the trainer the prutilized in a process	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. - —————	where it is gene
10.5 Written	the trainer the prutilized in a process	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
10.5 Written	the trainer the prutilized in a process	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	ocess, what essing sequ 	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	
10.5 Written	the trainer the prutilized in a process Test – Module 8	rocess, what essing sequ ——— Trainee	t it may be read ence. Pass/Fai Trainer	l. Date	

11.0 Module 9: Processing Technique – 1, 8-Diazafluoren-9-One (DFO), 1, 2 – Indanedione, and 1, 2 Indanedione TP

11.1 Background and Theory

1, 8-Diazafluoren-9-one (DFO) was originally prepared in 1950, but its reaction with amino acids was not explored until 1990, when it was first used as a fingerprint development reagent. It was observed that the application of DFO resulted in pink fingerprints that fluoresced. Fluorescence occurs when energy is supplied by an external source (in this case, an ALS) and is absorbed by a fluorescent chemical, creating an excited electronic state. In an effort to return to its ground state, the chemical emits energy that can be visualized as fluorescence. DFO fluoresces when illuminated between 450nm-570nm. The reagent is now widely used to develop friction ridge impressions composed of amino acids on porous surfaces.

The fingerprint developing qualities of 1, 2-Indanedione were first reviewed after a related compound, (6-methyl-thio-1, 2-indanedione) was found to produce fluorescent fingerprints. 1, 2-Indanedione was found to produce fingerprints similar to DFO. Prints treated with this chemical fluoresce when exposed to wavelengths of 450-570nm. As with DFO, 1, 2-Indanedione reacts with the amino acids present in fingerprints and is utilized on porous surfaces.

Special formulations of 1, 2-Indanedione have been created that allow for use on thermal papers. These formulations do not utilize an external heat source, decreasing the darkening of the substrate.

11.2 Objectives, Principles, and Knowledge

- Basic knowledge of the chemicals, the latent print matrices with which they react, potential safety hazards, and appropriate substrates for use.
- Demonstrate proper chemical application and preservation of developed prints.
- 11.2.3 Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

11.3 Health and Safety Hazards

DFO has not been fully investigated for potential health hazards, but is thought to be similar to ninhydrin, which may act as an irritant. Gloves, lab coats, and safety glasses should be worn when mixing and using DFO. The application of the DFO working solution should be performed in a fume hood, well-ventilated

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 31 of 109 Issu All printed copies are uncontrolled

- area, or while wearing an air-purifying respirator equipped with an organic vapor cartridge.
- Glacial acetic acid is *corrosive* and extremely irritating to the eyes and respiratory system. Avoid breathing the vapors and use in a fume hood or with adequate ventilation. Glacial acetic acid will cause burns if it comes in contact with skin.
- 11.3.3 Methanol needs to be handled carefully with gloves during mixing and use. Methanol is toxic in quantities as small as 30 ml and should not be allowed to come in contact with the skin, eyes, or mouth. It is possible for methanol to be absorbed through the skin. If methanol comes into contact with the eyes or mouth, the area should be flushed with generous amounts of water and a doctor may be consulted. Inhalation of methanol vapors should be kept at a minimum.
- 11.3.4 1,2 Indanedione may be harmful by: inhalation, ingestion and skin absorption.

 May cause skin and eye irritation.
- 11.3.5 Ethyl Acetate is hazardous by ingestion or inhalation and slightly hazardous in case of contact with skin or eyes. The substance is toxic to mucous membranes and the upper respiratory tract. Repeated or prolonged exposure to the substance can produce blood, kidneys, liver, or the central nervous system (CNS) damage.
- 11.3.6 HFE-7100 may be harmful if inhaled, swallowed or absorbed through skin. May cause skin, eye, and respiratory tract irritation. HFE-7100 is not considered a Hazardous chemical as defined by the OSHA Hazard Communication Standard, 29 CFR1910.1200.
- 11.3.7 Zinc chloride is hazardous. Avoid contact with skin and eyes. It is a known irritant, a permeator and is corrosive. It is classified as a possible human mutagen.
- Dichloromethane (Methylene Chloride) is hazardous. Avoid contact with skin and eyes. It is a known irritant, permeator and corrosive. Inflammation of the eye is characterized by redness, watering, and itching. It is classified as a possible human carcinogen.

11.4 Reading and Practical Exercises 11.4.1 Complete Module 9 Reading List Trainee Trainer Date 11.4.2 Practical Exercise I – continue adding to the game or other activity you developed in Module 1. Incorporate at least three topics from your reading that you feel are particularly notable. Pass/Fail. Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 32 of 109

11.4.5	Indanedione, 1,2		•		
		————— Trainee	Trainer	————— Date	
11.4.4	Practical Exercis				Э.
		Trainee	Trainer	Date	
11.4.5	Practical Exerci	se IV – Trainer	led lesson on th	e mixing of 1, 2	- Indanedione.
		Trainee	Trainer	Date	
11.4.6	Practical Exerci	se V – Trainer l	ed lesson on the	e mixing of 1, 2	Indanedione TP
	(Thermal Paper).			c es
		Trainee	Trainer	Date	
11.4.7	Practical Exerci	se VI – Trainer	led demonstrat	ion on the appli	ication and
	preservation of	DFO followed l	oy hands-on pro	cessing/preser	vation by the
	Trainee, utilizin	g training sam	ples. The traine	e will be able to	explain to the
	trainer the proc	ess, what it ma	y be reacting wi	ith, and where i	t is generally
	utilized in a pro	cessing sequen	ce. Pass/Fail.		
		Trainee	Trainer	Date	
11.4.8	Practical Exerci				lication and
	preservation of				
	processing/pres	- (•
	trainee will be a	-		=	=
	reacting with, an	nd where it is g	enerally utilized	d in a processin	g sequence.
	Pass/Fail.	~ ? ? ~ ~ (3~		
11.4.9	Practical Exerci	Trainee	Trainer	Date	dication and
11.4.5	preservation of			= =	
	followed by han	•	•	. , .	•
		-	o	-	iner the process
				-	d in a processing
X .	sequence. Pass	•	ind where it is g	enerally dulized	u iii a processiiig
	sequence. Tass,	ran.			
,0 ⁰ 0		Trainee	Trainer	Date	
Written T	est – Module 9				
		Trainee	Trainer	Date	Grade

12.0 Module 10: Processing Technique – Dye Stains – Rhodamine 6G and RAM

12.1 Background and Theory

Dye stains are chemicals that are used to help visualize or enhance latent prints developed with other methods. They do not develop prints on their own and are generally applied to non-porous surfaces after fuming with cyanoacrylate ester.

Rhodamine 6G is an extremely efficient and highly fluorescent dye stain. Rhodamine must be visualized using an alternate light source and fluoresces between 450nm and 525nm.

RAM is a dye stain consisting of Rhodamine 6G, Ardrox and MBD (7-(P-Methoxybenzlamino-4Notrobenz-2-Oxa-1, 3-Diazile). This combination allows the stain to fluoresce across a broad spectrum of wavelengths. Since it can be observed under various wavelengths, problematic backgrounds can be tuned out by using a wavelength that only fluoresces the friction ridge impression and not the background. As with rhodamine 6G, the print needs to have been previously developed by cyanoacrylate fuming before using the RAM stain.

12.2 Objectives, Principles, and Knowledge

- 12.2.1 Basic knowledge of the chemicals, the latent print matrices with which they react, potential safety hazards, and appropriate substrates for use.
- 12.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 12.2.3 Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

12.3 Health and Safety Hazards

- Rhodamine 6G, Ardrox P133D and MBD are classified as suspected animal carcinogens, but sufficient evidence of human carcinogenicity has not been established. Rhodamine 6G and RAM are thought to be relatively safe when exposure is at low levels. They should never be inhaled or allowed to get into the eyes or mouth, as they are irritants. If this should occur, the eyes or mouth should be flushed with a generous amount of water and a doctor may be consulted.
- 12.3.2 Methanol, isopropanol, and petroleum ether are highly *flammable*. All three chemicals need to be handled carefully with gloves during mixing and use of the stain. Methanol and isopropanol are toxic in quantities as small as 30 ml and should not be allowed to come in contact with the skin, eyes or mouth. It

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 34 of 109 Issuing Authority: Quality Manager

All printed copies are uncontrolled

is possible for methanol and isopropanol to be absorbed through the skin. If methanol, isopropanol or petroleum ether come into contact with the eyes or the mouth, the area should be flushed with generous amounts of water and a doctor may be consulted. Inhalation of vapors should be kept at a minimum and the stain should be used in a fume hood or a well-ventilated area.

- 12.3.3 Eye protection, a lab coat, and gloves should be worn. All mixing and application of chemicals should be done inside a ventilated laboratory fume hood. Excess reagent shall be collected and placed in the hazardous waste container located in the fume hood.
- 12.3.4 Acetonitrile may be fatal if swallowed, inhaled or absorbed through skin. It affects cardiovascular system, central nervous system, liver and kidneys and may cause irritation to skin, eyes and respiratory tract. It is also a flammable liquid and vapor.

12.4	Reading a	nd Practical Exe	ercises		CO)
	12.4.1	Complete Modul	le 10 Reading L	ist	. (1
					5
			Trainee	Trainer	Date
	12.4.2				me or other activity you
		-	-		ree topics from your reading
		that you feel are	particularly no	otable. Pass/Fai	il.
				<u> </u>	
			Trainee	Trainer	Date
	12.4.3				ata Sheet – Rhodamine 6G,
		Ardrox, MBD an	d carrier solver	its.	
			<u> x0</u>		
			Trainee	Trainer	Date
	12.4.4			led lesson on th	e mixing of Rhodamine 6G
		(methanol base)).		
		7.0.	Trainee	Trainer	Date
	12.4.5	X ·	se IV – Trainer l	led lesson on th	e mixing of Rhodamine 6G
		(water base).			
	/x				
	12.46	,	Trainee	Trainer	Date
	12.4.6	Practical Exercis	se V – Trainer le	ed lesson on the	e mixing of RAM.
0	Oz		Trainee	Trainer	————— Date
V	12.4.7	Practical Exercis			ion on the application and
*	12,				on processing/preservation by
		=	=	=	inee will be able to explain to
				=	g with, and where it is generally
		-		-	g with, and where it is generally
		utilized in a prod	cessing sequen	te. Pass/Fall.	
			Trainee	Trainer	Date
ent D	rint Evamine	er Training Manua	ı	Rovi	sion 11

Late

Issue Date: 12/15/2023

			
Trainee	Trainer	Date	Grade

Property of Idaho State Police Forensic Services

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 36 of 109 Issuing Author

13.0 Module 11: Processing Technique – Gentian Violet/Crystal **Violet**

13.1 Background and Theory

Gentian Violet or Crystal Violet is a biological stain used to dye epithelial cells and fatty components of latent print residues an intense purple color. This reagent is a toxic carcinogen and should only be used in small quantities. It can be used on the sticky side of tape (duct tape, clear plastic tape, packaging tape, black electrical tape) and items that are greasy or oily, to enhance prints.

13.2 Objectives, Principles, and Knowledge

- Basic knowledge of the chemical, the latent print matrices with which it reacts, 13.2.1 potential safety hazards, and appropriate substrates for use.
- 13.2.2 Demonstrate proper chemical application and preservation of developed prints.
- Demonstrate proper mixing, use of controls, documentation, storage, and 13.2.3 disposal.

13.3 Health and Safety Hazards

- Gentian Violet/Crystal Violet is a suspected human carcinogen. It is known to 13.3.1 affect the kidney, ureter, bladder, and thyroid of animals. It can be harmful if inhaled and is irritating to the eyes and skin.
- 13.3.2 Gentian Violet should not be used in large amounts.
- 13.3.3 A dust mask or respirator with dust filter should be used when working with the dry form. Gentian Violet should be prepared and used in a fume hood or well-ventilated area. The examiner should wear a lab coat, heavy-duty (nondisposable) gloves and safety glasses.

13.4 Reading	and Practical Exe	ercises			
13.4.1	Complete Modul	le 11 Reading	g List		
8		Trainee	Trainer	Date	
13.4.2	Practical Exercis	se I – continu	e adding to	the game or other acti	ivity you
	developed in Mo	dule 1. Inco	rporate at le	east three topics from	your reading
X	that you feel are	particularly	notable. Pa	ıss/Fail.	
		Trainee	Trainer	Date	
13.4.3	Practical Exercis	se II – locate	and read Sa	fety Data Sheet - Gent	ian Violet.
		Trainee	Trainer	Date	
Latent Print Examir	ner Training Manua		7 of 109	Revision 11 Issue Date: 12/15/20 Issuing Authority: Qu	
		i age 37	01 103	issuing Authority. Q	Junty Wallager

All printed copies are uncontrolled

13.4.4	Practical Exercis	se III – Trainer	· led lesson on th	ne mixing of G	entian Violet.
13.4.5	by the Trainee, ι	Gentian Violet utilizing traini rocess, what i	followed by haing samples. The timay be reacting	nds-on process trainee will b	olication and sing/preservation e able to explain to here it is generally
		Trainee	Trainer	Date	
13.5 Written	Γest – Module 11				S
		Trainee	Trainer	Date	Grade
Properti	3 of Idaho	Sinjer			

14.0 Module 12: Processing Technique – Iodine

14.1 Background and Theory

Iodine fuming is one of the oldest latent print methods still used today. It was advocated by Pierre Aubert in Paris in 1876. Iodine fuming exposes the evidentiary item to iodine fumes to develop latent prints. Iodine sublimates at low temperatures and the vapors are absorbed by the fats and oils in the latent print to turn it a yellow/brown color. Due to the sublimation of the iodine crystals, the print does not remain the yellow/brown color for very long. It is essential to photograph the print as quickly as possible after it is developed. It is considered a non-destructive technique.

14.2 Objectives, Principles, and Knowledge

- 14.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 14.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 14.2.3 Demonstrate proper use of controls, documentation, storage, and disposal.

14.3 Health and Safety Hazards

- 14.3.1 Safety is a serious concern when using the iodine fuming method. Iodine is toxic in any form. ALWAYS AVOID INHALING IODINE FUMES.
- Iodine fumes may irritate the skin and damage the respiratory tract.
 Headaches that can last for several days may result from exposure to iodine.
 Long-term effects to the thyroid gland may result from exposure.
- 14.3.3 Adequate ventilation when using the method is mandatory as the fumes are corrosive to metals and may discolor other surfaces that they come in contact with.
- 14.3.4 Iodine shall be purchased in glass ampoules. The ampoules shall stay sealed until use.

14.4 Reading and Practical Exercises 14.4.1 Complete Module 12 Reading List Trainee Trainer Date 14.4.2 Practical Exercise I – continue adding to the game or other activity you developed in Module 1. Incorporate at least three topics from your reading that you feel are particularly notable. Pass/Fail. Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 39 of 109

14.4.3	Practical Exerc	rise II – locate	and read Safety	y Data Sheet – I	Iodine.
	D 15	Trainee	Trainer	Date	
14.4.4	Trainee, utilizi trainer the pro	f Iodine folloving training sa	ved by hands-o mples. The trai	n processing/p nee will be abl g with, and whe	pplication and breservation by the eto explain to the ere it is generally
		Trainee	Trainer	Date	
	Test – Module 1				:1083
	Sofidaho	Trainee	Trainer	Date	Grade
			K0	0	
		Re	Olivos,		
	% C	Stater			
	70/19:01				
oper!	.7				

15.0 Module 13: Processing Technique – Leuco Crystal Violet (LCV)

15.1 Background and Theory

Leuco Crystal Violet (LCV) is a biological stain that reacts to the heme group in blood to cause the impression residues to turn an intense purple color. It should only be applied to thoroughly dried blood impressions. LCV gives an almost instantaneous visualization of latent prints in existing ambient light. Resulting prints should be photographed as soon as possible to avoid over development of the background.

15.2 Objectives, Principles, and Knowledge

- 15.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 15.2.2 Demonstrate proper chemical application and preservation of developed prints.
- Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

15.3 Health and Safety Hazards

- 15.3.1 Leuco Crystal Violet may be harmful by inhalation, ingestion or skin adsorption; may cause skin and eye irritation; may cause irritation to mucous membranes and upper respiratory tract.
- 15.3.2 Leuco Crystal Violet should not be used in large amounts.
- 15.3.3 A respirator should be used when working with the dry form. Leuco Crystal Violet should be prepared and used in a fume hood or well-ventilated area. The examiner should wear a lab coat, gloves and safety glasses.
- 15.3.4 In addition, examiners must be aware of the biological hazards associated with blood and other body fluids and take extra precautions to protect themselves.

15.4 Reading	and Practical Ex	ercises			
15.4.1	Complete Modu	le 13 Reading	List		
15.4.2		odule 1. Inco	rporate at l	the game or other activity your reast three topics from your rease/Fail.	
		Trainee	Trainer	Date	
15.4.3	Practical Exercise II – locate and read Safety Data Sheet – Leuco Crystal V				
	and carrier solv	ents.			
		Trainee	 Trainer	 Date	
Latent Print Exami	ner Training Manua		11411101	Revision 11	
				Issue Date: 12/15/2023	
		Page 41	of 109	Issuing Authority: Quality N	√lanager
	۰ ۱۱۸	arinted coniac	are unconti	collod	

All printed copies are uncontrolled

	ise III – Train	er led lesson or	n the mixing of L	euco Crystal Violet
	Trainee	Trainer	Date	
1 of mala B	Trainee	Trainer	Date	
preservation of processing/pretrainee will be reacting with, a	f Leuco Crystoneservation by able to explain	al Violet followe the Trainee, ut in to the trainer	ed by hands-on ilizing training s the process, wh	amples. The nat it may be
Pass/raii.	Trainee	Trainer	Date	5
Test – Module 1	3		•	dices
	Trainee	Trainer	Date Date	Grade
3 of loans	State	olice Kol	ensit	
	Formula A Formula B Practical Exerce preservation of processing/prestrainee will be reacting with, a Pass/Fail. Test – Module 1	Formula A Formula B Practical Exercise IV – Train preservation of Leuco Cryst processing/preservation by trainee will be able to explain reacting with, and where it is Pass/Fail. Trainee Test – Module 13 Trainee	Formula A Formula B Trainee Trainer Trainer Practical Exercise IV – Trainer led demonst preservation of Leuco Crystal Violet follower processing/preservation by the Trainee, uttrainee will be able to explain to the trainer reacting with, and where it is generally utility Pass/Fail. Trainee Trainee Trainer Trainer Trainer	Formula B Trainee Trainer Trainer Trainer Trainer Trainer Trainer Trainer Trainer Trainer Date Practical Exercise IV – Trainer led demonstration on the ap preservation of Leuco Crystal Violet followed by hands-on processing/preservation by the Trainee, utilizing training s trainee will be able to explain to the trainer the process, whereacting with, and where it is generally utilized in a process Pass/Fail. Trainee Trainer Trainer Date Trainer Date

16.0 Module 14: Processing Technique – Ninhydrin and ThermaNin

16.1Background and Theory

Ninhydrin (triketohydrindene hydrate) was first used in 1910 when Siegfried Ruhemann mistakenly prepared the compound. Ruhemann observed that the new compound reacted with amino acids to produce an intense purple color. Following Ruhemann's discovery, ninhydrin's use spread to analytical chemistry and biochemical applications. As early as 1916, the reaction with amino acids was used as an important test for the presence of protein in biological samples.

The technique is now one of the most popular methods for friction ridge detection on paper and other porous surfaces. The combination of heat and humidity accelerates the reaction of the proteins and amino acids with the ninhydrin.

Special formulations have been created that allow for use on thermal papers. These formulations do not utilize an external heat/humidity source, decreasing the darkening of the substrate.

16.2 Objectives, Principles, and Knowledge

- Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- Demonstrate proper chemical application and preservation of developed prints.
- Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

16.3 Health and Safety Hazards

- 16.3.1 Gloves, lab coat, and eye protection shall be worn when using or mixing ninhydrin or ThermaNin. Precautions should also be taken to avoid inhalation of the fumes.
- Hexane, is *extremely flammable* and the solution is to be used or mixed in a fume hood or in another well-ventilated area. Ensure that ninhydrin treated items are completely dry prior to exposing to the heat source.
- 16.3.3 Glacial acetic acid is *corrosive* and extremely irritating to the eyes and respiratory system. Avoid breathing the vapors and use in a fume hood or with adequate ventilation. Glacial acetic acid will cause burns if it comes in contact with skin.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 43 of 109 Issu

can be harmful if inhaled. Avoid breathing the vapors and use in a fume hood or with adequate ventilation. ThermaNin is combustible. It forms explosive mixtures with air on intense 16.3.5 heating in dry form. In event of a fire, ThermaNin will develop hazardous combustion gases or vapors. 16.3.6 Ethyl Acetate is hazardous if ingested or inhaled and slightly hazardous in case of contact with skin or eyes. The substance is toxic to mucous membranes and the upper respiratory tract. Repeated or prolonged exposure to the substance can damage the blood, kidneys, liver, or central nervous system (CNS). HFE-7100 may be harmful if inhaled, swallowed or absorbed through skin. 16.3.7 May cause skin, eye, and respiratory tract irritation. HFE-7100 is not considered a Hazardous chemical as defined by the OSHA Hazard Communication Standard, 29 CFR1910.1200. 16.4 Reading and Practical Exercises 16.4.1 Complete Module 14 Reading List Trainee Trainer 16.4.2 Practical Exercise I – continue adding to the game or other activity you developed in Module 1. Incorporate at least three topics from your reading that you feel are particularly notable. Pass/Fail. Trainer Practical Exercise II – locate and read Safety Data Sheet – Ninhydrin, 16.4.3 ThermaNin, and carrier solvents. Trainer Date Practical Exercise III - Trainer led lesson on the mixing of Ninhydrin stock and 16.4.4 working solutions. Trainee Trainer Date 16.4.5 Practical Exercise IV – Trainer led lesson on the mixing of ThermaNin. Trainee Trainer Date Practical Exercise V - Trainer led demonstration on the application and preservation of Ninhydrin followed by hands-on processing/preservation by the Trainee, utilizing training samples. The trainee will be able to explain to the trainer the process, what it may be reacting with, and where it is generally utilized in a processing sequence. Pass/Fail. Trainee Trainer Date 16.4.7 Practical Exercise VI - Trainer led demonstration on the application and preservation of ThermaNin developed latent prints followed by hands-on processing/preservation by the Trainee, utilizing training samples. The Latent Print Examiner Training Manual Revision 11 Issue Date: 12/15/2023 Page 44 of 109 Issuing Authority: Quality Manager

2-Propanol, also known as Isopropyl Alcohol, is *flammable*. It is an irritant and

16.3.4

All printed copies are uncontrolled

reacting with, and where it is generally utilized in a processing sequence.

Pass/Fail.

Trainee

Trainer

Date

16.5 Written Test - Module 14

Trainee

Trainer

Trainer

Date

Grade

trainee will be able to explain to the trainer the process, what it may be

Property of Idaho State Police Forensic Services

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 45 of 109 Issu All printed copies are uncontrolled

17.0 Module 15: Processing Technique – Powder Development of Latent Prints

17.1 Background and Theory

The development of latent prints using powder involves the application of fine particles that physically adhere to the aqueous or oily components in latent print residue. Powder is one of the most common methods of latent print development utilized on non-porous and some semi-porous surfaces. It is also one of the oldest dating back to 1891. At that time, available substances including charcoal, lead powder, soot, and cigar ashes, were used for latent print development.

Most commercial powders use two essential elements to provide adhesion to latent print residue: pigment and binder. The pigment in the powder provides effective visualization, giving contrast against the background surface. The binder provides for maximum and preferential adhesion to latent print residue. There are many different kinds of powders including, black powder, magnetic powder, white powder, fluorescent powder, and various colored powders. No powder is universally applicable to all types of evidence.

There are several different types and sizes of brushes that can be used when applying fingerprint powders. Types include fiberglass, feather and animal hairbrushes as well as magnetic wands. Certain types of brushes are used in conjunction with certain types of powders.

17.2 Objectives, Principles, and Knowledge

- 17.2.1 Understand the basic types of powders and brushes.
- 17.2.2 Knowledge of surfaces and environmental factors determining brush type, powder type, and color selection.
- 17.2.3 Understand the proper procedures for using different types of hair, fiberglass, and magnetic brushes.
- 17.2.4 Knowledge of equipment maintenance and safety procedures relative to powder development of latent prints.
- 17.2.5 Knowledge of lifting tape, gel lifters, hinge lifters, etc.

17.3 Health and Safety Hazards

- 17.3.1 Examiners are required to use the hoods or exhaust vents positioned at each workstation when performing powdering and lifting in the laboratory.
- 17.3.2 When fingerprint powders are to be used for an extended period of time, a dust mask or half face respirator with dust filters should be worn to minimize the inhalation of the powder particles.

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 46 of 109 Issu All printed copies are uncontrolled

		fingerprint powo	lers.			
17.4	Reading a	nd Practical Exe	rcises			
	17.4.1	Complete Modul	e 15 Reading Li	st		
					Data	
	17.4.2	Practical Exercis	Trainee e I – continue a	Trainer dding to the gar	Date ne or other activ	vitv vou
		developed in Mo				
		that you feel are	particularly no	table. Pass/Fail		
			Trainee	Trainer	 Date	9
	17.4.3	Practical Exercis		11011101		ssing to include
		standard, magne				0.8
	47.4.4	D 15	Trainee	Trainer	Date	
	17.4.4	Practical Exercis				=
		various tapes (cl and lifts (gel & h		s-MJ, casting me	eaiums (Mikrosi	1 & Accutrans),
		and mits (get & m	iligej.	(0)	•	
			Trainee	Trainer	Date	
	17.4.5	Practical Exercis	e IV – hands-on	powder and lif	ting exercises by	y the Trainee
		utilizing training	samples. The t	rainee will be al	ole to explain to	the trainer the
		process, what it	may be reacting	with, and when	re it is generally	utilized in a
		processing seque	ence. Pass/Fail.	()		
			- XO >			
			Trainee	Trainer	Date	
17.5	Written Te	est – Module 15				
		0				
			Trainee	Trainer	Date	Grade
		100				
		A.				
	121					
	OPERIN					
	Oz.					
V)						

Persons using fingerprint powders should monitor reactions (if any) to the

17.3.3

18.0 Module 16: Processing Technique – Physical Developer (PD)

18.1Background and Theory

Physical developer is a technique to detect friction ridge impressions on wet or dry porous items, including papers, tapes, and cardboard. The process involves an oxidation–reduction (redox) reaction whereby a solution of an iron salt reduces aqueous silver nitrate to finely divided metallic silver. The technique derives its name from the photographic developer used during film processing that undergoes a similar redox reaction. The physical developer develops the impressions as dark gray or black due to the adhesion of metallic silver particles on the fatty acid and lipid components of sweat residue. Prior to the introduction of physical developer in the 1970s, there was no reliable method for recovering prints from water-soaked documents.

18.2 Objectives, Principles, and Knowledge

- 18.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 18.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 18.2.3 Demonstrate proper mixing, documentation, storage, and disposal.

18.3 Health and Safety Hazards

- 18.3.1 Physical developer should only be used in a fume hood or well-ventilated area, as it is irritating to the respiratory tract.
- 18.3.2 Lab coats, gloves and safety glasses should be worn.
- 18.3.3 Standard laboratory protocol is followed for chemical handling.

18.4 Reading and Practical Exercises

reduing c	ma i racticai D.	ACT CISCS			
18.4.1	Complete Mod	ule 16 Readin	g List		
18.4.2	developed in N	Module 1. Inco	Ü	Date game or other act three topics from Fail.	5 5
18.4.3	Practical Exerc		Trainer and read Safety	Date y Data Sheet for p	hysical develope
		Trainee	Trainer	Date	

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 48 of 109

18.4.4	Practical Exercis Prewash.	se III – Trainei	eled lesson or	the mixing of PD and Maleic Acid
		Trainee	Trainer	Date
18.4.5	preservation of large	PD followed b g training sam ess, what it ma	y hands-on pr ples. The train ay be reacting	ration on the application and rocessing/preservation by the nee will be able to explain to the with, and where it is generally
		Trainee	Trainer	Date
18.5 Written	Test – Module 16			ices
		Trainee	Trainer	Date Grade
oroperti	3 of Idaho	State	Silce For	Date Grade

19.0 Module 17: Processing Technique - RECOVER LFT

19.1Background and Theory

RECOVER LFT is a technique that has been validated to detect friction ridge impressions on copper-based/brass items. The process involves converting tetrasulfur tetranitride (S_4N_4) to disulfur dinitride (S_2N_2) crystals by thermal cracking. The S_2N_2 crystals are then polymerized to form a solid state ($SN)_x$. It is believed that fingerprint residue reacts with the nitride during polymerization of S_2N_2 . RECOVER LFT developed friction ridges may present as a blue/black print on lighter background or as a light-colored print with blue/black background discoloration. While the interaction between S_2N_2 and the latent print matrix has yet to be fully explained, disulfur dinitride is an effective treatment on copper-based metal surfaces including those exposed to adverse environments including washing with water and detergents, bleach, and acetone.

19.2 Objectives, Principles, and Knowledge

- 19.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 19.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 19.2.3 Demonstrate proper documentation, storage, and disposal.

19.3 Health and Safety Hazards

- 19.3.1 RECOVER DEVELOP precursor contains copper bromide tetrathiatetrazocine. Heating (in air) may cause a fire. Keep RECOVER DEVELOP precursor out of extreme heat. The precursor may be harmful if contact with skin, eyes, or swallowing occurs.
- 19.3.2 Respiratory hazards may occur with inadequate or obstructed ventilation. The RECOVER instrument must be operated with in an active fume hood. Caution should be taken to avoid breathing fumes
- 19.3.3 Thermal hazards may occur with inadequate or obstructed ventilation. Caution should be taken to avoid contact with hot surfaces.
- 19.3.4 Lab coats, gloves and safety glasses should be worn.

19.4 F	keading ai	nd Practical Exe	rcises		
1	9.4.1	Complete Module	e 17 Reading Li	st	
			Trainee	Trainer	Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 50 of 109

Trainee Exercise II – locate	 Trainer	'Fail.	om your read
	Hainei	Date	
Trainee			RECOVER DE
Exercise III - Traindition of RECOVER Lition of RECOVER Litione, utilizing trainee, utilizing trainer the process, what a processing sequent	FT followed by l ning samples. T t it may be reac	hands-on proce he trainee will ting with, and v	essing/preserv be able to exp
Trainee	Trainer	Date	
dule 17		:(0	
Trainee	Trainer	Date	Grade
Trainee	olicopy		

20.0 Module 18: Processing Technique – Small Particle Reagent (SPR)

20.1 Background & Theory

Small particle reagent (SPR) is a technique used to develop latent friction ridge impressions on moist, non-porous surfaces. Two types of SPR are available: the conventional formula consisting of molybdenum (IV) disulfide and commercially available white SPR. This technique relies on the adherence of fine particles, within a suspension solution, to the fatty components of latent print residue. This is the same approach as fingerprint powder. This technique was originally discovered by J.R. Morris in 1981.

20.2 Objectives, Principles, and Knowledge

- 20.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 20.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 20.2.3 Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

20.3 Health and Safety Hazards

- There does not appear to be any health hazards associated with small particle reagent, but the process should be monitored to see if there are any allergies.
- 20.3.2 Lab coats, gloves and safety glasses should be worn.
- 20.3.3 Standard laboratory protocol is followed for chemical handling.

20.4 Reading and Practical Exercises

20.4.1	Complete M	lodule 18 Readir	ig List						
X	3	Trainee	Trainer	Date					
20.4.2	Practical Ex	ercise I – contin	ue adding to the	e game or other ac	tivity you				
~0	developed i	developed in Module 1. Incorporate at least three topics from your reading							
0,	that you fee	that you feel are particularly notable. Pass/Fail.							
		Trainee	Trainer	Date					
20.4.3	Practical Ex	Practical Exercise II - locate and read Safety Data Sheet - traditional and white							
	SPR.								
		Trainee	Trainer	Date					

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 52 of 109

20.4.4	Practical Exercis	se III – Trainer	led lesson on th	ne mixing of tra	ıditional SPR.
		Trainee	Trainer	Date	
20.4.5	by the Trainee,	traditional SPR utilizing trainin process, what it	followed by ha g samples. The may be reactin	nds-on process trainee will be	lication and sing/preservation able to explain to ere it is generally
		Trainee	Trainer	Date	
20.4.6	the Trainee, util	white SPR follo lizing training s process, what it	wed by hands-camples. The tra may be reactin	on processing/ inee will be ab	preservation by
		Trainee	Trainer	Date	
	Γest – Module 18		.0	nsi	
		Trainee	Trainer	Date	Grade
	of Idaho	Statern	il COPY		

21.0 Module 19: Processing Technique – Sticky Side Powder/Sticky Side Powder Equivalent

21.1 Background and Theory

Sticky-side powder is a liquid fingerprint powder method that develops latent prints on adhesive surfaces. Sticky-side powder detects epithelial cells and fatty/oily components of latent print residue left when handling adhesive surfaces. Sticky side powder can be used on almost any tape but works especially well on duct and electrical tape. Sticky side powder was developed in the mid-1990's when researchers at the National Identification Centre, Tokyo Metropolitan Police, were investigating methods for developing latent impressions on the adhesive side of tapes.

21.2 Objectives, Principles, and Knowledge

- 21.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- 21.2.2 Demonstrate proper chemical application and preservation of developed prints.
- 21.2.3 Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

21.3 Health and Safety Hazards

- 21.3.1 When using the powder in the dry form, precautions should be taken to prevent the powder from becoming airborne and possibly inhaled.
- 21.3.2 Lab coats, gloves, and safety glasses should be worn.
- 21.3.3 Standard laboratory protocol is followed for chemical handling.

21.4 Reading and Practical Exercises

21.4.1	Complete Mod	ule 19 Readin	ıg List		
21.4,2	developed in M	Iodule 1. Inco	•	Date e game or other act three topics fron	
21.4.3	Practical Exerc	Trainee rise II – locate Trainee	Trainer and read Safety Trainer	Date y Data Sheet – Stice Date	ky Side Powder.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 54 of 109

21.4.4	Practical Exerc	cise III – Train	ier led lesson or	n the mixing of	Sticky Side Powder.
21.4.5	preservation o processing/pro trainee will be	f Sticky Side I eservation by able to expla	Trainer Trainer Trainer Trainee, uting to the trainer Trainer Trainer	d by hands-on ilizing training the process, w	samples. The
21.5 Written	Test – Module 1				S
		Trainee	Trainer	Date	Grade
		× 210	Trainer		
Properti	3 of loans	Shire			

22.0 Module 20: Processing Technique – Sudan Black

22.1 Background and Theory

Sudan Black was originally used in laboratories for biological testing or chemical screening for fatty components. Sudan black was initially reported to detect the oily/fatty components of fingerprint residue by Misui, Katho, Shimada, and Wakasugi of the Criminal Science Laboratory in Nagoya-shi, Japan in 1980. It is a dye stain that produces a blue-black product and is used to develop latent friction ridge impressions on non-porous waxy substrates and surfaces contaminated with grease, dried beverages, and food residue. Sudan black will also enhance latent impressions developed by cyanoacrylate fuming.

22.2 Objectives, Principles, and Knowledge

- 22.2.1 Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- Demonstrate proper chemical application and preservation of developed prints.
- Demonstrate proper mixing, use of controls, documentation, storage, and disposal.

22.3 Health and Safety Hazards

- The Sudan Black working solution contains methanol. Methanol is toxic in quantities as small as 30 ml and should not be allowed to come in contact with the skin, eyes or mouth. It is possible for methanol to be absorbed through the skin. If methanol comes into contact with the eyes or mouth, the area should be flushed with generous amounts of water and a doctor may be consulted. Inhalation of methanol vapors should be kept at minimum.
- 22.3.2 Sudan Black should be used in a fume hood or well-ventilated area.
- 22.3.3 Lab coats, gloves and safety glasses should be worn.
- 22.3.4 Standard laboratory protocol is followed for chemical handling.

	, ,			
22.4 Reading	and Practical Exerc	cises		
22.4.1	Complete Module 2	20 Reading Li	st	
0	_			
	Tı	rainee	Trainer	Date
22.4.2			0	ne or other activity you ree topics from your reading
	that you feel are pa	articularly not	able. Pass/Fail	l.
		rainee	Trainer	————— Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 56 of 109

22.4.3	Practical Exercis carrier solvents.		d read Safety Da	ita Sheet – Suda	n Black and
		Trainee	Trainer	Date	
22.4.4	Practical Exercis	se III – Trainer	led lesson on the	e mixing of Suda	an Black.
		Trainee	Trainer	Date	
22.4.5	Practical Exercise preservation of State Trainee, utility the trainer the putilized in a production.	Sudan Black fol izing training sa rocess, what it	lowed by hands amples. The trai may be reacting	on processing, nee will be able	/preservation by to explain to
		Trainee	Trainer	Date	
22.5 Written 7	Sest – Module 20			Sell	
		Trainee	Trainer	Date	Grade

23.0 Module 21: Processing Technique – Cyanoacrylate Ester (Super Glue®)

23.1 Background and Theory

Cyanoacrylate ester (CAE), also known as "Super Glue®", is a technique used to develop latent friction ridge impressions on virtually all non-porous and some semi-porous surfaces, including glass, metal, coated papers, and all forms of plastics. This method is especially effective on rough or textured surfaces. CAE processing also prepares the evidence for the acceptance of powder and dyestains that may enable further visualization of the latent prints. Super Glue® was created in the 1950's by researchers who were trying to develop an acrylic polymer for the aircraft industry. In the late 1970's, researchers discovered its latent fingerprint development use, using the fumes of the glue. Shortly thereafter, the Bureau of Alcohol, Tobacco, and Firearms introduced this technique to North America and it quickly gained acceptance worldwide.

CAE fuming works by quickly bonding the CAE monomers to the latent print residues. The monomer on the fingerprint residue reacts with another CAE monomer in the vapor phase to form a dimer on the print. This reacts with another monomer to eventually form a polymer of CAE molecules. The overall development time is fast, especially when volatilization of the glue is accelerated (via heating or pretreatment).

23.2 Objectives, Principles, and Knowledge

- Basic knowledge of the chemical, the latent print matrices with which it reacts, potential safety hazards, and appropriate substrates for use.
- Demonstrate ability to properly utilize the CAE fuming chambers, wands, and vacuum chambers.
- 23.2.3 Demonstrate proper preservation of developed prints.
- 23.2.4 Demonstrate proper use of controls, documentation, storage, and disposal.

23.3 Health and Safety Hazards

- 23.3.1 CAE fuming should only be conducted in a filtered chamber or well-ventilated area. Precautions should be taken to avoid inhaling or allowing the vapors to contact the eyes, as the vapors can be irritating to the eyes, nose, and throat. Persons wearing contact lenses should not open CAE chambers without proper precautions. Non-vented goggles should be worn.
- 23.3.2 Precautions include properly sealed CAE chambers and evacuating the fumes from the chambers prior to removal of the questioned and test surfaces.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 58 of 109

skin so that it can be gently pulled apart. 23.4 Reading and Practical Exercises 23.4.1 Complete Module 21 Reading List Trainee Trainer Date 23.4.2 Practical Exercise I – continue adding to the game or other activity you developed in Module 1. Incorporate at least three topics from your reading that you feel are particularly notable. Pass/Fail. Trainee Trainer Date 23.4.3 Practical Exercise II – locate and read Safety Data Sheet – CAE Date Trainee 23.4.4 Practical Exercise III – Trainer led demonstration on the application of CAE using the fuming chambers followed by hands-on processing by the Trainee, utilizing training samples. Pass/Fail. Trainee Trainer Date 23.4.5 Practical Exercise IV – Trainer led demonstration on the application of CAE using the fuming wand followed by hands-on processing by the Trainee, utilizing training samples. Pass/Fail. Date Practical Exercise V - Trainer led demonstration on the application of CAE 23.4.6 using the vacuum chambers followed by hands-on processing by the Trainee, utilizing training samples. Pass/Fail. Trainer Practical Exercise VI – Demonstrate to the Trainer your ability to preserve 23.4.7 CAE developed prints. The trainee will be able to explain to the trainer the process, what it may be reacting with, and when it may be more or less advantageous to use the fuming chamber, fuming wand, or vacuum chamber. Pass/Fail. Trainee Trainer Date 23.5 Written Test - Module 21 Trainee Trainer Date Grade

Gloves should be worn to prevent the cyanoacrylate from contacting the skin. If liquid glue is allowed to contact the skin, adhesion may result. If the skin sticks together, immerse affected areas in warm water. This will loosen the

Latent Print Examiner Training Manual

23.3.3

Revision 11

Issue Date: 12/15/2023

24.0 Module 22: Digital Imaging

24.1 Background and Theory

Latent print images are frequently captured, processed and stored using digital devices. All of the techniques used in digital image processing have their roots in traditional photography and mathematics. The use of digital image processing can yield information not readily apparent in the original image and can assist in drawing a conclusion that might not have been reached otherwise. Image processing provides for higher image clarity and contrast.

24.2 Objectives, Principles, and Knowledge

- 24.2.1 Understand the capabilities and limitations of specific technologies that relate to digital imaging and storage of latent and inked prints.
- 24.2.2 Understand digital processing techniques using Adobe Photoshop to improve the visualization of latent print images.
- 24.2.3 Proficiency in the use of processing techniques to include, but not limited to: color reversal, position reversal, layers, contrast, image calibration/resolution, digital filters, and creation of enlargements.
- 24.2.4 Proficiency in the use of the current digital imaging system.

24.3 Health and Safety Hazards

As with other electrical appliances, guard against electrical shock. This can be accomplished by ensuring that all connections are proper and that no loose, damaged, or frayed wires exist.

24.4	Reading a	nd Practical Exercises
	24.4.1	Complete Module 22 Reading List

	19,0				
	6	Trainee	Trainer	Date	
24.4.2	The Trainee sh	ould attend a I	Digital Imagin	g course. (20 hour minim	ium -
X	attach copy of	certificate).			
~@,		Trainee	Trainer	Date	
24.4.3	Practical Exerc	cise I – Trainer	led lesson on	digital image processing	to include
	a demonstratio	on of commonly	y utilized tech	niques. The Trainee shall	l practice
•	processing tec	hniques on the	training imag	es. Processed images wil	l be
	evaluated by tl	he Trainer and	the Trainee w	ill discuss with the Train	er the
	reasons they c				
		Trainee	Trainer	Date	

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 60 of 109

24.4.4	developed in Mo Objectives, Princ	dule 1. Incorp	oorate at least t owledge Sectio	game or other activity you hree of the terms located in the n above as well as any others rly notable. Pass/Fail.
24.4.5		= =		 Date ll independently capture, digital imaging system, ten latent
	prints as assigne	ed by the Traii	ner. Pass/Fail.	
		Trainee	Trainer	Date
24.5 Written T	'est – Module 22			:(08)
		Trainee	Trainer	Date Grade
Property	ofloaho	Saleri	Sicoly Sicoly	Date Grade

25.0 Module 23: Biology and Physiology of Friction Ridge Skin

25.1 Background and Theory

A thorough understanding of the anatomy and physiology of friction ridge skin allows examiners to correctly analyze latent print impressions. Elements of biology and physiology explain why friction ridge skin is unique, why features of the skin persist, how the features of the skin age, how the skin responds to injury and why scars that form are unique. Understanding the pliability of friction ridge skin and how the skin reacts when it contacts a surface also provides valuable assistance during the examination of friction ridge impressions.

25.2	Obj	ectiv	es,	Pr	'inci	ples,	and	Kı	nov	vled	gf	e

- 25.2.1 Understand the biology and physiology of friction ridge skin.
- Understand the basic foundations of the science of friction ridge identification (persistence and discriminability).
- 25.2.3 Understand the basic anatomy and terminology of the hands and feet.
- Understand the general chemical composition of human perspiration as a means of understanding the composition of latent print residue.
- 25.2.5 Knowledge of genetic abnormalities of friction ridge skin (e.g. dysplasia, dissociated ridges).
- 25.2.6 Knowledge of alteration and mutilation of friction ridge skin.

25 2	Health	1	$C \cdot C \cdot$	TT	1
/5 4	Health	วทศ	Vatet	Hazar	กต
40.0	HCalui	anu	Jaicty	Hazar	us

25.3.1 N/A

25.4.1

25.4 Reading and Practical Exercises

	10,				
		Trainee	Trainer	Date	
25.4.2	Practical Exe	rcise I– continue	e adding to the	game or other acti	vity

Complete Module 23 Reading List

Practical Exercise I– continue adding to the game or other activity you developed in Module 1. Incorporate at least one question for each line of the Objectives, Principles, and Knowledge Section above as well as any others from your reading that you feel are particularly notable. Pass/Fail.

		Trainee	Trainer	Date	
25.4.3	Practical Exer	cise II – Find ar	nd read two ar	ticles (published	within the past 10
	years) on the	biology and phy	ysiology of fric	ction ridge skin. P	resent a synopsis
	of the papers	to the latent pri	int section. Pa	ss/Fail.	

Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 62 of 109

Trainee	Trainer	Date	Grade

Property of Idaho State Rolice Forensic Services

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 63 of 109

26.0 Module 24: Recording Inked Fingerprints, Palm Prints, and Footprints

26.1 Background and Theory

Recording inked fingerprints, palm prints and footprints is necessary for latent print examinations. These impressions can be made using a number of techniques, including traditional ink, Livescan, and powder/adhesive lift methods. Care and determination in recording the prints should always be exercised in order to obtain the best quality recordings for classification and/or comparison.

26.2 Objectives, Principles, and Knowledge

- Understand the various methods for recording known friction ridges for criminal history or personal identification including knowledge of chemical (inkless) systems, printer's ink, the black powder/adhesive lift (Handiprint®) method and electronic capture systems (Livescan).
- 26.2.2 Understand the quality of friction ridge detail produced by each method.
- 26.2.3 Understand the benefits associated with obtaining victim/elimination prints and complete friction ridge exemplars (major case prints).
- 26.2.4 Understand the proper method of completing fingerprint and palm print card information, sequence for recording fingers, and method of printing plain impressions.
- Demonstrate ability to properly use ink and brayer to record fingerprints, palm prints, and footprints (including equipment maintenance).
- Demonstrate ability to properly record complete friction ridge exemplars (major case prints).

26.3 Health and Safety Hazards

26.3.1 N/A

26.4 Reading and Practical Exercises

26.4.1 Complete Module 24 Reading List

26.4.2					
.08		Trainee	Trainer	Date	
26.4.2	Practical Exer	cise I – contini	ue adding to the	game or other	activity you
	1 1 1 1		1		

developed in Module 1. Incorporate at least three topics from your reading that you feel are particularly notable. Pass/Fail.

Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 64 of 109

26.4.3		east three ind	lividuals. Exem		Trainer followed by valuated by and
		Trainee	Trainer	—————— Date	
26.4.4	Practical Exerc Instruction by ' evaluated by a	Trainer follov	ved by hands-o	n application.	footprints) - Exemplars will be
				D. 1 -	
26.4.5	Practical Evero	Trainee ise IV - Black	Trainer Powder Adhes	Date ive Lift Method	d - Instruction by
20.1.3		ed by hands-o	n application. E		be evaluated by and
		Trainee	Trainer	Date	
26.4.6	Practical Exerc		an Terminal Fa	miliarity – Ove	erview led by
	Livescan termi	nal operator.		i,C	
		Trainee	Trainer	Date	
26.5 Written	Test – Module 2		Trainer	 Date	 Grade
Propert	3 of loans	Trainee	oli Copy	Bate	Grade

27.0 Module 25: Friction Ridge Pattern Recognition and Interpretation

27.1 Background and Theory

Friction ridge identification and classification has a long history rooted in scientific research and empirical observations. Various classification systems including Henry, Vucetich, and National Crime Information Center (NCIC), have been successfully used over the past 100 years. Today's classification systems rely mainly upon computers to digitize, categorize, recall, and identify matching tenprint cards. NCIC became operational in 1967.

While the use of computers has modernized fingerprint classification within the criminal justice system and forensic science, it is important that latent print examiners be able to recognize and articulate the various patterns and subpatterns, their use in analysis and comparison, as well as the history behind them.

27.2 Objectives, Principles, and Knowledge

- 27.2.1 Understand common terminology and definitions associated with friction ridge pattern recognition (arch, loop, and whorl).
- 27.2.2 Know frequency rates for each major fingerprint pattern type and which patterns are most likely to occur on which fingers.
- 27.2.3 Ability to differentiate between pattern types.
- 27.2.4 Awareness and understanding of the Henry Classification System to include: origin, FBI extensions, pattern interpretation, & parts of classification.
- 27.2.5 Awareness and understanding of other classification systems (NCIC Classification System, American System, and the Vucetich System)
- 27.2.6 Understand friction ridge formations as they relate to recognition, interpretation, and identification.

27.3 Health and Safety Hazards

27.3.1 N/A

27.4 Reading and Practical Exercises

Trainee	Trainer	Date

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 66 of 109 Issu All printed copies are uncontrolled

	on. Where w information f	ould you searc for these 5 prin	elieve these pat h first? Enter the ts into the follo .hugin.com/exa	he pattern, ridg wing website a	ge count, trac and evaluate
	Pass/Fail.			·	
27.4.3	Practical Eve	Trainee	Trainer nue adding to th	Date	ar activity vo
27.4.3			orporate at leas	_	= =
	that you feel	are particularly	y notable. Pass	/Fail.	•
		Trainee	Trainer	Date	S
27.4.4		_	erprint Classific	=	
	cards for both 80%.	h Primary Heni	ry and individu	al pattern type	s. Passing so
	00,0	Trainee	Trainer	Date	Grade
	Sofildan	State	olice K		
oroperti	A Of I gray,				

28.0 Module 26: Introduction to Latent Prints and the State of the Science

28.1 Background and Theory

Forensic scientists are entrusted with a great amount of responsibility. The public and the criminal justice system expect that forensic scientists be unbiased, intelligent, and thorough. In order to do so, scientists must take their responsibility seriously and uphold the ethics and values required for their position. Over the past decade, the news has been filled with stories of incompetence and out right misconduct. Crime labs in nearly every state have been affected and, in turn, the field of forensic science is facing more and more challenges. We are seeing them on multiple fronts from both the court system, in the form of Daubert hearings, to legislation requiring accreditation. Many resources are being put into exploring the state of the science and what the path forward should look like. From the 2009 NAS report on Strengthening Forensic Science in the United States to the formation of the Organization of Scientific Area Committees (OSACs), the field is rapidly changing.

- 28.2 Objectives, Principles, and Knowledge
 - 28.2.1 Knowledge of the professional duties, moral obligations, and code of ethics for Latent Print Examiners.
 - 28.2.2 Knowledge of the various professional organizations and certifications.
 - 28.2.3 Be familiar with the NAS and PCAST reports and the impact they are having on the field.
 - 28.2.4 Be familiar with the Friction Ridge OSAC and its activities.
- 28.3 Health and Safety Hazards

28.3.1 N/A

28.4	Read	ing and	l Practical	Exercises
------	------	---------	-------------	-----------

28.4.1 Complete Module 26 Reading List

Trainee Trainer Date

Practical Exercise I– continue adding to the game or other activity you developed in Module 1. Incorporate at least one question for each line of the Objectives, Principles, and Knowledge Section above as well as any others from your reading that you feel are particularly notable. Pass/Fail.

Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 68 of 109

28.4.3			atches exercise until all answer		e is 100% - exerci
		Trainee	Trainer	Date	Grade
28.4.4			te and read the tent print exam		and Standards of hed by the IAI.
	1101000101101				
28.4.5	Practical Eve	Trainee	Trainer application to	Date	DNWD-IAI
20.1.5	Tractical Exc				TIWD III.
		Trainee	Trainer	Date	
28.4.6			ttps://www.nis		
			forensic-science		
		a nve minute p hem. Pass/Fail.	resentation to t	ne iatent print	section on a topi
	refevant to the	ileili. 1 assy 1 all.		S	5
		Trainee	Trainer	Date	
20 E Writton	Tost Module	. 26		SI	
26.5 Willen	Test – Module			0	
		Trainee	Trainer	Date	Grade
Propert	S of loan	ostate P	Trainer		
Brox					

29.0 Module 27: Human Factors

29.1Background and Theory

The term "human factors" as it applies to forensic science, is the scientific discipline concerned with the understanding of interactions among humans and other elements of the forensic system including products, decisions, procedures, workspaces, and the overall environment encountered at work. It advances an understanding of the nature of errors in complex work settings and attempts to mitigate them by applying theory, principles, data, and method design to optimize overall performance and improve cognitive abilities with respect to judgment and decision making. Human factors research has its roots in post-World War I aviation psychology and was first applied to forensic science, and latent print examination in particular, in the mid 2000's. By 2008, the National Institute of Justice (NIJ) Office of Investigative and Forensic Sciences (OIFS) and the National Institute of Standards and Technology's (NIST's) Law Enforcement Standards Office (OLES) had put together an Expert Working Group on Human Factors in Latent Print Analysis. The Organization of Scientific Area Committees (OSAC) currently has a Human Factors Committee established to provide advice and guidance on human factors issues in forensics.

- 29.2 Objectives, Principles, and Knowledge
 - 29.2.1 Develop an understanding of the nature of errors in latent print examination.
 - 29.2.2 Identify the various human factors that lead to errors.
 - 29.2.3 Understand the role of human factors and their contributions to errors in latent print analysis.
 - 29.2.4 Understand how environmental conditions affect the quality of latent print examinations.
 - 29.2.5 Ability to define the different types of bias: cognitive bias, confirmation bias, contextual bias, etc.
 - 29.2.6 Ability to define the different types of errors: false positive, false negative, etc.

	29.3	Health	and	Safety	Hazards
--	------	--------	-----	--------	---------

29.3.1 N/A

29.4	Reading	and	Practical	Exercises
------	---------	-----	-----------	-----------

29.4.1 Complete Module 27 Reading List

Trainee Trainer Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

developed in Module 1. Incorporate at least one question for each line of the Objectives, Principles, and Knowledge Section above as well as any others from your reading that you feel are particularly notable. Pass/Fail. Trainee Trainer Date 29.5 Written Test - Module 27 Grade

Grade

Grade

Forensic Services

Property of Idaho State Police P Trainee Trainer Date Grade

Practical Exercise I- continue adding to the game or other activity you

29.4.2

30.0 Module 28: Analysis, Comparison, Evaluation, and Verification (ACE-V)

30.1 Background and Theory

The scientific method is a method of research in which a problem is identified, relevant data is gathered, and a hypothesis is formulated from the data and then tested. In forensic science, it is imperative to have a scientific technique for examination. Doing so ensures that evidence is treated equally, and conclusions are reliable and unbiased. The latent print section utilizes ACE-V as part of the examination methodology. ACE-V is an acronym that stands for analysis (A), comparison (C), evaluation (E) and verification (V). It is the process that latent print examiners utilize to reach a conclusion about a comparison examination.

Huber initially discussed ACE-V in 1959 when describing the philosophy of science and the correct use of the scientific method. Huber described hypothesis testing as analyzing, comparing, and evaluating and noted that verification was needed. In 1979, David Ashbaugh noted the applicability of the methodology to the latent print comparison process. In 1998, during the first Daubert hearing on fingerprint evidence, the members of the fingerprint community recognized the need to better articulate how they came to their conclusions. ACE-V was determined to be one such way to do so. Today, ACE-V has gained widespread recognition and implementation within the field.

30.2 Objectives, Principles, and Knowledge

- 30.2.1 Understand the scientific methodology and its application to friction ridge examination.
- 30.2.2 Understand the individual friction ridge structure (e.g., continuity, texture, pore, and edge definition) for determining the existence of individualizing details.
- Understand friction ridge characteristics (dots, ridge endings, and bifurcations), the varying definitions/interpretations assigned to combinations of those three ridge characteristics, and how they may be utilized in effecting identification.
- Understand the value of incipient ridge characteristics for use in latent print comparison/identification.
- 30.2.5 Understand the importance of elimination prints and the utility of completing comparisons of known individuals (e.g. victims) before searching a print in the MBIS system.
- Ability to recognize and utilize ridge flow configurations (size, pattern, focal points, etc.), scars, creases, and other friction ridge characteristics.

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 72 of 109

30.2.7	Ability to recognize, and if possible, determine the area from which the latent
30.2.8	fingerprints, palm prints, and foot/toe prints originated. Understand the nature of color reversals (entire print) and changes (within
30.2.0	the same print) and the ability to properly analyze these occurrences when
	they are encountered in latent print comparisons.
30.2.9	Understand the effects of pressure distortion, slippage, overlays, pre- and
	post- deposit artifacts (surface scratches, soil, brush strokes, etc.), and the
30.2.10	ability to properly analyze such disturbances/distortion. Understand the different policies and standards that exist regarding what
30.2.10	constitutes friction ridge identification in the U.S. and other countries and why
	no minimum number of ridge characteristics can be defined to effect an
	identification (i.e., positive opinion based on personal empirical experience in
	examining and comparing latent prints).
30.2.11	Knowledge of simultaneous or adjacent impressions and their value for identification.
30.2.12	Ability to analyze fragmentized friction ridge detail to determine its value
20.040	(comparison/identification, value/no value).
30.2.13	Knowledge of various methods used to record known friction ridge impressions and the ability to properly evaluate ridge structure based on each
	method.
30.2.14	Ability to properly conduct a comparison.
30.2.15	Ability to render an accurate source conclusion (identification, inconclusive,
	exclusion).
30.2.16	Understand the necessity for verification by another qualified latent print
30.2.17	examiner.
30.2.17	Understand the role of quality assurance measures in friction ridge examination.
30.2.18	Awareness of the impacts resulting from an erroneous conclusion.
30.2.19	Awareness of basic statistical models and the potential for their integration
	into the current friction ridge identification procedures.
30.3 Health ar	nd Safety Hazards
30.3.1	N/A
30.4 Reading	and Practical Exercises
30.4.1	Complete Module 28 Reading List
X	
30.4.2	Trainee Trainer Date The Trainee should attend an approved Latent Print Comparison Techniques
50.4.2	training course (36 hour minimum - attach certificate when completed).
	Trainee Trainer Date

Revision 11

Issue Date: 12/15/2023

30.4.3				_	logy or Complex
	Comparison co	ourse. (36 hour	minimum - at	tach certificate	e when completed).
		Trainee	Trainer	Date	
30.4.4	The Trainee sh	nould attend an	approved Pali	m Print trainir	ng course. (20 hour
	minimum - att	ach certificate	when complete	ed).	
		Trainee	 Trainer	 Date	
30.4.5	Practical Ever	cise I – continu			r activity you
30.4.3			=	_	from your reading
	=	re particularly	=	-	on our reading
	,				S
		Trainee	Trainer	Date	.0
30.4.6					O prints as to value
	=	dge detail "IRD			
	• ,	nd finger patte	rn/area of orig	gin. Passing so	ore is 90% for these
	two columns.			·.C	
		_			l of clarity (1, 2, 3),
	complexity, sh	ape clues, red f	flags/distortion	n, and orientat	ion are to be used by
	the Trainer to	assess analysis	skill and addi	tional training	needs.
22.4.7	B 15	Trainee	Trainer	Date	Grade
30.4.7		cise III – watch			•
					ridian. Pick one
		ch of the four fo	olders to view a	and discuss yo	ur observations with
	your trainer.	XO.	0		
		Trainee	Trainer	Date	Video #
	3/1/2	Trainee	Trainer	Date	Video #
	19.0	Trainee	Trainer	Date	Video #
	0)	Trainee	Trainer	Date	Video #
30.4.8	\ Practical Exerc	cise IV – Make t	wo new video:	s of your own	with a matrix of your
Cert	choosing and a	add them to the	e Trainee Video	o folder.	
604		Trainee	Trainer	Date	Video #
*		Trainee	Trainer	Date	Video #

	•)% - exercises w	vill be returned	l to the Trainee un
all answers ar	e correct.			
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date S	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
	Trainee	Trainer	Date	Packet #
30.5 Written Test – Module	28	iice 4		
	Trainee	Trainer	Date	Grade
30.6 Comparison Competen	. ()	- V)	-	_
mock case. Prints may		_	_	_
and non-matching prin		=		
such, Trainee will need	7	all appropriate	documentati	on and
attachments and issue	a report.			
O'S	Trainee	Trainer	Date	
30.7 Supervised Cases – Con	nplete 20 Sup	ervised Compa	arison Cases. '	Trainee shall
record all case number		-		
examiner.		•	J	. 0
10 ·				

Practical Exercise V - Complete comparison packets 1-10 as assigned by the

30.4.9

31.0 Module 29: Case Management and Reporting for Comparison and/or MBIS

31.1 Background and Theory

Forensic scientists are responsible for documenting the activities, methods, and results of their examinations in the case record. Written case records are recorded contemporaneously in ILIMS. All case documentation should be such that another qualified Latent Print Examiner could read the documentation and replicate the work. MBIS searches are also documented in ILIMS with supporting documentation attached so that they too may be evaluated by another qualified examiner.

31.2 Objectives, Principles, and Knowledge

- 31.2.1 Knowledge of and the ability to demonstrate proper procedures for maintaining chain of custody (documentation and physical control).
- Ability to navigate and query the various databases for location of criminal history records, fingerprint and palm print cards.
- 31.2.3 Ability to navigate and query ILIMS for latent print comparison and/or MBIS cases.
- Ability to demonstrate proper procedures for documentation of comparison casework. Documentation shall be such that another qualified Latent Print Examiner could evaluate what was done and replicate any comparisons.
- 31.2.5 Knowledge of and the ability to demonstrate proper procedures for reporting latent print comparison and MBIS examination findings in an accurate, concise, and clear manner.

31.3 Health and Safety Hazards

31.3.1 N/A

31.4 Reading and Practical Exercises

31.4.1 Complete Module 29 Reading List

0,	P				
·07		Trainee	Trainer	Date	
31.4.2	The Trainee so	hould attend a	Basic ILETS co	urse (attach certificate when	
		Trainee	Trainer	Date	
31.4.3		cise I – Obtain and obtaining k	J	d participate in Trainer led less ers.	son
		Trainee	Trainer	Date	

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 76 of 109

51.4.4	Trainer led dis	_	-	mparison repo	I IS III ILIMS –
		Trainee	 Trainer	 Date	
31.4.5	Practical Exerc				three comparison
011110	case reports. Pa		o onan maepen	acing product	o em co companioon
		Trainee	Trainer	Date	
31.4.6	Practical Exerc	ise IV – Techni	cal review train	ing for compa	rison cases -
	Trainer led dis	cussion and/or	demonstration	1.	
		Trainee	Trainer	Date	S
31.4.7	Practical Exerc	ise V – Trainee	shall perform a	administrative	and technical
			=	-	their Trainer. The
	Trainer will be			timately respo	nsible for the
	review on thes	e cases. Pass/F	ail.	5	
		Traine	Tuning	÷.C	
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
		Trainee	Trainer	Date	Case
	Silve	Trainee	Trainer	Date	Case
	Flor	Trainee	Trainer	Date	Case
	10	Trainee	Trainer	Date	Case
31.4.8				ing for MBIS c	ases - Trainer led
00,	discussion and	or demonstra	tion.		
40%		Trainee	Trainer	Date	
		Transce	Trainer	Dute	
5 Written	Test – Module 2	9			
		Trainee	Trainer	Date	Grade

32.0 Module 30: Court Procedures, Related Laws, Expert Testimony, Criminal and Civil Procedures Applicable to Latent Prints

32.1 Background and Theory

One of the most important parts of a forensic scientist's job is ensuring that the evidence that has been processed and evaluated is acceptable to the court system. ISPFS has numerous procedures to help ensure that evidence is handled and processed in an acceptable manner. It is also important to ensure that examiners are properly trained and prepared to testify as an expert witness within the field.

There are a number of important statutes and legal decisions that impact friction ridge testimony and the admission of evidence. It is important for latent print examiners to be familiar with some of the Federal Rules of Evidence, including Rules 701, 702, 703, and Rule 16. Important court cases include People v. Jennings, Frye v. United States, Daubert v. Merrel Dow Pharmaceuticals, US v. Byron Mitchell, US v Llera Plaza, and Mayfield v United States.

- 32.2 Objectives, Principles, and Knowledge.
 - 32.2.1 Understand the role of expert witness testimony.
 - 32.2.2 Knowledge of factors regarding the admissibility of evidence.
 - 32.2.3 Knowledge of relevant court cases and case histories.
 - 32.2.4 Understand the rules of discovery and evidence.
 - 32.2.5 Knowledge of applicable legal challenges to admissibility.
 - 32.2.6 Understand critical challenges to the discipline.
 - 32.2.7 Understand the advantages and disadvantages of different court chart types/methods (points, area bubbles, power point).
 - 32.2.8 Select appropriate prints and individual ridge characteristics for charting, create court charts, and utilize the digital imaging system to create court charts/exhibits.
 - 32.2.9 Ability to verbally articulate the friction ridge examination process and any resulting conclusions.
- 32.3 Health and Safety Hazards
 - 32.3.1 N/A
- 32.4 Reading and Practical Exercises
 - 32.4.1 Complete ISP FS core training court module

Trainee	Trainer	Date

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 78 of 109 Issuing Authority: Quality Manager

All printed copies are uncontrolled

32.4.2	Complete Mod	lule 30 Reading	List			
		Trainee	 Trainer	 Dat	 re	
32.4.3	The Trainee sl				training class when	
		ch copy of certif				
		Traine	e T	 'rainer	 Date	
32.4.4	Practical Exer	cise I – Write a t	hree-to-five	e-page pa	per on a recent court	
	development	as it relates to fi	ngerprints.	Pass/Fail		
		Trainee	Trainer	— — Dat	e e	
32.4.5	Practical Exer	cise II – Write oi	ne to two pa	aragraphs	outlining the argume	ents,
	decisions, and	impact of each	on the Scien	nce of Fric	ction Ridge Analysis fo	or eacl
	of the followin	g court cases: D	aubert v. M	errel Dov	v Pharmaceuticals, US	v.
	Byron Mitchel	l, US v Llera Plaz	za, and May	field v Un	ited States. Pass/Fail.	
22.4.6	Doortical Faces	Trainee	Trainer	Dat		
32.4.6			e your curri	cuium vii	ae utilizing the appro	priate
	template. Pass	5/Fa11.		40)		
		Trainee	Trainer	Dat	e	
32.4.7	Practical Exer	cise IV – Prepare	e a list of co	urt qualif	ying questions related	d to
		-		-	sample answers to th	
	= =	could be prese		7	-	
	1), ~ 0/	, ,	,	
	Processing	Trainee	Trainer	— — Dat	e	
		<u>x</u> O				
	Comparison	Trainee	Trainer	Dat	е	
32.4.8	Practical Exer	cise V –Discuss o	court dress	code, den	neanor, and etiquette	with
	your trainer. E	Explain to your t	rainer how	to procee	ed if there is an object	ion
	and what the p	protocol is for re	eferring to y	our note:	s. Discuss with your t	rainer
	what a "leadin	g question is" aı	nd why it is	importar	nt to remain accurate a	and
	transparent. I	Have your traine	er go throug	gh your qu	ıalifying questions. Pr	actice
	giving answer	s in a formalized	d manner. l	Pass/Fail		
X	4					
	_	Trainee	Trainer	Dat		
32.4.9			-	=	with the set of direct	
$\langle O \rangle$	=				t (NOTE: there are sep	
	=	_	=	-	Devise answers to thes	
	-	-		•	eting with your trainer	
		_	_	_	given to you. Practice	9
	giving answer	s in a formalized	d manner. l	Pass/Fail		
	Processing	Trainee	Trainer	— — Dat	e	
	Comparison	Trainee	 Trainer	— — Dat	 re	
Print Examir	ner Training Man	ual		Revision Issue Da	11 te: 12/15/2023	

Latent

32.4.10	Practical Exercise processing and of questions, directions, directions include case spetheir own making	comparison are t examination, c cific testimony.	separate). Mocl cross-examination. The trainee wil	c court will inclu on, and re-direc l utilize a compa	nde qualifying t. It will also arison chart of
	Processing	Trainee	Trainer	Date	
	Comparison	Trainee	Trainer	Date	
32.4.11	Practical Exercis	se VIII – Traine	e shall generate	a list of MBIS re	lated court
	qualifying quest	ions and provid	le sample answe	ers to those que	stions that
	could be present	ted in a court of	f law. Pass/Fail.	i	Co
		Trainee	Trainer	Date	*
32.4.12	Practical Exercis	se IX – Trainee s	shall attend the	testimony of tw	o forensic
	scientists. Pass/	/Fail.		ci ^C	
		Case #	Testifying Scientist	Date	
22.4.42	D 10	Case #	Testifying Scientist		
32.4.13	Practical Exercis				
	video from You				
	did well and wh				y will then
	discuss their cri	tique with their	trainer. Pass/F	all.	
		Trial/Video	Date	Trainer Discussion	Date
32.5 Written T	est – Module 30	~ ? ? ? ? ? ? ? ? ? ? ? ? ? ? ? ? ? ? ?	Date	Trainer Discussion	Date
	(19/3/10)	Trainee	Trainer	Date	Grade

33.0 Module 31: Multimodal Biometric Identification System (MBIS)

33.1 Background and Theory

Fingerprints are used as the foundation for criminal history records throughout the world. In 2016, the FBI's database was estimated to contain over 100 million fingerprint cards with the Idaho database having a little over half a million persons on file. Databases on all levels continue to grow with tens of thousands of individuals added to these repositories daily. These sophisticated computer file repositories are referred to as an Automated Fingerprint Identification System (AFIS) or Multimodal Biometric Identification System (MBIS). AFIS/MBIS is essentially a two part system: the ten-print system and the latent print system. The ten-print system is tasked with identifying sets of inked or Livescan fingerprints for criminal identification or employment purposes. The latent system is tasked with solving crimes through friction ridge impressions recovered from crime scenes or from items of evidence.

Idaho is a member of the Western Identification Network, Inc. (WIN). WIN was formed in 1988 to create a multi-state AFIS network. The members of WIN are Alaska, Montana, Oregon, Washington, Nevada, Utah, Wyoming, California and Idaho. WIN offers access to 20 million fingerprint records held within the western United States.

33.2 Objectives, Principles, and Knowledge

- 33.2.1 Understand automation technology and theory of operation to include:
 The history of the development of friction ridge automation technology;
 Theory of the operation of friction ridge automation technology, to include an understanding of distortion that may occur when three-dimensional friction ridge skin is captured as a two-dimensional image.
- Understand the function and use of image capture to include:

Types of friction ridge recordings (e.g. rolled, flat, simultaneous, palm);

Methods of friction ridge capture (e.g. ink, livescan);

Types of capture devices (e.g. livescan, flatbed, camera);

Point of capture variables (e.g. condition of fingers, condition of platen, rolling speed, movement);

Control measures needed to achieve quality friction ridge images (e.g. scan resolution, compression rate, equipment maintenance, calibration); Procedures for addressing amputations, temporary injuries, skin conditions, and rescans.

Understand the function and use of Multimodal Biometric Identification Systems (MBIS) to include:

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Issuing Authority: Quality Manager

Page 81 of 109 Issu All printed copies are uncontrolled

MBIS process related to acquisition, classification, searching, storage, retrieval, identification, and final reporting of friction ridge records;

Friction ridge search criteria (e.g. designated finger search, how many fingers, palm areas);

Importance of quality assurance on maintaining the integrity of friction ridge data;

Quality controls that ensure completeness, image quality, and data integrity.

Gain a working knowledge of the NEC Multimodal Biometric Identification System (MBIS) and the Integrated Automated Fingerprint Identification System (IAFIS) to include:

Who handles component maintenance and calibration;

System requirements and limitations including text data fields, fingerprint and palm print quality, finger sequence and image replacement, image rotation, and toleration for pattern interpretation;

Minutia recognition, placement, rotation, ridge counts, and other minutiae factors related to searching and matching;

Limitations of system interoperability;

Integration of friction ridge image, mug shot, scars, marks, tattoos, minutiae, other biometrics, as well as personal descriptors, and criminal history information;

Search parameters, pattern classification and referencing, minutiae extraction, search algorithms, significance in the range of candidate scores, threshold scoring, and candidate list comparisons, matching;

AFIS search capabilities in regards to latent print vs. ten print, ten print vs. latent print, latent print vs. latent print, ten print vs. ten print, and palm print vs. palm print;

"Lights out" processing of searches and mobile search capabilities;

Logical search progression (i.e. state, regional, national);

Filtering criteria used to establish logical candidates (e.g. finger position, sex, classification, race, offense, geographic location);

Search result contents (e.g. ranked order, unique identifier, finger or palm position);

Differences between AFIS digital images and original friction ridge impressions (e.g. potential loss of quality due to compression of image, monitor resolution, capture resolution);

Printer technology limitations vs. examinations from original friction ridge documents (e.g. paper quality, inked fingerprint cards);

AFIS processes related to latent print searches;

Various search options among databases within the system (e.g. image, feature);

Manual and automatic encoding of minutiae;

File penetration benefits and liabilities of partial vs. full data base searches; Record authentication processes (e.g. correct association of name, unique identifier, friction ridge images, and criminal history record).

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 82 of 109

3.4 Reading	and Practical E	xercises			
33.4.1	Complete Mod		g List		
		Trainee	 Trainer	Date	
33.4.2	The Trainee w	ill review the	AFIS Course Bi	nder and pass th	e associated te
		Trainee	Trainer	Grade	
33.4.3		-		rches through I	
	Working as "the Quality/Proce			fined by the ISP	HS 63
	Quarrey/11000				MIC
22.4.4	MDIC Communication	Trainee	Trainer	Date	E and become
33.4.4	-	=	-	endently search	
				lentification Syst	•
	= =	-	=	ninutia prints, d	=
	non-matching	prints. This c	ompetency test	will be entered	into ILIMS, as
	_	- 1. 1		1 107	0.1
	Trainee will ne		ent searches, at	tach proper MBI	S documentati
	_		ent searches, at	tach proper MBI	S documentati
	Trainee will ne	oort.	~		S documentati
	Trainee will ne		ent searches, at Trainer	tach proper MBI	S documentati
33.5 Written	Trainee will no and issue a rep	Trainee	~		S documentati
33.5 Written	Trainee will ne	Trainee	~		S documentati
33.5 Written	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer		S documentati Grade
33.5 Written	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
33.5 Written	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	
33.5 Written	Trainee will no and issue a rep Test – Module 3	Trainee	Trainer	Date	

34.0 Module 32: DNA Database Fingerprint Comparison

34.1 Background and Theory

Friction ridge identification and classification has a long history rooted in scientific research and empirical observations.

Various classification systems have been used over the past 100 years. Today's classification systems rely mainly upon computers to digitize, categorize, recall, and identify matching 10-print cards.

Examiners must be able to recognize and articulate the various patterns and subpatterns and understand their use in analysis and comparison.

The scientific method is a method of research in which a problem is identified, relevant data is gathered, and a hypothesis is formulated from the data and then tested. In forensic science, it is imperative to have a scientific technique for examination. Doing so ensures that evidence is treated equally, and conclusions are reliable and unbiased. The latent print section utilizes ACE-V as part of the examination methodology. ACE-V is an acronym that stands for analysis (A), comparison (C), evaluation (E) and verification (V). It is the process that latent print examiners utilize to reach a conclusion about a comparison examination.

34.2 Objectives, Principles, and Knowledge

- 34.2.1 Understand the basic biology and physiology of friction ridge skin.
- Understand the basic foundations of the science of friction ridge identification (persistence and discriminability).
- 34.2.3 Understand common terminology and definitions associated with friction ridge pattern recognition (arch, loop, and whorl).
- 34.2.4 Ability to differentiate between pattern types.
- 34.2.5 Understand friction ridge characteristics (dots, ridge endings, and bifurcations) the varying definitions/interpretations assigned to combinations of those three ridge characteristics, and how they may be utilized in effecting identification.
- Ability to successfully analyze and compare known fingerprint cards to plain inked fingerprint impressions.
- 34.2.7 Ability to render an accurate conclusion (identification, inconclusive, exclusion).
- 34.2.8 Understand the necessity for verification by another qualified latent print examiner.

Latent Print Examiner Training Manual

Revision 11 Issue Date: 12/15/2023

Page 84 of 109

34.3	34.3.1	d Safety Hazards N/A				
34.4	Reading a	nd Practical Exe	rcises			
	34.4.1	Complete Module	e 32 Reading Lis	t		
		_		Trainer	Date	
	34.4.2	Practical Exercise	e I – Trainer led	lesson on com	parison.	
				Trainer	Date	
	34.4.3	Practical Exercise documentation.	e II – Trainer lec	l lesson on DN	A database card	ces
			Trainee	Trainer	Date	
	34.4.4	Practical Exercise is 80%.	e III – Pattern re	ecognition – "10	00 fingerprints'	. Passing score
			Trainee	Trainer	Date	Grade
	34.4.5	Practical Exercise 100% of identific quality, there ma completed – this	cations effected a y be comparison	are correct. Duns that were at	e to examiner s	kill level or card
34.5	Written T	est – Module 32	Trainee	Trainer	Date	
		C	Trainee	Trainer	Date	Grade
	10 DNA D	on Competency T atabase Card Sar cation.		•		•
O [*]	operty		Trainee	Trainer	Date	

Appendix I – Reading Lists

Module 1 Reading List: History and Background of Friction Ridge Identification

Fingerprint Techniques - Andre Moenssens		
Chapter 1 - The History of Fingerprinting		
Advances in Fingerprint Technology, 2^{nd} edition - Lee, Gaensslen Chapter 1 - History and Development of Fingerprinting.	— jice	5
The Fingerprint Sourcebook – Scientific Working Group on Friction I Technology (SWGFAST), et al Chapter 1 - History	Ridge Analysis, Stud	dy and
Quantitative-Qualitative Friction Ridge Analysis - David R. Ashbaugh Chapter 2 - History of Fiction Ridge Identification	l. 	
Module 1 Supplemental Information Packet		
Module 2 Reading List: Other Scientific Personal Identification	Methods	
Biometrics Overview pdf Iris Recognition pdf		
Iris Recognition pdf		
Face Recognition pdf		
Vascular Pattern Recognition pdf		
Hand Geometry pdf		
Criminalistics, 9th edition Richard Saferstein		
Chapter 13, "DNA" Pages 380-418		
Chapter 16, "Document and Voice Examination" Pages 496-521		
Death Investigator's Handbook by Louis N. Eliopulos,		
Chapter 67 "Forensic Odontology Pages 679 – 693		

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 86 of 109 Issuing Authority: Quality Manager

"Handwriting and Handprinting Identifications." Pages 710-717
Module 3 Reading List: Safety Training
Latent Print Section Quality Manual sections: Safety
Chemicals, Supplies, and Reagent Preparation Equipment, Calibration, Maintenance and Repair
Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation Manual. 2014. Center for Applied Science and Technology. Section 3.2 – Working Safely
Module 4 Reading List: Case Management and Reporting for Processing
ISO/IEC 17025:2017 Section 7.8 Reporting of results
ISPFS Quality/Procedure Manual Section on "Technical records" Section on "Facilities and Environmental Conditions" Section on "Reporting of Results"
Latent Print Section Quality Manual - Documentation and Report Writing
Guideline - SWGFAST Document 5 Standard for Reporting Friction Ridge Examinations (Latent/Tenprint) or the OSAC successor document
Paper – "Integrating DNA Collection into the Latent Print Section."
Paper – "Evaluation of the impact of Different Visualization Techniques on DNA in Fingerprints."
ISPFS BEAST ILIMS General User Guide
Module 5 Reading List: Digital Preservation of Latent Prints
User's manual for the Nikon D810
Latent Print Evaminer Training Manual Revision 11

Page 87 of 109

User's manual for the Cannon EOS 6D	
User's manual for the Epson V700/V800/V850 pro	
Latent Print Section AM Section - Digital Imaging Procedure	
Foray Adams v6 User Manual and Adams Web Help files	
Guideline - SWGFAST Document 6 Standard for Friction Ridge Di successor document	igital Imaging or the OSAC
ASTM E3235-21 Standard Practice for Latent Print Evidence Ima	nging Resolution
Guideline - SWGIT Section 19 Issues Relating to Digital Image Co OSAC successor document	mpression and File Formats or the
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Ed Chapter 16, Digital Imaging – Sections 16.1-16.3	lition - Ramotowski
A Short Course in Photography, Digital – London & Stone Chapter 1 - Camera Chapter 2 - Lens Chapter 3 - Light and Exposure Home Office Center for Applied Science and Technology (CAST).	Fingerprint Visualisation
Manual. 2014. Center for Applied Science and Technology. Section 5.VE – Visual Examination National Centre for Forensic Studies - Fingermark Detection & En Lennard, Chapter 6 - Digital Imaging	nhancement 6 th Edition– Stoilovic &
Crime Scene Photography, 2 nd Edition – Robinson Chapter 1 – History of Forensic Imaging Chapter 2 - Composition and Cardinal Rules Chapter 3 - Basic Exposure (non-flash) Concepts Chapter 4 – Focus, Depth of Field, and Lenses Chapter 6 - Crime Scene Photography – "Close up Photographs" 3	336-341
Chapter 7 – Ultraviolet, Infrared and Fluorescence Chapter 10 - Digital Imaging Technologies	
The Fingerprint Sourcebook by Scientific Working Group on Fric Technology (SWGFAST), et al. Chapter 8 - The preservation of Friction Ridges.	etion Ridge Analysis, Study and
8	Revision 11 ssue Date: 12/15/2023

Fingerprints and other Ridge Skin Impressions, 2 nd Edition - Champod et al Section 3.5 – Photography
Module 6 Reading List: General Latent Print Processing
Registry Proposed - OSAC-2022-N-0033 Standard for Processing Evidence for the Detection of Friction Ridge Impressions
Latent Print Section AM - General Latent procedure to include Quick Reference Sequential Processing Guide
Latent Print Section Quality Manual – Evidence Control and Handling
Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation Manual. 2014. Center for Applied Science and Technology. Section 2.1 – An Introduction to Forensic Evidence Recovery Section 2.2 – Understanding Fingermarks Section 2.3 – Fingermark Visualisation Processes Section 2.5 – Using and Understanding Fingermark Evidence
The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al. Chapter 7 - Latent Print Development Chapter 11 - Equipment Fingerprint Detection with Lasers - Menzel Chapter 7 - Sections 7.1 & 7.2
Fingerprints and other Ridge Skin Impressions, 2 nd Edition - Champod et al Chapter 4 - Fingerprint Detection Techniques
Module 7 Reading List: Processing Technique – Alternate Light Sources
Latent Print Section AM - Alternate Light Source
Applicable ALS User Manuals
Fingerprints and other Ridge Skin Impressions, 2 nd Edition - Champod et al Sections 3.3-3.4 – Light theory

Revision 11

Issue Date: 12/15/2023

Page 89 of 109 Issuing Authority: Quality Manager

An Introduction to Lasers, Forensic Lights, and Fluorescent Fing Roland Menzel.	gerprint Detection Techniques, by A.
Fingerprint Detection with Lasers – Menzel Chapter 9 – Excitation Optimization and Filters	
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd E Chapter 7, Miscellaneous Methods and Challenging Surfaces - So	
National Centre for Forensic Studies - Fingermark Detection & Lennard	Enhancement 6th Edition- Stoilovic
Chapter 2 - General Nature of Light	
Chapter 3 - Optical Filters	
Chapter 4 - Optical Examination Techniques	
Chapter 5 - Forensic Light Sources	
Home Office Center for Applied Science and Technology (CAST) <i>Manual.</i> 2014. Center for Applied Science and Technology. Section 5.FE – Fluorescence Examination Note: additional readings for this section were covered in Modu	\
Module 8 Reading List: Processing Technique - Amido Blac	k
Latent Print Section AM - Amido Black	
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd E Chapter 6 Blood Reagents, Section 6.1 & 6.2 (pgs. 129-134 & 14	
Home Office Center for Applied Science and Technology (CAST) <i>Manual.</i> 2014. Center for Applied Science and Technology. Section 5.AD – Acid Dyes	. Fingerprint Visualisation
Paper – "Summary of Experiments Investigating the Impact of Fingerprint Reagents on PCR-based DNA Typing Profiles."	Fingerprint Processing and ————————————————————————————————————
Paper – "Chemical Enhancement of Fingerprints in Blood: An Education DNA, and Assessment of Chemical Hazards."	valuation of Methods, Effects on
Paper – "The Effect of Common Fingerprint Detection Techniqu	es on the DNA Typing of
Fingerprints Deposited on Different Surfaces. JFI, Vol. 54, No. 1,	-1 -
Latent Print Examiner Training Manual	Revision 11 Issue Date: 12/15/2023
Page 90 of 109	Issuing Authority: Quality Manager

Paper – Presumptive Testing for Blood on a Patent Print Developed with Amido Black."
Paper – "Deposition of Bloody Friction Ridge Impressions." JFI, Vol. 58, No. 3, 2008
Paper – "Developing Fingerprints in Blood: A Comparison of Several Chemical Techniques." JFI, Vol 57, No. 1, 2007
Note: additional readings for this section were covered in Module 6
Module 9 Reading List: Processing Technique – 1, 8-Diazafluoren-9-One (DFO), 1, 2 – Indanedione, and 1, 2 – Indanedione TP
Latent Print Section AM - DFO Latent Print Section AM 1, 2 – Indanedione Latent Print Section AM 1, 2 – Indanedione TP Fingerprint Detection with Lasers – Menzel Chapter 8 - Sections 8.3, 8.5, & 8.6
Fingerprint Detection with Lasers – Menzel Chapter 8 - Sections 8.3, 8.5, & 8.6
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowski Chapter 2 Amino Acid Reagents Sections: 2.4 & 2.5
Home Office Center for Applied Science and Technology (CAST). <i>Fingerprint Visualisation Manual</i> . 2014. Center for Applied Science and Technology. Section 5.DFO – DFO
Section 6.1.25-6.1.31 – Indandione
Paper – "Spectral Variations for Reaction Products Formed between Different Amino Acids and Latent Finger mark Detection Reagents on a Range of Cellulose-Based Substrates. JFI, Vol. 59, No. 3, 2009
Paper – "The Effectiveness of 1, 2-Indandione-Zinc Formulations and Comparison with HFE-Based 1, 8-diazafluoren-9-one for Fingerprint Development." JFI Vol. 59, No. 6, 2009
Paper – "DFO, Its Usage and Results," Masters, Morgan & Shipp
Paper – "1, 2-Indandiones: New Reagents for Visualizing the Amino Acid Components of Latent Prints." JFS Vol. 43, No. 4. 1998, pp. 744 – 747.

Revision 11

Issue Date: 12/15/2023

Paper – "Optimisation and Evaluation of 1, 2-ind Application to Real Samples." Forensic Science In	anedione For Use as a Fingermark Reagent and Its nternational. Vol. 168. 2007, pp. 14 – 26.
Paper – "Thermal Paper: Latent Friction Ridge D pp. 265-271	evelopment via 1, 2 Indanedione. JFI, Vol.53 (3),
Note: additional readings for this section were co	overed in Module 6
	S
Module 10 Reading List: Processing Techniqu	ie – Dye Stains – Rhodamine 6G and RAM
Latent Print Section AM - Rhodamine 6G	
Latent Print Section AM - RAM	
Butcht Time Section That Turns	:.0
Home Office Center for Applied Science and Tech	anology (CAST) Fingerprint Visualisation
Manual. 2014. Center for Applied Science and Te	
Section 5.SFDS – Superglue Fluorescent Dye Stai	
1 0	
Lee and Gaensslen's Advances in Fingerprint Ted	chnology, 3rd Edition - Ramotowski
Chapter 5 Vapor/Fuming Methods, Section 5.1.6	• []
onapoor o vapor/ rummig roomous, soomen errie	
Fingerprint Detection with Lasers – Menzel	
Chapter 7 – Section 7.3	
diapter / Section / is	
Note: additional readings for this section were co	overed in Module 6
Module 11 Reading List: Processing Techniqu	ie – Gentian Violet/Crystal Violet
Latent Print Section AM - Gentian Violet	
Latent Finit Section Aim - Gentian Violet	
Home Office Contact for Applied Science and Took	analogy (CAST) Eingermint Vigualization
Home Office Center for Applied Science and Tech Manual. 2014. Center for Applied Science and Te	
Section 5.BV3 – Basic Violet 3	chilology.
Section 3.DV3 - Dasic Violet 3	
Lee and Gaensslen's Advances in Fingerprint Tec	phnology 2rd Edition Domotovyski
<u> </u>	
Chapter 1, Vapor/Fuming Methods Section 5.1 (p	Jgs. 113-114)
Paper – "Development of Latent Fingerprints on Brightening."	Sticky Surfaces by Dye Staining or Fluorescent
Latent Print Examiner Training Manual	Revision 11

Issue Date: 12/15/2023

Note: additional readings for this section were covered in Module $\,6\,$

Module 12 Reading List: Processing Technique - Iodine	
Latent Print Section AM - Iodine	
Home Office Center for Applied Science and Technology (CAST). <i>Manual</i> . 2014. Center for Applied Science and Technology. Section 6.1.32-6.1.40 – Iodine Fuming	. Fingerprint Visualisation
The Science of Fingerprints - FBI. "Iodine Method." Pages 175-1	77
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Ed Chapter 5 Vapor/Fuming Methods, Section 5.2	dition - Ramotowski
Note: additional readings for this section were covered in Modu	le 6
Module 13 Reading List: Processing Technique - Leuco Crys	stal Violet (LCV)
Latent Print Section AM - Leuco Crystal Violet	
Home Office Center for Applied Science and Technology (CAST). <i>Manual</i> . 2014. Center for Applied Science and Technology. Section 6.1.42 – Leuco Crystal Violet	
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Ed Chapter 6 Blood Reagents, Section 6.1 & 6.2 (pgs. 148-149)	uitioii - Kaiiiotowski ————————————————————————————————————
Paper – "Lueco Crystal Violet: A Simple, Effective Blood Enhance	ement Reagent."
Note: additional readings for this section were covered in Modu	le 6
Module 14 Reading List: Processing Technique - Ninhydrin	and ThermaNin
Latent Print Section AM - Ninhydrin Latent Print Section AM ThermaNin	
CARON Fingerprint Development Chamber Operations Manual	
8	Revision 11 Issue Date: 12/15/2023

Lee and Gaensslen's Advances in Fingerprint Technology	
Chapter 2, Amino Acid Reagents Sections: 2.1, 2.4, & Chapter 7, Challenging Surfaces, Sections 7.2 (pgs. 1)	
chapter 7, changing surfaces, sections 7.2 (pgs. 1)	05-105)
The Science of Fingerprints - FBI. "Ninhydrin Metho	d." Pages 177-179
Home Office Center for Applied Science and Technol Manual. 2014. Center for Applied Science and Technol	
Section 5.Nin – Ninhydrin	
Paper – "Procedure to Develop Latent Prints on The	rmal Paper"
Paper – "Latent Fingerprints by a Superior Ninhydri	n Method"
Paper – "Ninhydrin Processing by Pat A. Wertheim"	
Paper – "Determining the Length of Time Required f No. 4, 2017	for Ninhydrin Development," JFI, 2017, Vol. 67,
Paper - "The Effectiveness of Ninhydrin Latent Print Regards to Climatic Conditions at the Time of Depos	
Paper – "Improved Results in the Development of La 58, No. 4, 2008	tent Fingerprints on Thermal Paper." JFI, Vol.
Paper - "A Limited Validation and Comparison of 1, 2 Development on Thermal Paper." JFI, Vol. 66(3), pp.	
Paper – "Thermal & Carbonless Papers: A Fundamer Development." JFI, Vol. 53(2), pp. 185-197	ntal Understanding for Latent Friction Ridge
Paper – "Chemical Fuming: A Practical Method for Fi	ngerprint Development on Thermal Paper." JFI
Note: additional readings for this section were cover	red in Module 6
Module 15 Reading List: Processing Technique -	Powder Development of Latent Prints
Latent Print Section AM - Powder Detection Method Latent Print Section AM - Lifting Methods	S
Latent Print Examiner Training Manual	Revision 11 Issue Date: 12/15/2023

Page 94 of 109

Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowski Chapter 1, Powder Methods Section 1.1 (pgs. 1-5)
The Science of Fingerprinting - FBI. Chapter 14, "Powdering and Lifting Latent Impressions." Pages 173-174
Fingerprint Techniques, by Andre A. Moenssens, Chapter 4, "Latent Prints," Pages 106-114
Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation Manual. 2014. Center for Applied Science and Technology.
Section 5.Lif – Lifting
Section 5.Pow – Powders
Section 6.2.12 – Powders (Fluorescent)
Paper – "Evaluation of Fingerprint Powders." JFI, Vol. 56, No. 2, 2006
Paper – Beware of the Possibility of Fingerprint Techniques Transferring DNA," Journal of Forensic
Science, Vol.50, No.6, 2005
Module 15 Supplemental Information Packet
Note: additional readings for this section were covered in Module 6
Module 16 Reading List: Processing Technique – Physical Developer (PD)
Latent Print Section AM - PD
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowski Chapter 3 Metal Deposition Methods: Section 3.2
Home Office Center for Applied Science and Technology (CAST). <i>Fingerprint Visualisation Manual</i> . 2014. Center for Applied Science and Technology. Section 5.PD – Physical Developer
Paper – "Physical Developer" - David Burow
Paper – "Physical Developer: A Practical and Productive Latent Print Developer"
Paper – "PD, Maleic Acid and Synperonic N"

Revision 11

Issue Date: 12/15/2023

Paper – "The Efficacy of Commercial vs. N Sequential Enhancement of Friction Ridge 2010	•	
Paper – "Physical developer method for do of Forensic Sciences	etection of latent fin	gerprints: A review." Egyptian Journal
Note: additional readings for this section	were covered in Mod	lule 6
		65
Module 17 Reading List: Processing Te	chnique – RECOVEI	RLFT
Latent Print Section AM – RECOVER LFT		
Foster+Freeman, "RECOVER Latent Finger	rprint Technology U	ser Manual," 2019
Idaho State Police Forensic Services, Later Validation, 2022	nt Print Section, Fost	ter+Freeman RECOVER LFT
Paper "Recovery of Fingermarks from Fire Devices using S_2N_2 – A proof of Concept St		
Module 18 Reading List: Processing Te	chnique – Small Par	rticle Reagent (SPR)
	100	
Latent Print Section AM - SPR Home Office Center for Applied Science ar Manual. 2014. Center for Applied Science Section 5.SPR - Small Particle Reagent		T). Fingerprint Visualisation
I and a second of the second o	era maraka alam 2001	Elitica Bassas all'
Lee and Gaensslen's Advances in Fingerpr Chapter 1, Powder Methods Section 1.2.1	int Technology, 3rd	
Paper – "Development of Latent Prints Us White (SPR-W) on Adhesives." JFI, Vol. 55	•	e (TiO2) in Small Particle Reagent,
Paper - "Report of Validation Testing" Sirc	chie SPR-W by Albuq	uerque Police
Paper – "Small Particle Reagent" by Pat A.	Wertheim	
Paper – "Lightning Powder Co. Technical I	Note Small Particle R	deagent"
Note: additional readings for this section v	were covered in Mod	lule 6
Latent Print Examiner Training Manual	Page 96 of 109	Revision 11 Issue Date: 12/15/2023 Issuing Authority: Quality Manager

All printed copies are uncontrolled

Equivalent Latent Print Section AM - Sticky Side Powder Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation Manual. 2014. Center for Applied Science and Technology. Section 5.PS – Powder Suspension Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowsk Chapter 1, Powder Methods Section 1.2.2 & 1.2.3 Paper – "Homemade Solution for Processing Latent Prints on the Adhesive Side of Tape." Paper - "A Black Powder method to Process Adhesive Tapes." Paper - "Anomalous Results with Sticky Side Powder." Paper - "A New Approach to Unraveling Tangled Adhesive Tape or Potential Detection of Latent Prints and Recovery of Trace Evidence Paper - "Does CA Fuming Interfere with Powder Suspension Processing?" JFI, Vol. 59, No. 2, 2009 Paper – "The Effects of Cyanoacrylate Fuming and Rhodamine 6G on the Adhesive side of Tape when Processing with Adhesive-side Powders" JFI, Vol. 70, No. 1, 2020 Note: additional readings for this section were covered in Module 6 Module 20 Reading List: Processing Technique - Sudan Black Latent Print Section AM - Sudan Black Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation *Manual.* 2014. Center for Applied Science and Technology. Section 5.SB3 – Solvent Black 3 Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowski Chapter 4 Lipid Reagents, Section 4.1

Module 19 Reading List: Processing Technique – Sticky Side Powder/Sticky Side Powder

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Friction Ridge Skin, by James F. Cowger, "Locating, Developing, Preserving, and Collecting Evidence Prints." Page 104
Note: additional readings for this section were covered in Module 6
Module 21 Reading List: Processing Technique - Cyanoacrylate Ester (Super Glue)
Latent Print Section AM - Cyanoacrylate Ester
MEGAfume User Manual
SAFEFUME Cyanoacrylate Fuming Chamber Operating Manual
"Fast Vac" – Operating Instructions – CAE Vacuum chambers
"AMETEK" – Use and Installation of Pressure Gauges – CAE Vacuum chambers ————————————————————————————————————
Home Office Center for Applied Science and Technology (CAST). Fingerprint Visualisation Manual. 2014. Center for Applied Science and Technology. Section 5.SF – Superglue Fuming
Lee and Gaensslen's Advances in Fingerprint Technology, 3rd Edition - Ramotowski Chapter 5 Vapor/Fuming Methods, Section 5.1 (pgs. 98-105 & 115-116)
Chapter 11 Cyanoacrylate Fuming Method
Paper – "A Modified Cyanoacrylate Technique Utilizing Treated Neutral Filter Paper for Developing Latent Fingerprints"
Paper - "Fivis by 3M – Instructions and Notes"
Paper - "Effects of Cyanoacrylate Processing on Cocaine HCL Trace Analysis"
Note: additional readings for this section were covered in Module 6
Module 22 Reading List: Digital Imaging
Latent Print Section AM - Digital Imaging Procedure
FORAY Adams V6 User Mannual and review AdamsWeb Help File

Revision 11

Issue Date: 12/15/2023

Review Current Adobe Photoshop user manual (highlighted sections on Photoshop manual reading document)
Techniques of Crime Scene Investigation, Sixth Edition - Barry A. J. Fisher Page 113
Crime Scene Photography, 2 nd Edition – Robinson Chapter 11 - Digital Imaging Processing of Evidentiary Photography
A Short Course in Photography, Digital – London & Stone Chapter 4 - Digital Workplace Basics Chapter 5 - Image Editing
Criminalistics 9 th edition, An Introduction to Forensic Science - Richard Saferstein. Pages 452-454, 509-510
Advances in Fingerprint Technology, 2^{nd} edition - Lee & Gaensslen. Page 267
Guideline - SWGFAST Document 6 Standard for Friction Ridge Impression Digital Imaging (Latent/Tenprint) or the OSAC successor document
ASTM Standard Terminology for Digital and Multimedia Evidence Examination E2916-19 ^{e1} (See Trainer or Discipline Lead)
Guideline - SWGIT Section 5 Guidelines for Image Processing or the OSAC successor document
Guideline - SWGIT Section 11 Best Practices for Documenting Image Enhancement or the OSAC successor document
Paper – "Digital Enhancement of Latent Prints using Adobe Photoshop Black & White Adjustments." JFI, Vol. 59, No. 4, 2009
Paper – "Image Enhancement and Adobe Photoshop: Using Calculations to Extract Image Detail." JFI, Vol. 57, No. 4, 2007
Paper – "Techniques for Digital Enhancement of Latent Prints Obscured by Disruptive Backgrounds." JFI, Vol. 54, No. 2, 2004
Paper – "Computer Fingerprint Enhancement: The Joy of Lab Color." JFI, Vol. 62, No. 5, 2012
Paper – "Adapting Narrow Bandpass Filters to Photography." JFI, Vol. 62, No. 3, 2012
Paper – "Improved Multiple Exposure and Panoramic Photography of Latent Fingerprints." JFI, Vol. 63, No. 1, 2013

Revision 11 Issue Date: 12/15/2023

Page 99 of 109

Module 23 Reading List: Biology and Physiology of Friction Ridge Skin

The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al.
Chapter 2 - Anatomy and Physiology of Adult Friction Ridge Skin Chapter 3 - Embryology and Morphology of Friction Ridge Skin ———————————————————————————————————
Scott's Fingerprint Mechanics - Robert D. Olsen Sr., Pages 114-125
Fingerprint Techniques – Andre Moenssens Chapter 2 - The Nature of Friction Skin Chapter 11, Pages294-297
Finger Prints, Palms and Soles - Harold Cummins and Charlie Midlo Chapter 10 - Embryology Chapter 12 - Inheritance
Advances in Fingerprint Technology, 2 nd Edition - Lee & Gaensslen, Chapter 3 - Composition of Latent Print Residue
Quantitative-Qualitative Friction Ridge Analysis - David R. Ashbaugh. Chapter 3 - Friction Ridge Medium
Fingerprints and Other Friction Ridge Skin Impression - Christophe Champod et. al. Chapter 1 - Friction Ridge Skin
Paper – "The Critical Stage of Friction Ridge Skin and Pattern Formation - Kasey Wertheim and Alice Maceo
Paper – "Qualitative Assessment of Skin Deformation: A Pilot Study." JFI, Vol. 59, No. 4, 2009
Paper – "Discriminability of Fingerprints of Twins." JFI, Vol. 58, No. 1, 2008
Paper – "Fingerprint Patterns: A Study on the Finger and Ethnicity Prioritized Order of Occurrence. JFI, Vol. 55, No. 4, 2005
Paper – "Permanent Intentional Fingerprint Mutilation" - Kasey Wertheim
Paper – "An Extreme Case of Fingerprint Mutilation." JFI, Vol. 48, No. 4, 1998
Paper – "Fingerprint Formation," Kucken, Journal of Theoretical Biology, Vol. 235, No. 1, 2005

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

Page 100 of 109

Module 24 Reading List: Recording Inked Fingerprints, Palm Prints, and Footprints Latent Print Section AM Section – Taking Known Exemplars Scott's Fingerprint Mechanics - Robert D. Olsen Sr. Chapter 2 - Taking Finger, Palm, and Footprints Fingerprint Techniques - Andre A. Moenssens Chapter 5, "Recording Prints." Pages 137-145. The Science of Fingerprints - FBI Chapter 9, "Techniques for Taking Good Fingerprints." Pages 111-115_ Chapter 10, "Problems in Taking Inked Fingerprints." Pages 116-128 Finger Prints, Palm and Soles - Harold Cummins, Charles Midlo Chapter 3, "Methods of Printing." Pages 45-55 Friction Ridge Skin - James F. Cowger Chapter 2, "Taking Inked Prints." Pages 9-33 The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al. Chapter 4, "Recording Living and Postmortem Friction Ridge Skin Exemplars," sections 4.1-4.3 Module 25 Reading List: Friction Ridge Pattern Recognition and Interpretation Criminalistics, 9th edition - Richard Saferstein Chapter 14 "Classification of Fingerprints." Pages 435-436 Scott's Fingerprint Mechanics - Robert D. Olsen Sr. Chapter 1 Sections 7 Fingerprint Classification," 8 "Space Value on Fingerprint Cards," and 9 "Fingerprint Patterns are Complex Yet Simple." Pages 17-21 Friction Ridge Skin, by James F. Cowger Chapter 3 - Classification Fingerprint Techniques - Andre A. Moenssens Chapter 3 - Pattern Interpretation Chapter 6 - Fingerprint Classification in the United States Fingerprints and the Law - Andre Moenssens Chapter 2, "Fingerprint Principles and Techniques." Pages 10-23

Latent Print Examiner Training Manual

Revision 11

Issue Date: 12/15/2023

The Science of Fingerprints - The FBI. Chapters - 2-8. Pages 5-110	
The Fingerprint Sourcebook - Scientific Working G Technology (SWGFAST), et al. Chapter 5 - Systems of Fingerprint Classification	roup on Friction Ridge Analysis, Study and
Module 26 Reading List: Introduction to Latent	Prints and the State of the Science
The Fingerprint Sourcebook by Scientific Working Technology (SWGFAST), et al.	Group on Friction Ridge Analysis, Study and
Chapter 14 - Scientific Research Supporting the Fo	undations of Friction Ridge Examinations
Strengthening Forensic Science in the United State Identifying the Needs of the Forensic Sciences Com Executive Summary Friction Ridge Analysis pages 136-145	
PCAST Forensic Science in Criminal Courts: Ensuri Methods. Executive Summary Section 5.4 Latent Fingerprint Analysis 2016 Addendum to the PCAST report of Forensic S	
NIST/NIJ Latent Print Examination and Human Fac Approach, the Report of the Expert Working Group Chapter 1	
International Association of Identification "Positio and Probability Modeling" & "Resolution 2016-4"	n Statement on Conclusions, Qualified Opinions,
Module 27 Reading List: Human Factors	
The Fingerprint Sourcebook by Scientific Working Technology (SWGFAST), et al. Chapter 15: Special Abilities and Vulnerabilities in	
Latent Print Examiner Training Manual	Revision 11 Issue Date: 12/15/2023

Page 102 of 109 Issu All printed copies are uncontrolled

Latent Print Examination and Human Factors: Improving the Practice through a Systems Approach The Report of the Expert Working Group on Human Factors in Latent Print Analysis 2012 Chapters 2 – Human Factors and Errors Chapter 3 - Interpreting Latent Prints Chapter 7 – A Systems Approach to the Work Environment Chapter 8 – Training and Education ———————————————————————————————————
Paper "The Authority of Fingerprint Experts: Is it Based on Belief or Science?" JFI, Vol. 59, No. 6, 2009
Paper – "Why Experts Make Errors." JFI Vol. 56, No. 4, 2006
Paper – "A Report of Latent Print Examiner Accuracy During Comparison Training Exercises." JFI, Vol. 56, No. 1, 2006
Paper – "Subjective- The Misused Word." William Leo. JFI Vol. 58, No. 1, 2008
Paper - "Accuracy and Reliability of Forensic Latent Fingerprint Decisions." Ulery et al. PNAS, Vol. 61, No. 4, 2011
Paper - "Latent Fingerprint Quality: A Survey of Examiners." Hicklin et al. JFI, Vol. 61, No. 4, 2011
Paper - "Measuring what Latent Fingerprint Examiners Consider Sufficient Information for Individualization Determinations." Ulery et al. PLoS ONE, Vol. 9, No. 11, 2014
Paper - "Understanding the sufficiency of information for fingerprint value determinations." Ulery et al. Forensic Science International, Vol. 226, No. 1, 2013
Paper - "Inter-examiner variation of minutia markup on latent fingerprints." Ulery et al. Forensic Science International, Vol. 264, March, 2016
Paper - "Repeatability and Reproducibility of Decisions by Latent Fingerprint Examiners." Ulery et al. PLoS ONE, Vol. 7, No. 3, 2012
Paper - "Changes in latent fingerprint examiner' markup between Analysis and Comparison." Ulery et al. Forensic Science International, Vol. 247, 2014
Paper - "The forensic confirmation bias: Problems, perspectives and proposed solutions." Kassin et al. Journal of Applied Research in Memory and Cognition, Vol. 2, 2013
Paper – "Confirmation Bias, Ethics and Mistakes in Forensics," JFI,Vol. 56, No. 4, 2006

Revision 11 Issue Date: 12/15/2023

Page 103 of 109

Paper – Contextual bias and cross-contamination in the forensi investigations, plea bargains, trials and appeals." Law, Probabil	-
Module 28 Reading List: Analysis, Comparison, Evaluation,	and Verification (ACE-V)
ISPFS Latent Print Section AM – Friction Ridge Examination Me	ethodology
Guideline - SWGFAST Document 10 Standards for Examining For Resulting Conclusions (Latent/Tenprint) or the OSAC successor	
Registry Proposed - OSAC-2022-N-0038 Standard for Feature S Examination	Selectin in Friction Ridge
Registry Proposed - OSAC-2021-N-0020 Best Practice Recomm Examinations	endation for Limited Friction Ridge
ANSI/ASB Standard Recommendation 145, First Edition 2023 - Friction Ridge Examination	Standard for Consultation during
ANSI/ASB Best Practice Recommendation 144, First Edition 20 Recommendations for the Verification Component in Friction R	
ANSI/ASB Best Practice Recommendation 142, First Edition 20 Recommendations for the Resolution of Conflicts in Friction Ric	
Friction Ridge Skin - James F. Cowger	
Chapter 6 - The Basis for Comparison"	
Chapter 7 - Comparing Prints	
Chapter 8 - Some Comparisons of Evidence Prints	
Scott's Fingerprint Mechanics - Robert D. Olsen Sr. Pages 5-46,	171-175
Fingerprint Techniques - Andre A. Moenssens, Chapter 10 - Comparison of Fingerprints	
Palm Print Comparison Techniques course packet - Ron Smith	
Advances in Fingerprint Technology, 2^{nd} Edition - Lee & Gaenss Chapter 2 - Identification of Latent Prints	slen.
Latent Print Examiner Training Manual Page 104 of 109	Revision 11 Issue Date: 12/15/2023 Issuing Authority: Quality Manager

The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al.
Chapter 9 - Examination Process
Chapter 12 - Quality Assurance
Quantitative-Qualitative Friction Ridge Analysis - David R. Ashbaugh Chapters 4 - The Identification Process Chapter 5 - Poroscopy and Edgeoscopy
Analysis of Distortion in Latent Prints course packet – Alice Maceo
Fingerprints and Other Ridge Skin Impressions - Champod, et. al., Second Edition Chapter 2 – The Friction Ridge Identification Process
Paper - "Detection of Forged and Fabricated Latent Prints" Pat A. Wertheim, JFI Vol. 44, No. 6. 1994
Paper- "Fingerprints What They Can & Cannot Do!" Allan McRoberts The Print Vol. 10(6), June 1994 Pages 1-3
Paper - "The Ability Equation" Pat A. Wertheim
Paper - "Forensic Individualization of Images Using Quality and Quantity of Information." John Vanderkolk, JFI, Vol. 49. No. 3, 1999
Paper - "ACE-V and the Scientific Method." JFI Vol. 60 No.1, 2010
Paper – "The Investigation of the Reproducibility of Third-Level Characteristics," JFI Vol. 61, No.2, 2011.
Paper - "Scientific Comparison and Identification of Fingerprint Evidence." Pat. Wertheim. Fingerprint Whorld Vol. 26, No. 101, July 2000
Paper - "Distortion Versus Dissimilarity in Friction Skin Identification." William Leo. JFI, Vol. 48, No. 2, 1998
Paper - "A Performance Study of the ACE-V Process: A Pilot Study to Measure the Accuracy, Precision, Reproducibility, Repeatability, and Biasability of Conclusions Resulting from the ACE-V Process." JFI, Vol. 59, No. 2, 2009
Paper - "Incipient Ridges and the Clarity Spectrum" David R. Ashbaugh. JFI Vol.42. No. 2 1992
Paper - "Level 3 Details and Their Role in Fingerprint Identification: A Survey among Practitioners." JFI, Vol.58. No. 5, 2008

Revision 11

Issue Date: 12/15/2023

Page 105 of 109

Paper - "The Etiology of ACE-V and its Proper Use: An Exploration of the Relationship between ACE-V and the Scientific Method of Hypothesis Testing." JFI, Vol. 56 No. 3, 2006
Paper – "Palmar Flexion Crease Identification" David R. Ashbaugh Identification Canada Jan/Feb/March 1992
Paper – "Coins in the Pocket: A Simple Explanation of Quantitative – Qualitative Friction Ridge Analysis." JFI, Vol. 55, No. 3, 2005
Paper – "Assessing the Clarity of Friction Ridge Impressions." Forensic Science International, Vol.226, No. 1, 2012
Paper – "Fingerprint analysis for the determination of hand origin (right/left) using the axis slant in whorl patterns." Forensic Science Research, Vol.7, No. 2, 2022
Paper – "Testing the accuracy and reliability of palmar friction ridge comparisons – A black box study." Forensic Science International, 318:110457, 2021.
Module 29 Reading List: Case Management and Reporting for Comparison and/or MBIS
Latent Print Section Quality Manual - Casework Documentation and Report Writing ISPFS Quality/Procedure Manual Section on "Technical records" ———————————————————————————————————
Section on "Reporting of results" ASCLD/LAB-International Supplemental Requirements for the Accreditation of Forensic Science Testing Laboratories Appendix C- Latent Print Examination Records
Guideline - SWGFAST Document 8 Standard for the Documentation of Analysis, Comparison, Evaluation, and Verification (ACE-V) (Latent) or the OSAC successor document
Guideline - SWGFAST Document 5 Standard for Reporting Friction Ridge Examinations (Latent/Tenprint) or the OSAC successor document
The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al. Chapter 10 - Documentation of Friction Ridge Impressions from the Scene to the Conclusion
Latent Print Evaminer Training Manual Revision 11

Issue Date: 12/15/2023

Page 106 of 109

Latent Print Examination and Human Factors: Improving the The Report of the Expert Working Group on Human Factors i Chapter 5 - Reports and Documentation	9 - 1 - 1 -
Module 30 Reading List: Court Procedures, Related Laws, Civil Procedures Applicable to Latent Prints	, Expert Testimony, Criminal and
Guideline - SWGIT Section 17 Digital Imaging Technology Issusuccessor document	ues for the Courts or the OSAC
Friction Ridge Skin - James F. Cowger, Chapter 9 – Reporting and Testifying to Conclusions	Ser
Fingerprint Techniques - Andre A. Moenssens, Pages 270-280	
Fingerprints and the Law - Andre A. Moenssens Chapters 7 - Fingerprint Evidence in Criminal Cases Chapter 8 - Fingerprints in Non-Criminal Cases Chapter 9 - The Prosecutor's Approach to Fingerprint Evidence Chapter 10 - The Defense approach to Fingerprint Evidence Chapter 11 - The Fingerprint Witness in Court The Fingerprint Sourcebook by Scientific Working Group on I	X ·
Technology (SWGFAST), et al. Chapter 13 – Fingerprints and the Law	Friction Ridge Analysis, Study and
Law for the Expert Witness - Daniel A. Bronstein	
Advances in Fingerprint Technology, 2^{nd} Edition - Lee and Ga Chapter 10 – The Expert Fingerprint Witness	ensslen
Crime Scene Photography, 2 nd Edition – Robinson Chapter 12 – Legal Issues Related to Photographs and Digital	Images
Latent Print Examination and Human Factors: Improving the the Report of the Expert Working Group on Human Factors in Chapter 6 – Testimony	
National Commission of Forensic Science: Presentation of Exp Recommendations, 2012	pert Testimony Policy
Latent Print Examiner Training Manual Page 107 of 109	Revision 11 Issue Date: 12/15/2023 Issuing Authority: Quality Manager

All printed copies are uncontrolled

Paper – "Qualifying as an Expert Fingerprint Witness: Design Testimony." Pat A. Wertheim. JFI, Vol. 40, No. 2 1990	ing a Set of Questions to Assist in Court
Department of Justice Uniform Language for Testimony and Discipline 08/20. https://www.justice.gov/olp/page/file/10	=
Module 31 Reading List: Multimodal Biometric Identifica	ntion System (MBIS)
ISPFS Latent Print Section AM – MBIS	
The Fingerprint Sourcebook by Scientific Working Group on Technology (SWGFAST), et al.	Friction Ridge Analysis, Study and
Chapter 6 - Automated Fingerprint Identification System (AF	ris)
Criminalistics, 9 th edition - Richard Saferstein Chapter 14 - "AFIS" Pages 436-440	Keilis.
Advances in Fingerprint Technology $2^{\rm nd}$ edition Lee, Gaensslo Chapter 8 – Automated Fingerprint Identification and Imagin	
NEC – Integra-ID Integrated Biometric Workstation Latent U	ser Guide (current version)
NEC – IBW Latent Quick Reference (current version)	
NEC- IBW Application Keyboard Shortcuts (current version)	
NEC –Archive Quick Reference (current version)	
NEC WIN Best Practices for Latent Examiners V 1.0, 11/6/20	20
Universal Latent Workstation Training Guide ULW Version 6	.6.7, July 2020
Universal Latent Workstation (ULW) Version 6.6.7 Supplement successor document	ental Instructions October 2017 or its
Latent Print Examination and Human Factors: Improving the The Report of the Expert Working Group on Human Factors in Chapters 4 – Looking Ahead to Emerging and Improving Tech	in Latent Print Analysis 2012
PowerPoint "ULW-WEB"	
Latent Print Examiner Training Manual Page 108 of 109	Revision 11 Issue Date: 12/15/2023 Issuing Authority: Quality Manager

All printed copies are uncontrolled

Paper – "A Latent Print Examiner's Guide to IAFIS" JFI, Vol. 57, No. 4, 2007
Paper – "Determination of AFIS "sufficiency" in friction ridge examination" Forensic Science International, Vol. 263, 2016
Paper – "Fingerprint error rates on close non-matches" Journal of Forensic Sciences: 2021; 66(1): 129–134.
Module 32 DNA Database Comparison Training
Friction Ridge Skin, by James F. Cowger Pages 129-206.
Guideline - SWGFAST Document 10 Standards for Examining Friction Ridge Impressions and Resulting Conclusions (Latent/Tenprint) or the OSAC successor document
The Fingerprint Sourcebook by Scientific Working Group on Friction Ridge Analysis, Study and Technology (SWGFAST), et al. Chapter 9 - Examination Process
Forensic Pathways webpage/article – Confirmation Bias, Ethics, and Mistakes in Forensics "The eyes are not responsible when the mind does the seeing."
Latent Print Section AM - Friction Ridge Examination Methodology
Latent Print Section Quality Manual – Documentation and Report Writing Sections 9.8 and 9.9
Property of Idaho Inite