

**REVIEWED**

1/14/2019

By Britany Wylie at 1:46 pm, Jan 24, 2019

**Worklist: 2882**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
M2018-5719	2	135151	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6160	1	134751	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6171	1	134779	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6183	1	134810	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6185	1	134834	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6228	1	134926	AM 25 Blood Multi-Drug Screen by LC-QQQ	
M2018-6364	3	136184	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3423	1	133560	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3591	1	135464	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3604	2	137055	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3652	1	135713	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3653	1	135716	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3654	1	135763	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3655	1	135766	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2018-3656	1	135770	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0006	1	135807	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0007	1	135810	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0052	1	136169	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0053	1	136172	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0054	1	136187	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0060	1	136195	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0084	1	136298	AM 25 Blood Multi-Drug Screen by LC-QQQ	
P2019-0085	1	136301	AM 25 Blood Multi-Drug Screen by LC-QQQ	

B

**Worklist: 2882**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
P2019-0099	1	136337	AM 25 Blood Multi-Drug Screen by LC-QQQ
P2019-0103	1	136403	AM 25 Blood Multi-Drug Screen by LC-QQQ



B

# AM# 25: Multi-Drug Screen in Blood by LC-MS/MS

Extraction Date: 1/17/19  
Plate lot#: 0543908

Analyst: Sarah Pickle  
Plate Expiration: 11/28/19

**Mobile phase A:** 10mM Amm Form  
0.5M Ammonium Hydroxide

**Mobile phase B:** 0.1% Formic Acid in MeOH  
Ethyl Acetate LC Methanol

**Blank Blood Lot:** 445283-1

**Column:** Phenomenex Phenyl Hexyl (4.6x50mm, 2.6um)

**LCMS-QQQ ID:** 59740

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: Data Path Name: \_\_\_\_\_

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **250µL blood (calibrated pipette) Pipette ID: #16** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **250µL 0.5M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **300µL of blood+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate**.
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **900uL ethyl acetate**.
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% LC MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Open quantitation software and create a new quantitation batch.  
Batch Name: 011719 MDS SP  
Worklist path: MDS P1, MDS P2 final
- 2. Make necessary changes to integration limits
- 3. Evaluate samples, S/N of primary transition >5 and S/N of secondary transition >3 or evaluation of peak symmetry and resolution. Within +/- 2% or 0.1 min RT of administrative control. Calculated concentration >5.
- 4. Did all QCs pass for each analyte? Y / N Yes
- 5. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:



# Idaho State Police Forensic Services

P

---

## AM #25 Blood Multi-Drug Screen by LCMS-QQQ

---

### Methanol External Control Solution (Lot: 042618)

*100 ul of 1mg/mL stock was added of each drug to 9600 ul of LC MeOH.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	177145	
Oxycodone	Cerilliant	FE01181701	March 2022
Diazepam	Cerilliant	FE05201602	June 2021
Metoprolol	Cerilliant	FN06091510	July 2020
Meperidine	Cerilliant	FE01271602	June 2021
Prepared:	04/26/18		
Prepared By:	Tamara Salazar		
Expires:	04/26/19		

### Blood External Control Solution (Lot: WS122618)

*100 ul of methanol external control solution was added to 9900 ul of blood.*

*Approximately 50ng/mL of each compound.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-1
Methanol External Control Solution		042618
Prepared:	12/26/18	
Prepared by:	Tamara Salazar	
Expires:	04/26/19	

Compound Method		p1 cal 1a			p1 cal 1a Results			Qualifier 1 M...			Qualifier 1 Results			ISTD Method	
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	MI	Transition	MI	S/N	Area	Name	RT	Resp.	
6-MAM	328.2 -> 165...	1/17/2019 2:49 PM	4.732	2903	10.0000	52.79	328.2 -> 211...	345.28	1641	6-MAM-D6	4.712	101106			
7-Aminoclonazepam	286.1 -> 121...	1/17/2019 2:49 PM	5.451	97878	10.0000	381.31	286.1 -> 222...	2214.45	13891	7-Aminoclonazepam-D4	5.449	371678			
Benzoylcegonine	290.1 -> 168...	1/17/2019 2:49 PM	5.243	12916	10.0000	2048.53	290.1 -> 105...	3568.60	19683	Benzoylcegonine-D8	5.243	673248			
Metoprolol	268.2 -> 116...	1/17/2019 2:49 PM	5.290	67247	10.0000	16082.28	268.2 -> 56.3	236.04	71950	Noroxycodone-D3	4.714	366224			
a-hydroxyalprazolam	325.1 -> 297...	1/17/2019 2:49 PM	6.320	23161	10.0000	220.29	325.1 -> 215...	1054.25	6293	a-hydroxyalprazolam-D5	6.319	160668			
alpha-PVP	232.2 -> 91.0	1/17/2019 2:49 PM	5.336	145432	10.0000	80.47	232.2 -> 77.1	41.93	1015...	alpha-PVP-d8	5.334	1145303			
Alprazolam	309.1 -> 281...	1/17/2019 2:49 PM	6.425	254231	10.0000	26622.22	309.1 -> 205...	125.91	81536	Alprazolam-D5	6.404	1106673			
Amphetamine	136.1 -> 91.1	1/17/2019 2:49 PM	4.649	318079	10.0000	60.47	136.1 -> 119...	84.06	1483...	Amphetamine-D11	4.588	341050			
Cocaine	304.2 -> 182...	1/17/2019 2:49 PM	5.385	332060	10.0000	114859...	304.2 -> 105...	120.36	48132	Cocaine-D3	5.384	779394			
Bupropion	240.1 -> 184...	1/17/2019 2:49 PM	5.539	115044	10.0000	366.04	240.1 -> 131...	1011.85	2008...	Bupropion-D9	5.539	38711			
Carisoprodol	261.2 -> 176...	1/17/2019 2:49 PM	6.067	34938	10.0000	10659.87	261.2 -> 55.3	22.83	70574	Carisoprodol-D7	6.045	119721			
Citalopram	325.2 -> 109...	1/17/2019 2:49 PM	5.873	408602	10.0000	92830.80	325.2 -> 262...	10763...	52693	Citalopram-D6	5.872	526018			
Clonazepam	316.1 -> 269...	1/17/2019 2:49 PM	6.281	96876	10.0000	205.42	316.1 -> 213...	241.68	17704	Clonazepam-D4	6.260	93026			
Codeine	300.2 -> 214...	1/17/2019 2:49 PM	4.596	18986	10.0000	69883.75	300.2 -> 165...	18495...	21880	Codeine-D6	4.572	107773			
Cyclobenzaprine	276.2 -> 215...	1/17/2019 2:49 PM	6.104	113453	10.0000	15582.84	276.2 -> 231...	2096.99	57276	Cyclobenzaprine-D3	6.104	746304			
Dextromethorphan	272.2 -> 171...	1/17/2019 2:49 PM	5.842	158330	10.0000	1134.06	272.2 -> 147...	17452...	1357...	Dextromethorphan-D3	5.841	407449			
Dextroproprian	258.2 -> 157...	1/17/2019 2:49 PM	5.231	196948	10.0000	24961.23	258.2 -> 133...	351.54	97884	Dextroproprian-D3	5.230	584230			
Diazepam	285.1 -> 193...	1/17/2019 2:49 PM	6.674	136800	10.0000	75.94	285.1 -> 153...	176.65	1860...	Diazepam-D5	6.652	612482			
Dihydrocodeine	302.2 -> 198...	1/17/2019 2:49 PM	4.534	62031	10.0000	145.27	302.2 -> 128...	162.02	52837	Dihydrocodeine-D6	4.531	310883			
Diphenhydramine	256.2 -> 167...	1/17/2019 2:49 PM	5.804	168127	10.0000	26.05	256.2 -> 152...	Infinity	1597...	Diphenhydramine-D3	5.824	1464455			
Doxylamine	271.2 -> 182...	1/17/2019 2:49 PM	5.432	210324	10.0000	85.94	271.2 -> 167...	87.71	3456...	Doxylamine-D5	5.452	2226142			
EDDP	279.2 -> 235...	1/17/2019 2:49 PM	5.860	749332	10.0000	133598...	279.2 -> 250...	935.74	3300...	EDDP-D3	5.840	4252572			
Fentanyl	337.2 -> 105...	1/17/2019 2:49 PM	5.790	23840	10.0000	28.29	337.2 -> 188...	1990.05	8120	Fentanyl-D5	5.789	2095143			
Fluoxetine	310.1 -> 44.3	1/17/2019 2:49 PM	6.099	859249	10.0000	313.21	310.1 -> 148...	85.62	10059	Fluoxetine-D6	6.098	2938287			
Hydrocodone	300.2 -> 198...	1/17/2019 2:49 PM	4.798	66817	10.0000	356.72	300.2 -> 128...	Infinity	53290	Hydrocodone-D6	4.795	479533			
Hydromorphone	286.2 -> 184...	1/17/2019 2:49 PM	4.175	55866	10.0000	Infinity	286.2 -> 157...	Infinity	59688	Hydromorphone-D6	4.174	236495			
Ketamine	238.1 -> 125...	1/17/2019 2:49 PM	5.193	311419	10.0000	106.70	238.1 -> 220...	19.69	49120	Ketamine-D4	5.191	1193047			
Mirtazapine	266.2 -> 195...	1/17/2019 2:49 PM	5.536	361012	10.0000	1129.78	266.2 -> 72.3	95.41	1786...	Tramadol-13C-D3	5.248	3118507			
Mirtagryline	399.2 -> 238...	1/17/2019 2:49 PM	5.910	15553	10.0000	2447.75	399.2 -> 174...	190.15	1604...	Tramadol-13C-D3	5.248	3118507			
Meprobamate	219.1 -> 158...	1/17/2019 2:49 PM	5.541	30079	10.0000	41.75	219.1 -> 97.1	40.15	17759	Meprobamate-D7	5.541	6597			
Methadone	310.2 -> 265...	1/17/2019 2:49 PM	6.119	368995	10.0000	48913.37	310.2 -> 105...	2187.00	3279...	Methadone-D9	6.118	1167840			
Methamphetamine	150.1 -> 91.1	1/17/2019 2:49 PM	4.789	328335	10.0000	34.35	150.1 -> 119...	5.69	1053...	Methamphetamine-D11	4.748	785367			
Morphine	286.2 -> 201...	1/17/2019 2:49 PM	3.916	12844	10.0000	Infinity	286.2 -> 165...	19.81	9887	Morphine-D6	3.895	24798			
Norbuprenorphine	414.3 -> 101...	1/17/2019 2:49 PM	5.665	13137	10.0000	88.93	414.3 -> 57.3	40.15	11997	Norbuprenorphine-D3	5.644	23162			
Nordiazepam	271.1 -> 140...	1/17/2019 2:49 PM	6.534	187126	10.0000	64.11	271.1 -> 208...	91.92	46592	Nordiazepam-D5	6.512	426938			
Norfentanyl	233.2 -> 84.2	1/17/2019 2:49 PM	5.173	340935	10.0000	970.61	233.2 -> 55.3	23367...	1414...	Norfentanyl-D5	5.152	1480218			
Norhydrocodone	286.2 -> 199...	1/17/2019 2:49 PM	4.760	11030	10.0000	2000.65	286.2 -> 171...	1383.95	6117	Norhydrocodone-D3	4.760	127597			
Noroxycodone	302.1 -> 284...	1/17/2019 2:49 PM	4.716	72642	10.0000	29.46	302.1 -> 186...	72.94	18541	Noroxycodone-D3	4.714	366224			
O-desmethyl-tramadol	250.2 -> 58.3	1/17/2019 2:49 PM	4.745	549330	10.0000	79.70	250.2 -> 42.3	43.91	18695	O-desmethyl-tramadol...	4.766	1334362			
Oxazepam	287.1 -> 240...	1/17/2019 2:49 PM	6.326	96599	10.0000	Infinity	287.1 -> 268...	6.24	76524	Oxazepam-D5	6.326	549976			
Oxycodone	316.2 -> 298...	1/17/2019 2:49 PM	4.733	197935	10.0000	134.96	316.2 -> 241...	94.14	46939	Oxycodone-D6	4.712	494956			
Oxymorphone	302.1 -> 284...	1/17/2019 2:49 PM	4.716	72642	10.0000	29.46	302.1 -> 227...	30.57	9531	Oxymorphone-D3	4.714	366224			
Phentermine	150.1 -> 65.2	1/17/2019 2:49 PM	4.952	55075	10.0000	Infinity	150.1 -> 133...	39.93	5010	Phentermine-D5	4.951	795320			
Promethazine	285.1 -> 86.2	1/17/2019 2:49 PM	6.062	449242	10.0000	912.55	285.1 -> 71.3	214.52	93323	Promethazine-D3	6.061	1259711			
Quetiapine	384.2 -> 253...	1/17/2019 2:49 PM	6.053	303605	10.0000	112433...	384.2 -> 221...	33810...	1632...	Quetiapine-D8	6.011	369284			
Setraline	306.1 -> 158...	1/17/2019 2:49 PM	6.304	329894	10.0000	494.88	306.1 -> 275...	29.73	1020...	Setraline-D3	6.303	872047			
Temazepam	301.1 -> 255...	1/17/2019 2:49 PM	6.468	273561	10.0000	1910.40	301.1 -> 283...	Infinity	83560	Temazepam-D5	6.467	1697586			
Tramadol	264.2 -> 58.3	1/17/2019 2:49 PM	5.271	626610	10.0000	169.39	264.2 -> 43.3	6.82	7687	Tramadol-13C-D3	5.248	3118507			

Compound Name	Transition	Acq. Date/Time	RT	Resp.	Calc. Conc.	S/N	MI	Qualifier 1 M...	Qualifier 1 Results	ISTD Method	ISTD Results
Trazodone	372.2 -> 176...	1/17/2019 2:49 PM	5.972	520650	10.0000	86032.46		372.2 -> 148...	2506.63	Trazodone-D6	5.952 1108511
Venlafaxine	278.2 -> 58.3	1/17/2019 2:49 PM	5.616	669439	10.0000	115.82		278.2 -> 260...	1421.74	Venlafaxine-D6	5.615 868241
Zolpidem	308.2 -> 235...	1/17/2019 2:49 PM	5.712	453328	10.0000	660.76		308.2 -> 263...	12678...	Zolpidem-D6	5.711 2826537

12

Compound Method		p1 external control			p1 external control Results			Qualifier 1 M...			Qualifier 1 Results			ISTD Method	
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	M/I	Transition	M/I	S/N	Area	Name	RT	Resp.	
6-MAM	328.2 -> 165...	1/17/2019 3:28 PM						328.2 -> 211...				6-MAM-D6	4.712	89012	
7-Aminoclonazepam	286.1 -> 121...	1/17/2019 3:28 PM						286.1 -> 222...				7-Aminoclonazepam-D4	5.449	332144	
Benzoylclonazepam	290.1 -> 168...	1/17/2019 3:28 PM						290.1 -> 105...				Benzoylclonazepam-D8	5.243	46702	
Metoprolol	268.2 -> 116...	1/17/2019 3:28 PM	5.269	697753	103.5293	7687.27		268.2 -> 56.3		931.42	664109	Noroxycodone-D3	4.693	314435	
a-hydroxyalprazolam	325.1 -> 297...	1/17/2019 3:28 PM						325.1 -> 215...				a-hydroxyalprazolam-D5	6.319	133515	
alpha-PVP	232.2 -> 91.0	1/17/2019 3:28 PM						232.2 -> 77.1				alpha-PVP-d8	5.334	976080	
Alprazolam	309.1 -> 281...	1/17/2019 3:28 PM						309.1 -> 205...				Alprazolam-D5	6.404	913048	
Amphetamine	136.1 -> 91.1	1/17/2019 3:28 PM						136.1 -> 119...				Amphetamine-D11	4.588	323318	
Cocaine	304.2 -> 182...	1/17/2019 3:28 PM						304.2 -> 105...				Cocaine-D3	5.384	652783	
Bupropion	240.1 -> 184...	1/17/2019 3:28 PM						240.1 -> 131...				Bupropion-D9	5.539	31186	
Carisoprodol	261.2 -> 176...	1/17/2019 3:28 PM						261.2 -> 55.3				Carisoprodol-D7	6.045	108517	
Citalopram	325.2 -> 109...	1/17/2019 3:28 PM						325.2 -> 262...				Citalopram-D6	5.872	446734	
Clonazepam	316.1 -> 269...	1/17/2019 3:28 PM						316.1 -> 213...				Clonazepam-D4	6.260	87160	
Codeine	300.2 -> 215...	1/17/2019 3:28 PM						300.2 -> 165...				Codeine-D6	4.572	103486	
Cyclobenzaprine	276.2 -> 214...	1/17/2019 3:28 PM						276.2 -> 231...				Cyclobenzaprine-D3	6.104	615789	
Dextromethorphan	272.2 -> 171...	1/17/2019 3:28 PM						272.2 -> 147...				Dextromethorphan-D3	5.841	382166	
Dextroproprian	258.2 -> 157...	1/17/2019 3:28 PM						258.2 -> 133...				Dextroproprian-D3	5.210	493673	
Diazepam	285.1 -> 193...	1/17/2019 3:28 PM	6.674	12643...	116.3790	1016625...		285.1 -> 153...		165294...	17854...	Diazepam-D5	6.673	486423	
Dihydrocodeine	302.2 -> 198...	1/17/2019 3:28 PM						302.2 -> 128...				Dihydrocodeine-D6	4.531	260986	
Diphenhydramine	256.2 -> 167...	1/17/2019 3:28 PM						256.2 -> 152...				Diphenhydramine-D3	5.803	12743...	
Doxylamine	271.2 -> 182...	1/17/2019 3:28 PM						271.2 -> 167...				Doxylamine-D5	5.432	19016...	
EDDP	279.2 -> 235...	1/17/2019 3:28 PM						279.2 -> 250...				EDDP-D3	5.840	37703...	
Fentanyl	337.2 -> 105...	1/17/2019 3:28 PM						337.2 -> 188...				Fentanyl-D5	5.789	17886...	
Fluoxetine	310.1 -> 44.3	1/17/2019 3:28 PM						310.1 -> 148...				Fluoxetine-D6	6.098	26680...	
Hydrocodone	300.2 -> 198...	1/17/2019 3:28 PM						300.2 -> 128...				Hydrocodone-D6	4.774	421464	
Hydromorphone	286.2 -> 184...	1/17/2019 3:28 PM						286.2 -> 157...				Hydromorphone-D6	4.153	187934	
Ketamine	238.1 -> 125...	1/17/2019 3:28 PM						238.1 -> 220...				Ketamine-D4	5.191	10942...	
Mirtazapine	266.2 -> 195...	1/17/2019 3:28 PM						266.2 -> 72.3				Tramadol-13C-D3	5.248	27781...	
Mirtazapine	399.2 -> 238...	1/17/2019 3:28 PM						399.2 -> 174...				Tramadol-13C-D3	5.248	27781...	
Meperbamate	219.1 -> 158...	1/17/2019 3:28 PM						219.1 -> 97.1				Meperbamate-D7	5.561	2563	
Methadone	310.2 -> 265...	1/17/2019 3:28 PM						310.2 -> 105...				Methadone-D9	6.118	943137	
Methamphetamine	150.1 -> 91.1	1/17/2019 3:28 PM						150.1 -> 119...				Methamphetamine-D11	4.748	771855	
Morphine	286.2 -> 201...	1/17/2019 3:28 PM						286.2 -> 165...				Morphine-D6	3.875	19075	
Norbuprenorphine	414.3 -> 101...	1/17/2019 3:28 PM						414.3 -> 57.3				Norbuprenorphine-D3	5.624	23635	
Nordiazepam	271.1 -> 140...	1/17/2019 3:28 PM						271.1 -> 208...				Nordiazepam-D5	6.492	253971	
Nortentanyl	233.2 -> 84.2	1/17/2019 3:28 PM						233.2 -> 55.3				Nortentanyl-D5	5.152	12361...	
Norhydrocodone	286.2 -> 199...	1/17/2019 3:28 PM						286.2 -> 171...				Norhydrocodone-D3	4.760	119201	
Noroxycodone	302.1 -> 284...	1/17/2019 3:28 PM						302.1 -> 186...				Noroxycodone-D3	4.693	314435	
O-desmethyl-tramadol	250.2 -> 58.3	1/17/2019 3:28 PM						250.2 -> 42.3				O-desmethyl-tramadol...	4.745	10898...	
Oxazepam	287.1 -> 240...	1/17/2019 3:28 PM						287.1 -> 268...				Oxazepam-D5	6.326	523121	
Oxycodone	316.2 -> 298...	1/17/2019 3:28 PM	4.733	18934...	111.4650	334.55		316.2 -> 241...		236.82	459873	Oxycodone-D6	4.712	424786	
Oxymorphone	302.1 -> 284...	1/17/2019 3:28 PM						302.1 -> 227...				Oxymorphone-D3	4.693	314435	
Phentermine	150.1 -> 65.2	1/17/2019 3:28 PM						150.1 -> 133...				Phentermine-D5	4.951	755627	
Promethazine	285.1 -> 86.2	1/17/2019 3:28 PM						285.1 -> 71.3				Promethazine-D3	6.061	10522...	
Quetiapine	384.2 -> 253...	1/17/2019 3:28 PM						384.2 -> 221...				Quetiapine-D8	6.032	341135	
Seftalaine	306.1 -> 158...	1/17/2019 3:28 PM						306.1 -> 275...				Seftalaine-D3	6.303	699700	
Temazepam	301.1 -> 255...	1/17/2019 3:28 PM						301.1 -> 283...				Temazepam-D5	6.467	12409...	
Tramadol	264.2 -> 58.3	1/17/2019 3:28 PM						264.2 -> 43.3				Tramadol-13C-D3	5.248	27781...	



Compound Method		p1 negative		p1 negative Results			ISTD Method					
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	MI	Qualifier 1 M...	Qualifier 1 R...	Name	RT	Resp.
6-MAM	328.2 -> 165...	1/17/2019 3:09 PM						328.2 -> 211...		6-MAM-D6	4.712	87653
7-Aminoclonazepam	286.1 -> 121...	1/17/2019 3:09 PM						286.1 -> 222...		7-Aminoclonazepam-D4	5.449	332294
Benzoylcegonine	290.1 -> 168...	1/17/2019 3:09 PM						290.1 -> 105...		Benzoylcegonine-D8	5.243	51462
Metoprolol	268.2 -> 116...	1/17/2019 3:09 PM						268.2 -> 56.3		Noroxycodone-D3	4.693	337871
a-hydroxyalprazolam	325.1 -> 297...	1/17/2019 3:09 PM						325.1 -> 215...		a-hydroxyalprazolam-D5	6.299	148208
alpha-PVP	232.2 -> 91.0	1/17/2019 3:09 PM						232.2 -> 77.1		alpha-PVP-d8	5.334	10652...
Alprazolam	309.1 -> 281...	1/17/2019 3:09 PM						309.1 -> 205...		Alprazolam-D5	6.404	912386
Amphetamine	136.1 -> 91.1	1/17/2019 3:09 PM						136.1 -> 119...		Amphetamine-D11	4.588	336379
Cocaine	304.2 -> 182...	1/17/2019 3:09 PM						304.2 -> 105...		Cocaine-D3	5.384	741987
Bupropion	240.1 -> 184...	1/17/2019 3:09 PM						240.1 -> 131...		Bupropion-D9	5.518	34114
Carisoprodol	261.2 -> 176...	1/17/2019 3:09 PM						261.2 -> 55.3		Carisoprodol-D7	6.045	113144
Citalopram	325.2 -> 109...	1/17/2019 3:09 PM						325.2 -> 262...		Citalopram-D6	5.851	463398
Clonazepam	316.1 -> 269...	1/17/2019 3:09 PM						316.1 -> 213...		Clonazepam-D4	6.281	93638
Codeine	300.2 -> 215...	1/17/2019 3:09 PM						300.2 -> 165...		Codeine-D6	4.572	96932
Cyclobenzaprine	276.2 -> 214...	1/17/2019 3:09 PM						276.2 -> 231...		Cyclobenzaprine-D3	6.104	691725
Dextromethorphan	272.2 -> 171...	1/17/2019 3:09 PM						272.2 -> 147...		Dextromethorphan-D3	5.841	383670
Dextroproprian	258.2 -> 157...	1/17/2019 3:09 PM						258.2 -> 133...		Dextroproprian-D3	5.230	512787
Diazepam	285.1 -> 193...	1/17/2019 3:09 PM						285.1 -> 153...		Diazepam-D5	6.652	507095
Dihydrocodeine	302.2 -> 198...	1/17/2019 3:09 PM						302.2 -> 128...		Dihydrocodeine-D6	4.531	272693
Diphenhydramine	256.2 -> 167...	1/17/2019 3:09 PM						256.2 -> 152...		Diphenhydramine-D3	5.803	12919...
Doxylamine	271.2 -> 182...	1/17/2019 3:09 PM						271.2 -> 167...		Doxylamine-D5	5.432	20038...
EDDP	279.2 -> 235...	1/17/2019 3:09 PM						279.2 -> 250...		EDDP-D3	5.840	39048...
Fentanyl	337.2 -> 105...	1/17/2019 3:09 PM						337.2 -> 188...		Fentanyl-D5	5.789	18933...
Fluoxetine	310.1 -> 44.3	1/17/2019 3:09 PM						310.1 -> 148...		Fluoxetine-D6	6.098	27782...
Hydrocodone	300.2 -> 198...	1/17/2019 3:09 PM						300.2 -> 128...		Hydrocodone-D6	4.774	439860
Hydromorphone	286.2 -> 184...	1/17/2019 3:09 PM						286.2 -> 157...		Hydromorphone-D6	4.174	205453
Ketamine	238.1 -> 125...	1/17/2019 3:09 PM						238.1 -> 220...		Ketamine-D4	5.191	11469...
Mirtazapine	266.2 -> 195...	1/17/2019 3:09 PM						266.2 -> 72.3		Tramadol-13C-D3	5.248	28968...
Mirtazapine	399.2 -> 238...	1/17/2019 3:09 PM						399.2 -> 174...		Tramadol-13C-D3	5.248	28968...
Meprobamate	219.1 -> 158...	1/17/2019 3:09 PM						219.1 -> 97.1		Meprbamate-D7	5.541	4345
Methadone	310.2 -> 265...	1/17/2019 3:09 PM						310.2 -> 105...		Methadone-D9	6.118	10116...
Methamphetamine	150.1 -> 91.1	1/17/2019 3:09 PM						150.1 -> 119...		Methamphetamine-D11	4.748	772915
Morphine	286.2 -> 201...	1/17/2019 3:09 PM						286.2 -> 165...		Morphine-D6	3.895	20407
Norbuprenorphine	414.3 -> 101...	1/17/2019 3:09 PM						414.3 -> 57.3		Norbuprenorphine-D3	5.644	27920
Nordiazepam	271.1 -> 140...	1/17/2019 3:09 PM						271.1 -> 208...		Nordiazepam-D5	6.512	258691
Norfentanyl	233.2 -> 84.2	1/17/2019 3:09 PM						233.2 -> 55.3		Norfentanyl-D5	5.152	13216...
Norhydrocodone	286.2 -> 199...	1/17/2019 3:09 PM						286.2 -> 171...		Norhydrocodone-D3	4.760	123484
Noroxycodone	302.1 -> 284...	1/17/2019 3:09 PM						302.1 -> 186...		Noroxycodone-D3	4.693	337871
O-desmethyl-tramadol	250.2 -> 58.3	1/17/2019 3:09 PM						250.2 -> 42.3		O-desmethyl-tramadol...	4.726	11632...
Oxazepam	287.1 -> 240...	1/17/2019 3:09 PM						287.1 -> 268...		Oxazepam-D5	6.326	507886
Oxycodone	316.2 -> 298...	1/17/2019 3:09 PM						316.2 -> 241...		Oxycodone-D6	4.712	457125
Oxymorphone	302.1 -> 284...	1/17/2019 3:09 PM						302.1 -> 227...		Oxymorphone-D3	4.693	337871
Phentermine	150.1 -> 65.2	1/17/2019 3:09 PM						150.1 -> 133...		Phentermine-D5	4.931	779350
Promethazine	285.1 -> 86.2	1/17/2019 3:09 PM						285.1 -> 71.3		Promethazine-D3	6.061	11456...
Quetiapine	384.2 -> 253...	1/17/2019 3:09 PM						384.2 -> 221...		Quetiapine-D8	6.011	364687
Sepraline	306.1 -> 158...	1/17/2019 3:09 PM						306.1 -> 275...		Sepraline-D3	6.303	782260
Temazepam	301.1 -> 255...	1/17/2019 3:09 PM						301.1 -> 283...		Temazepam-D5	6.467	12578...
Tramadol	264.2 -> 58.3	1/17/2019 3:09 PM						264.2 -> 43.3		Tramadol-13C-D3	5.248	28968...

Compound Method	p1 negative		p1 negative Results				Qualifier 1 M...	Qualifier 1 R...	ISTD Method	ISTD Results				
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	MI	Transition	MI	S/N	Area	Name	RT	Resp.
Trazodone	372.2 -> 176...	1/17/2019 3:09 PM						372.2 -> 148...				Trazodone-D6	5.952	10234...
Venlafaxine	278.2 -> 58.3	1/17/2019 3:09 PM						278.2 -> 260...				Venlafaxine-D6	5.615	815894
Zolpidem	308.2 -> 235...	1/17/2019 3:09 PM						308.2 -> 263...				Zolpidem-D6	5.690	25296...

*P*

Compound Method		p2 cal 1b		p2 cal 1b Results		Qualifier 1 M...		Qualifier 1 Results		ISTD Method		ISTD Results		
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	MI	Transition	MI	S/N	Area	Name	RT	Resp
Acetylfentanyl	323.2 -> 188...	1/17/2019 2:59 PM	5.578	11723	10.0000	33.94	323.2 -> 105...	63.73	29315	Acetylfentanyl-D5	5.577	2044140		
alpha-hydroxymidazola...	342.1 -> 324...	1/17/2019 2:59 PM	6.367	99315	10.0000	30441.25	342.1 -> 203...	16607.37	68556	alpha-hydroxymidazola...	6.366	932797		
Amiripivline	278.2 -> 91.2	1/17/2019 2:59 PM	6.173	478966	10.0000	Infinity	278.2 -> 105...	Infinity	1656...	Amiripivline-D3	6.172	1101605		
Acetylornofentanyl	219.1 -> 84.2	1/17/2019 2:59 PM	4.694	15126	10.0000	17024.44	219.1 -> 56.3	36.33	6895	Acetylornofentanyl-D5	4.693	873818		
Carbamazepine	237.1 -> 193...	1/17/2019 2:59 PM	6.059	68389	10.0000	639.45	237.1 -> 194...	153620...	4361...	Carbamazepine-13C6	6.057	2714705		
Chlordiazepoxide	300.1 -> 227...	1/17/2019 2:59 PM	6.514	79613	10.0000	79.00	300.1 -> 282...	Infinity	1218	Chlordiazepoxide-D5	6.494	1367690		
Sufentanil	387.2 -> 238...	1/17/2019 2:59 PM	6.041	29792	10.0000	611.06	387.2 -> 111...	35.50	17250	Sufentanil-D5	6.019	2137108		
Desipramine	267.2 -> 72.3	1/17/2019 2:59 PM	6.134	609894	10.0000	384.11	267.2 -> 44.3	58302.27	2843...	Desipramine-D3	6.134	3654371		
Doxepin	280.2 -> 107...	1/17/2019 2:59 PM	5.929	132026	10.0000	87.35	280.2 -> 77.2	Infinity	81392	Doxepin-D3	5.928	423349		
Estazolam	295.1 -> 267...	1/17/2019 2:59 PM	6.333	201971	10.0000	294.17	295.1 -> 205...	29547.20	1075...	Estazolam-D5	6.331	1367476		
Flunitrazepam	314.1 -> 268...	1/17/2019 2:59 PM	6.390	189078	10.0000	Infinity	314.1 -> 239...	20096.03	62960	Flunitrazepam-D7	6.370	80350		
Flurazepam	388.2 -> 315...	1/17/2019 2:59 PM	5.878	365505	10.0000	946.27	388.2 -> 317...	19.43	48016	Flurazepam-D4	5.614	325265		
Imipramine	281.2 -> 86.2	1/17/2019 2:59 PM	6.131	638286	10.0000	167363...	281.2 -> 58.3	Infinity	3742...	Imipramine-D3	6.130	2413283		
Lepamisole	205.1 -> 91.2	1/17/2019 2:59 PM	4.719	110837	10.0000	211.46	205.1 -> 178...	35.57	93466	Sufentanil-D5	6.019	2137108		
Mepricitline	278.2 -> 91.2	1/17/2019 2:59 PM	6.173	478966	10.0000	Infinity	278.2 -> 117...	28697.10	1054...	Mepricitline-D3	6.172	1101605		
MDA	180.1 -> 163...	1/17/2019 2:59 PM	4.802	128714	10.0000	Infinity	180.1 -> 105...	58.58	47154	MDA-D5	4.802	433328		
MDEA	208.1 -> 163...	1/17/2019 2:59 PM	5.063	241593	10.0000	Infinity	208.1 -> 105...	41.28	1334...	MDEA-D6	5.063	949068		
MDMA	194.1 -> 163...	1/17/2019 2:59 PM	4.902	273949	10.0000	Infinity	194.1 -> 105...	148.08	1157...	MDMA-D6	4.902	83308		
Meperidine	248.2 -> 220...	1/17/2019 2:59 PM	5.402	143280	10.0000	288.68	248.2 -> 174...	307.48	95295	Meperidine-D4	5.401	670032		
Methylphenidate	234.2 -> 84.2	1/17/2019 2:59 PM	5.325	589301	10.0000	119.93	234.2 -> 56.3	Infinity	1500...	Meperidine-D4	5.401	670032		
Midazolam	326.1 -> 223...	1/17/2019 2:59 PM	6.409	70122	10.0000	339.79	326.1 -> 249...	484.74	88471	Midazolam-D4	6.408	958460		
Normeperidine	234.1 -> 42.3	1/17/2019 2:59 PM	5.427	146590	10.0000	Infinity	234.1 -> 160...	31895.02	85469	Normeperidine-D4	5.425	651206		
7-Aminoflunitrazepam	284.1 -> 135...	1/17/2019 2:59 PM	5.643	214889	10.0000	66356.29	284.1 -> 93.2	8280.95	22314	7-aminoflunitrazepam...	5.622	1078915		
Nortriptyline	264.2 -> 91.2	1/17/2019 2:59 PM	6.175	198232	10.0000	47448.46	264.2 -> 233...	159.74	96569	Nortriptyline-D3	6.175	799954		
Phenazepam	349.0 -> 206...	1/17/2019 2:59 PM	6.446	81545	10.0000	25646.38	349.0 -> 183...	19992.42	66970	Phenazepam-D4	6.446	438004		
Phencyclidine	244.2 -> 86.2	1/17/2019 2:59 PM	5.710	155881	10.0000	37591.41	244.2 -> 91.1	Infinity	2183...	Phencyclidine-D5	5.689	1133508		
Phenylethylamine	253.1 -> 104...	1/17/2019 2:59 PM	5.973	6422	10.0000	1575.83	253.1 -> 182...	679.92	1845	Phenylethylamine-D5	5.973	30600		
Pseudoephedrine	166.1 -> 148...	1/17/2019 2:59 PM	4.422	4562202	10.0000	Infinity	166.1 -> 133...	Infinity	9673...	Pseudoephedrine-D3	4.421	205583...		
Tapentadol	222.2 -> 107...	1/17/2019 2:59 PM	5.286	281624	10.0000	237.28	222.2 -> 121...	88.15	1245...	Tapentadol-D3	5.285	1083468		
Zaleplon	306.1 -> 236...	1/17/2019 2:59 PM	6.168	172698	10.0000	36431.83	306.1 -> 264...	209.59	1754...	Zaleplon-D4	6.168	522461		
Zopiclone	389.1 -> 244...	1/17/2019 2:59 PM	5.615	65826	10.0000	442.69	389.1 -> 216...	57.16	25591	Zopiclone-D4	5.614	325265		

Compound Method	p2 external control	p2 external control Results					Qualifier 1 M...	Qualifier 1 Results	ISTD Method	ISTD Results				
Name	Transition	Acq. Date/Time	RT	Resp	Calc. Conc.	S/N	MI	Transition	MI	S/N	Area	Name	RT	Resp.
Acetylferfentanyl	323.2 -> 188...	1/17/2019 3:38 PM						323.2 -> 105...				Acetylferfentanyl-D5	5.577	1651125
alpha-hydroxymidazolam	342.1 -> 324...	1/17/2019 3:38 PM						342.1 -> 203...				alpha-hydroxymidazolam-...	6.366	793689
Amiripryline	278.2 -> 91.2	1/17/2019 3:38 PM						278.2 -> 105...				Amiripryline-D3	6.172	946935
Acetyl norfentanyl	219.1 -> 84.2	1/17/2019 3:38 PM						219.1 -> 56.3				Acetyl norfentanyl-D5	4.693	755901
Carbamazepine	237.1 -> 193...	1/17/2019 3:38 PM						237.1 -> 194...				Carbamazepine-13C6	6.078	2388273
Chlordiazepoxide	300.1 -> 227...	1/17/2019 3:38 PM						300.1 -> 282...				Chlordiazepoxide-D5	6.494	1231680
Sufentanil	387.2 -> 238...	1/17/2019 3:38 PM						387.2 -> 111...				Sufentanil-D5	6.019	1779412
Desipramine	267.2 -> 72.3	1/17/2019 3:38 PM						267.2 -> 44.3				Desipramine-D3	6.134	3229835
Doxepin	280.2 -> 107...	1/17/2019 3:38 PM						280.2 -> 77.2				Doxepin-D3	5.928	345936
Estazolam	295.1 -> 267...	1/17/2019 3:38 PM						295.1 -> 205...				Estazolam-D5	6.331	1234028
Flunitrazepam	314.1 -> 268...	1/17/2019 3:38 PM						314.1 -> 239...				Flunitrazepam-D7	6.370	77612
Flurazepam	388.2 -> 315...	1/17/2019 3:38 PM						388.2 -> 317...				Zopiclone-D4	5.614	190584
Imipramine	281.2 -> 86.2	1/17/2019 3:38 PM						281.2 -> 58.3				Imipramine-D3	6.130	2098365
Levamisole	205.1 -> 91.2	1/17/2019 3:38 PM						205.1 -> 178...				Sufentanil-D5	6.019	1779412
Maprotiline	278.2 -> 91.2	1/17/2019 3:38 PM						278.2 -> 117...				Amiripryline-D3	6.172	946935
MDA	180.1 -> 163...	1/17/2019 3:38 PM						180.1 -> 105...				MDA-D5	4.802	412413
MDEA	208.1 -> 163...	1/17/2019 3:38 PM						208.1 -> 105...				MDEA-D6	5.063	788378
MDMA	194.1 -> 163...	1/17/2019 3:38 PM						194.1 -> 105...				MDMA-D6	4.882	77284
Meperidine	248.2 -> 220...	1/17/2019 3:38 PM	5.402	1313149	113.2907	293.48		248.2 -> 174...		Infini...	9664...	Meperidine-D4	5.401	542038
Methylphenidate	234.2 -> 84.2	1/17/2019 3:38 PM						234.2 -> 56.3				Meperidine-D4	5.401	542038
Midazolam	326.1 -> 223...	1/17/2019 3:38 PM						326.1 -> 249...				Midazolam-D4	6.408	825735
Normeperidine	234.1 -> 42.3	1/17/2019 3:38 PM						234.1 -> 160...				Normeperidine-D4	5.425	566616
7-Aminoflunitrazepam	284.1 -> 135...	1/17/2019 3:38 PM						284.1 -> 93.2				7-Aminoflunitrazepam-D7	5.622	893238
Notriptyline	264.2 -> 91.2	1/17/2019 3:38 PM						264.2 -> 233...				Notriptyline-D3	6.175	655883
Phenazepam	349.0 -> 206...	1/17/2019 3:38 PM						349.0 -> 183...				Phenazepam-D4	6.446	403975
Phencyclidine	244.2 -> 86.2	1/17/2019 3:38 PM						244.2 -> 91.1				Phencyclidine-D5	5.689	958420
Phenytoln	253.1 -> 104...	1/17/2019 3:38 PM						253.1 -> 182...				Phenytoln-D10	5.952	24061
Pseudoephedrine	166.1 -> 148...	1/17/2019 3:38 PM						166.1 -> 133...				Pseudoephedrine-D3	4.421	188172...
Tapentadol	222.2 -> 107...	1/17/2019 3:38 PM						222.2 -> 121...				Tapentadol-D3	5.285	917000
Zaleplon	306.1 -> 236...	1/17/2019 3:38 PM						306.1 -> 264...				Zaleplon-D4	6.168	485359
Zopiclone	389.1 -> 244...	1/17/2019 3:38 PM						389.1 -> 216...				Zopiclone-D4	5.614	190584

1

Compound Method	p2 negative			p2 negative Results			ISTD Method			ISTD Results				
Name	Transition	Acq. Date-Time	RT	Resp.	Calc. Conc.	S/N	MI	Qualifier 1 M...	MI	S/N	Area	Name	RT	Resp.
Acetylfentanyl	323.2 -> 188...	1/17/2019 3:18 PM						323.2 -> 105...				Acetylfentanyl-D5	5.577	1682114
alpha-hydroxymidazolam	342.1 -> 324...	1/17/2019 3:18 PM						342.1 -> 203...				alpha-hydroxymidazolam-...	6.366	920649
Amtripyline	278.2 -> 91.2	1/17/2019 3:18 PM						278.2 -> 105...				Amtripyline-D3	6.172	1003074
Acetylnorfentanyl	219.1 -> 84.2	1/17/2019 3:18 PM						219.1 -> 56.3				Acetylnorfentanyl-D5	4.693	804886
Carbamazepine	237.1 -> 193...	1/17/2019 3:18 PM						237.1 -> 194...				Carbamazepine-13C6	6.078	2548851
Chlordiazepoxide	300.1 -> 227...	1/17/2019 3:18 PM						300.1 -> 282...				Chlordiazepoxide-D5	6.494	1239275
Sufentanil	387.2 -> 238...	1/17/2019 3:18 PM						387.2 -> 111...				Sufentanil-D5	6.019	1867520
Desipramine	267.2 -> 72.3	1/17/2019 3:18 PM						267.2 -> 44.3				Desipramine-D3	6.134	3411387
Doxepin	280.2 -> 107...	1/17/2019 3:18 PM						280.2 -> 77.2				Doxepin-D3	5.928	360874
Estazolam	295.1 -> 267...	1/17/2019 3:18 PM						295.1 -> 205...				Estazolam-D5	6.331	1416696
Flunitrazepam	314.1 -> 268...	1/17/2019 3:18 PM						314.1 -> 239...				Flunitrazepam-D7	6.370	81979
Furazepam	388.2 -> 315...	1/17/2019 3:18 PM						388.2 -> 317...				Zopiclone-D4	5.614	290032
Imipramine	281.2 -> 86.2	1/17/2019 3:18 PM						281.2 -> 58.3				Imipramine-D3	6.130	2170756
Levamisole	205.1 -> 91.2	1/17/2019 3:18 PM						205.1 -> 178...				Sufentanil-D5	6.019	1867520
Maprotiline	278.2 -> 91.2	1/17/2019 3:18 PM						278.2 -> 117...				Amtripyline-D3	6.172	1003074
MDA	180.1 -> 163...	1/17/2019 3:18 PM						180.1 -> 105...				MDA-D5	4.802	425327
MDEA	208.1 -> 163...	1/17/2019 3:18 PM						208.1 -> 105...				MDEA-D6	5.063	812292
MDEA	194.1 -> 163...	1/17/2019 3:18 PM						194.1 -> 105...				MDEA-D6	4.902	70983
Meperidine	248.2 -> 220...	1/17/2019 3:18 PM						248.2 -> 174...				Meperidine-D4	5.401	602543
Methylphenidate	234.2 -> 84.2	1/17/2019 3:18 PM						234.2 -> 56.3				Meperidine-D4	5.401	602543
Midazolam	326.1 -> 223...	1/17/2019 3:18 PM						326.1 -> 249...				Midazolam-D4	6.408	915870
Normeperidine	234.1 -> 42.3	1/17/2019 3:18 PM						234.1 -> 160...				Normeperidine-D4	5.425	599698
7-Aminoflunitrazepam	284.1 -> 135...	1/17/2019 3:18 PM						284.1 -> 93.2				7-Aminoflunitrazepam-D7	5.643	957426
Nortriptyline	264.2 -> 91.2	1/17/2019 3:18 PM						264.2 -> 233...				Nortriptyline-D3	6.175	727229
Phenazepam	349.0 -> 206...	1/17/2019 3:18 PM						349.0 -> 183...				Phenazepam-D4	6.446	456009
Phencyclidine	244.2 -> 86.2	1/17/2019 3:18 PM						244.2 -> 91.1				Phencyclidine-D5	5.689	1018817
Phenylethylamine	253.1 -> 104...	1/17/2019 3:18 PM						253.1 -> 182...				Phenylethylamine-D10	5.973	27880
Pseudoephedrine	166.1 -> 148...	1/17/2019 3:18 PM						166.1 -> 133...				Pseudoephedrine-D3	4.421	183341...
Tapentadol	222.2 -> 107...	1/17/2019 3:18 PM						222.2 -> 121...				Tapentadol-D3	5.285	996733
Zaleplon	306.1 -> 236...	1/17/2019 3:18 PM						306.1 -> 264...				Zaleplon-D4	6.168	502085
Zopiclone	389.1 -> 244...	1/17/2019 3:18 PM						389.1 -> 216...				Zopiclone-D4	5.614	290032